

1-StepTM Slow TMB-ELISA

34024

Number Description

1-Step Slow TMB-ELISA, 250mL

Storage: Upon receipt store product at 4°C. Product is shipped at ambient temperature.

Introduction

The Thermo Scientific 1-Step Slow TMB-ELISA contains soluble TMB (3,3′,5,5′-tetramentylbenzidine) that is ideal for ELISA-based applications. TMB is a chromogenic substrate for horseradish peroxidase. Once oxidized by the enzyme, this substrate yields a blue product that absorbs at 370nm and 652nm. Adding sulfuric or phosphoric acid to stop the reaction changes the product to yellow that absorbs at 450nm and provides greater sensitivity than the blue product. The 1-Step Slow TMB-ELISA is supplied ready to use with no mixing required.

Procedure for using the 1-Step Slow TMB-ELISA

- 1. Equilibrate the TMB to room temperature. Add 100-150μL of the 1-Step Slow TMB-ELISA to each microplate well.
- 2. Incubate plate at room temperature for 5-30 minutes or until the desired color develops.
- 3. Stop reaction by adding 100-150µL of 2M sulfuric acid to each well.
- 4. Measure the absorbance of each well at 450nm.

Troubleshooting

Problem	Possible Cause	Solution
Some wells are blue while others are blue/green	Common TMB appearance	Reduce the HRP concentration, use a less sensitive substrate or reduce
A precipitate forms in the wells	The HRP concentration was too high	incubation time
High background	or incubation time was too long	

Related Thermo Scientific Products

15075	Reagent Reservoirs, 200/pkg
15082	Microtube Racked System, 960 tubes
15036	Sealing Tape for 96-Well Plates, 100/pkg
21124	Streptavidin, Horseradish Peroxidase Conjugated, 2mg
31002	NeutrAvidin [®] Protein, Horseradish Peroxidase Conjugated, 2mg
37538	StartingBlock TM (PBS) Blocking Buffer, 1L
37542	StartingBlock (TBS) Blocking Buffer, 1L



This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.