ė	inv	vitro	gen™
8	1117	'IU U	qen

CellSensor[®] NFAT CHO-K1 DA Cell-based Assay

Cat. nos. K1447

Shipping: Dry Ice

Storage: Liquid Nitrogen


Protocol part no. K1078.pps

Rev. date: 6 November 2010

Table	e of Contents Page
1.	Description1
2.	Overview of GeneBLAzer [®] GPCR Cell-based Assays2
3.	Overview of Division-Arrested Cells
4.	Materials Supplied
5.	Materials Required . 4 5.1 Optional Equipment and Materials . 4
6.	Thawing DA Cells5
7.	Media Requirements
8.	Assay Procedure.68.1Quick Assay Reference Guides.68.2Detailed Assay Protocol78.3Substrate Preparation, Loading and Incubation .88.4Detection8
9.	Data Analysis
10.	References10
11.	Purchaser Notification

1. Description

CellSensor[®] NFAT CHO-K1 DA (Division-arrested) cells are a division-arrested version of the CellSensor[®] NFAT-bla CHO-K1 cell line. CellSensor[®] NFAT-bla CHO-K1 contains a beta-lactamase reporter gene under control of a NFAT response element stably integrated into CHO-K1 cells. NFAT CHO-K1 DA cells have been functionally validated for Z'-Factor and EC₅₀ concentrations of Thapsigargin.



	Division- arrested cells
EC ₅₀	2.50 nM
Z'-factor at EC ₁₀₀	0.87

Dose response of NFAT CHO-K1 DA cells to Thapsigargin.

2. Overview of GeneBLAzer[®] GPCR Cell-based Assays

The GeneBLAzer[®] GPCR Cell-based Assay provides a highly accurate, sensitive, and easy-to-use method of monitoring cellular response to drug candidates or other stimuli (1). The core of the GeneBLAzer[®] assay technology is a beta-lactamase (*bla*) fluorescence resonance energy transfer (FRET) substrate that generates a ratiometric reporter response with minimal experimental noise. In addition to the dual-color (blue/green) readout of stimulated and unstimulated cells, this ratiometric method reduces the absolute and relative errors that can mask the underlying biological response of interest. Such errors include variations in cell number, transfection efficiency, substrate concentration, excitation path length, fluorescence detectors, and volume changes. The GeneBLAzer[®] GPCR assay technology has been proven in high-throughput screening (HTS) campaigns for a range of target classes, including G-protein coupled receptors (GPCRs) (2, 3), nuclear receptors (4-6) and kinase signaling pathways (7). The utility of division-arrested cells in HTS has also been demonstrated (8-11).

The GeneBLAzer[®] assay technology uses a mammalian-optimized *bla* reporter gene combined with a FRET-enabled substrate to provide reliable and sensitive detection in cells (1) (Figure 1). Cells are loaded with an engineered fluorescent substrate containing two fluorophores: coumarin and fluorescein. In the absence of *bla* expression, the substrate molecule remains intact. In this state, excitation of the coumarin results in fluorescence resonance energy transfer to the fluorescein moiety and emission of green fluorescent light. However, in the presence of *bla* expression, the substrate is cleaved separating the fluorophores and disrupting energy transfer. Excitation of the coumarin in the presence of *bla* enzyme activity results in a blue fluorescence signal.




Figure 1.

Fluorescent detection of beta-lactamase reporter gene response using a FRET-enabled substrate. After substrate loading, in the absence of beta-lactamase expression, cells appear green fluorescent. In the presence of beta-lactamase expression, the substrate is cleaved and cells appear blue fluorescent.

3. Overview of Division-Arrested Cells

Many of Invitrogen's cell lines are available in two forms: dividing or division arrested. Invitrogen's division-arrest technology allows the use of frozen cells, made from the exact same cell line sold in its dividing form, as ordinary, cost-effective assay reagents for screening. Division-arrested (DA) cells exhibit response profiles, similar to dividing cells thus ensuring you obtain the correct pharmacological profile.

DA cells may be plated and assayed within 24 hours of thawing. Cell numbers for DA cells increase only marginally after plating, thereby removing the variability caused by cell division during the course of an assay and providing more consistent results.

4. Materials Supplied

Product:	Name	Size	Catalog #
(One of these)	 CellSensor[®] NFAT CHO-K1 DA Assay Kit Each system contains sufficient division-arrested cells and substrate to assay one 384-well plate. Includes: CellSensor[®] NFAT CHO-K1 DA cells (K1447A) LiveBLAzer[™]-FRET B/G Loading Kit, 70 µg Solution D, 1 ml Protocol Certificate of Analysis 		K1447
Shipping Condition:	Dry ice		
Storage Condition of Cells:	Liquid nitrogen. Immediately upon receipt, cells must be stored in liquid nitrogen or thawed for immediate use. Cells stored at -80°C can quickly lose viability.		hawed for
Cell Phenotype:	Epithelial		
Mycoplasma Testing:	Negative		
BioSafety Level:	1		

5. Materials Required

Use the table below to determine the additional media and reagents required for use with each kit:

Media/Reagents	Recommended	Part #	Required Separately?	
	Source		NFAT CHO-K1 DA Assay Kit (K1447)	
LiveBLAzer [™] -FRET B/G Loading Kit: LiveBLAzer [™] -FRET B/G Substrate (CCF4-AM) DMSO for Solution A Solution B	Invitrogen	K1427 (70 μg) K1095 (200 μg) K1096 (1 mg)	No (included in kit)	
Solution C		K1030 (5 mg)		
Solution D	Invitrogen	K1156 (1 ml) K1157 (25 ml)	No (included in kit)	
DMEM (high-glucose), with GlutaMAX [™]	Invitrogen	10569-010	Yes	
DMSO	Fluka	41647	Yes	
Fetal bovine serum (FBS), dialyzed, (DO NOT SUBSTITUTE !)	Invitrogen	26400-036	Yes	
Non-essential amino acids (NEAA)	Invitrogen	11140-050	Yes	
Penicillin/Streptomycin (antibiotics)	Invitrogen	15140-122	Yes	
HEPES (1 M, pH 7.3)	Invitrogen	15630-080	Yes	
Thapsigargin	Calbiochem	586005	Yes	

The following table lists additional items required for use with all kits:

Consumables	Recommended Source	Part #
Black-wall, clear-bottom, 384-well assay plates (with low fluorescence background)	Corning	3712
Compressed air	Various	
Equipment	Recommended Source	
Fluorescence plate reader with bottom-read capabilities	Various	
Filters if required for plate reader (see Section 8.4.1)	Chroma Technologies	

5.1 Optional Equipment and Materials

- Epifluorescence- or fluorescence-equipped microscope, with appropriate filters
- Microplate centrifuge

6. Thawing DA Cells

- *Note:* Refer to Section 7, Media Requirements for specific media recipes.
- *Note:* Once cells are thawed per the instructions below, cells must be counted and the density adjusted to the appropriate level as specified in **Section 8.2.2**, **Assay Procedure**, prior to analysis.
- 1. Rapidly thaw the vial of cells by placing at 37°C in a water bath with gentle agitation for 1–2 minutes. Do not submerge vial in water.
- 2. Decontaminate the vial by wiping with 70% ethanol before opening in a Class II biological safety cabinet.
- 3. Transfer the vial contents drop-wise into 10 ml of Assay Medium in a sterile 15-ml conical tube.
- 4. Centrifuge cells at $200 \times g$ for 5 minutes.
- 5. Aspirate supernatant and resuspend the cell pellet in 1 ml fresh Assay Medium. (For the 1 x 384-well plate size of cells, dilute in 1 ml.)
- 6. Count the cells.
- 7. Adjust the cell density with Assay Medium to the appropriate cell density as specified in Section 8.2.2. Proceed to **Section 8, Assay Procedure**, for guidance on using cells in an assay.

7. Media Requirements

- *Note:* Unless otherwise stated, have all media and solutions at least at room temperature (we recommend 37°C for optimal performance) before adding to cells.
- *Note:* Make **NO MEDIA SUBSTITUTIONS**, as these cell lines have been specifically validated for optimal assay performance with these media.

Component	Assay Medium
DMEM	90%
Dialyzed FBS (Do not substitute!)	10%
NEAA	0.1 mM
HEPES (pH 7.3)	25 mM
Penicillin (antibiotic)	100 U/ml
Streptomycin (antibiotic)	100 µg/ml

8. Assay Procedure

The following instructions outline the recommended procedure for determining activity of compounds as modulators of using LiveBLAzerTM-FRET B/G Substrate as the readout. If alternative substrates are used (*e.g.*, ToxBLAzerTM DualScreen or LyticBLAzerTM Loading kits), follow the loading protocol provided with the product.

8.1 Quick Assay Reference Guides

For a more detailed assay protocol, see Section 8.2.

Agonist Assay Quick Reference Guide

	Unstimulated Wells	Stimulated Wells	Cell-free Wells	Test Compound Wells
Step 1 Plate cells, incubate	32 μl cells in Assay Medium (5,000 cells/well)	32 µl cells in Assay Medium (5,000 cells/well)	32 μl Assay Medium (no cells)	32 μl cells in Assay Medium (5,000 cells/well)
		Incubate cells for 16-2	4 hrs. at 37°C/ 5%CO ₂	
Step 2 Add Agonist or Test Compounds	8 μl Assay Medium with 0.5% DMSO	8 μl 5X agonist in Assay Medium with 0.5% DMSO	8 μl Assay Medium with 0.5% DMSO	8 μl 5X Test Compounds in 0.5% DMSO
Step 3 Incubate cells	Incubate in a humidified 37°C/5% CO ₂ incubator for 5 hours			
Step 4 Prepare 6X Substrate Mix	6 µl of 1 mM LiveBLAzer [™] -FRET B/G (CCF4-AM) Substrate + 60 µl of solution B, mix. Add 904 µl of Solution C, mix. Add 30 µl of Solution D, mix.			
Step 5 Add Substrate Mixture	8 μl per well			
Step 6 Incubate Substrate Mix. + cells	2 hours at room temperature in the dark			
Step 7 Detect activity	See Section 8.4			
Step 8 Analyze data	See Section 9			

Antagonist Assay Quick Reference Guide

	Unstimulated Wells	Stimulated Wells	Antagonist Control Wells	Cell-free Wells	Test Compound Wells
Step 1 Plate cells, incubate	32 μl cells in Assay Medium (5,000 cells/well)	32 μl cells in Assay Medium (5,000 cells/well)	32 μl cells in Assay Medium (5,000 cells/well)	32 μl Assay Medium (no cells)	32 μl cells in Assay Medium (5,000 cells/well)
		Incubate ce	ells for 16-24 hrs. at 37°C	C/ 5%CO ₂	
Step 2 Add Antagonist or Test Compounds, incubate	4 μl Assay Medium with 0.5% DMSO	4 μl Assay Medium with 0.5% DMSO	4 μl 10X antagonist in Assay Medium with 0.5% DMSO	4 μl Assay Medium with 0.5% DMSO	4 μl 10X Test Compounds in Assay Medium with 0.5% DMSO
	Incubate plate with Antagonist for 30 minutes before proceeding				
Step 3 Add Agonist	4 μl Assay Medium with 0.5% DMSO	4 μl 10X agonist in Assay Medium with 0.5% DMSO		4 μl 10X agonist in Assay Medium with 0.5% DMSO	4 μl 10X agonist in Assay Medium with 0.5% DMSO
Step 4 Incubate cells	Incubate in a humidified 37°C/5% CO ₂ incubator for 5 hours				
Step 5 Prepare 6X Substrate Mix	6 μl of 1 mM LiveBLAzer [™] -FRET B/G (CCF4-AM) Substrate + 60 μl of solution B, mix. Add 904 μl of Solution C, mix. Add 30 μl of Solution D, mix.				
Step 6 Add Substrate Mixture	8 μl per well				
Step 7 Incubate Mixture	2 hours at room temperature in the dark				
Step 8 Detect activity	See Section 8.4				
Step 9 Analyze data	See Section 9				

8.2 Detailed Assay Protocol

Plate layouts and experimental outlines will vary; in screening mode, we recommend using at least three wells for each control: Unstimulated Control, Stimulated Control, and Cell-free Control.

Note: Some solvents may affect assay performance. Assess the effects of solvent before screening. The cell stimulation procedure described below is carried out in the presence of 0.5% DMSO to simulate the effect that a Test Compound's solvent might have on the assay. If you use other solvents and/or solvent concentrations, optimize the following assay conditions appropriately.

8.2.1 Precautions

- Work on a dust-free, clean surface. Always handle the 384-well, black-wall, clear-bottom assay plate by the sides; do not touch the clear bottom of the assay plate.
- If pipetting manually, you may need to centrifuge the plate briefly at room temperature (for 1 minute at $14 \times g$) after additions to ensure all assay components are on the bottom of the wells.

8.2.2 Plating Cells

- 1. Thaw DA cells and resuspend in Assay Medium to a density of 156,250 cells/ml.
- Add 32 µl per well of the Assay Medium to the Cell-free Control wells. Add 32 µl per well of the cell suspension to the Test Compound wells, the Unstimulated Control wells, and Stimulated Control wells. Incubate cells at 37°C/ 5% CO₂ for 16-24 hours. Proceed to Section 8.2.3 for an Agonist assay or Section 8.2.4 for an Antagonist assay.

8.2.3 Agonist Assay Plate Setup

- *Note:* This subsection provides directions for performing an Agonist assay. See **Section 8.2.4** for directions for performing an Antagonist assay.
- 1. Prepare a stock solution of 0.5% DMSO in Assay Medium.
- 2. Prepare a 5X stock of Test Compounds in Assay Medium with 0.5% DMSO.
- 3. Prepare a 5X stock of agonist in Assay Medium with 0.5% DMSO. We recommend running a dose response curve to determine the optimal concentration of the agonist solution.
- 4. Add 8 μl of the stock solution of 0.5% DMSO in Assay Medium to the Unstimulated Control and Cell-free Control wells.
- 5. Add 8 µl of the 5X stock solution of agonist to the Stimulated Control wells.
- 6. Add 8 µl of the 5X stock of Test Compounds to the Test Compound wells.
- 7. Incubate the Agonist assay plate in a humidified 37°C/5% CO₂ incubator for 5 hours. Then proceed to **Section 8.3** for Substrate Loading and Incubation.

8.2.4 Antagonist Assay Plate Setup

- *Note:* This subsection provides directions for performing an Antagonist assay. See **Section 8.2.3** for directions for performing an Agonist assay.
- 1. Prepare a stock solution of 0.5% DMSO in Assay Medium.
- 2. Prepare a 10X stock of Test Compounds in Assay Medium with 0.5% DMSO.
- 3. Prepare a 10X stock of agonist in Assay Medium with 0.5% DMSO. We recommend running a dose response curve to determine the optimal agonist concentration. For antagonist assays, we recommend stimulating cells with an agonist concentration in the EC₅₀-EC₈₀ range.
- 4. Prepare a 10X stock of antagonist in Assay Medium with 0.5% DMSO. We recommend running a dose response curve to determine the optimal inhibition concentration for the Antagonist solution.
- 5. Add 4 µl of the 10X stock of Test Compounds to the Test Compound wells.
- 6. Add 4 μl of the stock solution of 0.5% DMSO to the Stimulated Control wells, the Unstimulated Control wells, and the Cell-free Control wells.
- Add 4 µl of the 10X stock of antagonist in Assay Medium with 0.5% DMSO to the Antagonist Control wells.
- 8. If desired, incubate the Test Compounds with the cells humidified 37°C/5% CO₂ incubator before proceeding. Typically, a 30-minute incubation is sufficient.
- 9. Add 4 µl of the 10X stock solution of agonist to the Test Compound wells, the Stimulated Control wells, and the Antagonist Control wells.
- 10. Add 4 µl of Assay Medium with 0.5% DMSO to the Unstimulated Control and Cell-free Control wells.

11. Incubate the Antagonist assay plate in a humidified 37°C/5% CO₂ incubator for 5 hours. Then proceed to **Section 8.3** for Substrate Loading and Incubation.

8.3 Substrate Preparation, Loading and Incubation

This protocol is designed for loading cells with LiveBLAzer[™]-FRET B/G Substrate Mixture (CCF4-AM) Substrate Mixture. If you use alternative substrates, follow the loading protocol provided with the substrate.

Prepare LiveBLAzer[™]-FRET B/G Substrate Mixture (CCF4-AM) Substrate Mixture and load cells in the absence of direct strong lighting. Turn off the light in the hood.

- Prepare Solution A: 1 mM LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM) Substrate Mixture in dry DMSO by adding 912 µl of DMSO per mg of dry substrate. Store the aliquots of the stock solution at -20°C until use. The molecular weight of the LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM) is 1096 g/mol.
- 2. Prepare 6X Loading Solution:
 - a. Add 6 µl of Solution A to 60 µl of Solution B and vortex.
 - b. Add 904 µl of Solution C to the above solution and vortex.
 - c. Add 30 µl of Solution D to the above solution and vortex.
- 3. Remove assay plate from the humidified $37^{\circ}C/5\%$ CO₂ incubator.
 - *Note:* Handle the plate gently and do not touch the bottom.
- 4. Add 8 µl of the 6X Substrate Mixture to each well.
- 5. Cover the plate to protect it from light and evaporation.
- 6. Incubate at room temperature for 2 hours.

8.4 Detection

Make measurements at room temperature from the bottom of the wells, preferably in 384-well, black-wall, clearbottom assay plates with low fluorescence background. Before reading the plate, remove dust from the bottom with compressed air.

8.4.1 Instrumentation, Filters, and Plates

- Fluorescence plate reader with bottom reading capabilities.
- Recommended filters for fluorescence plate reader:

Excitation filter:	409/20 nm
Emission filter:	460/40 nm
Emission filter:	530/30 nm

8.4.2 Reading an Assay Plate

- 1. Set the fluorescence plate reader to bottom-read mode with optimal gain and 5 reads.
- 2. Allow the lamp in the fluorescence plate reader to warm up for at least 10 minutes before making measurements.
- 3. Use the following filter selections:

	Scan 1	Scan 2
Purpose:	Measure fluorescence in the Blue channel	Measure FRET signal in the Green channel
Excitation filter:	409/20 nm	409/20 nm
Emission filter:	460/40 nm	530/30 nm

9. Data Analysis

9.1 Background Subtraction and Ratio Calculation

We recommend that you subtract the background for both emission channels (460 nm and 530 nm).

- 1. Use the assay plate layout to identify the location of the Cell-free Control wells. These Control wells are used for background subtraction.
- 2. Determine the average emission from the Cell-free Control wells at both 460 nm (Average Blue Background) and 530 nm (Average Green Background).
- 3. Subtract the Average Blue background from all of the Blue emission data.
- 4. Subtract the Average Green background from all of the Green emission data.
- 5. Calculate the Blue/Green Emission Ratio for each well, by dividing the background-subtracted Blue emission values by the background-subtracted Green emission values.

9.2 Visual Observation of Intracellular Beta-lactamase Activity Using LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM)

Note: Microscopic visualization of cells will cause photobleaching. Always read the assay plate in the fluorescence plate reader before performing microscopic visualization.

An inverted microscope equipped for epifluorescence and with either a xenon or mercury excitation lamp may be used to view the LiveBLAzer[™]-FRET B/G Substrate (CCF4-AM) signal in cells. To visually inspect the cells, you will need a long-pass filter passing blue and green fluorescence light, so that your eye can visually identify whether the cells are fluorescing green or blue.

Recommended filter sets for observing beta-lactamase activity are described below and are available from Chroma Technologies (800-824-7662, www.chroma.com).

Chroma Set # 41031

Excitation filter:	HQ405/20x (405 ± 10)
Dichroic mirror:	425 DCXR
Emission filter:	HQ435LP (435 long-pass)

Filter sizes vary for specific microscopes and need to be specified when the filters are ordered. For epifluorescence microscopes, a long-pass dichroic mirror is needed to separate excitation and emission light and should be matched to the excitation filter (to maximally block the excitation light around 405 nm, yet allow good transmission of the emitted light).

10. References

- 1. Zlokarnik, G., *et al*, Quantitation of Transcription and Clonal Selection of Single Living Cells with Beta-Lactamase as Reporter, (1998) *Science*; 279: p84-88.
- Kunapuli P, Ransom R, Murphy K, Pettibone D, Kerby J, Grimwood S, Zuck P, Hodder P, Lacson R, Hoffman I, Inglese J, Strulovici B, Development of an Intact Cell Reporter Gene Beta-lactamase Assay for G Protein-coupled Receptors, (2003) *Analytical Biochem.*; 314: p16-29.
- 3. Xing, H., Pollok, B., et al, A Fluorescent Reporter Assay For The Detection of Ligands Acting Through G1 Proteincoupled Receptors, (2000) J. Receptor & Signal Transduction Research; 20: p189-210.
- 4. Qureshi, S., *et al*, A One-Arm Homologous Recombination Approach for Developing Nuclear Receptor Assays in Somatic Cells, (2003) *Assay and Drug Dev. Tech*; **1**: p755-766.
- Peekhaus, N. et al, A Beta-Lactamase-Dependent Gal4-Estrogen Receptor Transactivation Assay for the Ultra-High Throughput Screening of Estrogen Receptor Agonists in a 3,456-Well Format, (2003) Assay and Drug Dev Tech, 1: p789-800.
- 6. Chin, J., et al, Miniaturization of Cell-Based, Beta-Lactamase-Dependent FRET Assays to Ultra-High Throughput Formats to Identify Agonists of Human Liver X Receptors, (2003) Assay and Drug Dev. Tech.; 1: p777-787.
- Whitney M, Rockenstein E, Cantin G, Knapp T, Zlokarnik G, Sanders P, Durick K, Craig FF, Negulescu PA., A Genome-wide Functional Assay of Signal Transduction in Living Mammalian Cells, (1998) Nat. Biotechnol.;16: p1329-1333.
- 8. Fursov, et al, (2005) Improving Consistency of cell-based assays by using division-arrested cells. Assay & Drug Development Technologies, Vol 3, pp 7-15
- 9. Kunapuli, *et al*, (2005) Application of division arrest technology to cell-based HTS: comparison with frozen and fresh cells. *Assay & Drug Development Technologies*, Vol 3, pp17-26.
- 10. Digan, *et al*, (2005) Evaluation of division-arrested cells for cell-based high-throughput screening and profiling. J. Biomolecular Screening. Vol. 10, pp 615 623.
- 11. Vasudevan, *et al*, (2005) **Improving high-content-screening assay performance by using division-arrested cells**. *Assay* & Drug Development Technologies, **Vol 3**, p.p.515.

11. Purchaser Notification

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008; Phone (760) 603-7200 or email: outlicensing@lifetech.com.

Limited Use Label License No. 150: GeneBLAzer[®] Technology

This product and/or its use is the subject of one or more of U.S. Patent Nos. 5,741,657, 5,955,604, 6,291,162, and 6,472,205 and foreign equivalents licensed to Life Technologies Corporation. The right to use this product for internal research specifically excludes the right to use this product to identify, discover, and profile compounds that act as a flavor, fragrance or taste-enhancers and modify a target identified in taste, olfaction, or pheromone detection, which compound does not require FDA approval of an NDA for claims of safety and efficacy. The right to use methods claimed in the foregoing patents with this product for research purposes can only be acquired by the use of GeneBLAzer® substrates purchased from Life Technologies Corporation or its authorized distributors.

Limited Use Label License No. 238: NFAT Technology

This product is the subject of one or more of US Patents 5,837,840, 6,096,515, 6,197,925, 6,150,099, pending applications and foreign equivalents licensed by Life Technologies Corporation.

Limited Use Label License No. 315: Division Arrested Cell Products

The purchaser of this product may not remove or obtain components from cells sold under this label license or use these cells in a manner inconsistent with protocols provided by Life Technologies Corporation for use of division arrested cells. The purchaser of this product may not attempt to expand division arrested cells and may only use the purchased amount of this product.

Use of Genetically Modified Organisms (GMO)

Information for European Customers The CHO-K1 DA cell lines are genetically modified to contain p3XNFATBlaX-Zeo and As a condition of sale, use of this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

Zeocin[™] is a trademark of CAYLA Corporation.

© 2008, 2010 Invitrogen Corporation. All rights reserved. Reproduction forbidden without permission