



00-2201      6 x 1 ml  
PNPP (100X)-Substrate  
for Alk. Phosphatase  
Lot No.

**PNPP (p-nitrophenyl phosphate)**  
**SUBSTRATE FOR ALKALINE PHOSPHATASE**

**FORM:**

Powder (100 mg). Dissolve the whole content in 1 ml of deionized water and keep it frozen as a stock solution.

**RECOMMENDED WORKING DILUTION:**

To use as a soluble substrate for alkaline phosphatase in enzyme immunoassay, dilute the stock PNPP solution 1:100 with 0.1 M 2-amino-2-methyl-1,3-propanediol, pH 10.3. Discard the unused portion of the diluted PNPP. Read optical density at 405 nm.

**STORAGE:**

Store at 4°C.

**REMARKS:**

For research use only. Not for human or drug use. No preservative added.

[www.invitrogen.com](http://www.invitrogen.com)

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)

PI002201

(Rev 10/08) DCC-08-1089

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, [www.invitrogen.com](http://www.invitrogen.com)). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.