PNPP - phosphatase substrate

Catalog Numbers 34045, 34047, 34064, and 37620

Doc. Part No. 2160338 Pub. No. MAN0011219 Rev. B.0

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Product description

Thermo Scientific[™] PNPP (*p*-nitrophenyl phosphate disodium salt) is a widely used substrate for detecting alkaline phosphatase in ELISA-based applications. The reaction produces a water-soluble yellow product that absorbs light at 405 nm.

Contents and storage

Product	Cat. No.	Amount	Storage
Phosphatase Substrate Kit:	37620		
• Diethanolamine Substrate Buffer (5X) (5.1 M diethanolamine, pH 9.8)		Buffer: 225 mL	Buffer: 4°C
• PNPP Tablets: Each tablet contains 5 mg of <i>p</i> -nitrophenyl phosphate disodium salt with a balance of inert filler		Tablets: 105	Tablets: -20°C
Diethanolamine Substrate Buffer (5X) (5.1 M diethanolamine, pH 9.8)	34064	225 mL	4°C
PNPP Tablets: Each tablet contains 5 mg of p -nitrophenyl phosphate disodium salt with a balance of inert filler	34047	105 tablets	-20°C
PNPP	34045	25 g	4°C

The tablets are shipped on wet ice, and the buffer is shipped at ambient temperature.

Use PNPP phosphatase substrate in an ELISA

CAUTION! Avoid PNPP contact with eyes, skin, and clothing. For handling information, see the MSDS.

- 1. Equilibrate PNPP to room temperature before opening.
- 2. Dilute the Diethanolamine Substrate Buffer from 5X to 1X.
- Dissolve 10 mg of PNPP or two PNPP Tablets in 10 mL of 1X Diethanolamine Substrate Buffer.
 Note: Diethanolamine is effective at concentrations from 10 mM to 1 M.
- 4. Add 100 μ L of the PNPP solution to each microplate well.
- 5. Incubate plate at room temperature for 30 minutes, or until sufficient color develops.
- 6. Add 50 μ L of 2 N NaOH to each well to stop the reaction.
- 7. Measure the absorbance at 405 nm. Alternatively, monitor reaction in a kinetic ELISA reader.

Related products

Product	Cat. No.
EZ-Link [™] Maleimide Activated Alkaline Phosphatase Kit	31493
Reagent Reservoirs, 200/pkg	15075
Microtube Racked System, 960 tubes	15082
Sealing Tape for 96-Well Plates, 100/pkg	15036
1-Step™ PNPP, 100 mL	37621
Streptavidin, Alkaline Phosphatase Conjugated, 3 mg	21323
NeutrAvidin [™] Protein, Alkaline Phosphatase Conjugated, 2 mg	31002

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Pierce Biotechnology, Inc. | Thermo Fisher Scientific | 3747 N. Meridian Road | Rockford, Illinois 61101 USA

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.