

Insulin-Transferrin-Selenium Supplement (100X) Insulin-Transferrin-Selenium-A Supplement (100X) Insulin-Transferrin-Selenium-X Supplement (100X)

Description

Insulin-Transferrin-Selenium (ITS) supplementation of many conventional synthetic nutrient media permits substantial reduction in the Fetal Bovine Serum (FBS) requirement for routine maintenance and low density attachment of many adherent cell types. ITS Supplement contains Insulin, Transferrin, and Sodium Selenite, prepared in Earle's Balanced Salt Solution (EBSS) without Phenol Red. In addition, ITS-A and ITS-X contain Sodium Pyruvate and Ethanolamine, respectively. ITS and ITS-A are designed to supplement RPMI-1640 and Earle's Minimal Essential Medium while ITS-X is intended for use as a supplement for F-12 Nutrient Mixture, Dulbecco's Modified Eagle Medium/F-12, DMEM, and Earle's Minimal Essential Medium. All ITS supplements will enhance the growth of various adherent cell types at FBS concentrations less than 4%.

Product	Catalog no.	Amount	Storage	Shelf life*
Insulin-Transferrin-Selenium Supplement (100X)	41400-045			
Insulin-Transferrin-Selenium-A Supplement (100X)	51300-044	10 mL	2°C to 8°C	18 months
Insulin-Transferrin-Selenium-X Supplement (100X)	51500-056			

^{*} Shelf life duration is determined from Date of Manufacture.

Product use

For Research Use Only. Not for use in diagnostic procedures.

Important information

Insulin, transferrin, and selenium are components shown to promote optimal performance of serum-free media and reduce the need for serum supplementation of conventional media used in cell culture

- Insulin has pleiotropic anabolic effects on mammalian cells, promoting glucose and amino acid uptake, lipogenesis, monovalent cation and phosphate transport, protein, and nucleic acid synthesis.
- Transferrin serves as a carrier for iron and may also help to reduce toxic levels of oxygen radicals and peroxide.
- Selenite is a co-factor for glutathione peroxidase and other proteins and is used as an anti-oxidant in media.
- Pyruvate is an important precursor to amino acid, fatty acid, and cholesterol biosynthetic pathways. Pyruvate can be utilized in the Krebs' cycle, and in gluconeogenesis.
- Ethanolamine can serve as a precursor for the biosynthesis of phosphoglycerides, which are essential to the structure of the plasma membrane and cellular organelles.

Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Caution: Human origin materials are non-reactive (donor level) for anti-HIV 1 & 2, anti-HCV, and HB_sAg . Handle in accordance with established bio-safety practices.

Use

- Each 10 mL vial of Insulin-Transferrin-Selenium 100X
 Supplement is sufficient for up to one liter of medium.
- In general, it is necessary to add 2–4% FBS to achieve optimal growth, although some adherent cultures may require less serum supplementation following initial adaptation.
- Store Insulin-Transferrin-Selenium supplemented medium in the dark at 2°C to 8°C.

Formulation*

Component	Concentration (g/L)				
	ITS	ITS-A	ITS-X		
Insulin	1.00	1.00	1.00		
Transferrin	0.55	0.55	0.55		
Sodium Selenite	0.00067	0.00067	0.00067		
Sodium Pyruvate	1	11.00	_		
Ethanolamine	_	_	0.20		

^{*} Prepared in EBSS without Phenol Red.

Explanation of symbols and warnings

The symbols present on the product label are explained below:

***	LOT	MM-1111	REF
Manufacturer	Batch code	Use By:	Catalog number
<u></u>	i	*	STERILE A
Caution, consult accompanying documents	Consult instructions for use	Temperature Limitation	Sterilized using aseptic processing techniques

Publication Number MAN0007340 Rev. 1.00

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Important licensing information

These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit www.lifetechnologies.com/support For further assistance, email techsupport@lifetech.com

©2014 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

DISCLAIMER - LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

life technologies™