PRODUCT INFORMATION Exonuclease III (Exo III)

Thermo

Pub. No. MAN0011997 Rev. Date 15.July.2016 (Rev. B.00)

#_ Lot: _

Expiry Date: _

Components	#EN0191		
Exonuclease III, 200 U/µL	4000 U		
10X Reaction Buffer	0.2 mL		

Store at -20 °C

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

Description

Exonuclease III (Exo III) exhibits four different catalytic activities (1):

• 3'→5' exodeoxyribonuclease activity specific for doublestranded DNA:

Exo III degrades dsDNA from blunt ends, 5'-overhangs or nicks, releasing 5'-mononucleotides from the 3'-ends of DNA strands and producing stretches of single-stranded DNA. It is not active on 3'-overhang ends of DNA that are at least four bases long and do not carry a 3'-terminal C-residue (2); on single-stranded DNA, or on phosphorothioate-linked nucleotides.

- 3'-phosphatase activity: Exo III removes the 3'-terminal phosphate and generates a 3'-OH group.
- RNase H activity: Exo III exonucleolytically degrades the RNA strand in DNA-RNA hybrids.
- Apurinic/apyrimidinic- endonuclease activity: Exo III cleaves phosphodiester bonds at apurinic or apyrimidinic sites to produce 5'-termini that are base-free deoxyribose 5'-phosphate residues.

Applications

- Creation of unidirectional deletions in DNA fragments in conjunction with S1 Nuclease (2, 3), see protocol on back page.
- Generation of a single-stranded template for dideoxy sequencing of DNA (4).
- Site-directed mutagenesis (5).
- Cloning of PCR products (6).
- Preparation of strand-specific probes.

Source

E.coli cells carrying a cloned *E.coli xth* gene.

Molecular Weight

31 kDa monomer.

Definition of Activity Unit

One unit of the enzyme catalyzes the release of 1 nmol of acid soluble reaction products from *E.coli* [³H]-DNA in 30 min at 37 °C.

Enzyme activity is assayed in the following mixture: 50 mM Tris-HCI (pH 8.0), 5 mM MgCl₂, 1 mM DTT and 0.05 mM sonicated *E.coli* [³H]-DNA.

Storage Buffer

The enzyme is supplied in: 50 mM Tris-HCl (pH 8.0), 50 mM KCl, 1 mM DTT and 50% (v/v) glycerol.

10X Reaction Buffer

660 mM Tris-HCl (pH 8.0 at 30 °C), 6.6 mM MgCl₂. Inhibition and Inactivation

- Inhibitors: metal chelators, p-chloromercuri benzoate (50-90% inhibitory at 0.1 mM) (7).
- Inactivated by heating at 70 °C for 10 min.

Note

- The rate of DNA digestion by Exo III depends upon temperature, salt concentration and the molar ratio of DNA to enzyme in a reaction mixture (4, 8). Optimal reaction conditions should be determined experimentally.
- Activity in Thermo Scientific REase Buffers, % (in comparison to activity in assay buffer)

Buffers	Activity, %
Thermo Scientific FastDigest, FastDigest™	
Green, O, R, 2X Thermo Scientific Tango,	
BamHI, EcoRI	0-25
G, 1X Tango™	25-50
B, Ecl136II, Kpnl, Pacl, Sacl	100

CERTIFICATE OF ANALYSIS

Endodeoxyribonuclease Assay

No detectable degradation was observed after incubation of supercoiled plasmid DNA with Exonuclease III.

Quality authorized by:

Jurgita Zilinskiene

(continued on back page)

Protocol for Generation of Unidirectional Deletions in DNA Fragments (for 25 time points)

- Digest 5-10 µg of DNA with a pair of restriction enzymes; one which generates a blunt or 5'-overhanging end adjacent to the target sequence, and another which produces a resistant 4 base 3'-overhang close to the priming site.
- Extract the reaction mixture with 1 volume of phenol:chloroform (1:1) and then with 1 volume of chloroform:isoamyl alcohol (24:1). Precipitate DNA by adding 0.1 volume of NaCl/glycogen solution (1.1 M NaCl, 0.25 mg/mL glycogen) and 2 volumes of 95% ethanol. Mix and centrifuge at 10,000 rpm for 10 minutes. Wash the pellet with 1 mL of 70% ethanol. Dry the pellet.
- 3. Dissolve DNA in 50 μL of 1X Exo III reaction buffer.
- 4. Prepare mixture (S1 Nuclease mix) by adding the following:

5X Reaction Buffer for S1 Nuclease*	40 µL
S1 Nuclease (#EN0321)	0.5 µL (50 U)
Water, nuclease-free (#R0581)	to 200 μL

200 mM sodium acetate (pH 4.5 at 25 °C), 1.5 M NaCl and 10 mM ZnSO₄)

 Add 7.5 µL aliquots of the S1 Nuclease mix into each of 25 numbered tubes, place on ice. 6. <u>Select incubation temperature and time as follows:</u>

Temperature, °C	25	30	37	45	
Digestion rate, bp/min	80	210	420	600	

Mononucleotides are released at base-dependent rates in the order: C>A=T>G. DNA termini of different are degraded at different rates, so these guidelines should only be used as an approximation.

- Warm the DNA solution to the digestion temperature. Add 500 units Exo III and mix rapidly. Remove 2 µL aliquots at the chosen time intervals and add to the tubes containing the cold S1 Nuclease mix.
- 8. After all samples have been taken, move the tubes to room temperature for 30 minutes. Add 1 μ L S1 "STOP" solution (300 mM Tris base, 50 mM EDTA) and heat at 70 °C for 10 minutes to inactivate S1 Nuclease.
- 9. Use 7 μL of each sample for gel electrophoresis.

References

- 1. Rogers, S.G. and Weiss, B., Exonuclease III of *Escherichia coli* K-12, an AP endonuclease, Methods Enzymol., 65, 201-211, 1980.
- 2. Hoheisel, J.D., On the Activities of *Escherichia coli* Exonuclease III, Anal.Biochem., 209, 238-246, 1993.
- 3. Henikoff, S., Unidirectional digestion with exonuclease III creates targeted breakpoints for DNA sequencing, Gene, 28, 351-359, 1984.
- 4. Guo, Li-He., Wu, R., New rapid methods for DNA sequencing based on exonuclease III digestion followed by repair synthesis, Nucleic Acids Res., 10, 2065-2084, 1982.
- 5. Vandeyar, M.A., A simple and rapid method for the selection of oligodeoxynucleotide-directed mutants, Gene, 65, 129-133, 1988.
- Li, Ch., Evans, R.M., Ligation independed cloning irrespective of restriction site compatibility, Nucleic Acids Res., 25, 4165-4166, 1997.
- 7. Eun, H.-M. Enzymology Primer for RecombinantDNA Technology, Academic Press. Inc., 1996.
- 8. Tomb, J.F., Barcak, G.J., Regulating the 3'-5' activity of exonuclease III by varying the sodium chloride concentration, BioTechniques, 7, 932-933, 1989.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, nontransferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only.* The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to <u>www.thermofisher.com</u> for Material Safety Data Sheet of the product.

 $\ensuremath{\mathbb{C}}$ 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.