

Dynamis[™] Medium

Description

DynamisTM Medium is specifically designed to offer the highest batch and fed-batch culture performance and yield with recombinant CHO cells in a chemically defined environment. The medium is formulated without hypoxanthine and thymidine for use in dihydrofolate reductase (DHFR)-amplified systems, without L-glutamine for use in glutamine synthetase systems, and without phenol red to minimize estrogen-like effects of phenol red. Furthermore, the glucose concentration is formulated to minimize potential lactic acid accumulation under typical culture conditions. The chemically defined, protein-free, animal origin component-free Dynamis TM Medium provides the power to achieve high titers, start process development faster, and streamline or simplify transfer to manufacturing scale.

Product	Catalog no.	Amount	Storage	Shelf life
Dynamis [™] Medium	A26615-01	1000 mL	2°C to 8°C; Protect from light	12 months

^{*} Shelf life duration is determined from Date of Manufacture.

Product use

Caution: For manufacturing, processing, or repacking.

Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Culture conditions

Media: Dynamis[™] Medium

Cell line: CHO cells
Culture type: Suspension

Temperature range: 36°C to 38°C

Incubator atmosphere: Humidified atmosphere of 8% CO₂ in air. Ensure that proper gas exchange is achieved in culture vessels and minimize exposure of cultures to light.

Culture vessels: Shake flasks, spinner bottles (rpm may vary with shaker platform/impeller design and should be empirically determined for optimal cell growth), or bioreactor.

Prepare complete medium

 $Dynamis^{^{\text{\tiny{TM}}}} Medium\ requires\ a septic\ supplementation\ with\ L-glutamine\ or\ GlutaMAX^{^{\text{\tiny{TM}}}}\text{-}I\ prior\ to\ use.$

- Add GlutaMAX[™]-I or L-glutamine at 2–8 mM final concentration to the medium before use.
- 2. Add 10 mL/L of HT Supplement for use in applications not requiring DHFR amplification.
- 3. Glucose supplementation may be required for terminal batch cultures and should be determined empirically.
- 4. Add Anti-Clumping Agent (1 mL/L) to the medium after transfection to reduce cell aggregation, if required.

Note: Consider reducing L-glutamine concentration for fed batch or perfusion protocols, or if the cell line in use is sensitive to ammonia. Addition of a surfactant such as Pluronic® F-68 is not required.

Recover frozen cells

- 1. Rapidly thaw (<1 minute) frozen cells in a 37°C water bath.
- Transfer the entire contents of the cryovial into a 125-mL shake flask containing 30 mL pre-warmed complete Dynamis™ Medium.
- 3. Incubate at 37° C in a humidified atmosphere of 8% CO₂ in air on an orbital shaker platform rotating at 115-135 rpm.

4. Maintain a cell density of 0.5×10^6 – 1×10^6 viable cells/mL for the first two passages following recovery; thereafter, return to your normal maintenance schedule.

Note: Do not centrifuge the cells after thawing as they are extremely fragile upon recovery from cryopreservation.

Subculture cells

- Determine viable cell density using a Countess® Automated Cell Counter (alternate automated or manual methods may also be used).
- 2. Ensure that the cell density is $\ge 1 \times 10^6$ viable cells/mL, viability is $\ge 90\%$, and the growth rate is in mid-logarithmic phase prior to subculturing.
- 3. Calculate the volume of cell culture and medium necessary to seed a flask at 2 × 10⁵–3 × 10⁵ viable cells/mL in a total volume of 30 mL fresh Dynamis[™] Medium per 125-mL shake

Note: If cell density does not reach 1×10^6 viable cells/mL within 5 days of recovery, centrifuge cells at $100 \times g$ for 5 minutes and resuspend the cell pellet in 20–30 mL of fresh complete DynamisTM Medium.

- 4. Incubate at 37° C in a humidified atmosphere of 8% CO₂ in air on an orbital shaker platform rotating at 115-135 rpm.
- For optimal performance and cell growth, dilute the cells at a seeding density of 3 × 10⁵ viable cells/mL every 3–4 days with fresh Dynamis[™] Medium.

Note: We recommend thawing a fresh, low-passage vial of cells for use every 3 months or 30 passages.

Adaptation of CHO cells to Dynamis™ Medium

We recommend adapting CHO cells to DynamisTM Medium using sequential adaptation. However, some CHO cell lines will adapt directly from other serum-free medium. It is critical that cell viability be $\geq 90\%$ and the growth rate be in mid-logarithmic phase prior to initiating adaptation procedures.

Direct adaptation

- For direct adaptation of CHO cells grown in other serum-free medium into Dynamis[™] Medium, dilute cells into 100%
 Dynamis[™] Medium using a seeding density of 3 × 10⁵–4 × 10⁵ viable cells/mL when subculturing (see Subculture cells).
- 2. Continue to subculture cells at $3\times10^5-4\times10^5$ viable cells/mL (every 3–4 days) until consistent growth is achieved. Once cell growth has been demonstrated, the seeding density may be reduced to $2\times10^5-3\times10^5$ viable cells/mL during the final stages of adaptation.

Publication Number MAN0013436 Revision 1.0

After several passages in Dynamis[™] Medium, the viable cell count should reach at least 2 × 10⁶ cells/mL with ≥85% viability within 3–4 days of seeding culture. At this stage, the culture is considered to be adapted to Dynamis[™] Medium.

Note: If suboptimal performance is achieved using the direct adaptation method, use the sequential adaptation method.

Sequential adaptation

- 1. During sequential adaptation of CHO cells grown in conventional 5–10% serum-supplemented medium or other serum-free medium, use a seeding density of 3×10^5 – 4×10^5 viable cells/mL.
- 2. Monitor cell growth using Countess® Automated Cell Counter until viable cell density reaches $\ge 1 \times 10^6$ cells/mL.
- 3. Dilute cells using a 25:75 ratio of complete Dynamis[™] Medium to the original medium. We recommend maintaining backup cultures in the original ratio medium until success with the new ratio medium is achieved. At each subsequent passage, dilute cells with stepwise increasing ratios of complete Dynamis[™] Medium to original medium (25:75, 50:50, 75:25, and 90:10, followed by 100% Dynamis[™] Medium). Multiple passages at each step may be needed to achieve consistent growth.
- 4. After several passages in 100% Dynamis™ Medium, the viable cell count should reach at least 2 × 106 cells/mL with ≥85% viability within 3–4 days of seeding culture. At this stage, the culture is considered to be adapted to Dynamis™ Medium.

Cryopreservation

- 1. Prepare the desired quantity of cells, harvesting in mid-log phase of growth with >90% viability. Reserve the conditioned medium to prepare cryopreservation medium.
- 2. Determine the viable cell density and calculate the required volume of cryopreservation medium to give a final cell density of $>1 \times 10^7$ cells/mL.
- 3. Prepare the required volume of cryopreservation medium of 92.5% Dynamis™ Medium (50:50 ratio of fresh-complete to conditioned media) + 7.5% DMSO and store at 4°C until use. Important: Prepare cryopreservation medium on the day of use.
- 4. Harvest cells by centrifugation at $100 \times g$ for 5–10 minutes. Resuspend the pellet in the pre-determined volume of 4°C cryopreservation medium.
- 5. Dispense aliquots of this suspension into cryovials according to the manufacturer's specifications.
- Achieve cryopreservation in an automated or manual controlled rate freezing apparatus following standard procedures (1°C decrease per minute).
- 7. Transfer frozen cells to liquid nitrogen (vapor phase) storage at -200°C to -125°C.

Note: Check viability of cryopreserved cells 24 hours after storage of vials in liquid nitrogen (see **Recover frozen cells**).

Related Products

Product	Cat. no.
L-glutamine, 200mM (100X), liquid	25030
GlutaMAX [™] -I, (100X), liquid	35050
HT Supplement, (100X), liquid	11067
Anti-Clumping Agent	0010057
Freedom® CHO-S® Kit	A13696-01
CHO-S® Cells (cGMP banked) and Media Kit	A11557-01
250X Cholesterol Lipid Concentrate	12531
FreeStyle [™] MAX Reagent	16447
EfficientFeed™ C+ AGT™ Supplement	A25031
Trypan Blue Stain	15250
Countess® Automated Cell Counter	C10227

Explanation of Symbols and Warnings

\triangle	\mathcal{X}	誉	\boxtimes	\prod i
aution, consult accompanying documents	Temperature Limitation	Keep away from light	Use By:	Consult instructions for use
LOT	REF	444	STERILE A	Read SDS
<u> </u>			5 LNIL	Keau 3D3

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Limited Use Label License No. 545: Chemically Defined Cell Culture Media

Notice to Purchaser: This product is covered by US and ex-US patents owned by Life Technologies, corresponding to chemically defined cell culture media and uses thereof. The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product to (a) perform internal research for the sole benefit of the purchaser; (b) manufacture protein (or other biological material) for resale; and (c) perform research or manufacturing services conducted by the purchaser on a fee for service or contract basis for or on behalf of third parties. No additional rights are granted. By purchasing this product, the purchaser agrees not to: (1) transfer or resell the product in any form; (2) use the product as a therapeutic agent or diagnostics test component; (3) reverse engineer the product or cause the product to be reverse engineered; or (4) use the product to perform tests other than what is indicated in this Limited Use Label License on a contract or fee per test basis for or on behalf of third parties. The purchaser is responsible for obtaining all regulatory approvals necessary for any therapeutic or diagnostic use of the protein (or biological material) manufactured using this product. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit www.lifetechnologies.com/support. For further assistance, email techsupport@lifetech.com.

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Pluronic is a registered trademark of BASF Corporation.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

