

# Champion<sup>™</sup> pET302/NT-His and pET303/CT-His Vectors

#### Vectors for high-level, inducible expression of N- and C-terminal 6x Histagged protein in *E. coli*

Catalog no. K6302-03

**Rev. Date: 7 July 2010** Manual part no. 25-0957

MAN0000576

**User Manual** 

# **Table of Contents**

Important Information	iv
Accessory Products	v
Overview	1
Methods	3
Cloning	3
Expression and Analysis	7
Appendix	10
Recipes	10
Map of pET302/NT-His	11
Map of pET303/CT-His	12
Features of pET302/NT-His and pET303/CT-His	13
Map of pET303/CT-Rac Kinase	14
Technical Service	15
Purchaser Notification	17
References	18

# Important Information

Shipping and Storage	Champion <sup>™</sup> pET302/NT-His, pET303/CT-His and pET303/CT-Rac Kinase vectors are shipped at room temperature. Upon receipt, store lyophilized vectors at -20°C.
Contents	The Champion <sup>™</sup> pET302/NT-His, pET303/CT-His vector kit contains the N- and C-terminal His vectors and an expression control plasmid as listed below:

Item	Concentration	Amount
pET302/NT-His	lyophilized in TE, pH 8.0	6 µg
pET303/CT-His	lyophilized in TE, pH 8.0	6 µg
pET303/CT-Rac Kinase	lyophilized in TE, pH 8.0	10 µg

# **Accessory Products**

#### Additional Products

Additional products that may be used with pET302/NT-His and pET303/CT-His are available from Invitrogen. Ordering information is provided below.

Product	Amount	Catalog no.
BL21(DE3) Chem. Competent Cells	20 x 50 µl	C6000-03
One Shot <sup>®</sup> BL21 Star (DE3) Chem. Competent Cells	20 x 50 µl	C6010-03
One Shot <sup>®</sup> BL21(DE3) pLysS Chem. Competent Cells	20 x 50 µl	C6060-03
MagicMedia <sup>™</sup> E. coli Expression	1 L SoluPouch™	K6801
Medium	5 x 1 L SoluPouch™	K6802
	1 L liquid	K6803
Ampicillin	5 g	Q100-16
Carbenicillin	5 g	10177-012

#### Detection of Recombinant Fusion Protein

Purification and If your gene of interest is in frame with the N- or C-terminal polyhistidine (6x His) tag, you may detect your fusion protein with an antibody to the polyhistidine tag. You may also purify your recombinant fusion protein using a metal chelating system.

Product	Quantity	Catalog no.
Mouse anti-His Tag monoclonal antibody	100 µg	37-2900
ProBond <sup>™</sup> Nickel-chelating Resin	50 ml	R801-01
	150 ml	R801-15
ProBond <sup>™</sup> Purification System	6 purifications	K850-01
Ni-NTA Agarose	10 ml	R901-15
	25 ml	R901-25
	100 ml	R901-10
Purification Columns	50	R640-50
(10 ml polypropylene columns)		

# Overview

Description	Champion <sup>TM</sup> pET302/NT-His (5.7 kb) and pET303/CT-His (5.3 kb) vectors allow you to clone your gene of interest using restriction enzyme digestion and ligation. Both Champion <sup>TM</sup> pET302/NT-His and pET303/CT-His are included, allowing you to choose the best configuration ( <i>i.e.</i> N- or C-terminal polyhistidine tag) and to optimize protein expression levels. The vectors are designed to allow high-level, inducible expression of recombinant fusion proteins in <i>E. coli</i> using the pET system. A control expression plasmid, pET303/CT-Rac Kinase, is included to optimize protein expression.
The pET Expression System	The pET system was originally developed by Studier and colleagues and takes advantage of the high activity and specificity of the bacteriophage T7 RNA polymerase to allow regulated expression of heterologous genes in <i>E. coli</i> from the T7 promoter (Rosenberg <i>et al.</i> , 1987; Studier & Moffatt, 1986; Studier <i>et al.</i> , 1990). For more information about T7-regulated expression, see the next page.
Features	Champion <sup>™</sup> pET302/NT-His and pET303/CT-His contain the following elements:
	• T7 <i>lac</i> promoter for high-level expression of the gene of interest in <i>E. coli</i> (see next page for more information)
	• Multiple cloning site for restriction enzyme digestion and ligation of gene of interest
	<ul> <li>N- or C-terminal 6x His tag for detection and purification of protein</li> </ul>
	• Ampicillin resistance gene for selection in <i>E. coli</i>
	• pBR322 origin for low-copy replication and maintenance of the plasmid in <i>E. coli</i>
	<ul> <li><i>lacI</i> gene encoding the lac repressor to reduce basal transcription from the T7<i>lac</i> promoter</li> </ul>
	For maps of pET302/NT-His and pET303/CT-His, see pages 11-12.
	Continued on next page

# **Overview**, continued

T7-Regulated Expression	pET302/NT-His and pET303/CT-His contain elements from bacteriophage T7 to control expression of heterologous genes in <i>E. coli</i> . In the vector, expression of the gene of interest is controlled by a strong bacteriophage T7 promoter that has been modified to contain a <i>lac</i> operator sequence (see below). In bacteriophage T7, the T7 promoter drives expression of gene 10. T7 RNA polymerase specifically recognizes this promoter. To express the gene of interest, it is necessary to deliver T7 RNA polymerase to the cells by inducing expression of the polymerase.
T7 <i>lac</i> Promoter	pET302/NT-His and pET303/CT-His contain the T7 <i>lac</i> promoter to drive expression of the gene of interest. The T7 <i>lac</i> promoter consists of a <i>lac</i> operator sequence placed downstream of the T7 promoter. The <i>lac</i> operator serves as a binding site for the lac repressor (encoded by the <i>lac</i> I gene) and functions to further repress T7 RNA polymerase-induced basal transcription of the gene of interest in BL21 strains.
BL21 Strains	The BL21(DE3) <i>E. coli</i> strain is specifically designed for expression of genes regulated by the T7 promoter. Basal level expression of T7 polymerase, particularly in BL21(DE3) cells, may lead to plasmid instability if your gene of interest is toxic to <i>E. coli</i> . You may also use BL21 Star <sup>™</sup> (DE3) and BL21 Star <sup>™</sup> (DE3)pLysS strains if your protein is toxic to <i>E. coli</i> . See page v for ordering information.

# Methods

Cloning	
Introduction	The following information is provided to help you clone your gene of interest into pET302/NT-His and pET303/CT- His. For basic information on DNA ligations, <i>E. coli</i> transformations, restriction analysis, DNA sequencing and DNA biochemistry, see <i>Current Protocols in Molecular Biology</i> (Ausubel <i>et al.</i> , 1994).
Resuspending and Propagating the Vectors	Resuspend pET302/NT-His and pET303/CT-His to 150 ng/µl in sterile water. Store at -20°C. To propagate the vectors, use this stock solution to transform a <i>recA</i> , <i>endA E</i> . <i>coli</i> strain like DH5 $\alpha$ , TOP10F' or equivalent. Transformants are selected on LB plates containing 50-100 µg/ml ampicillin. Be sure to prepare a glycerol stock of a transformant containing plasmid for long-term storage (see page 10).

### **Cloning**, continued

# Cloning into pET302/NT-His

pET302/NT-His is an N-terminal fusion vector and contains an ATG initiation codon and a Shine-Dalgarno ribosome binding site (RBS) with optimal spacing for proper translation. Your gene of interest must:

- Be in frame with the N-terminal tag (you may need to add additional nucleotides between your gene of interest and the restriction site)
- Contain a stop codon if you are cloning with *Bam*HI

The multiple cloning region of pET302/NT-His is shown below:

				Laco	operator	
	_	T7 promoter/pr	iming site			
11	CCGCGAAATT GGCGCTTTAA	AATACGACTC TTATGCTGAG	ACTATAGGGG TGATATCCCC	AATTGTGAGC TTAACACTCG	GGATAACAAT CCTATTGTTA	TCCCCTCTAG AGGGGAGATC
		R	ibosome bindir	ng site	►6X His	Tag Pml
71	АААТААТТТТ ТТТАТТАААА	GTTTAAACTT CAAATTTGAA	TAAGAAGGAG ATTCTTCCTC	ATATACATAT TATATGTATA	GCATCATCAT CGTAGTAGTA	CATCATCACG GTAGTAGTGC
	EcoRI Xho	I		Avrll stop code	BamHI	T7 terminator
131	TGAATTCGCT ACTTAAGCGA	CGAGATCGAT GCTCTAGCTA	GATATTCGAG CTATAAGCTC	CCTAGGTATA GGATCCATAT	ATCGGATCCG TAGCCTAGGC	GCTGCTAACA CGACGATTGT
				T7 reverse p	riming site	7
191	AAGCCCGAAA TTCGGGCTTT	GGAAGCTGAG CCTTCGACTC	TTGGCTGCTG AACCGACGAC	CCACCGCTGA GGTGGCGACT	GCAATAACTA CGTTATTGAT	GCATAACCCC CGTATTGGGG
251	TTGGGGCCTC AACCCCGGAG	TAAACGGGTC ATTTGCCCAG	TTGAGGGGTT AACTCCCCAA	TTTTGCTGAA AAAACGACTT	AGGAGGAACT TCCTCCTTGA	ATATCCGGAT TATAGGCCTA
311	ATCCCGCAAG TAGGGCGTTC	AGGCCCGGCA TCCGGGCCGT	GTACCGGCAT CATGGCCGTA			

Continued on next page

### **Cloning**, continued

# Cloning into pET303/CT-His is an C-terminal fusion vector and contains a Shine-Dalgarno ribosome binding site (RBS) with optimal spacing for proper translation. Your gene of interest must include:

- An ATG initiation codon
- If you wish to include the 6x His tag, your gene **should not** contain a stop codon.
- If you do not wish to fuse your gene of interest to the 6xHis tag, your gene **should** contain a stop codon The multiple cloning region of pET303/CT-His is shown below:

11	CCGCGAAATT GGCGCTTTAA	<b>17 promoter/pri</b> AATACGACTC TTATGCTGAG	Ming site ACTATAGGGG TGATATCCCC	Lac o AATTGTGAGC TTAACACTCG	ggataacaat CCTATTGTTA	TCCCCTCTAG AGGGGAGATC	ТААТААТТТТ АТТАТТАААА
81	Ribosc GTTTAACTTT CAAATTGAAA	me binding site		NsilXhol	6X His Ta CCACCACCAC GGTGGTGGTG	g CACCACTGAG GTGGTGACTC	ATCCGGCTGC TAGGCCGACG
	<b>T7</b> (			T7		4.0	
	T7 termi				erse priming si		
151	TAACAAAGCC ATTGTTTCGG	CGAAAGGAAG GCTTTCCTTC	CTGAGTTGGC GACTCAACCG	TGCTGCCACC	GCTGAGCAAT CGACTCGTTA	AACTAGCATA TTGATCGTAT	ACCCCTTGGG TGGGGAACCC
151	TAACAAAGCC	CGAAAGGAAG		TGCTGCCACC	GCTGAGCAAT	AACTAGCATA	
151 221	TAACAAAGCC	CGAAAGGAAG		TGCTGCCACC	GCTGAGCAAT	AACTAGCATA	

Continued on next page

# Cloning, continued

Ligation	vector gene o	you have determined a cloning strategy, digest the with the selected restriction enzymes. Ligate your f interest into the vector using standard molecular y techniques.			
Transformation	pET303 compe DNA f	igating your gene of interest into pET302/NT-His or 3/CT-His, transform the ligation mixture into tent <i>E. coli</i> . Select 10-20 clones and analyze plasmid for the presence and orientation of your insert by noing or appropriate restriction enzyme digestion.			
Sequencing	or C-te constru	firm that your gene of interest is in frame with the N- erminal 6X His tag, you may sequence your expression uct, if desired. Refer to the diagrams on pages 4-5 for ation of the primer binding sites.			
	For your convenience, Invitrogen offers a custom primer synthesis service. For more information, go to www.invitrogen.com or contact Technical Service (page 15).				
Primer		Sequence			

Primer	Sequence
T7 Promoter	5'-TAATACGACTCACTATAGGG-3'
T7 Reverse	5'-TAGTTATTGCTCAGCGGTGG-3'

The Next Step	Once you have generated your expression clone, you will
	need to transform it into a BL21 E. coli strain for expression
	studies. See page 15 for recommended BL21 host strains and
	media, and proceed to Expression and Analysis, next
	section.

# **Expression and Analysis**

Introduction	This section provides general guidelines for expressing and analyzing your protein of interest. For detailed information on transforming your BL21 strain, inducing expression, and analyzing samples, refer to your specific BL21 <i>E. coli</i> strain manual.
Basic Strategy	The basic steps needed to induce expression of your gene in a BL21 <i>E. coli</i> strain are outlined below.
	<ol> <li>Isolate plasmid DNA using standard procedures and transform your construct into BL21 cells. Use pET303/CT- GW/Rac Kinase included with the kit as a positive control for transformation and protein expression (see next page).</li> </ol>
	2. Grow the transformants and induce expression over several hours. Take several time points to determine the optimal time of expression. Alternatively, you can grow <i>E. coli</i> using MagicMedia <sup>™</sup> (see Recommendation, below)
- UNMENO	MagicMedia <i>E. coli</i> Expression Medium allows high yield of T7-regulated heterologous protein expression without time- consuming steps such as monitoring O.D. or adding inducing agents such as IPTG. MagicMedia <sup>™</sup> is available separately from Invitrogen, see page v for ordering information or go to www.invitrogen.com for more details.
Plasmid Preparation	You may prepare plasmid DNA using your method of choice. We recommend using the PureLink <sup>™</sup> HiPure Plasmid Midiprep Kit for isolation of pure plasmid DNA. Note that since you are purifying a low-copy number plasmid, you may need to increase the amount of bacterial culture that you use to prepare your plasmid construct.
	Continued on next page

# **Expression and Analysis, continued**

Choosing a Selection Agent	trar find cark ence into the rapi con	most purposes, ampicillin works well for selection of asformants and expression experiments. However, if you I that your expression level is low, you may want to use benicillin instead. The resistance gene for ampicillin odes a protein called $\beta$ -lactamase. This protein is secreted the medium where it hydrolyzes ampicillin, inactivating antibiotic. Since $\beta$ -lactamase is catalytic, ampicillin is idly removed from the medium, resulting in non-selective ditions. If your plasmid is unstable, this may result in the sof plasmid and low expression levels.
Using Carbenicillin	stuc amj pre cark exp <b>Not</b> con	benicillin is generally more stable than ampicillin, and dies have shown that using carbenicillin in place of picillin may help to increase expression levels by venting loss of the expression plasmid. If you wish to use penicillin, perform your transformation and expression eriments in LB containing 50 $\mu$ g/ml carbenicillin. <b>a:</b> If your gene of interest is highly toxic, increasing the centration of carbenicillin used from 50 $\mu$ g/ml to $\mu$ g/ml may help to increase expression levels.
Expression Control Vector pET303/CT- Rac Kinase	con This tagg deta	C303/CT-Rac Kinase is provided for use as a positive trol vector for protein synthesis in a suitable <i>E. coli</i> host. s vector allows expression of a 6x His C-terminally ged fusion protein of 57.7 kDa (Jones <i>et al.</i> , 1991). For ails about the vector, see page 14. To propagate and ntain the plasmid:
	1.	Resuspend the vector in 10 $\mu l$ of sterile water to prepare a 1 $\mu g/\mu l$ stock solution.
	2.	Use the stock solution to transform a <i>recA</i> , <i>endA E</i> . <i>coli</i> strain like TOP10, DH5 $\alpha^{\text{TM}}$ -T1 <sup>R</sup> , or equivalent. Use 10 ng of plasmid for transformation.
	3.	Select transformants on LB agar plates containing 50- 100 µg/ml ampicillin.
	4.	Prepare a glycerol stock of a transformant containing plasmid for long-term storage (see page 10).

# **Expression and Analysis, continued**

#### Detection of Recombinant Fusion Proteins

To detect expression of your recombinant fusion protein by western blot analysis, you may use antibodies against the polyhistidine tag (such as mouse anti-His monoclonal antibody, available separately from Invitrogen) or an antibody to your protein of interest. For more information, see page v, go to www.invitrogen.com or contact Technical Service (page 15).



The N-terminal peptide containing the 6x His tag will add approximately 1.5 kDa to your protein. The C-terminal peptide containing the 6x His tag will add approximately 0.9 kDa to your protein.

#### Purification of Recombinant Fusion Proteins

The presence of the N- or C-terminal 6x His tag in your recombinant fusion protein allows use of a metal-chelating resin such as ProBond<sup>™</sup> to purify your fusion protein. The ProBond<sup>™</sup> Purification System and bulk ProBond<sup>™</sup> resin are available from Invitrogen (see page v for ordering information). Invitrogen also offers Ni-NTA Agarose for purification of proteins containing a 6x His tag. **Note:** Other metal-chelating resins and purification methods are suitable.

# Appendix

# Recipes

Making Glycerol Stocks


- 1. Grow 1-2 ml of the *E. coli* strain to be frozen in SOB medium overnight with antibiotic selection when appropriate.
- 2. Combine 0.85 ml of the overnight culture with 0.15 ml of sterile glycerol (sterilized by autoclaving).
- 3. Mix well by vortexing.
- 4. Transfer to an appropriate freezing vial (preferably a screw cap, air-tight gasket).

Freeze in an ethanol-dry ice bath or liquid nitrogen and then transfer to -80°C for long-term storage.

# Map of pET302/NT-His

# Map of pET302/NT-His

The map below shows the elements of pET302/NT-His. The complete sequence of pET302/NT-His is available from www.invitrogen.com or by contacting Technical Service (page 15.)




(c) = complementary strand

# Map of pET303/CT-His

# Map of pET303/CT-His

The map below shows the elements of pET303/CT-His. The complete sequence of pET303/CT-His is available from www.invitrogen.com or by contacting Technical Service (page 15).



*bla* promoter: bases 775-879 Ampicillin (*bla*) resistance gene: bases 874-1734 pBR322 origin: bases 1945-2678 (c) *ROP* ORF: bases 2920-3011 (c) *lacl* ORF: bases 3914-5032 (c)

(c) = complementary strand


# Features of pET302/NT-His and pET303/CT-His

Features of the<br/>VectorsThe pET302/NT-His and pET303/CT-His vectors contain<br/>the following elements. Features have been functionally<br/>tested.

testea.		
Feature	Benefit	
T7 promoter	Allows high-level, inducible expression of your recombinant protein in <i>E. coli</i> strains expressing the T7 RNA polymerase	
T7 primer binding site	Allows sequencing of the insert	
<i>lac</i> operator ( <i>lac</i> O)	Binding site for lac repressor that serves to reduce basal expression of the recombinant protein	
Ribosome binding site	Optimally spaced from the initiation site for efficient translation of insert	
N-terminal or C-terminal 6x His tag	Allows purification of the recombinant protein on metal-chelating resin such as ProBond <sup>™</sup>	
T7 reverse primer binding site	Allows sequencing of the insert	
T7 transcription termination region	Sequence from bacteriophage T7 which allows efficient transcription termination.	
bla promoter	Allows expression of the ampicillin resistance gene	
Ampicillin resistance gene	Allows selection of transformants in <i>E. coli</i>	
pBR322 origin	Allows replication and maintenance in E. coli	
ROP ORF	Interacts with the pBR322 origin to facilitate low-copy replication in <i>E. coli</i> .	
lacI ORF	Encodes lac repressor which binds to the T7 <i>lac</i> promoter to block basal transcription of the gene of interest. Also binds the <i>lac</i> UV5 promoter in BL21 strains containing the $\lambda$ DE3 lysogen to repress transcription of T7 RNA polymerase	

### Map of pET303/CT-Rac Kinase

Map of pET303/CT-Rac Kinase The map below shows the elements of pET303/CT-Rac Kinase. The complete sequence of pET303/CT-Rac Kinase is available from www.invitrogen.com or by contacting Technical Service (page 15).



(c) = complementary strand

# **Technical Service**

Web Resources	<ul> <li>Visit the Invitrogen website at www.invitrogen.com for:</li> <li>Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.</li> <li>Complete technical support contact information</li> <li>Access to the Invitrogen Online Catalog</li> <li>Additional product information and special offers</li> </ul>
Contact Us	For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our website ( <u>www.invitrogen.com</u> ).

Corporate Headquarters:	Japanese Headquarters:	European Headquarters:
5791 Van Allen Way	LOOP-X Bldg. 6F	Inchinnan Business Park
Carlsbad, CA 92008 USA	3-9-15, Kaigan	3 Fountain Drive
Tel: 1 760 603 7200	Minato-ku, Tokyo 108-0022	Paisley PA4 9RF, UK
Tel (Toll Free): 1 800 955 6288	Tel: 81 3 5730 6509	Tel: +44 (0) 141 814 6100
Fax: 1 760 602 6500	Fax: 81 3 5730 6519	Tech Fax: +44 (0) 141 814 6117
E-mail:	E-mail:	E-mail:
tech_support@invitrogen.com	jpinfo@invitrogen.com	eurotech@invitrogen.com
···· · · · · · · · · · · · · · · · · ·	<u> </u>	

SDS

Safety Data Sheets (SDSs) are available at www.invitrogen.com/sds.

**Certificate of Analysis** The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to <u>www.invitrogen.com/support</u> and search for the Certificate of Analysis by product lot number, which is printed on the box.

continued on next page

# **Technical Service**, Continued

Limited Warranty	Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives. All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. <u>This warranty limits the Company's liability to only the price of the product</u> . No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order. Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives. Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above
	consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

### **Purchaser Notification**

Limited Use Label License No. 22: Vectors and Clones Encoding Histidine Hexamer This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

Limited Use Label License No. 30: T7 Expression System The composition and/or use of this product may be claimed in U.S. Patent No. 5,693,489 licensed to Life Technologies Corporation by Brookhaven Science Associates, LLC. The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy, and is the subject of patents and patent applications assigned to Brookhaven Science Associates, LLC (BSA,). By provisions of the Distribution License Agreement granted to Life Technologies covering said patents and patent applications, Life Technologies grants you a nonexclusive sub-license under patents assigned to BSA for the use of this technology, including the enclosed materials, based upon the following conditions: 1 - these materials are to be used for noncommercial research purposes only. A separate license under patents owned by BSA is required for any commercial use, including the use of these materials for research purposes or production purposes by any commercial entity. Information about commercial license may be obtained from The Office of Technology Transfer, Brookhaven National Laboratory, Bldg. 475D, P.O. Box 5000, Upton, New York 11973-5000. Phone (516) 344-7134. 2 - No materials that contain the cloned copy of the T7 gene 1, the gene for T7 RNA polymerase, may be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of this sub-license and agrees to be bound by its terms. This limitation applies to strains BL21(DE3), BL21(DE3)pLysS and BL21(DE3)pLysE, CE6, BL21-SI Competent Cells and any derivatives that are made of them. You may refuse this sub-license by returning this product unused in which case Life Technologies accept return of the product with a full refund. By keeping or using this product, you agree to be bound by the terms of this license.

### References

- Jones, P. F., Jakubowicz, T., Pitossi, F. J., Maurer, F., and Hemmings, B. A. (1991) Molecular cloning and identification of a serine/threonine protein kinase of the second messenger subfamily. PNAS 88, 4171-4175
- Kozak, M. (1987) An Analysis of 5´-Noncoding Sequences from 699 Vertebrate Messenger RNAs. Nucleic Acids Res. 15, 8125-8148
- Kozak, M. (1990) Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. Proc. Natl. Acad. Sci. USA 87, 8301-8305
- Kozak, M. (1991) An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. J. Cell Biology 115, 887-903
- Landy, A. (1989) Dynamic, Structural, and Regulatory Aspects of Lambda Sitespecific Recombination. Ann. Rev. Biochem. 58, 913-949
- Rosenberg, A. H., Lade, B. N., Chui, D.-S., Lin, S.-W., Dunn, J. J., and Studier, F. W. (1987) Vectors for Selective Expression of Cloned DNAs by T7 RNA Polymerase. Gene 56, 125-135
- Studier, F. W., and Moffatt, B. A. (1986) Use of Bacteriophage T7 RNA Polymerase to Direct Selective High-Level Expression of Cloned Genes. J. Mol. Biol. 189, 113-130
- Studier, F. W., Rosenberg, A. H., Dunn, J. J., and Dubendorff, J. W. (1990) Use of T7 RNA Polymerase to Direct Expression of Cloned Genes. Meth. Enzymol. 185, 60-89

©2010 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

### Notes

### Notes



Corporate Headquarters 5791 Van Allen Way Carlsbad, CA 92008 T: 1 760 603 7200 F: 1 760 602 6500 E: tech\_support@invitrogen.com

For country-specific contact information visit our web site at www.invitrogen.com