

PRODUCT INFORMATION

T4 Polynucleotide Kinase (T4 PNK)

Pub. No. MAN0011986

Rev. Date 15.July.2016 (Rev. B.00)

#_

Lot: _ Expiry Date: _

Components	#EK0031	#EK0032
T4 Polynucleotide Kinase (T4 PNK), 10 U/ μ L	500 U	2500 U
10X Reaction Buffer A <i>(for forward reaction)</i>	0.4 mL	2 \times 1 mL
10X Reaction Buffer B <i>(for exchange reaction)</i>	0.2 mL	1 mL
24% PEG 6000 Solution	0.2 mL	1 mL

Store at -20 °C

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

Description

T4 Polynucleotide Kinase (T4 PNK) catalyzes the transfer of the γ -phosphate from ATP to the 5'-OH group of single- and double-stranded DNAs and RNAs, oligonucleotides or nucleoside 3'-monophosphates (forward reaction). The reaction is reversible. In the presence of ADP T4 Polynucleotide Kinase exhibits 5'-phosphatase activity and catalyzes the exchange of phosphate group between 5'-P-oligo-/polynucleotides and ATP (exchange reaction) (1). The enzyme is also a 3'-phosphatase (2).

Applications

- Labeling 5'-termini of nucleic acids (3, 4) (see protocols on back page) to be used as:
 - probes for hybridization,
 - probes for transcript mapping,
 - markers for gel-electrophoresis,
 - primers for DNA sequencing,
 - primers for PCR.
- 5'-phosphorylation of oligonucleotides, PCR products, other DNA or RNA prior to ligation.
- Phosphorylation of PCR primers.
- Detection of DNA modification by the [³²P]-postlabeling assay (5, 6).
- Removal of 3'-phosphate groups (2).

Source

E.coli cells with a cloned *pseT* gene of bacteriophage T4.

Molecular Weight

The enzyme is a homotetramer. It consists of four identical subunits of 28.9 kDa.

Definition of Activity Unit

One unit of the enzyme transfers 1 nmol of γ -phosphate from ATP to 5'-OH DNA in 30 min at 37 °C.

Enzyme activity is assayed in the following mixture:

100 mM Tris-HCl (pH 8.0), 10 mM MgCl₂, 5 mM DTT, 0.5 mM 5'-OH DNA, 0.05 mM ATP and 0.1 MBq/mL [γ -³³P]-ATP.

Storage Buffer

The enzyme is supplied in: 20 mM Tris-HCl (pH 7.5), 25 mM KCl, 0.1 mM EDTA, 2 mM DTT and 50% (v/v) glycerol.

10X Reaction Buffer A (for forward reaction)

500 mM Tris-HCl (pH 7.6 at 25 °C), 100 mM MgCl₂, 50 mM DTT, 1 mM spermidine.

10X Reaction Buffer B (for exchange reaction)

500 mM imidazole-HCl (pH 6.4 at 25 °C), 180 mM MgCl₂, 50 mM DTT, 1 mM spermidine and 1 mM ADP.

Inhibition and Inactivation

- Inhibitors: metal chelators, phosphate and ammonium ions, KCl and NaCl at a concentration higher than 50 mM.
- Inactivated by heating at 75 °C for 10 min or by addition of EDTA.

Note

- 5'-termini of nucleic acids can be labeled by either the forward or the exchange reaction (1).
- Polyethylene glycol (PEG) and spermidine improve the rate and efficiency of the phosphorylation reaction (7). PEG is used in the exchange reaction.
- As T4 Polynucleotide Kinase is inhibited by ammonium ions, use sodium acetate to precipitate DNA prior to phosphorylation (1, 2).
- Activity in Thermo Scientific Buffers, % (in comparison to activity in buffer A):

Buffers	Activity, %
for restriction enzymes: Thermo Scientific FastDigest, FastDigest™ Green, O, G, 1X Thermo Scientific Tango, 2X Tango™, BamHI, EcoRI B, R, KpnI Ecl136II, PaeI, SacI	100 75-100 50-75
for <i>Taq</i> buffer with KCl	100
for RT buffers	100

(continued on back page)

Protocol for DNA/RNA 5'-end labeling by T4 PNK in the forward reaction

1. Prepare the following reaction mixture:

Dephosphorylated DNA <i>or</i> Oligonucleotide	1-20 pmol of 5'-termini 10-50 pmol
10X reaction buffer A	2 μ L
$[\gamma\text{-}^{32}\text{P}$ or $\gamma\text{-}^{33}\text{P}]$ -ATP	20 pmol
T4 Polynucleotide Kinase	1 μ L (10 U)
Water, nuclease-free (#R0581)	to 20 μ L
Total volume	20 μ L

2. Incubate at 37 °C for 30 min.

3. Add 1 μ L 0.5 M EDTA (pH 8.0) and incubate at 75 °C for 10 min.

4. Separate labeled DNA from unincorporated label by gel filtration on Sephadex G-50.

Protocol for DNA 5'-end labeling by T4 PNK in the exchange reaction

1. Prepare the following reaction mixture:

Linear DNA	1-20 pmol of 5'-termini
10X reaction buffer B	2 μ L
$[\gamma\text{-}^{32}\text{P}$ or $\gamma\text{-}^{33}\text{P}]$ -ATP	40 pmol
24% (w/v) PEG 6000 solution	4 μ L
T4 Polynucleotide Kinase	1 μ L (10 U)
Water, nuclease-free (#R0581)	to 20 μ L
Total volume	20 μ L

2. Incubate at 37 °C for 30 min.

3. Add 1 μ L 0.5M EDTA (pH 8.0) and incubate at 75 °C for 10 min.

4. Separate labeled DNA from unincorporated label by gel filtration on Sephadex G-50.

Note

- If ethanol solution of $[\gamma\text{-}^{32}\text{P}$ or $\gamma\text{-}^{33}\text{P}]$ -ATP is used, dry the required amount of ATP under vacuum and dissolve in water, nuclease-free.
- The ATP concentration should be at least 1 μ M in the forward reaction and at least 2 μ M in the exchange reaction (3, 4).

Protocol for Phosphorylation of DNA

1. Prepare the following reaction mixture:

Linear ds DNA <i>or</i> Oligonucleotide	1-20 pmol of 5'-termini 10-50 pmol
10X reaction buffer A for T4 Polynucleotide Kinase	2 μ L
ATP, 10 mM*	2 μ L
T4 Polynucleotide Kinase	1 μ L (10 U)
Water, nuclease-free (#R0581)	to 20 μ L
Total volume	20 μL

* Prepare 10 mM ATP solution by combining 10 μ L of 100 mM ATP solution (#R0441) and 90 μ L of Water, nuclease-free.

2. Mix thoroughly, spin briefly and incubate at 37 °C for 20 min.

3. Heat at 75 °C for 10 min.

Note

Visit www.thermofisher.com for molar calculations.


CERTIFICATE OF ANALYSIS

Endodeoxyribonuclease Assay

No detectable degradation was observed after incubation of supercoiled plasmid DNA with T4 Polynucleotide Kinase.

Ribonuclease Assay

No detectable degradation was observed after incubation of [3H]-RNA with T4 Polynucleotide Kinase.

Quality authorized by:  Jurgita Zilinskiene

References

1. Berkner, K.L., Folk, W.R., Polynucleotide kinase exchange reaction, *J. Biol. Chem.*, 252, 3176-3184, 1977.
2. Richardson, C.C., Bacteriophage T4 polynucleotide kinase, *The Enzymes* (Boyer, P.D., ed.), 14, 299-314, Academic Press, San Diego, 1981.
3. Sambrook, J., Russell, D.W., *Molecular Cloning: A Laboratory Manual*, the third edition, Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York, 2001.
4. *Current Protocols in Molecular Biology*, vol. 1 (Ausubel, F.M., et al., ed.), John Wiley & Sons, Inc., Brooklyn, New York, 3.10.2-3.10.5, 1994-2004.
5. Phillips, D.H., Detection of DNA modifications by the ³²P-postlabelling assay, *Mutation Res.*, 378, 1-12, 1997.
6. Keith, G., Dirheimer, G., Postlabeling: a sensitive method for studying DNA adducts and their role in carcinogenesis, *Curr. Opin. Biotechnol.*, 6, 3-11, 1995.
7. Harrison, B., Zimmerman, S.B., T4 polynucleotide kinase: macromolecular crowding increases the efficiency of reaction at DNA termini, *Anal. Biochem.*, 158, 307-315, 1986.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermofisher.com for Material Safety Data Sheet of the product.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.