SuperScript[™] III First-Strand Synthesis SuperMix

Catalog Number 18080-400

Doc. Part No. 18080400.pps Pub. No. MAN0001347 Rev. A.0

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **thermofisher.com/support**.

Product description

The Invitrogen^{\mathbb{N}} SuperScript^{\mathbb{N}} III First–Strand Synthesis SuperMix is an optimized SuperMix formulation for first-strand cDNA synthesis from purified poly(A)+ or total RNA. RNA targets from 100 bp to >12 kb can be detected with this system. The amount of starting material can vary from 0.1 pg to 5 µg of total RNA.

The kit includes SuperScript[™] III/RNaseOUT[™] Enzyme Mix, 2X First-Strand Reaction Mix, and Annealing Buffer. SuperScript[™] III Reverse Transcriptase, included in the Enzyme Mix, is a version of M–MLV RT that has been engineered to reduce RNase H activity and provide increased thermal stability. The enzyme can be used to synthesize cDNA at a temperature range of 45–55°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases. Because SuperScript[™] III Reverse Transcriptase is not significantly inhibited by ribosomal and transfer RNA, it can be used to synthesize cDNA from total RNA.

RNaseOUT[™] Recombinant RNase Inhibitor is included in the Enzyme Mix to safeguard against degradation of target RNA due to ribonuclease contamination.

The 2X First-Strand Reaction Mix includes 10 mM MgCl₂ and 1 mM of each dNTP in a buffer formulation that has been optimized for first-strand synthesis of cDNA. The Annealing Buffer is used in the initial template-primer annealing step. Separate tubes of $oligo(dT)_{20}$ and random hexamers are also provided.

cDNA synthesis can be performed using either total RNA or poly(A)+-selected RNA primed with oligo(dT), random primers, or a gene-specific primer. Recommendations for PCR amplification of the first-strand cDNA are provided in the "Guidelines for PCR" on page 2.

Contents and storage

Contents	Amount (50 reactions)	Storage
SuperScript [™] III/RNaseOUT [™] Enzyme Mix	100 µL	
2X First-Strand Reaction Mix (contains 10 mM MgCl ₂ and 1 mM each dNTP)	500 μL	
Annealing Buffer	50 µL	Store at -20°C
Oligo(dT) ₂₀ (50 μM)	50 µL	
Random hexamers (50 ng/µL)	50 μL	

Required materials not supplied

- 0.2-mL, nuclease-free, thin-walled PCR tubes
- Thermal cycler, preheated to 65°C
- Ice
- RNase/DNase-free water
- Microcentrifuge
- Vortex mixer

Procedural guidelines

Guidelines for RNA

- High-quality, intact RNA is essential for full-length, high-quality cDNA synthesis. This kit is designed for use with 0.1 pg to 5 μg of total RNA or 0.1 pg to 500 ng of poly(A)+ RNA.
- RNaseOUT[™] Recombinant RNase Inhibitor is included in the Enzyme Mix to safeguard against degradation of target RNA due to ribonuclease contamination of the RNA preparation.
- To isolate total RNA, we recommend TRIzol[™] Reagent (Cat. No. 15596-018) or the PureLink[™] RNA Mini Kit (Cat. No. 12183018A). Isolation of poly(A)+ RNA is typically not necessary, although it may improve the yield of specific cDNAs; we recommend Dynabeads[™] magnetic beads with DynaMag[™] magnets.
- Small amounts of genomic DNA in the RNA preparation may be amplified along with the target cDNA. If your application requires removal of all genomic DNA from your RNA preparation, we recommend using Invitrogen[™] ezDNase (Cat. No. 11766-051). ezDNase is ribonuclease-free and allows fast, efficient, and specific removal of genomic DNA from RNA samples without compromising the quality or quantity of RNA.



Guidelines for cDNA synthesis reactions

- For difficult or high GC-content templates, use a 55°C cDNA synthesis temperature.
- Keep all components, reaction mixes, and samples on ice. After preparation of the annealing mixture, transfer the mixture to the preheated thermal cycler.

Guidelines for primers

First-strand cDNA synthesis can be primed using random hexamers, oligo(dT), or gene-specific primers (GSPs):

• Random hexamers are the most nonspecific priming method, and are typically used for difficult or high GC-content mRNA. Using random hexamers, all RNAs in a population are templates for first-strand cDNA synthesis, and PCR primers confer specificity during PCR.

Note: For most RT-PCR applications, 50 ng of random hexamers per 5 µg of total RNA is adequate. Increasing hexamers to 250 ng per 5 µg of RNA may increase yield of small PCR products (<500 bp), but may decrease the yield of longer PCR products and full-length transcripts.

 Oligo(dT), a more specific priming method, is used to hybridize to 3' poly(A) tails, which are found in the vast majority of eukaryotic mRNAs. Since poly(A)+ RNA constitutes approximately 1–2% of total RNA, the amount and complexity of cDNA is considerably less than with random hexamers. We recommend using the oligo(dT)₂₀ provided in the kit.

Note: Oligo(dT) is recommended over random hexamers or GSPs for new mRNA targets. Oligo(dT) produces an RT-PCR product more consistently than random hexamers or GSPs.

• The most specific priming method uses a gene-specific primer (GSP) for the sequence of interest. First-strand synthesis can be primed with the PCR primer that hybridizes nearest to the 3' terminus of the mRNA. Note that some GSPs fail to prime cDNA synthesis even though they work in PCR on DNA templates. If gene-specific priming fails in RT-PCR, repeat first-strand synthesis using oligo(dT) as the primer.

Guidelines for PCR

The first-strand cDNA from the synthesis reaction may be amplified directly using PCR. We recommend using 10% of the first-strand reaction (2 μ L) for PCR. However, for some targets, increasing the amount of first-strand reaction to up to 10 μ L may result in increased product yield.

We recommend the following DNA polymerases (see "Ordering information" on page 4):

- Platinum[™] Taq DNA Polymerase provides automatic hot-start conditions for increased specificity and sensitivity. It is recommended for targets up to 4 kb.
- Platinum[™] Taq High Fidelity DNA Polymerase provides increased fidelity and higher yields for targets up to 15 kb.
- Platinum[™] *Pfx* DNA Polymerase possesses a proofreading 3′ to 5′ exonuclease activity and provides maximum fidelity for PCR. It is recommended for targets up to 12 kb.

Consult the product documentation provided with each DNA polymerase for recommended protocols and optimization guidelines.

Documentation is also available on our website at thermofisher.com/support.

Synthesize first-strand cDNA

The following procedure is designed to convert 0.1 pg to 5 µg of total RNA or 0.1 pg to 500 ng of poly(A)+ RNA into first-strand cDNA:

- 1. Mix and briefly centrifuge each component before use. Preheat the thermal cycler to 65°C.
- 2. Combine the following in a 0.2–mL, thin-walled PCR tube on ice:

Component	Amount
Up to 5 µg total RNA	nμL
Primer (50 μM oligo(dT) ₂₀ , <i>or</i> 2 μM gene-specific primer, <i>or</i> 50 ng/μL random hexamers)	1 µL
Annealing Buffer	1 µL
RNase/DNase-free water	to 8 μL

3. Incubate in a thermal cycler at 65°C for 5 minutes, then immediately place on ice for at least 1 minute. Collect the contents of the tube by brief centrifugation.

4. Add the following to the tube on ice:

Component	Amount
2X First-Strand Reaction Mix	10 µL
SuperScript™ III/RNaseOUT™ Enzyme Mix	2 μL

5. Vortex the sample briefly to mix, and collect by brief centrifugation. Incubate as follows:

- Oligo(dT)₂₀ or GSP primed: 50 minutes at 50°C
- Random hexamer primed: 5–10 minutes at 25°C, followed by 50 minutes at 50°C
- 6. Terminate the reactions at 85°C for 5 minutes. Chill on ice.
- 7. Store the cDNA synthesis reaction at –20°C, or proceed directly to PCR.

Troubleshooting

Observation	Possible cause	Recommended action
Unexpected bands after electrophoretic analysis		Pretreat RNA with ezDNase (Cat. No. 11766-051), as described in the ezDNase documentation.
		Design primers that anneal to sequence in exons on both sides of an intron or at the exon/exon boundary of the mRNA to differentiate between amplified cDNA and potential contaminating genomic DNA.
	Nonspecific annealing of primers	Vary the annealing conditions. Use Platinum [™] <i>Taq</i> DNA Polymerase for automatic hot-start PCR.
		Optimize magnesium concentration for each template and primer combination.
	Primers formed dimers	Design primers without complementary sequences at the 3' ends.
No or faint bands after analysis of amplified products	Procedural error in first-strand cDNA synthesis	Repeat the procedure, being careful to follow each step. Be careful to include the Annealing Buffer when adding primers and template for optimal yield.
	RNase contamination	Maintain aseptic conditions to prevent RNase contamination. RNaseOUT™ is included in the Enzyme Mix to inhibit RNases.
	Polysaccharide coprecipitation of RNA	Precipitate RNA with lithium chloride to remove polysaccharides.
	Target mRNA contains strong transcriptional pauses	Use random hexamers instead of oligo(dT) in the first-strand reaction.
		Increase the reaction temperature to 55°C.
		Prewarm the 2X First-Strand Reaction Mix to 55°C and hold the primer- template mixture at 55°C, then add the Enzyme Mix. Perform cDNA synthesis at 55°C.
		Use PCR primers closer to the 3' terminus of the target cDNA.
	Too little first-strand product was used in PCR	Use up to 10 μL of the first-strand reaction.
	GSP was used for first-strand synthesis	Try another GSP or switch to oligo(dT). Make sure the GSP is the antisense sequence.
	Inhibitors of RT present	Remove inhibitors by ethanol precipitation of mRNA preparation before the first-strand reaction. Include a 70% (v/v) ethanol wash of the mRNA pellet.
		Note: Inhibitors of RT include sodium dodecyl sulfate (SDS), EDTA, guanidinium salts, formamide, sodium pyrophosphate, and spermidine.

Ordering information

The following products are also available. Unless otherwise indicated, all materials are available through thermofisher.com.

Item	Amount	Source	
Platinum [™] <i>Taq</i> DNA Polymerase	100 units	10966-018	
	250 units	10966-026	
	500 units	10966-034	
Platinum [™] <i>Taq</i> High Fidelity DNA Polymerase	100 units	11304-011	
	500 units	11304-029	
Platinum [™] <i>Pfx</i> DNA Polymerase	100 units	11708-013	
	250 units	11708-021	
	500 units	11708-039	
PureLink™ RNA Mini Kit	50 reactions	12183018A	
TRIzol [™] Reagent	100 mL	15596-026	
	200 mL	15596-018	
ezDNase	100 units	11766-051	
Dynabeads [™] magnetic beads with DynaMag [™] magnets	To order, visit thermofisher.co thermofisher.com/magnets.	To order, visit thermofisher.com/dynabeads and thermofisher.com/magnets .	
Custom Primers	To order, visit thermofisher.co	To order, visit thermofisher.com .	

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

The information in this guide is subject to change without notice.

DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0001347

Revision	Date	Description
A.0	24 May 2016	Updated products (ezDNase and Dynabeads [™] magnetic beads)
		Format, style, and legal updates
_	20 May 2010	Baseline for this revision history

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies Corporation | Carlsbad, CA 92008 USA | Toll Free in USA 1 800 955 6288

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TRIzol is a trademark of Molecular Research Center, Inc.

