


# T7 RNA Polymerase-Plus™ Enzyme Mix

Catalog Number AM2716, AM2718

Doc. Part no. Pub. No. 4393878 Rev. C

Contents	Quantity	Storage conditions
T7 RNA Polymerase-Plus™ Enzyme Mix	AM2716: 5000 Units (20 U/μL) AM2718: 30,000 Units (200 U/μL)	Store at -20°C. <i>Do not store in a frost-free freezer.</i>
10X Transcription Buffer: 400 mM Tris pH 7.8, 200 mM NaCl, 60 mM MgCl <sub>2</sub> , 20 mM Spermidine HCl, 100 mM DTT	500 μL	

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

## Product description

T7 RNA Polymerase-Plus™ Enzyme Mix contains T7 RNA polymerase and RNase inhibitor for broad-spectrum protection against enzymatic degradation of the synthesized RNA.

**Source:** An *E. coli* strain harboring a plasmid that overexpresses T7 RNA Polymerase.

**Unit (U) definition:** One unit is the amount of T7 RNA Polymerase-Plus™ Enzyme Mix required to catalyze the incorporation of 1 nmol of nucleoside triphosphate into acid-insoluble material in 60 minutes at 37°C.

**Storage buffer (not included):** 20 mM KPO<sub>4</sub>, pH 7.7, 100 mM NaCl, 1 mM EDTA, 10 mM DTT, 50% glycerol (v/v), and detergent.

## Using T7 RNA Polymerase-Plus™ Enzyme Mix

T7 RNA Polymerase-Plus™ Enzyme Mix is ideal for the synthesis of RNA for blot and *in situ* hybridization probes, RPA analysis, *in vitro* translation, and antisense RNA. T7 RNA polymerase is highly specific for its own promoter, a conserved 23 bp sequence that is not efficiently recognized by SP6 or T3 RNA polymerases. It will transcribe large amounts of RNA from DNA sequences (for example, plasmids, polymerase chain reaction (PCR) fragments, or hybridized oligonucleotides) downstream of its promoter, without cross-talk from nearby SP6 or T3 promoters. Thus, RNA molecules transcribed from a linear template will be of a defined length.

Using circular plasmid DNA as a template will result in heterogeneous transcripts of multiple lengths.

## Labeled transcription reactions

The yield and amount of full-length RNA transcript obtained depend on the ratio of template DNA to the concentration of the limiting ribonucleoside triphosphate (rNTP) in the transcription reaction. Typically, three nucleotides are present at 500 μM and the nucleotide used for labeling is at various concentrations, depending on the desired specific activity of the probe (Tabor and Richardson, 1985). The limiting nucleotide should generally be present at a minimum of 3 μM to maximize synthesis of full-length RNA transcripts. Under standard assay conditions, >50% of the label is incorporated in 30 minutes into RNA.

We find that temperature is not a critical variable, although 37°C is frequently recommended as the appropriate incubation temperature. In fact, lower temperatures seem to favor the synthesis of full-length transcripts under conditions of limiting nucleotide concentration. It may be convenient to run these reactions at room temperature.

The following reaction conditions will yield labeled RNA probe suitable for use with Northern blots containing a moderately abundant mRNA.

### Typical labeled transcription reaction conditions (20-μL reactions):

- 1 μg template DNA
- 2 μL 10X Transcription Buffer (included)
- 500 μM (final) rNTPs (A, G, C)
- 50 μCi [ $\alpha$ -<sup>32</sup>P]UTP (800 mCi/mmol, 10 mCi/mL)
- 20 U T7 RNA Polymerase-Plus™ Enzyme Mix
- Nuclease-free water to a final volume of 20 μL

Incubate 30 minutes at 37°C.

Treat the reaction with DNase I to remove DNA template (optional; see below) or simply stop the transcription reaction by adding 2  $\mu$ L of 0.2 M EDTA and/or heating to 65°C.

## Unlabeled transcription reactions

In this reaction, rNTP levels are not limiting, and large amounts of RNA are synthesized throughout the incubation period. Frequently, more than 4  $\mu$ g RNA may be synthesized per  $\mu$ g of input DNA.

**Note:** This protocol may be altered to include nonisotopically-labeled rNTPs. See **Technical Bulletin 173, Methods for Nonisotopic Labeling** for a protocol and sources of these nucleotides.

### Typical unlabeled transcription reaction conditions (20- $\mu$ L reaction):

- 1  $\mu$ g template DNA
- 2  $\mu$ L 10X Transcription Buffer (included)
- 500  $\mu$ M (final) rNTPs (A, C, G, U)
- 20 U T7 RNA Polymerase-Plus™ Enzyme Mix
- Nuclease-free water to a final volume of 20  $\mu$ L

Incubate 60 minutes at 37°C.

---

## References

Tabor, S. and Richardson, C.C. (1985) *Proc. Natl. Acad. Sci. USA* 82, 1074–1078.

Treat the reaction with DNase I to remove DNA template (optional; see below) or simply stop the transcription reaction by adding 2  $\mu$ L of 0.2 M EDTA and/or heating to 65°C.

## (Optional) Removal of DNA template

Remove the DNA template by digestion with 2 units of DNase I (RNase-free; Cat. no. AM2222, AM2224) or TURBO™ DNase (Cat. no. AM2238, AM2239) for 15 minutes at 37°C. Inactivate the DNase by adding 2  $\mu$ L 0.2 M EDTA and heating at 70°C for 10 minutes, or by phenol/chloroform extraction.

For information about other post-reaction options, such as removal of unincorporated nucleotides, see the user guide for the MAXIscript® Kit (Cat. no. AM1312), available at [www.lifetechnologies.com](http://www.lifetechnologies.com).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.lifetechnologies.com/termsandconditions](http://www.lifetechnologies.com/termsandconditions). If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

### For Research Use Only. Not for use in diagnostic procedures.

The information in this guide is subject to change without notice.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

NOTICE TO PURCHASER: LIMITED USE LABEL LICENSE: Research Use Only : The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

© 2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit [lifetechnologies.com/support](http://lifetechnologies.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[lifetechnologies.com](http://lifetechnologies.com)