


# CountBright™ and CountBright™ Plus Absolute Counting Beads

Catalog Numbers C36950 and C36995

Pub. No. MAN0018850 Rev. A.0

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](http://thermofisher.com/support).

## Product description

Flow cytometry provides a rapid method to quantify cell characteristics, however most flow cytometers cannot directly provide the concentration or absolute count of cells in a sample. Absolute cell counts have been widely used in quantifying cell populations and disease progression, including studies of stem cells, HIV/AIDS, and leukemia.

CountBright™ and CountBright™ Plus Absolute Counting Beads are microparticles that can be used with most cell types, including no-wash/lysed whole blood. The microspheres are compatible with UV (350 nm), violet (405 nm), blue (488 nm), yellow (532 nm), green (561 nm), red (633 nm), and NIR (810 nm, CountBright™ Plus only ) excitation sources and emit fluorescence ranging from 385 nm to 860 nm.

## Contents and storage

Product	Catalog No.	Amount	Concentration	Size	Excitation	Emission	Storage
CountBright™ Absolute Counting Beads	C36950	5 mL	~0.5 × 10 <sup>5</sup> beads/50 µL Refer to bottle for lot-specific concentration.	7 µm	350-635 nm	385-800 nm	Store at 2-8°C protected from light.
CountBright™ Plus Absolute Counting Beads	C36995			4 µm	350-810 nm	385-860 nm	

## Experimental protocol

Cell count accuracy based on CountBright™ and CountBright™ Plus Absolute Counting Beads depends on sample handling and precise delivery of the volume of beads. The beads must be mixed well to ensure a uniform suspension of microspheres. This can be achieved by vortexing for 30 seconds before removing an aliquot. The microsphere suspension can be pipetted by standard techniques, but more viscous solutions, such as blood, require reverse pipetting for accurate volume delivery. Cell suspensions may be diluted, but should be assayed without wash steps after the addition of the beads.

- Process cells as needed before counting. This may include:
  - Staining cells with fluorochrome-conjugated antibodies
  - Fixation and/or permeabilization of cells
  - Lysis of red blood cells
- Allow the absolute counting beads to warm to room temperature. Vortex the bead suspension for 30 seconds to completely resuspend.
- Immediately after vortexing, add 50 µL of the counting bead suspension to each sample to be counted and vortex (the volume of each sample needs to be noted prior to addition of beads).

**Note:** Any buffer appropriate for cells (e.g., flow cytometry staining buffer or other buffered saline solutions, lysis buffer, fixation buffer, permeabilization buffer) can be used to bring the volume up to at least 300 µL. This volume is sufficient to dilute the small amount of detergent present in the beads storage buffer.

- Run the sample on the flow cytometer making sure to set the forward scatter (FSC) threshold low enough to include the beads on the linear-forward scatter vs. the linear-side scatter plot (SSC).
- Adjust fluorescence detector voltages such that stained cells and the absolute counting beads are appropriately on scale.
- Use normal gating strategies to identify the cell population to be enumerated.
- Draw a gate on the absolute counting beads and collect at least 1,000 bead events to ensure a statistically significant determination of sample volume (e.g., Fig. 1).

**Note:** A fluorescence threshold may also be used to analyze cells and beads.

- Using the count statistics from these 2 gates, the concentration of the original cell sample may be determined by the equations listed below.




Figure 1 Normal human peripheral blood was stained with Anti-Human CD3 Alexa Fluor™ 488 and Anti-Human CD19 PE-Cy7, then lysed with 1-Step™ Fix/Lyse Solution. CountBright™ Plus Absolute Counting Beads (boxed, purple) were added to the sample prior to collection on a flow cytometer.

## Calculations

**Note:** Only the ratio of the measured volumes of cells and absolute counting beads added to a sample is accounted for in the counting equations. Additional volumes of staining buffer, lysis buffer, etc., do not need to be taken into account when calculating absolute counts. Absolute count is defined as the concentration of cells in the original sample added to the counting tube.

- When cells and absolute counting beads are mixed in a 1:1 ratio of volumes, as in the preceding protocol, the following equation applies:

$$\text{Absolute count} \left( \frac{\text{cells}}{\mu\text{L}} \right) = \frac{\text{Cell count}}{\text{Counting beads count}} \times \text{Counting beads concentration from bottle} \left( \frac{\text{beads}}{\mu\text{L}} \right)$$

- When the cell to bead ratio is not 1:1, use the following equation (Note: Cell volume is the volume noted in Step 3 of the protocol.):

$$\text{Absolute count} \left( \frac{\text{cells}}{\mu\text{L}} \right) = \frac{(\text{Cell count} \times \text{Counting beads volume})}{(\text{Counting beads count} \times \text{Cell volume})} \times \text{Counting beads concentration} \left( \frac{\text{beads}}{\mu\text{L}} \right)$$

Sample calculation:

A 1 mL volume of cells was stained. Afterwards, 50 μL of absolute counting beads was added at a bottle concentration of  $0.49\text{--}0.5 \times 10^5$  beads/50 μL. The number of cell events was determined to be 1,700; bead events numbered 1,030.

$$\frac{1,700 \text{ cells}}{1,030 \text{ beads}} \times \frac{49,500 \text{ beads}}{1,000 \mu\text{L}} = 81.7 \text{ cells}/\mu\text{L}$$

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).



Life Technologies Corporation | 29851 Willow Creek | Eugene, OR 97402

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](http://thermofisher.com/symbols-definition).

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

**Important Licensing Information:** These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.