

GeneBLAzer® Gateway® Fusion Kits for *In Vivo* or *In Vitro* Detection

Gateway®-adapted destination vectors for
fluorescence detection of β -lactamase reporter
activity in mammalian cells

Catalog numbers 12578-035, 12578-043, 12578-050, and
12578-068

Revision date 29 March 2012

Publication Part number 25-0664

MAN0000390

For Research Use Only. Not intended for any animal or human
therapeutic or diagnostic use.

Contents

Kit Contents and Storage	v
Introduction	1
Product Overview	1
The GeneBLAzer™ Technology	3
Methods	4
Generating an Entry Clone	4
Creating an Expression Clone	7
Performing the LR Recombination Reaction	10
Transfecting Cells	13
Creating Stable Cell Lines	15
Detecting Recombinant Fusion Proteins	17
Appendix	19
Recipes	19
Blasticidin	20
Map and Features of pcDNA™ 6.2/cGeneBLAzer™-DEST	21
Map and Features of pcDNA™ 6.2/nGeneBLAzer™-DEST	23
Map of pcDNA™ 6.2/cGeneBLAzer™-GW/lacZ	25
Map of pcDNA™ 6.2/nGeneBLAzer™-GW/lacZ	26
Accessory Products	27
Technical Support	28
Purchaser Notification	29
Gateway® Clone Distribution Policy	30
References	31

Kit Contents and Storage

Types of Kits

This manual is supplied with the following products. All products are also supplied with the GeneBLAzer® Detection Kits manual.

Product	Catalog no.
GeneBLAzer® C-terminal Gateway® Fusion Kit for <i>In Vitro</i> Detection	12578-035
GeneBLAzer® C-terminal Gateway® Fusion Kit for <i>In Vivo</i> Detection	12578-043
GeneBLAzer® N-terminal Gateway® Fusion Kit for <i>In Vitro</i> Detection	12578-050
GeneBLAzer® N-terminal Gateway® Fusion Kit for <i>In Vivo</i> Detection	12578-068

Kit Components

The GeneBLAzer® Gateway® Fusion Kits include the following components. For a detailed description of the contents of each component, see page vi.

Component	Catalog no.			
	12578-035	12578-043	12578-050	12578-068
pcDNA™ 6.2/cGeneBLAzer™-DEST Gateway® Vector	√	√		
pcDNA™ 6.2/nGeneBLAzer™-DEST Gateway® Vector			√	√
GeneBLAzer® <i>In Vitro</i> Detection Kit	√		√	
GeneBLAzer® <i>In Vivo</i> Detection Kit		√		√

Shipping and Storage

The GeneBLAzer® Gateway® Fusion Kits are shipped as detailed below. Upon receipt, store as indicated.

Box	Shipping	Storage
GeneBLAzer® Gateway® Vector	Room temperature	Vectors: –20°C
GeneBLAzer® <i>In Vitro</i> Detection Kit	Dry Ice	CCF2-FA: –20°C, desiccated and protected from light
GeneBLAzer® <i>In Vivo</i> Detection Kit	Room temperature	CCF2-AM: –20°C, desiccated and protected from light Solutions: Room temperature, protected from light

Continued on next page

Kit Contents and Storage, Continued

GeneBLAzer® Gateway® Vectors

The vectors provided with the GeneBLAzer® Gateway® Fusion Kits are listed below. Store at –20°C.

Item	Concentration	Amount
Gateway® Destination Vector (pcDNA™ 6.2/cGeneBLAzer™-DEST or pcDNA™ 6.2/nGeneBLAzer™-DEST)	150 ng/μL in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	6 μg
Control Plasmid (pcDNA™ 6.2/cGeneBLAzer™-GW/ <i>lacZ</i> or pcDNA™ 6.2/nGeneBLAzer™-GW/ <i>lacZ</i>)	0.5 μg/μL in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	10 μg

GeneBLAzer® Detection Kits

The GeneBLAzer® Gateway® Fusion Kits include either the GeneBLAzer® *In Vitro* Detection Kit or the GeneBLAzer® *In Vivo* Detection Kit for fluorescence detection of β-lactamase reporter activity. Refer to the GeneBLAzer® Detection Kits manual for detailed information pertaining to each kit and a description of the reagents provided in each kit.

Product Use

For research use only. Not intended for any human or animal therapeutic or diagnostic use.

Introduction

Product Overview

Description of the System

The GeneBLAzer® Gateway® Fusion Kits contain Gateway®-adapted destination vectors designed for use with the GeneBLAzer® Technology. The pcDNA™ 6.2/GeneBLAzer®-DEST vector supplied with each kit facilitates *in vivo* or *in vitro* detection of β -lactamase reporter activity in mammalian cells using a unique fluorescent substrate. Use of the GeneBLAzer® Technology provides a highly sensitive and accurate method to quantitate gene expression in mammalian cells.

Features

The pcDNA™ 6.2/cGeneBLAzer™-DEST and pcDNA™ 6.2/nGeneBLAzer™-DEST vectors contain the following elements:

- Human cytomegalovirus immediate-early (CMV) promoter/enhancer for high-level expression in a wide range of mammalian cells
- β -lactamase *bla*(M) reporter gene for C-terminal (pcDNA™ 6.2/cGeneBLAzer™-DEST) or N-terminal (pcDNA™ 6.2/nGeneBLAzer™-DEST) fusion to the gene of interest
- Two recombination sites, *attR1* and *attR2*, downstream of the CMV promoter for recombinational cloning of the gene of interest from an entry clone
- Chloramphenicol resistance gene located between the two *attR* sites for counterselection
- The *ccdB* gene located between the two *attR* sites for negative selection
- The V5 epitope tag for detection using Anti-V5 antibodies (pcDNA™ 6.2/nGeneBLAzer™-DEST only)
- The Herpes Simplex Virus thymidine kinase polyadenylation signal for proper termination and processing of the recombinant transcript
- *f1* intergenic region for production of single-strand DNA in F plasmid-containing *E. coli*
- SV40 early promoter and origin for expression of the Blasticidin resistance gene and stable propagation of the plasmid in mammalian hosts expressing the SV40 large T antigen
- Blasticidin resistance gene for selection of stable cell lines
- The pUC origin for high copy replication and maintenance of the plasmid in *E. coli*
- The ampicillin resistance gene for selection in *E. coli*

For a map of pcDNA™ 6.2/cGeneBLAzer™-DEST or pcDNA™ 6.2/nGeneBLAzer™-DEST, refer to pages 21 and 23, respectively.

Continued on next page

Product Overview, Continued

The Gateway® Technology

The Gateway® Technology is a universal cloning method that takes advantage of the site-specific recombination properties of bacteriophage lambda (Landy, 1989) to provide a rapid and highly efficient way to move your gene of interest into multiple vector systems. To express your gene of interest in mammalian cells using Gateway® Technology, simply:

1. Clone your gene of interest into a Gateway® entry vector to create an entry clone.
2. Generate an expression clone by performing an LR recombination reaction between the entry clone and a Gateway® destination vector (e.g. pcDNA™ 6.2/cGeneBLAzer™-DEST or pcDNA™ 6.2/nGeneBLAzer™-DEST).
3. Transfect your expression clone into the cell line of choice for transient or stable expression of your gene of interest.

For more information on Gateway®, refer to the Gateway® Technology manual. This manual is available for downloading from www.lifetechnologies.com/support or by contacting Technical Support (page 28).

Advantages of the GeneBLAzer® Technology

Using the GeneBLAzer® Technology and the GeneBLAzer® Detection System as a reporter of gene expression in mammalian cells provides the following advantages:

- Suitable for use as a sensitive reporter of gene expression in living mammalian cells using fluorescence microscopy.
- Provides ratiometric readout to minimize differences due to variability in cell number, substrate concentration, fluorescence intensity, and emission sensitivity.
- Compatible with a wide variety of *in vivo* and *in vitro* applications including microplate-based transcriptional assays and flow cytometry.
- Provides a flexible and simple assay development platform for gene expression in mammalian cells.
- Using a non-toxic substrate allows continued cell culturing after quantitative analysis.

For more information on GeneBLAzer® Technology, see page 3.

The GeneBLAzer® Technology

Components of the GeneBLAzer® System

The GeneBLAzer® System facilitates fluorescence detection of β -lactamase reporter activity in mammalian cells, and consists of two major components:

- The β -lactamase reporter gene, *bla*(M), a truncated form of the *E. coli bla* gene. When fused to a gene of interest, the *bla*(M) gene can be used as a reporter of gene expression in mammalian cells. For more information about the *bla*(M) gene, see below.
 - A fluorescence resonance energy transfer (FRET)-enabled substrate, CCF2 to facilitate fluorescence detection of β -lactamase activity. In the absence or presence of β -lactamase reporter activity, cells loaded with the CCF2 substrate fluoresce green or blue, respectively. Comparing the ratio of blue to green fluorescence in a population of live cells or in a cell extract of your sample to a negative control provides a means to quantitate gene expression. For more information about the CCF2 substrate and how FRET works, refer to the GeneBLAzer® Detection Kits manual.
-

β -Lactamase (*bla*) Gene

β -lactamase is the product encoded by the ampicillin resistance gene (*bla*) and is the bacterial enzyme that hydrolyzes penicillins and cephalosporins. The *bla* gene is present in many cloning vectors and allows ampicillin selection in *E. coli*. β -lactamase enzyme activity is not found in mammalian cells.

bla(M) Gene

The GeneBLAzer® Technology uses a modified *bla* gene as a reporter in mammalian cells. This *bla* gene is derived from the *E. coli TEM-1* gene present in many cloning vectors (Zlokarnik *et al.*, 1998), and has been modified in the following ways:

- 72 nucleotides encoding the first 24 amino acids of β -lactamase were deleted from the N-terminal region of the gene. These 24 amino acids comprise the bacterial periplasmic signal sequence, and deleting this region allows cytoplasmic expression of β -lactamase in mammalian cells.
- The amino acid at position 24 was mutated from His to Asp to create an optimal Kozak sequence for optimal translation initiation.

This modified reporter gene is named *bla*(M).


Note: The *TEM-1* gene also contains 2 mutations (at nucleotide positions 452 and 753) that distinguish it from the *bla* gene in pBR322 (Sutcliffe, 1978).

Methods

Generating an Entry Clone

Introduction

To recombine your gene of interest into pcDNA™6.2/cGeneBLAzer™-DEST or pcDNA™6.2/nGeneBLAzer™-DEST, you will need an entry clone containing the gene of interest. Many entry vectors including pENTR™/D-TOPO® are available to facilitate generation of entry clones (see page 27 for ordering). For more information, refer to www.lifetechnologies.com or contact Technical Support (page 28). Refer to the manual for the specific entry vector you are using for detailed instructions to construct an entry clone.



Note

If you wish to express a human or mouse gene of interest, we recommend using an Ultimate™ Human ORF (hORF) or Ultimate™ Mouse ORF (mORF) Clone available for purchase. Each Ultimate™ ORF Clone is a fully-sequenced clone provided in a Gateway® entry vector that is ready-to-use in an LR recombination reaction with a pcDNA™6.2/GeneBLAzer®-DEST vector. In addition, each clone contains a **TAG** stop codon, making it fully compatible for use in the Tag-On-Demand™ System. For more information about the Ultimate™ ORF Clones available, refer to www.lifetechnologies.com or contact Technical Support (page 28).

Kozak Consensus Sequence

If you will be expressing your protein from pcDNA™6.2/cGeneBLAzer™-DEST, your insert in the entry clone should contain a Kozak translation initiation sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1990; Kozak, 1991). An example of a Kozak consensus sequence is provided below. The ATG initiation codon is shown underlined.

(G/A)NNATGG

Other sequences are possible, but the G or A at position -3 and the G at position +4 are the most critical for function (shown in bold).

Continued on next page

Generating an Entry Clone, Continued

Points to Consider for pcDNA™ 6.2/ cGeneBLAzer™ - DEST

pcDNA™ 6.2/cGeneBLAzer™-DEST allows expression of recombinant proteins containing a C-terminal β -lactamase reporter; however, you may use this vector to express a native protein. You may also use this vector in the Tag-On-Demand™ System (see page **Error! Bookmark not defined.**). Consider the following when generating your entry clone.

If you wish to...	Then your insert...
include the β -lactamase reporter	<ul style="list-style-type: none">• should contain a Kozak initiation sequence (see page 4)• should not contain a stop codon• should be in frame with the <i>bla</i>(M) reporter gene after recombination (see page 8 for a diagram)
include the β -lactamase reporter for use in the Tag-On-Demand™ System	<ul style="list-style-type: none">• should contain a Kozak initiation sequence (see previous page)• should contain a TAG stop codon• should be in frame with the <i>bla</i>(M) reporter gene after recombination (see page 8 for a diagram)
not include the β -lactamase reporter	<ul style="list-style-type: none">• should contain a Kozak initiation sequence (see page 4)• should contain a stop codon

Continued on next page

Generating an Entry Clone, Continued

Points to Consider for pcDNA™ 6.2/ nGeneBLAzer™ - DEST

pcDNA™ 6.2/nGeneBLAzer™-DEST allows expression of recombinant proteins containing an N-terminal β -lactamase reporter and a C-terminal V5 epitope tag, if desired, and contains an ATG initiation codon within the context of a Kozak consensus sequence (see page 4). You may use this vector in the Tag-On-Demand™ System (see page **Error! Bookmark not defined.**). Consider the following when generating your entry clone.

If you wish to...	Then your insert...
include the β -lactamase reporter	<ul style="list-style-type: none">• should not contain a Kozak initiation sequence• should be in frame with the <i>bla</i>(M) reporter gene after recombination (see page 9 for a diagram)
include the V5 epitope tag	<ul style="list-style-type: none">• should not contain a stop codon• should be in frame with the V5 epitope after recombination (see page 9 for a diagram)
include the V5 epitope for use in the Tag-On-Demand™ System	<ul style="list-style-type: none">• should contain a TAG stop codon• should be in frame with the V5 epitope after recombination (see page 9 for a diagram)
not include the V5 epitope tag	<ul style="list-style-type: none">• should contain a stop codon

Creating an Expression Clone

Introduction

After you have generated an entry clone, perform the LR recombination reaction to transfer the gene of interest into your pcDNA™ 6.2/GeneBLAzer®-DEST vector to create your expression clone. To ensure that you obtain the best possible results, we recommend that you read this section and the next section entitled **Performing the LR Recombination Reaction** (pages 10–12) before beginning.

Experimental Outline

To generate an expression clone:

1. Perform an LR recombination reaction using the *attL*-containing entry clone and the *attR*-containing pcDNA™ 6.2/GeneBLAzer®-DEST vector.
 2. Transform the reaction mixture into a suitable *E. coli* host.
 3. Select for expression clones (refer to pages 8–9 for diagrams of the recombination region of the resulting expression clones).
-

Propagating the pcDNA™ 6.2/GeneBLAzer®-DEST Vectors

If you wish to propagate and maintain the pcDNA™ 6.2/GeneBLAzer®-DEST vectors, we recommend using Library Efficiency® DB3.1™ Competent Cells for transformation (see page 27 for ordering). The DB3.1™ *E. coli* strain is resistant to CcdB effects and can support the propagation of plasmids containing the *ccdB* gene.

Note: DO NOT use general *E. coli* cloning strains including TOP10 or DH5α™ for propagation and maintenance of the pcDNA™ 6.2/GeneBLAzer®-DEST vectors as these strains are sensitive to CcdB effects.

Continued on next page

Creating an Expression Clone, Continued


Recombination Region for pcDNA™ 6.2/cGeneBLAzer™ -DEST

The recombination region of the expression clone resulting from pcDNA™ 6.2/cGeneBLAzer™-DEST × entry clone is shown below.

Note: If you are using pcDNA™ 6.2/cGeneBLAzer™-DEST in the Tag-On-Demand™ System, your gene of interest must contain a TAG stop codon (see page **Error! Bookmark not defined.**).

Features of the Recombination Region:

- Shaded regions correspond to DNA sequences transferred from the entry clone into pcDNA™ 6.2/cGeneBLAzer™-DEST by recombination. Non-shaded regions are derived from the pcDNA™ 6.2/cGeneBLAzer™-DEST vector.
- Bases 922 and 2605 of the pcDNA™ 6.2/cGeneBLAzer™-DEST vector sequence are marked.



Continued on next page

Creating an Expression Clone, Continued

Recombination Region for pcDNA™ 6.2/nGeneBLAzer™ -DEST

The recombination region of the expression clone resulting from pcDNA™ 6.2/nGeneBLAzer™-DEST × entry clone is shown below.

Note: If you are using pcDNA™ 6.2/nGeneBLAzer™-DEST in the Tag-On-Demand™ System, your gene of interest must contain a TAG stop codon (see page [Error! Bookmark not defined.](#)).

Features of the Recombination Region:

- Shaded regions correspond to DNA sequences transferred from the entry clone into pcDNA™ 6.2/nGeneBLAzer™-DEST by recombination. Non-shaded regions are derived from the pcDNA™ 6.2/nGeneBLAzer™-DEST vector.
- Bases 1719 and 3402 of the pcDNA™ 6.2/nGeneBLAzer™-DEST vector sequence are marked.

```

      CAAT                                TATA      3' end of CMV promoter      Putative
771  CAAATGGGCG GTAGGCGTGT ACGGTGGGAG GTCTATATAA GCAGAGCTCT CTGGCTAACT      transcriptional start
                                     T7 promoter
831  AGAGAACCCA CTGCTTACTG GCTTATCGAA ATTAATACGA CTCACTATAG GGAGACCCAA
                                     β-lactamase bla(M) reporter
891  GCTGGCTAGT TAAGCTGAGC ATG GAC CCA GAA ACG CTG GTG AAA GTA AAA GAT
      Met Asp Pro Glu Thr Leu Val Lys Val Lys Asp
      1694
944  GCT GAA GAT ... // ... AAG CAT TGG CTG TTA TCA ACA AGT TTC TAC AAA
      CGA CTT CTA ... // ... TTC GTA ACC GAC AAT AGT TGT TCA AAC ATG TTT
      Ala Glu Asp ... // ... Lys His Trp Leu Leu Ser Thr Ser Leu Tyr Lys
      1719 attB1
1727 AAA GCA GGC TNN --- --- NAC CCA GCT TTC TTG TAC AAA GTG GTT GAT AAC
      TTT CGT CCG ANN GENE NTG GGT CGA AAG AAC ATG TTT CAC CAA CTA TTG
      Lys Ala Gly ... ... Pro Ala Phe Leu Tyr Lys Val Val Asp Asn
      3402 attB2
3424 GGG AAG CCT ATC CCT AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT
      Gly Lys Pro Ile Pro Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly
      V5 epitope
      TK polyA Reverse priming site
3475 TAG TAA TGA GTTTAAACGG GGGAGGCTAA CTGAAACACG GAAGGAGACA ATACCGGAAG
      *** *** ***
```


Performing the LR Recombination Reaction

Introduction

Once you have obtained an entry clone containing your gene of interest, you may perform an LR recombination reaction between the entry clone and your pcDNA™6.2/GeneBLAzer®-DEST vector, and transform the reaction mixture into a suitable *E. coli* host (see below) to select for an expression clone. We recommend including the pENTR™-gus positive control supplied with the LR Clonase® enzyme mix in your experiments to help you evaluate your results.

E. coli Host

You may use any *recA*, *endA* *E. coli* strain including TOP10, DH5α™, or equivalent for transformation. **Do not** transform the LR reaction mixture into *E. coli* strains that contain the F' episome (e.g. TOP10F'). These strains contain the *ccdA* gene and will prevent negative selection with the *ccdB* gene.



Note

The presence of the EM7 promoter and the Blasticidin resistance gene in the pcDNA™6.2/GeneBLAzer®-DEST vectors allows for selection of *E. coli* transformants using Blasticidin. For selection, use Low Salt LB agar plates containing 100 µg/mL Blasticidin (see the **Appendix**, page 19 for a recipe). For Blasticidin to be active, the salt concentration of the medium must remain low (<90 mM) and the pH must be 7.0.

Blasticidin is available separately for purchase (see page 27 for ordering information). See the **Appendix**, page 20 for instructions on how to prepare and store Blasticidin.

Materials Needed

- Purified plasmid DNA of your entry clone (50–150 ng/ µL in TE, pH 8.0)
 - pcDNA™6.2/cGeneBLAzer™-DEST or pcDNA™6.2/nGeneBLAzer™-DEST vector (150 ng/µL in TE, pH 8.0)
 - LR Clonase® enzyme mix (keep at –80°C until immediately before use; see page 27 for ordering)
 - 5X LR Clonase® Reaction Buffer (supplied with the LR Clonase® enzyme mix)
 - pENTR™-gus positive control, optional (50 ng/µL in TE, pH 8.0; supplied with the LR Clonase® enzyme mix)
 - TE Buffer, pH 8.0 (10 mM Tris-HCl, pH 8.0, 1 mM EDTA)
 - 2 µg/µL Proteinase K solution (supplied with the LR Clonase® enzyme mix; thaw and keep on ice until use)
 - Appropriate competent *E. coli* host and growth media for expression
 - S.O.C. Medium
 - LB agar plates containing the appropriate antibiotic to select for expression clones
-

Continued on next page

Performing the LR Recombination Reaction, Continued

Setting Up the LR Recombination Reaction

1. Add the following components to 1.5 mL microcentrifuge tubes at room temperature and mix.

Note: To include a negative control, set up a second sample reaction and substitute TE Buffer, pH 8.0 for the LR Clonase® enzyme mix (see Step 4).

Component	Sample	Positive Control
Entry clone (100–300 ng/reaction)	1–10 µL	–
Destination vector (150 ng/µL)	2 µL	2 µL
pENTR™-gus (50 ng/µL)	–	2 µL
5X LR Clonase® Reaction Buffer	4 µL	4 µL
TE Buffer, pH 8.0	to 16 µL	8 µL

2. Remove the LR Clonase® enzyme mix from –80°C and thaw on ice (~ 2 minutes).
3. Vortex the LR Clonase® enzyme mix briefly twice (2 seconds each time).
4. To each sample above, add 4 µL of LR Clonase® enzyme mix. Mix well by pipetting up and down.

Reminder: Return LR Clonase® enzyme mix to –80°C immediately after use.

5. Incubate reactions at 25°C for 1 hour.
Note: Extending the incubation time to 18 hours typically yields more colonies.
6. Add 2 µL of the Proteinase K solution to each reaction. Incubate for 10 minutes at 37°C.
7. Transform 1 µL of the LR recombination reaction into a suitable *E. coli* host (follow the manufacturer's instructions) and select for expression clones.
Note: You may store the LR reaction at –20°C for up to 1 week before transformation, if desired.

What You Should See

If you use *E. coli* cells with a transformation efficiency of 1×10^8 cfu/µg, the LR reaction should give >5000 colonies if the entire LR reaction is transformed and plated.

Confirming the Expression Clone

The *ccdB* gene mutates at a very low frequency, resulting in a very low number of false positives. True expression clones will be ampicillin-resistant and chloramphenicol-sensitive. Transformants containing a plasmid with a mutated *ccdB* gene will be both ampicillin- and chloramphenicol-resistant. To check your putative expression clone, test for growth on LB plates containing 30 µg/mL chloramphenicol. A true expression clone will not grow in the presence of chloramphenicol.

Continued on next page

Performing the LR Recombination Reaction, Continued

Sequencing Primers for pcDNA™ 6.2/cGene BLAzer™ -DEST

To confirm that your gene of interest is in frame with the *bla*(M) reporter gene, you may sequence your expression construct, if desired. Keep the following in mind when designing your sequencing primers:

- Use a forward primer which hybridizes within the 3' end of your gene of interest to sequence through the *att*B2 site and the 5' region of the *bla*(M) reporter gene (see page 8 for a diagram).
 - Do not use a reverse primer that hybridizes within the *bla*(M) reporter gene. Any primer that hybridizes within the *bla*(M) reporter gene will also hybridize within the ampicillin resistance gene, contaminating your results.
Note: Because you will not be using a reverse primer, you will only be able to sequence the sense strand of your expression construct.
 - Use the T7 Promoter primer to sequence through the *att*B1 site and into the 5' region of your gene of interest. Refer to the diagram on page 8 for the location of the T7 Promoter primer binding site.
-

Sequencing Primers for pcDNA™ 6.2/nGene BLAzer™ -DEST

To confirm that your gene of interest is in frame with the *bla*(M) reporter gene or the V5 epitope tag, you may sequence your expression construct, if desired. Keep the following in mind when designing your sequencing primers:

- Use a reverse primer which hybridizes within the 5' end of your gene of interest to sequence through the *att*B1 site and the 3' region of the *bla*(M) reporter gene (see page 9 for a diagram).
 - Do not use a forward primer that hybridizes within the *bla*(M) reporter gene. Any primer that hybridizes within the β -lactamase reporter gene will also hybridize within the ampicillin resistance gene, contaminating your results.
Note: Because you will not be using a forward primer, you will only be able to sequence the anti-sense strand of your expression construct.
 - Use the TK polyA Reverse primer to sequence through the *att*B2 site and into the V5 epitope. Refer to the diagram on page 9 for the location of the TK polyA Reverse primer binding site.
-

Transfecting Cells

Introduction

This section provides general information for transfecting your expression clone into the mammalian cell line of choice. We recommend that you include a positive control vector (pcDNA™ 6.2/cGeneBLAzer™-GW/*lacZ* or pcDNA™ 6.2/nGeneBLAzer™-GW/*lacZ*) and a mock transfection (negative control) in your experiments to evaluate your results.



If you plan to detect β -lactamase reporter activity *in vivo* using the GeneBLAzer® *In Vivo* Detection Kit (supplied with Catalog nos. 12578-043 and 12578-068 only), note that a number of factors including cell type and cell density can influence the degree of the fluorescence signal detected. We recommend taking these factors into account when designing your transfection experiment. For more information, refer to the section entitled **General Guidelines to Use the GeneBLAzer® *In Vivo* Detection Kit** in the GeneBLAzer® Detection Kits manual.

Plasmid Preparation

Once you have generated your expression clone, you must isolate plasmid DNA for transfection. Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink® HQ Mini Plasmid Purification Kit (see page 27 for ordering) or CsCl gradient centrifugation.

Positive Control

pcDNA™ 6.2/cGeneBLAzer™-GW/*lacZ* or pcDNA™ 6.2/nGeneBLAzer™-GW/*lacZ* is provided as a positive control vector for mammalian cell transfection and expression (see pages 25 and 26 for a map) and may be used to optimize recombinant protein expression levels in your cell line. These vectors allow expression of the β -galactosidase gene with either an N-terminal or C-terminal fusion to the β -lactamase reporter.

To propagate and maintain the plasmid:

1. Use the stock solution to transform a *recA*, *endA* *E. coli* strain like TOP10, DH5 α ™, or equivalent.
 2. Select transformants on LB agar plates containing 50–100 μ g/mL ampicillin.
 3. Prepare a glycerol stock of a transformant containing plasmid for long-term storage.
-

Continued on next page

Transfecting Cells, Continued

Methods of Transfection

For established cell lines (e.g. HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. We recommend that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

Methods for transfection include calcium phosphate (Chen & Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1989; Felgner & Ringold, 1989) and electroporation (Chu *et al.*, 1987; Shigekawa & Dower, 1988). For high efficiency transfection in a broad range of mammalian cell lines, we recommend using Lipofectamine® 2000 Reagent (see page 27 for ordering). For more information about Lipofectamine® 2000 and the other transfection reagents available for purchase, refer to www.lifetechnologies.com or contact Technical Support (page 28).

Creating Stable Cell Lines

Introduction

The pcDNA™ 6.2/GeneBLAzer®-DEST vectors contain the Blasticidin resistance gene to allow selection of stable cell lines. If you wish to create stable cell lines, transfect your construct into the mammalian cell line of choice and select for foci using Blasticidin. General information and guidelines are provided below.



To obtain stable transfectants, we recommend that you linearize your pcDNA™ 6.2/GeneBLAzer®-DEST construct before transfection. While linearizing the vector may not improve the efficiency of transfection, it increases the chances that the vector does not integrate in a way that disrupts elements necessary for expression in mammalian cells. To linearize your construct, cut at a unique site that is not located within a critical element or within your gene of interest.

Determining Blasticidin Sensitivity

To successfully generate a stable cell line expressing your protein of interest, you need to determine the minimum concentration of Blasticidin required to kill your untransfected host cell line by performing a kill curve experiment (see the procedure below). Typically, concentrations ranging from 2.5 to 10 µg/mL Blasticidin are sufficient to kill most untransfected mammalian cell lines. Blasticidin is available separately (see page 27 for ordering). Refer to page 20 for instructions on how to prepare and store Blasticidin.

1. Plate cells at approximately 25% confluence. Prepare a set of 6 plates.
 2. On the following day, replace the growth medium with fresh growth medium containing varying concentrations of Blasticidin (e.g. 0, 1, 3, 5, 7.5, and 10 µg/mL Blasticidin).
 3. Replenish the selective media every 3–4 days, and observe the percentage of surviving cells.
 4. Count the number of viable cells at regular intervals to determine the appropriate concentration of Blasticidin that prevents growth within 10–14 days after addition of Blasticidin.
-

Continued on next page

Creating Stable Cell Lines, Continued

Generating Stable Cell Lines

Once you have determined the appropriate Blasticidin concentration to use for selection, you can generate a stable cell line expressing your pcDNA[™]6.2/GeneBLAzer[®]-DEST construct.

1. Transfect the mammalian cell line of interest with the pcDNA[™]6.2/cGeneBLAzer[™]-DEST or pcDNA[™]6.2/nGeneBLAzer[™]-DEST expression construct using your transfection method of choice.
 2. 24 hours after transfection, wash the cells and add fresh growth medium.
 3. 48 hours after transfection, split the cells into fresh growth medium such that they are no more than 25% confluent. If the cells are too dense, the antibiotic will not kill the cells. Antibiotics work best on actively dividing cells.
 4. Incubate the cells at 37°C for 2–3 hours until they have attached to the culture dish.
 5. Remove the growth medium and replace with fresh growth medium containing Blasticidin at the predetermined concentration required for your cell line.
 6. Feed the cells with selective media every 3–4 days until Blasticidin-resistant colonies can be identified.
 7. Pick at least 5 Blasticidin-resistant colonies and expand them to assay for recombinant protein expression.
-

Detecting Recombinant Fusion Proteins

Introduction

Depending on the kit you are using, you will assay for β -lactamase reporter activity through *in vivo* or *in vitro* detection methods. A brief description of each detection method is provided below. For detailed information, refer to the GeneBLAzer® Detection Kits manual. If you have generated a pcDNA™6.2/nGeneBLAzer™-DEST expression construct that contains your gene of interest fused to the V5 epitope tag, you may also detect your recombinant fusion protein by Western blot analysis using the Anti-V5 Antibodies.

In Vitro Detection

Using the GeneBLAzer® *In Vitro* Detection Kit allows you to quantitate the amount of intracellular β -lactamase in cells based on the β -lactamase activity in lysates.

To detect β -lactamase activity in mammalian cell lysates, use the CCF2-FA substrate. CCF2-FA is the non-esterified, free acid form of CCF2, and is recommended for *in vitro* use because it is readily soluble in aqueous solution and may be added directly to pre-made cell lysates. Once added to cell lysates, you may quantitate the CCF2-FA fluorescence signal using a fluorescence plate reader or a fluorometer.

To prepare cell lysates from mammalian cells containing the *bla*(M) reporter gene, you **must** use a method that will preserve the activity of the β -lactamase enzyme. Refer to the GeneBLAzer® Detection Kits manual for detailed guidelines and protocols to prepare CCF2-FA solution, prepare cell lysates and samples, and detect CCF2 signal.

In Vivo Detection

Using the GeneBLAzer® *In Vivo* Detection Kit allows you to measure β -lactamase reporter activity in live mammalian cells. Once β -lactamase reporter activity has been measured, cells may be cultured further for use in additional assays or other downstream applications.

To detect β -lactamase activity in live mammalian cells, use the CCF2-AM substrate. CCF2-AM is the membrane-permeable, esterified form of CCF2, and is recommended for *in vivo* use because it is non-toxic, lipophilic, and readily enters the cell. Once cells are “loaded” with CCF2-AM, you may quantitate the CCF2 fluorescence signal using a variety of methods.


Refer to the GeneBLAzer® Detection Kits manual for detailed guidelines and protocols to prepare CCF2-AM solution, load cells with CCF2-AM substrate, and detect CCF2 signal.

Continued on next page

Detecting Recombinant Fusion Proteins, Continued

Detecting the V5 Epitope Tag

If you are using pcDNA[™]6.2/nGeneBLAzer[™]-DEST and you have fused your gene of interest to the V5 epitope tag, you may detect expression of your recombinant fusion protein using the Anti-V5 Antibody, Anti-V5-HRP Antibody, or Anti-V5-AP Antibody available (see page 27 for ordering). In addition, the Positope[™] Control Protein is available for use as a positive control for detection of fusion proteins containing a V5 epitope (see page 27). The ready-to-use WesternBreeze[®] Chromogenic Kits and WesternBreeze[®] Chemiluminescent Kits are available to facilitate detection of antibodies by colorimetric or chemiluminescent methods. For more information, refer to www.lifetechnologies.com/support or contact Technical Support (page 28).



Note

Expression of your protein fused to the β -lactamase reporter and/or to the V5 epitope tag will increase the size of your recombinant protein. The table below lists the increase in the molecular weight of your recombinant protein that you should expect from a particular fusion. Note that the expected sizes take into account any additional amino acids between the gene of interest and the fusion peptide (see pages 8 and 9 for diagrams).

Vector	Fusion	Expected Size Increase (kDa)
pcDNA [™] 6.2/cGeneBLAzer [™] -DEST	β -lactamase (C-terminal)	30 kDa
pcDNA [™] 6.2/nGeneBLAzer [™] -DEST	β -lactamase (N-terminal)	30 kDa
	V5 (C-terminal)	3 kDa

Assay for β -Galactosidase

If you use pcDNA[™]6.2/cGeneBLAzer[™]-GW/*lacZ* or pcDNA[™]6.2/nGeneBLAzer[™]-GW/*lacZ* as a positive control vector, you may assay for β -galactosidase expression by Western blot analysis or activity assay (Miller, 1972). We offer the β -Gal Antiserum, the β -Gal Assay Kit, and the β -Gal Staining Kit for fast and easy detection of β -galactosidase expression (see page 27 for ordering information).

Appendix

Recipes

LB (Luria-Bertani) Medium and Plates

Composition:

1.0% Tryptone
0.5% Yeast Extract
1.0% NaCl
pH 7.0

1. For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 mL deionized water.
2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.
3. Autoclave on liquid cycle for 20 minutes at 15 psi. Allow solution to cool to 55°C and add antibiotic (100 µg/mL ampicillin) if needed.
4. Store at room temperature or at 4°C.

LB agar plates

1. Prepare LB medium as above, but add 15 g/L agar before autoclaving.
 2. Autoclave on liquid cycle for 20 minutes at 15 psi.
 3. After autoclaving, cool to ~55°C, add antibiotic (100 µg/mL of ampicillin), and pour into 10 cm plates.
 4. Let harden, then invert and store at 4°C.
-

Low Salt LB Medium with Blasticidin

Low Salt LB Medium:

10 g Tryptone
5 g NaCl
5 g Yeast Extract

1. Combine the dry reagents above and add deionized, distilled water to 950 mL. Adjust pH to 7.0 with 1 N NaOH. Bring the volume up to 1 liter. For plates, add 15 g/L agar before autoclaving.
 2. Autoclave on liquid cycle at 15 psi and 121°C for 20 minutes.
 3. Allow the medium to cool to at least 55°C before adding the Blasticidin to 100 µg/mL final concentration.
 4. Store plates at 4°C in the dark. Plates containing Blasticidin are stable for up to 2 weeks.
-


Blasticidin

Blasticidin

Blasticidin S HCl is a nucleoside antibiotic isolated from *Streptomyces griseochromogenes* which inhibits protein synthesis in both prokaryotic and eukaryotic cells (Takeuchi *et al.*, 1958; Yamaguchi *et al.*, 1965). Resistance is conferred by expression of either one of two Blasticidin S deaminase genes: *bsd* from *Aspergillus terreus* (Kimura *et al.*, 1994) or *bsr* from *Bacillus cereus* (Izumi *et al.*, 1991). These deaminases convert Blasticidin S to a non-toxic deaminohydroxy derivative (Izumi *et al.*, 1991).

Molecular Weight, Formula, and Structure

The formula for Blasticidin S is $C_{17}H_{26}N_8O_5 \cdot HCl$, and the molecular weight is 458.9. The diagram below shows the structure of Blasticidin.



Handling Blasticidin

Always wear gloves, mask, goggles, and protective clothing (e.g. a laboratory coat) when handling Blasticidin. Weigh out Blasticidin and prepare solutions in a hood.

Preparing and Storing Stock Solutions


Blasticidin may be obtained separately in 50 mg aliquots (see page 27 for ordering). Blasticidin is soluble in water. Sterile water is generally used to prepare stock solutions of 5–10 mg/mL.

- Dissolve Blasticidin in sterile water and filter-sterilize the solution.
- Aliquot in small volumes suitable for one time use (see next to last point below) and freeze at -20°C for long-term storage or store at 4°C for short-term storage.
- Aqueous stock solutions are stable for 1–2 weeks at 4°C and 6–8 weeks at -20°C .
- pH of the aqueous solution should be 7.0 to prevent inactivation of Blasticidin.
- Do not subject stock solutions to freeze/thaw cycles (**do not store in a frost-free freezer**).
- Upon thawing, use what you need and store the thawed stock solution at 4°C for up to 2 weeks.
- Medium containing Blasticidin may be stored at 4°C for up to 2 weeks.

Map and Features of pcDNA™ 6.2/cGeneBLAzer™ -DEST

Map

The map below shows the elements of pcDNA™ 6.2/cGeneBLAzer™ -DEST. DNA from the entry clone replaces the region between bases 918 and 2605. The nucleotide sequence of this vector is available for downloading from www.lifetechnologies.com or by contacting Technical Support (page 28).



Comments for pcDNA™ 6.2/cGeneBLAzer™ -DEST 7519 nucleotides

CMV promoter: bases 232-819

T7 promoter/priming site: bases 863-882

attR1 site: bases 915-1039

Chloramphenicol resistance gene: bases 1148-1807

ccdB gene: bases 2149-2454

attR2 site: bases 2495-2619

β-lactamase *bla*(M) reporter gene: bases 2630-3424

TK polyadenylation signal: bases 3447-3718

f1 origin: bases 3754-4182

SV40 early promoter and origin: bases 4209-4517

EM7 promoter: bases 4572-4638

Blasticidin resistance gene: bases 4639-5037

SV40 early polyadenylation signal: bases 5195-5325

pUC origin (c): bases 5708-6378

Ampicillin resistance gene (c): bases 6523-7383

Ampicillin promoter (c): bases 7384-7482

(c) = complementary strand

Continued on next page

Map and Features of pcDNA™ 6.2/cGeneBLAzer™ -DEST, Continued

Features


pcDNA™ 6.2/cGeneBLAzer™ -DEST (7519 bp) contains the following elements. All features have been functionally tested.

Feature	Benefit
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Allows efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter/priming site	Allows <i>in vitro</i> transcription in the sense orientation and sequencing through the insert
<i>attR1</i> and <i>attR2</i> sites	Allows recombinational cloning of the gene of interest from an entry clone
Chloramphenicol resistance gene	Allows counterselection of plasmid
<i>ccdB</i> gene	Allows negative selection of plasmid
β -lactamase <i>bla(M)</i> reporter gene	Allows fusion of the β -lactamase reporter to the C-terminus of your protein for use as a reporter of gene expression (Zlokarnik <i>et al.</i> , 1998)
Herpes Simplex Virus Thymidine Kinase (TK) polyadenylation signal	Allows efficient transcription termination and polyadenylation of mRNA (Cole & Stacy, 1985)
f1 origin	Allows rescue of single-stranded DNA
SV40 early promoter and origin	Allows efficient, high-level expression of the Blasticidin resistance gene and episomal replication in cells expressing the SV40 large T antigen
EM7 promoter	Allows expression of the Blasticidin resistance gene in <i>E. coli</i>
Blasticidin (<i>bsd</i>) resistance gene	Allows selection of stable transfectants in mammalian cells (Kimura <i>et al.</i> , 1994)
SV40 early polyadenylation signal	Allows efficient transcription termination and polyadenylation of mRNA
pUC origin	Allows high-copy number replication and growth in <i>E. coli</i>
Ampicillin resistance gene	Allows selection of transformants in <i>E. coli</i>

Map and Features of pcDNA™ 6.2/nGeneBLAzer™-DEST

Map

The map below shows the elements of pcDNA™ 6.2/nGeneBLAzer™-DEST. DNA from the entry clone replaces the region between bases 1719 and 3402. The nucleotide sequence of this vector is available for downloading from www.lifetechnologies.com or by contacting Technical Support (page 28).



Comments for pcDNA™ 6.2/nGeneBLAzer™-DEST 7564 nucleotides

CMV promoter: bases 232-819
T7 promoter: bases 863-882
β-lactamase *bla*(M) reporter gene: bases 911-1702
attR1 site: bases 1712-1836
Chloramphenicol resistance gene: bases 1945-2604
ccdB gene: bases 2946-3251
attR2 site: bases 3292-3416
V5 epitope: bases 3424-3465
TK polyadenylation signal: bases 3492-3763
TK polyA reverse priming site: bases 3499-3517
f1 origin: bases 3799-4227
SV40 early promoter and origin: bases 4254-4562
EM7 promoter: bases 4617-4683
Blasticidin resistance gene: bases 4684-5082
SV40 early polyadenylation signal: bases 5240-5370
pUC origin (c): bases 5753-6423
Ampicillin resistance gene (c): bases 6568-7428
Ampicillin promoter (c): bases 7429-7527

(c) = complementary strand

Continued on next page

Map and Features of pcDNA™ 6.2/nGeneBLAzer™ -DEST, Continued

Features


pcDNA™ 6.2/nGeneBLAzer™ -DEST (7564 bp) contain the following elements. All features have been functionally tested.

Feature	Benefit
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Allows efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter	Allows <i>in vitro</i> transcription in the sense orientation
β-lactamase reporter gene	Allows fusion of the β-lactamase reporter to the N-terminus of your protein for use as a reporter of gene expression (Zlokarnik <i>et al.</i> , 1998)
<i>attR1</i> and <i>attR2</i> sites	Allows recombinational cloning of the gene of interest from an entry clone
Chloramphenicol resistance gene	Allows counterselection of plasmid
<i>ccdB</i> gene	Allows negative selection of plasmid
V5 epitope	Allows detection of the recombinant fusion protein by the Anti-V5 antibodies (Southern <i>et al.</i> , 1991).
Herpes Simplex Virus Thymidine Kinase (TK) polyadenylation signal	Allows efficient transcription termination and polyadenylation of mRNA (Cole & Stacy, 1985)
TK polyA reverse priming site	Allows sequencing through the insert
f1 origin	Allows rescue of single-stranded DNA
SV40 early promoter and origin	Allows efficient, high-level expression of the Blasticidin resistance gene and episomal replication in cells expressing the SV40 large T antigen
EM7 promoter	Allows expression of the Blasticidin resistance gene in <i>E. coli</i>
Blasticidin (<i>bsd</i>) resistance gene	Allows selection of stable transfectants in mammalian cells (Kimura <i>et al.</i> , 1994)
SV40 early polyadenylation signal	Allows efficient transcription termination and polyadenylation of mRNA
pUC origin	Allows high-copy number replication and growth in <i>E. coli</i>
Ampicillin resistance gene	Allows selection of transformants in <i>E. coli</i>

Map of pcDNA™ 6.2/cGeneBLAzer™ -GW//lacZ

Description

pcDNA™ 6.2/cGeneBLAzer™-GW//lacZ (8981 bp) is a control vector containing the *lacZ* gene. pcDNA™ 6.2/cGeneBLAzer™-GW//lacZ was constructed using the Gateway® LR recombination reaction between an entry clone containing the *lacZ* gene and pcDNA™ 6.2/cGeneBLAzer™-DEST. *lacZ* is expressed as a fusion to the β-lactamase reporter protein. The molecular weight of the β-galactosidase fusion protein is approximately 147 kDa. The nucleotide sequence of this vector is available for downloading from www.lifetechnologies.com or by contacting Technical Support (page 28).



Comments for pcDNA™ 6.2/cGeneBLAzer™-GW//lacZ 8981 nucleotides


CMV promoter: bases 232-819
T7 promoter/priming site: bases 863-882
attB1 site: bases 915-939
lacZ gene: bases 981-4040
attB2 site: bases 4057-4081
β-lactamase *bla*(M) reporter gene: bases 4092-4886
TK polyadenylation signal: bases 4909-5180
f1 origin: bases 5216-5644
SV40 early promoter and origin: bases 5671-5979
EM7 promoter: bases 6034-6100
Blasticidin resistance gene: bases 6101-6499
SV40 early polyadenylation signal: bases 6657-6787
pUC origin (c): bases 7170-7840
Ampicillin resistance gene (c): bases 7985-8845
Ampicillin promoter (c): bases 8846-8944

(c) = complementary strand

Map of pcDNA™ 6.2/nGeneBLAzer™ -GW//lacZ

Description

pcDNA™ 6.2/nGeneBLAzer™ -GW//lacZ (9041 bp) is a control vector containing the *lacZ* gene. pcDNA™ 6.2/nGeneBLAzer™ -GW//lacZ was constructed using the Gateway® LR recombination reaction between an entry clone containing the *lacZ* gene and pcDNA™ 6.2/nGeneBLAzer™ -DEST. *lacZ* is expressed as a fusion to the β-lactamase reporter protein. Note that the *lacZ* gene contains a TAA stop codon and is **not** fused to the V5 tag. The molecular weight of the β-galactosidase fusion protein is approximately 148 kDa. The nucleotide sequence of this vector is available for downloading from www.lifetechnologies.com or by contacting Technical Support (28).



Comments for pcDNA™ 6.2/nGeneBLAzer™ -GW//lacZ 9041 nucleotides

CMV promoter: bases 232-819
T7 promoter: bases 863-882
β-lactamase *bla*(M) reporter gene: bases 911-1702
*attB*1 site: bases 1712-1736
lacZ gene: bases 1778-4852
*attB*2 site: bases 4869-4893
V5 epitope: bases 4901-4942
TK polyadenylation signal: bases 4969-5240
TK polyA reverse priming site: bases 4976-4994
f1 origin: bases 5276-5704
SV40 early promoter and origin: bases 5731-6039
EM7 promoter: bases 6094-6160
Blasticidin resistance gene: bases 6161-6559
SV40 early polyadenylation signal: bases 6717-6847
pUC origin (c): bases 7230-7900
Ampicillin resistance gene (c): bases 8045-8905
Ampicillin promoter (c): bases 8906-9004

(c) = complementary strand

Accessory Products

Additional Products

Additional products that may be used with the GeneBLAzer™ Gateway® Fusion Kits are available for purchase. Ordering information is provided below.

Product	Amount	Catalog no.
pENTR™ /D-TOPO® Cloning Kit	20 reactions	K2400-20
GeneBLAzer® <i>In Vitro</i> Detection Kit	100 µg	12578-126
GeneBLAzer® <i>In Vivo</i> Detection Kit	50 µg	12578-134
Gateway® LR Clonase® Enzyme Mix	20 reactions	11791-019
	100 reactions	11791-043
Library Efficiency® DB3.1™ Competent Cells	1 mL (5 × 0.2 mL)	11782-018
One Shot® TOP10 Chemically Competent Cells	10 reactions	C4040-10
	20 reactions	C4040-03
One Shot® TOP10 Electrocompetent Cells	10 reactions	C4040-50
	20 reactions	C4040-52
PureLink® HQ Mini Plasmid Purification Kit	100 reactions	K2100-01
Lipofectamine® 2000	1.5 mL	11668-019
	0.75 mL	11668-027
Blasticidin	50 mg	R210-01
β-Gal Assay Kit	100 reactions	K4155-01
β-Gal Staining Kit	1 kit	K1465-01
β-Gal Antiserum*	50 µL	R901-25

*The amount supplied is sufficient to perform 25 Western blots using 10 mL working solution per reaction.

Detecting Fusion Protein

A number of antibodies are available for purchase to detect expression of your fusion protein from the pcDNA™ 6.2/nGeneBLAzer™-DEST vector. Horseradish peroxidase (HRP) or alkaline phosphatase (AP)-conjugated antibodies allow one-step detection using colorimetric or chemiluminescent detection methods. The fluorescein isothiocyanate (FITC)-conjugated antibody allows one-step detection in immunofluorescence experiments.

Antibody	Epitope	Catalog no.
Anti-V5	Detects 14 amino acid epitope derived from the P and V proteins of the paramyxovirus, SV5 (Southern et al., 1991). GKPIPNNLLGLDST	R960-25
Anti-V5-HRP		R961-25
Anti-V5-AP Antibody		R962-25
Anti-V5-FITC Antibody		R963-25
Positope™ Control Protein		R900-50

Technical Support

Obtaining support

For the latest services and support information for all locations, go to www.lifetechnologies.com/support.

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
 - Search through frequently asked questions (FAQs)
 - Submit a question directly to Technical Support (techsupport@lifetech.com)
 - Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
 - Obtain information about customer training
 - Download software updates and patches
-

Safety Data Sheets (SDS)

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/support.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.

Limited warranty

Life Technologies Corporation is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about a Life Technologies product or service, contact our Technical Support Representatives.

All Life Technologies products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Life Technologies makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

**Limited Use Label
License No:
Research Use
Only**

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

**Limited Use Label
License:
ULB ccdB
Selection
Technology**

ccdB selection technology is described in Bernard et al., "Positive Selection Vectors Using the F Plasmid ccdB Killer Gene" Gene 148 (1994) 71-74. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). For licensing information for use in other than research, please contact: outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Gateway® Clone Distribution Policy

Introduction

The information supplied in this section is intended to provide clarity concerning Life Technologies' policy for the use and distribution of cloned nucleic acid fragments, including open reading frames, created using Life Technologies' commercially available Gateway® Technology.

Gateway® Entry Clones

Life Technologies understands that Gateway® entry clones, containing *attL1* and *attL2* sites, may be generated by academic and government researchers for the purpose of scientific research. Life Technologies agrees that such clones may be distributed for scientific research by non-profit organizations and by for-profit organizations without royalty payment to Life Technologies.

Gateway® Expression Clones

Life Technologies also understands that Gateway® expression clones, containing *attB1* and *attB2* sites, may be generated by academic and government researchers for the purpose of scientific research. Life Technologies agrees that such clones may be distributed for scientific research by academic and government organizations without royalty payment to Life Technologies. Organizations other than academia and government may also distribute such Gateway® expression clones for a nominal fee (\$10 per clone) payable to Life Technologies.

Additional Terms and Conditions

We would ask that such distributors of Gateway® entry and expression clones indicate that such clones may be used only for research purposes, that such clones incorporate the Gateway® Technology, and that the purchase of Gateway® Clonase® from Life Technologies is required for carrying out the Gateway® recombinational cloning reaction. This should allow researchers to readily identify Gateway® containing clones and facilitate their use of this powerful technology in their research. Use of Life Technologies' Gateway® Technology, including Gateway® clones, for purposes other than scientific research may require a license and questions concerning such commercial use should be directed to licensing department outlicensing@lifetech.com.

References

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989) Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* 264, 8222-8229
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994) *Current Protocols in Molecular Biology*, Greene Publishing Associates and Wiley-Interscience, New York
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985) A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* 41, 521-530
- Chen, C., and Okayama, H. (1987) High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* 7, 2745-2752
- Chu, G., Hayakawa, H., and Berg, P. (1987) Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nucleic Acids Res.* 15, 1311-1326
- Cole, C. N., and Stacy, T. P. (1985) Identification of Sequences in the Herpes Simplex Virus Thymidine Kinase Gene Required for Efficient Processing and Polyadenylation. *Mol. Cell. Biol.* 5, 2104-2113
- Felgner, P. L., Holm, M., and Chan, H. (1989) Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* 32, 115-121
- Felgner, P. L. a., and Ringold, G. M. (1989) Cationic Liposome-Mediated Transfection. *Nature* 337, 387-388
- Izumi, M., Miyazawa, H., Kamakura, T., Yamaguchi, I., Endo, T., and Hanaoka, F. (1991) Blasticidin S-Resistance Gene (*bsr*): A Novel Selectable Marker for Mammalian Cells. *Exp. Cell Res.* 197, 229-233
- Kimura, M., Takatsuki, A., and Yamaguchi, I. (1994) Blasticidin S Deaminase Gene from *Aspergillus terreus* (*BSD*): A New Drug Resistance Gene for Transfection of Mammalian Cells. *Biochim. Biophys. ACTA* 1219, 653-659
- Kozak, M. (1987) An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nucleic Acids Res.* 15, 8125-8148
- Kozak, M. (1990) Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* 87, 8301-8305
- Kozak, M. (1991) An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biology* 115, 887-903
- Landy, A. (1989) Dynamic, Structural, and Regulatory Aspects of Lambda Site-specific Recombination. *Ann. Rev. Biochem.* 58, 913-949
- Miller, J. H. (1972) *Experiments in Molecular Genetics*, Cold Spring Harbor Laboratory, Cold Spring Harbor, New York
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987) Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Molec. Cell. Biol.* 7, 4125-4129
- Shigekawa, K., and Dower, W. J. (1988) Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* 6, 742-751
- Southern, J. A., Young, D. F., Heaney, F., Baumgartner, W., and Randall, R. E. (1991) Identification of an Epitope on the P and V Proteins of Simian Virus 5 That Distinguishes Between Two Isolates with Different Biological Characteristics. *J. Gen. Virol.* 72, 1551-1557

Continued on next page

References, Continued

- Sutcliffe, J. G. (1978) Nucleotide Sequence of the Ampicillin Resistance Gene of *Escherichia coli* Plasmid pBR322. *Proc. Nat. Acad. Sci. USA* 75, 3737-3741
- Takeuchi, S., Hirayama, K., Ueda, K., Sakai, H., and Yonehara, H. (1958) Blasticidin S, A New Antibiotic. *The Journal of Antibiotics, Series A* 11, 1-5
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977) Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* 11, 223-232
- Yamaguchi, H., Yamamoto, C., and Tanaka, N. (1965) Inhibition of Protein Synthesis by Blasticidin S. I. Studies with Cell-free Systems from Bacterial and Mammalian Cells. *J. Biochem (Tokyo)* 57, 667-677
- Zlokarnik, G., Negulescu, P. A., Knapp, T. E., Mere, L., Burres, N., Feng, L., Whitney, M., Roemer, K., and Tsien, R. Y. (1998) Quantitation of Transcription and Clonal Selection of Single Living Cells with Beta-Lactamase as Reporter. *Science* 279, 84-88

©2012 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the properties of Life Technologies Corporation or their respective owners.



Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit www.invitrogen.com/support or email techsupport@invitrogen.com

www.lifetechnologies.com

