TaqMan® Rodent GAPDH Control Reagents

VIC™ Probe

Protocol



© Copyright 2000, Applied Biosystems

For Research Use Only. Not for use in diagnostic procedures.

All rights reserved.

Printed in the U.S.A.

NOTICE TO PURCHASER: DISCLAIMER OF LICENSE

This product is optimized for use in the Polymerase Chain Reaction ("PCR") and 5' nuclease detection methods covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd. No license under these patents to use the PCR process or 5' nuclease detection methods is conveyed expressly or by implication to the purchaser by the purchase of this product. A license to use the PCR process for certain research and development activities accompanies the purchase of certain Applied Biosystems reagents when used in conjunction with an authorized thermal cycler, or is available from PE Corporation. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc.,1145 Atlantic Avenue, Alameda, California 94501.

ABI PRISM and MicroAmp are registered trademarks of PE Corporation or its subsidiaries in the U.S. and certain other countries. Applied Biosystems, BioLIMS, and MultiScribe are trademarks of PE Corporation or its subsidiaries in the U.S. and certain other countries.

AmpErase and TaqMan are registered trademarks and AmpliTaq Gold is a trademark of Roche Molecular Systems, Inc.

All other trademarks are the sole property of their respective owners.

Contents

Introduction
Purpose of the Kit
Materials and Equipment
Module Components2
Core Kits Supplied by the User
Materials Supplied by User4
Amplicon Size5
Performance Specifications5
Storage and Stability5
Preventing Contamination
Introduction6
AmpErase UNG6
UNG in Two-Step Reactions6
UNG in One-Step Reactions
General Guidelines for PCR
Fluorescent Contaminants
Preparing Reaction Mix Components
Overview9
Preparation of Reagents
Reaction Mix Preparation
PCR Using the TaqMan PCR Core Reagents
Introduction10
PCR Reaction Mix10
Thermal Cycling Parameters for PCR11
One-Step RT-PCR Using The TaqMan EZ RT-PCR Reagents
Introduction12

Description of One-Step RT-PCR
EZ Reaction Mix
Reaction Preparation
Thermal Cycling for One-Step EZ RT-PCR
Run Types
One-Step RT-PCR Using The TaqMan Gold RT-PCR Kit
Introduction 1
Description of One-Step RT-PCR
One-Step Reaction Mix
Reaction Preparation
Thermal Cycling for One-Step Gold RT-PCR
Two-Step RT-PCR Using The TaqMan Gold RT-PCR Kit
Two-Step RT-PCR Reaction
RT-Reaction Mix Preparation
RT-Reaction Preparation
RT-Reaction Thermal Cycling Parameters
PCR Reaction Mix
PCR Reaction Preparation
PCR-Reaction Thermal Cycling Parameters
Interpretation of Results
Normalization
Multi-componenting
Rn and ΔRn Values
Real Time Detection
Plate Read Analysis
Determination of Threshold ΔRn Value
Troubleshooting
Technical Support
To Reach Us On the Web
Hours for Telephone Technical Support
To Reach Us by Telephone or Fax in North America
Documents-on-Demand
To Reach Us by E-Mail
Regional Offices Sales and Service

Appendix A.	References		32
-------------	------------	--	----

Introduction

Purpose of the Kit The TaqMan® Rodent GAPDH Control Reagents are designed to detect rat, mouse, and Chinese hamster glyceraldehyde-3-phosphate dehydrogenase (GAPDH) genes. Either DNA or RNA can be used as a target template.

> The TagMan Rodent GAPDH Control Reagents can be used with TaqMan® PCR Core Reagents or TaqMan® Universal PCR Master Mix to perform polymerase chain reaction (PCR). Both TaqMan® EZ RT-PCR Core Reagents and the TagMan® Gold RT-PCR Kit (one-step and two-step) can be used to perform reverse transcription-polymerase chain reaction (RT-PCR).

Materials and Equipment

Module The TaqMan Rodent GAPDH Control Reagents (P/N 4308313) contain Components the probe and primers sufficient to perform one thousand 50-µL reactions. The module provides enough Control RNA (Rodent) for 100 control reactions. The module includes the following components:

Component	Description
Rodent GAPDH Probe (VIC)	One tube containing 500 µL of 20 µM probe in TE buffer.
Rodent GAPDH Forward Primer	One tube containing 500 µL of 10 µM primer in TE buffer.
Rodent GAPDH Reverse Primer	One tube containing 500 µL of 10 µM primer in TE buffer.
Control RNA (Rodent)	One tube containing 100 µL of 50 ng/µL total mouse RNA in 100 mM Tris-HCl, pH 7.0, 100 mM NaCl, and 1 mM EDTA.

IMPORTANT The TagMan VIC dye must be configured as a Pure Dye on the ABI PRISM® 7700 Sequence Detection System for it to appear on the Reporter pull-down menu. See User Bulletin #4: Generating New Spectra Components (P/N 4306234) pages 6-7 to configure TaqMan VIC as a Pure Dye.

Core Kits Supplied One of the TaqMan core reagent kits listed in the following table is by the User required in addition to the reagents supplied in the TaqMan Rodent GAPDH Control Reagents.

Application	TaqMan Core Reagents ^a	Source
PCR	TaqMan PCR Core Reagents	Applied Biosystems (P/N N808-0228)
	TaqMan Universal PCR Master Mix	Applied Biosystems (P/N 4304437)
RT-PCR	TaqMan Gold RT-PCR Kit	Applied Biosystems (P/N N808-0232)
	TaqMan EZ RT-PCR Core Reagents	Applied Biosystems (P/N N808-0236)
RT	TaqMan [®] Reverse Transcription Reagents	Applied Biosystems (P/N N808-0234)

a. See your local Applied Biosystems representative for a current listing of available reagents.

Equipment Item	Source
ABI PRISM® 7700 Sequence Detection System	See your local Applied Biosystems representative for the instrument
ABI PRISM® 7200 Sequence Detection System/GeneAmp® PCR System 9600	best suited to meet your needs.
TaqMan® LS-50B PCR Detection System/GeneAmp PCR System 9600	

Consumable or Reagent Item	Source
Sequence Detection Primers	Applied Biosystems
Minimum 4000 pmol, unlabeled Minimum 40,000 pmol, unlabeled Minimum 130,000 pmol, unlabeled	(P/N 4304970) (P/N 4304971) (P/N 4304972)
Custom TaqMan Probes	Applied Biosystems
5000 pmol 15,000–25,000 pmol 50,000–100,000 pmol	(P/N 450025) (P/N 450024) (P/N 450003)
MicroAmp® Optical 96-Well Reaction Plate and Optical Caps	Applied Biosystems (P/N 403012)
MicroAmp Optical 96-Well Reaction Plate	Applied Biosystems (P/N N801-0560)
MicroAmp Optical Caps	Applied Biosystems (P/N N801-0935)
MicroAmp Optical Tubes	Applied Biosystems (P/N N801-0933)
96-Well Microplate (Portvair) ^a	Applied Biosystems (P/N L225-1692)
Microcentrifuge	Major laboratory suppliers (MLS)
Centrifuge with adapter for 96-well plate	MLS
Polypropylene tubes	MLS

Consumable or Reagent Item	Source
Pipettors, positive-displacement or air-displacement	MLS
Pipette tips, with filter plugs	MLS
Vortexer	MLS
Disposable gloves	MLS
Tris EDTA (TE) Buffer	MLS
RNase free water	MLS

a. The ABI PRISM 7700 and ABI PRISM 7200 Sequence Detectors use the MicroAmp Optical 96-Well Reaction Plate and MicroAmp Optical Caps. The TaqMan LS-50B PCR Detection System uses the 96-Well Microplate (Portvair). Do not use MicroAmp Optical Tubes with the ABI PRISM 7200 Sequence Detector.

Amplicon Size The amplicon generated from the GAPDH gene for both rodent RNA and DNA is 177 bp in length.

Performance The following performance specifications apply to the TagMan Rodent **Specifications** GAPDH Control Reagents.

Description	Performance Specification
Detection limit of total RNA	Two picograms of Control RNA (Rodent) per 50-µL reaction
Dynamic range of Rodent GAPDH on ABI PRISM 7700 Sequence Detection System	Five log magnitude on the ABI PRISM 7700 Sequence Detection System

Stability

Storage and Upon receipt, store the TaqMan Rodent GAPDH Control Reagents at -15 to -25 °C in a constant-temperature freezer. If stored under the recommended conditions, the product will maintain performance for one year from time of receipt.

Preventing Contamination

Introduction Due to the high throughput and repetitive nature of the 5´nuclease assay, special laboratory practices are necessary in order to avoid false-positive amplifications (Kwok and Higuchi, 1989).

AmpErase UNG AmpErase® UNG (uracil-N-glycosylase, UNG), is a pure nuclease-free, 26-kDa recombinant enzyme encoded by the Eschericia coli uracil-Nglycosylase gene. This gene has been inserted into an E. coli host to direct expression of the native form of the enzyme (Kwok and Higuchi, 1989).

> UNG acts on single- and double-stranded dU-containing DNA. It acts by hydrolyzing uracil-glycosidic bonds at dU-containing DNA sites. The enzyme causes the release of uracil, thereby creating an alkalisensitive apyrimidic site in the DNA. The enzyme has no activity on RNA or dT-containing DNA.

Reactions

UNG in Two-Step When two-step RT-PCR is performed with the TagMan Gold RT-PCR Kit, AmpErase UNG treatment can prevent the reamplification of carry-over PCR products. When dUTP replaces dTTP during PCR amplification, AmpErase UNG treatment can remove up to 200,000 copies of the carry-over GAPDH amplicon per 50-µL reaction.

UNG in One-Step Amperase UNG can be used to remove amplicon contamination in one-Reactions step RT-PCR when using the TaqMan EZ RT-PCR Core Reagents. The rTth DNA Polymerase contained in the kit is thermally stable and is used at temperatures at which Amperase UNG is inactive. Because one-step EZ RT-PCR utilizes dUTP, amplicons generated during this reaction contain uridine residues.

> AmpErase UNG cannot be used when one-step RT-PCR is performed using the TagMan Gold RT-PCR Kit. UNG is active at the temperatures for reverse transcription. The active UNG enzyme would remove uracil bases that are incorporated into the newly synthesized complementary DNA (cDNA) strand (AmpErase UNG Product Insert, 1993).

> If contamination is suspected from previous PCR runs, performing PCR with and without AmpErase UNG will help to identify the source of contamination. To do this, set up parallel No Template Control PCR reactions with and without UNG. A positive signal in the reaction without UNG indicates contamination of reaction components.

General Certain laboratory practices are necessary in order to avoid Guidelines for false-positive amplifications (Kwok and Higuchi, 1989). This is because PCR the PCR process is capable of amplifying single DNA molecules (Saiki et al., 1985; Mullis and Faloona, 1987).

- Wear a clean lab coat (one never worn while handling amplified PCR products or doing sample preparation) and clean gloves when preparing samples for PCR amplification.
- Change gloves whenever contamination is possible.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation
 - PCR setup
 - PCR amplification and detection
- Never bring amplified PCR products into the PCR setup area.
- Open and close all sample tubes carefully. Try not to splash or spray PCR samples.
- Keep reactions and components capped as much as possible.
- Clean lab benches and equipment regularly with 10% bleach solution.

Contaminants

Fluorescent Since fluorescent contaminants may interfere with this assay and give false-positive results, it may be necessary to include a No Amplification Control tube that contains sample and no enzyme. If the absolute fluorescence of the No Amplification Control is greater than that of the No Template Control after PCR, fluorescent contaminants may be present in the sample or in the heat block of the thermal cycler.

Preparing Reaction Mix Components

Overview This procedure is optimized for the TaqMan Rodent GAPDH Control Reagents. The number of reactions depends upon the plate set up by the user. Between 2 pg and 200 ng of total RNA may be used for a onestep RT-PCR reaction. The No Template Control reaction is the complete RT-PCR formulation without the target RNA.

Preparation of The kit is stored at -15 to -25 °C. Prior to use, thaw all reagents except Reagents the enzymes and the RNase Inhibitor. When the reagents are thawed, keep them on ice. Keep the enzymes in a freezer until immediately prior to use. Mix kit components, except enzymes and RNase Inhibitor, by vortexing and using a microcentrifuge to briefly spin down the tube contents. Protect the fluorescent dve labeled probe from excessive exposure to light. When finished with the kit, return it to the -15 to –25 °C freezer.

Reaction Mix Preparing a Reaction Mix of RT-PCR or PCR components is Preparation recommended in order to increase the accuracy of the results. The use of a Reaction Mix reduces the number of reagent transfers and minimizes volume loss due to pipetting.

Step	Action
1	Prepare Reaction Mix by combining all the nonenzymatic components listed in the appropriate table.
2	Mix the components by pipetting up and down.
3	Vortex briefly.
4	Add the enzymatic components (for example: AmpliTaq Gold TM DNA Polymerase, MultiScribe TM Reverse Transcriptase, RNase Inhibitor) listed for the appropriate reaction mix.
5	Mix the components by inverting the microcentrifuge tube.

PCR Using the TaqMan PCR Core Reagents

Introduction PCR can be performed using the TaqMan PCR Core Reagents (P/N N808-0228) by following the instructions described in the TaqMan PCR Core Reagents protocol (P/N 402823).

PCR Reaction Mix The ingredients of a 50 μ L reaction, PCR Reaction Mix are listed in the table below. To make the PCR Reaction Mix, follow the instructions described in "Reaction Mix Preparation" on page 9.

Component	Volume/Tube (μ L)	Final Concentration
RNase-free water	see below ^a	_
10X TaqMan Buffer A	5	1X
25 mM Magnesium chloride	11	5.5 mM
10 mM deoxyATP	1	200 μΜ
10 mM deoxyCTP	1	200 μΜ
10 mM deoxyGTP	1	200 μΜ
20 mM deoxyUTP	1	400 μM
10 μM Rodent GAPDH Forward Primer	0.5	100 nM
10 µM Rodent GAPDH Reverse Primer	0.5	100 nM
20 μM Rodent GAPDH Probe (VIC)	0.5	200 nM
AmpErase UNG	0.5	0.01 Units/µL
AmpliTaq Gold DNA Polymerase	0.25	0.025 Units/μL

a. The volume of RNase-free water will be (27.75 μ L – DNA sample volume).

PCR

Thermal Cycling Use the following thermal cycling parameters for PCR:

Parameters for

Step	UNG Digestion	AmpliTaq Gold Activation	PCR	
	HOLD	HOLD	CYCLE (4	0 cycles)
			Denature	Anneal/ Extend
Time	2 min	10 min	15 sec	1 min
Temp	50 °C	95 °C	95 °C	60 °C

One-Step RT-PCR Using The TaqMan EZ RT-PCR Reagents

Introduction RT-PCR can be performed with the TaqMan EZ RT-PCR Core Reagents (P/N N808-0236) by following the instructions described in the TaqMan EZ RT-PCR Core Reagents protocol (P/N 402877).

One-Step RT-PCR

Description of One-step RT-PCR is RT as well as PCR in a single buffer system using the rTth DNA polymerase enzyme for both steps. The reaction proceeds without the addition of reagents between the RT and PCR steps. This offers the convenience of a single-tube preparation for RT and PCR amplification. The carryover prevention enzyme, AmpErase UNG (uracil-N-glycosylase, UNG), can be used with one-step RT-PCR using the TaqMan EZ RT-PCR Core Reagents.

EZ Reaction Mix The ingredients used for a 50-μL RT-PCR reaction mix using the TaqMan EZ RT-PCR Core Reagents is listed in the table below.

Component	Volume/Tube (μL)	Final Concentration
RNase-free water	see below ^a	_
5X TaqMan EZ Buffer	10	1X
25 mM Manganese acetate	6	3 mM
10 mM deoxyATP	1.5	300 μΜ
10 mM deoxyCTP	1.5	300 μΜ
10 mM deoxyGTP	1.5	300 μΜ
20 mM deoxyUTP	1.5	600 μM
10 µM Rodent GAPDH Forward Primer	0.5	100 nM
10 µM Rodent GAPDH Reverse Primer	0.5	100 nM
20 μM Rodent GAPDH Probe (VIC)	0.5	200 nM
AmpErase UNG	0.5	0.01 U/μL
rTth DNA Polymerase	2	0.1 U/μL

a. The volume of RNase-free water will be (24 μ L – RNA sample volume).

To run a control sample, use 1 µL of the the Control RNA (Rodent) provided with the TagMan Rodent GAPDH Control Reagents.

Reaction For a description of how to prepare the reaction, refer to the TagMan EZ Preparation RT-PCR Core Reagents protocol (P/N 402877).

RT-PCR

Thermal Cycling The following thermal cycling parameters are optimized for the Rodent for One-Step EZ GAPDH system. See thermal cycler manuals for details on operation.

Step	UNG Activation	RT	UNG Deactivation	PC	R
	HOLD	HOLD	HOLD	CYCLE (4	0 cycles)
				Denature	Anneal/ Extend
Time	2 min	30 min	5 min	15 sec	1 min
Temp	50 °C	60 °C	95 °C	94 °C	60 °C

Run Types The TaqMan EZ RT-PCR Kit is designed for two run types: Real Time Detection and Plate Read Detection. Real Time Detection monitors fluorescence during each PCR cycle. Plate Read (endpoint) Detection collects one fluorescence scan per tube after PCR is completed. The features of these detection systems are listed below.

Run Type	Fluorescence Detection	Analysis Results
Real Time	At each PCR cycle	Quantitation of initial template amount
Plate Read Only	Only endpoint	Detection of final amplified products

One-Step RT-PCR Using The TaqMan Gold RT-PCR Kit

Introduction RT-PCR can be performed with the TaqMan Gold RT-PCR Kit (P/N N808-0232) by following the instructions in the TagMan Gold RT-PCR Kit protocol (P/N 402876).

> The TaqMan Gold RT-PCR Kit can be used to perform one-step or twostep RT-PCR.

One-Step RT-PCR

Description of One-step RT-PCR is RT as well as PCR in a single buffer system. The reaction proceeds without the addition of reagents between the RT and PCR steps. This offers the convenience of a single-tube preparation for RT and PCR amplification. However, the carryover prevention enzyme, AmpErase UNG cannot be used with one-step RT-PCR when using the TaqMan Gold RT-PCR Kit.

 $\begin{array}{ccc} \textbf{One-Step Reaction} & \textbf{The ingredients for a 50 } ~\mu \textbf{L RT-PCR reaction mix using the TaqMan} \\ \textbf{Mix} & \textbf{Gold RT-PCR Kit are listed in the table below.} \end{array}$

Component	Volume/Tube (μL)	Final Concentration
RNase-free water	see below ^a	_
10X TaqMan Buffer A	5	1X
25 mM Magnesium chloride	11	5.5 mM
10 mM deoxyATP	1.5	300 μM
10 mM deoxyCTP	1.5	300 μM
10 mM deoxyGTP	1.5	300 µM
20 mM deoxyUTP	1.5	600 µM
10 µM Rodent GAPDH Forward Primer	0.5	100 nM
10 µM Rodent GAPDH Reverse Primer	0.5	100 nM
20 μM Rodent GAPDH Probe (VIC)	0.5	200 nM
RNase Inhibitor	1	0.4 U/μL
MultiScribe Reverse Transcriptase	0.25	0.25 U/μL
AmpliTaq Gold DNA Polymerase	0.25	0.025 U/μL

a. The volume of RNase-free water will be (25 μL – RNA sample volume).

To run a control sample, use 1 μL of the the Control RNA (Rodent) provided with the TaqMan Rodent GAPDH Control Reagents.

Reaction For a description of how to prepare the reaction, refer to the TaqMan **Preparation** Gold RT-PCR Core Reagents protocol (P/N 402876).

Thermal Cycling The following thermal cycling parameters are optimized for the one-step for One-Step Gold RT-PCR for the Rodent GAPDH system. See thermal cycler manuals RT-PCR for details on operation.

Step	RT	AmpliTaq Gold activation	PCR	
	HOLD	HOLD	CYCLE (4	40 cycles)
			Denature	Anneal/ Extend
Time	30 min	10 min	15 sec	1 min
Temperature	48 °C	95 °C	95 °C	60 °C

Two-Step RT-PCR Using The TaqMan Gold RT-PCR Kit

Two-Step RT-PCR Two-step RT-PCR is performed in two separate reactions, reverse Reaction transcription (RT) of the RNA to the cDNA, followed by PCR amplification. This is useful when detecting multiple transcripts from a single cDNA reaction, or when storing a portion of the cDNA for later use. The advantage of the two-step RT-PCR amplification is that the UNG enzyme can be used to prevent carryover contamination.

RT-Reaction Mix Preparation

The ingredients of a 10-µL RT-reaction mix are listed below. A 10-µL reaction volume is sufficient to reverse transcribe up to 200 ng of total RNA. When more than 200 ng of total RNA is used for RT, the reaction volume must be increased. When changing the reaction volume, make sure that the final concentration is consistent with that described in the table below.

Component	Volume/Tube (μL)	Final Concentration
RNase-free water	see below ^a	_
10X TaqMan RT Buffer	1	1X
25 mM Magnesium chloride	2.2	5.5 mM
10 mM dNTP mixture	2.0	500 µM each
Random hexamers or Oligo d(T) ₁₆ ^b	0.5	2.5 μΜ
RNase Inhibitor	0.2	0.4 U/μL
MultiScribe Reverse Transcriptase	0.25	1.25 U/μL

a. The volume of RNase-free water will be $(3.85 \, \mu L - RNA \, sample \, volume)$.

b. If using a sequence-specific reverse primer for the first strand cDNA synthesis, the final concentration should be 200 nM

To run a control sample, use 1 µL of the the Control RNA (Rodent) provided with the TaqMan Rodent GAPDH Control Reagents.

RT-Reaction **Preparation**

For a description of how to prepare the RT-reaction, refer to the TagMan Gold RT-PCR Core Reagents protocol (P/N 402876). The RT reagents are available in a seperate module called TagMan Reverse Transcription Reagents (P/N N808-0234).

Thermal Cycling Parameters

RT-Reaction The cycling parameters for a 10-µL RT-reaction mix using the TaqMan Gold RT-PCR Kit are listed below.

Step	Incubationa	RT	RT Inactivation
	HOLD	HOLD	HOLD
Time	10 min	30 min	5 min
Temperature	25 °C	48 °C	95 °C

a. If using random hexamers or oligo $d(T)_{16}$ primers for first-strand cDNA synthesis, a primer incubation step is necessary at 25 °C for 10 minutes. If using a sequence specific reverse primer, the incubation step is not necessary.

PCR Reaction Mix The ingredients of a 50-µL PCR reaction are listed below.

Component	Volume/Tube (μ L)	Final Concentration
RNase-free water	see below ^a	_
10X TaqMan Buffer A	5	1X
25 mM Magnesium chloride	11	5.5 mM
10 mM deoxyATP	1	200 μM
10 mM deoxyCTP	1	200 μM
10 mM deoxyGTP	1	200 μM
20 mM deoxyUTP	1	400 μM
10 μM Rodent GAPDH Forward Primer	0.5	100 nM
10 µM Rodent GAPDH Reverse Primer	0.5	100 nM
20 μM Rodent GAPDH Probe (VIC)	0.5	200 nM
AmpErase UNG	0.5	0.01 U/μL
AmpliTaq Gold DNA Polymerase	0.25	0.025 U/μL

a. The volume of RNase-free water will be (27.75 μ L – cDNA sample volume).

To run a control sample, use 5 μL of the completed RT-reaction from the previous step.

PCR Reaction For a description of how to prepare the PCR reaction, refer to the Preparation TaqMan Gold RT-PCR Core Reagents protocol (P/N 402876).

Parameters

 $\begin{tabular}{ll} \textbf{PCR-Reaction} & \textbf{The cycling parameters for the PCR step (step two) of a 50-μL, two-step} \end{tabular}$ Thermal Cycling RT-PCR reaction using the TaqMan Gold RT-PCR Kit are listed below.

Step	UNG Incubation	AmpliTaq Gold Activation	PC	CR
	HOLD	HOLD	CYCLE (4	10 cycles)
			Denature	Anneal/ Extend
Time	2 min	10 min	15 sec	1 min
Temp	50 °C	95 °C	95 °C	60 °C

Interpretation of Results

Normalization The Passive Reference is a dye included in the 10X TaqMan Buffer A and in the 5X TagMan EZ Buffer and does not participate in the 5' nuclease assay. The Passive Reference provides an internal reference to which the reporter dye signal can be normalized during data analysis. Normalization is necessary to correct for fluorescent fluctuations due to changes in concentration or volume.

Multicomponenting

Multicomponenting is the term used to distinguish the contribution each individual dye makes to the fluorescent spectra. The overlapping spectra from the three pure dye components generate the composite spectrum. This spectrum represents one fluorescent reading from one well. For the TagMan Rodent GAPDH Control system, the three dyes used for Multicomponent Analysis are:

Reporter	VIC
Quencher	TAMRA
Passive Reference	ROX

 R_n and ΔR_n Values Normalization is accomplished by dividing the emission intensity of the reporter dye by the emission intensity of the Passive Reference to obtain a ratio defined as the R_n (normalized reporter) for a given reaction tube.

> R_n⁺ is the R_n value of a reaction containing all components including the template.

> R_n is the R_n value of an unreacted sample. This value may be obtained from the early cycles of a Real Time run, those cycles prior to a detectable increase in fluorescence. This value may also be obtained from a reaction not containing template.

> ΔR_n is the difference between the R_n^+ value and the R_n^- value. It reliably indicates the magnitude of the signal generated by the given set of PCR conditions.

The following equation expresses the relationship of these terms:

$$\Delta \mathsf{R}_\mathsf{n} = (\mathsf{R}_\mathsf{n}^+) - (\mathsf{R}_\mathsf{n}^-)$$

where:

$R_n^+ =$	Emission Intensity of Reporter	PCR with template
''n -	Emission Intensity of Passive Reference	1 Off Will template
R_n^- =	Emission Intensity of Reporter Emission Intensity of Passive Reference	PCR without template or early cycles of a Real Time reaction

Detection

Real Time The threshold cycle or C_T value is the cycle at which a statistically significant increase in ΔR_n is first detected. Threshold is defined as the average standard deviation of $R_{\rm n}$ for the early cycles, multiplied by an adjustable factor.

> On the graph of $\rm R_n$ versus cycle number (figure below), the threshold cycle occurs when the Sequence Detection Application begins to detect the increase in signal associated with an exponential growth of PCR product.




Plate Read Plate Read analysis relies on endpoint data to calculate ΔR_n . To ensure Analysis statistically high confidence levels using either the GAPDH probe or a custom probe, run the protocol with at least three No Template Controls per microplate. A positive result is defined as ΔR_n values greater than the threshold ΔR_n . The threshold ΔR_n is calculated by multiplying the standard deviation of three R_n-values by an appropriate multiplier. The value for the appropriate multiplier is selected from a table of t-distribution values and depends upon the desired confidence level (Beyer, 1984). For 99% confidence levels, the multiplier is 6.965.

> When more than three No Template Controls are run, the multiplier for the standard deviation decreases. Refer to a table of *t*-distribution values for the appropriate multiplier (Beyer, 1984).

To evaluate reproducibility, calculate the coefficient of variation on replicate samples. Inconsistent results (coefficients of variation exceeding 10%) may be caused by pipetting errors and incomplete mixing.

Threshold ∆R_n Value

Determination of To determine threshold ΔR_n value:

- On one microplate, measure both reporter and passive reference fluorescence, and determine the normalized reporter for each No Template Control tube (R_n^-) .
- Determine the mean and standard deviation of R_n⁻.
- When running three No Template Controls, multiply the R_n⁻ standard deviation by 6.965 to determine the threshold ΔR_n for the system.

For example, if the mean R_n- for a system is 0.5 with a standard deviation of 0.02, the threshold ΔR_n is 6.965 \times 0.02 or 0.14. Any $\Delta R_n > 0.14$ is a positive result and indicates that the sample contains target.

Troubleshooting

Observation	Possible Cause	Recommended Action
$\Delta R_n \le No template$	Inappropriate reaction conditions	Troubleshoot RT-PCR optimization.
Control ΔR_n , and no amplification plot	Incorrect dye components chosen	Check dye component prior to data analysis.
	Reaction component omitted	Check that all the correct reagents were added.
	Incorrect primer or probe sequence	Resynthesize with appropriate sequence.
	Degraded template or no template added	Repeat with fresh template.
	Reaction inhibitor present	Repeat with purified template.
$\Delta R_n \le No Template$	Amplicon contamination of reagents	Check technique and equipment to
Control ΔR_n , and both reactions show an amplification plot	Template contamination of reagents	confine contamination. Use fresh reagents.
Shifting R _n value during	Fluorescent emissions have not	Reset lower value of baseline range.
the early cycles of PCR (cycles 0–5)	stabilized to new buffer conditions of reaction mix. This does not affect PCR, or the final results.	Add probe to the buffer component and allow it to equilibrate at room temperature prior to Reagent Mix formulation.
Abnormal amplification plot:	C _T value <15, amplification signal detected in early cycles	Reset upper value of baseline range.
0.100 ΔRn -0.450 0 Cycle 40		Dilute the sample to increase the \mathbf{C}_{T} value.
Multicomponent signal for ROX is not linear	Pure dye component's spectra are incorrect	Rerun pure dye spectra.
	Incorrect dye components choosen	Choose correct dyes for data analysis.
Small ∆R _n	PCR efficiency is poor	Reoptimize reaction conditions.
	Low copy number of target	Increase starting copy number.

Observation	Possible Cause	Recommended Action
C _T value is higher than	Less template added than expected	Increase sample amount.
expected	Sample is degraded	Evaluate sample integrity.
C _T value is lower than	More sample added than expected	Reduce sample amount.
expected	Template or amplicon contamination	Review "General Guidelines for PCR" on page 8.
Standard deviation of C _T value >0.16	Inaccurate pipetting	Prepare a Reagent Mix. Refer to "Reaction Mix Preparation" on page 9.
		Use positive-displacement pipettors.

Technical Support

To Reach Us on the Web

To Reach Us on the Applied Biosystems web site address is:

http://www.appliedbiosystems.com/techsupport

We strongly encourage you to visit our web site for answers to frequently asked questions, and to learn more about our products. You can also order technical documents and/or an index of available documents and have them faxed or e-mailed to you through our site (see the "Documents on Demand" section below).

Hours for Telephone Technical Support

Hours for In the United States and Canada, technical support is available at the **Telephone** following times.

Product	Hours
Chemiluminescence	9:00 a.m. to 5:00 p.m. Eastern Time
LC/MS	9:00 a.m. to 5:00 p.m. Pacific Time
All Other Products	5:30 a.m. to 5:00 p.m. Pacific Time

See the "Regional Offices Sales and Service" section below for how to contact local service representatives outside of the United States and Canada.

To Reach Us by Telephone or Fax in North America

Call Technical Support at 1-800-831-6844, and select the appropriate option (below) for support on the product of your choice at any time during the call. (To open a service call for other support needs, or in case of an emergency, press **1** after dialing 1-800-831-6844.)

For Support On This Product	Dial 1-800-831-6844, and		
ABI PRISM® 3700 DNA Analyzer	Press	FAX	
Analyzei	8	650-638-5981	
ABI PRISM® 3100 Genetic Analyzer	Press	FAX	
7 111017 201	26	650-638-5891	
DNA Synthesis	Press	FAX	
	21	650-638-5981	

For Support On This Product	Dial 1-800-831-6844	, and
Fluorescent DNA Sequencing	Press	FAX
Sequencing	22	650-638-5891
Fluorescent Fragment Analysis (includes	Press	FAX
GeneScan® applications)	23	650-638-5891
Integrated Thermal Cyclers	Press	FAX
	24	650-638-5891
BioInformatics (includes BioLIMS™, BioMerge™, and	Press	FAX
SQL GT™ applications)	25	505-982-7690
PCR and Sequence	Press	FAX
Detection	5, or call 1-800-762-4001, and press 1 for PCR, or 2 for Sequence Detection	240-453-4613
FMAT	Telephone	FAX
	1-800-899-5858, and press 1, then press 6	508-383-7855
Peptide and Organic Synthesis	Press	FAX
Synthesis	31	650-638-5981
Protein Sequencing	Press	FAX
	32	650-638-5981
Chemiluminescence	Telephone	FAX
	1-800-542-2369 (U.S. only), or	781-275-8581 (Tropix)
	1-781-271-0045 (Tropix)	9:00 a.m. to 5:00 p.m. ET

For Support On This Product	Dial 1-800-831-6844	Dial 1-800-831-6844, and	
LC/MS	Telephone	FAX	
	1-800-952-4716	650-638-6223	
		9:00 a.m. to 5:00 p.m. PT	

Documents on Pree 24-hour access to Applied Biosystems technical documents, including MSDSs, is available by fax or e-mail.

You can access Documents on Demand through the internet or by telephone:

If you want to order	Then
through the	Use http://www.appliedbiosystems.com/techsupport
internet	You can search for documents to order using keywords.
	Up to five documents can be faxed or e-mailed to you by title.
by phone from the United States or	a. Call 1-800-487-6809 from a touch-tone phone. Have your fax number ready.
Canada	 b. Press 1 to order an index of available documents and have it faxed to you. Each document in the index has an ID number. (Use this as your order number in step "d" below.)
	c. Call 1-800-487-6809 from a touch-tone phone a second time.
	d. Press 2 to order up to five documents and have them faxed to you.
by phone from outside the	a. Dial your international access code, then 1-858-712-0317, from a touch-tone phone.
United States or Canada	Have your complete fax number and country code ready (011 precedes the country code).
	 b. Press 1 to order an index of available documents and have it faxed to you. Each document in the index has an ID number. (Use this as your order number in step "d" below.)
	c. Call 1-858-712-0317 from a touch-tone phone a second time.
	d. Press 2 to order up to five documents and have them faxed to you.

E-Mail areas.

To Reach Us by Contact technical support by e-mail for help in the following product

For this product area	Use this e-mail address
Chemiluminescence	info@appliedbiosystems.com
Genetic Analysis	galab@appliedbiosystems.com
LC/MS	apisupport@sciex.com
PCR and Sequence Detection	pcrlab@appliedbiosystems.com
Protein Sequencing, Peptide and DNA Synthesis	corelab@appliedbiosystems.com

Regional Offices If you are outside the United States and Canada, you should contact Sales and Service your local Applied Biosystems service representative.

The Americas			
United States Applied Biosystems 850 Lincoln Centre Drive		Latin America (Del.A. Obregon, Mexico)	
	y, California 94404	Tel:(305) 670-4350 Fax: (305) 670-4349	
,	650) 570-6667 800) 345-5224		
Fax: (650) 572-2743		

Europe				
Austria (Wien)		Hungary (Budapest)		
Tel:	43 (0)1 867 35 75 0	Tel: 36 (0)1 2		
Fax:	43 (0)1 867 35 75 11	Fax: 36 (0)1 2	70 8288	
Belgiu	m	Italy (Milano)		
	32 (0)2 712 5555 32 (0)2 712 5516	Tel: 39 (0)39 8 Fax: 39 (0)39 8		
Czech Republic and Slovakia (Praha)		The Netherlands (Nieuwerkerk a/d IJssel)		
_	420 2 61 222 164 420 2 61 222 168	Tel: 31 (0)180 Fax: 31 (0)180		
Denmark (Naerum)		Norway (Oslo)		
	45 45 58 60 00 45 45 58 60 01	Tel: 47 23 12 Fax: 47 23 12		

Furon	Europe			
Finland (Espoo) Tel: 358 (0)9 251 24 250		Poland, Lithuania, Latvia, and Estonia (Warszawa)		
Fax:	358 (0)9 251 24 243	Tel: Fax:	- ()	
France	(Paris)	Portug	al (Lisboa)	
	33 (0)1 69 59 85 85 33 (0)1 69 59 85 00	Tel: Fax:	(-)	
Germa	ny (Weiterstadt)	Russia	(Moskva)	
Tel: Fax:	49 (0) 6150 101 0 49 (0) 6150 101 101	Tel: Fax:	7 095 935 8888 7 095 564 8787	
Spain ((Tres Cantos)	South	South Africa (Johannesburg)	
Tel: Fax:	34 (0)91 806 1210 34 (0)91 806 1206	Tel: Fax:	27 11 478 0411 27 11 478 0349	
	Sweden (Stockholm) Tel: 46 (0)8 619 4400		United Kingdom (Warrington, Cheshire)	
	46 (0)8 619 4401	Tel: Fax:	44 (0)1925 825650 44 (0)1925 282502	
Switze	rland (Rotkreuz)	South East Europe (Zagreb, Croatia)		
	41 (0)41 799 7777 41 (0)41 790 0676	Tel: Fax:	385 1 34 91 927 385 1 34 91 840	
Middle Eastern Countries and North Africa (Monza, Italia)			(English Speaking) and West Fairlands, South Africa)	
	39 (0)39 8389 481 39 (0)39 8389 493	Tel: Fax:		
	er Countries Not Listed ngton, UK)			
	44 (0)1925 282481 44 (0)1925 282509			

Japan

Japan (Hatchobori, Chuo-Ku, Tokyo)

Tel: 81 3 5566 6100 Fax: 81 3 5566 6501

Eastern Asia, China, Oceania				
Australia (Scoresby, Victoria)		Malaysia (Petaling Jaya)		
-	61 3 9730 8600 61 3 9730 8799	_	60 3 758 8268 60 3 754 9043	
China (Beijing)		Singap	Singapore	
-	86 10 6238 1156 86 10 6238 1162	_	65 896 2168 65 896 2147	
Hong k	Kong	Taiwan (Taipei Hsien)		
1	852 2756 6928 852 2756 6968	Tel: Fax:	886 2 2698 3505 886 2 2698 3405	
Korea (Seoul)		Thailaı	nd (Bangkok)	
	82 2 593 6470/6471 82 2 593 6472	Tel: Fax:	66 2 719 6405 66 2 319 9788	

Appendix A. References

Beyer, W.H., Ed., 1984. *CRC Standard Mathematical Tables*. 27th ed. CRC Press, Boca Raton, FL.

Kwok, S., and Higuchi, R. 1989. Avoiding false positives with PCR. *Nature* 339: 237–238.

Mullis, K.B., and Faloona, F.A. 1987. Specific synthesis of DNA *in vitro* via a polymerase-catalyzed chain reaction. In: Wu, R., ed. *Methods in Enzymology*, Vol. 155, Academic Press, Inc., San Diego, CA, pp 335–350.

Saiki, R. K., Scharf, S. J., Faloona, F. A., Mullis, K. B., Horn, G. T., Erlich, H. A., and Amheim, N. 1985. Enzymatic amplification of β -globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science* 230: 1350–1354.

Worldwide Sales Offices

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches into 150 countries on six continents. For international office locations, please call our local office or refer to our web site at www.appliedbiosystems.com.

Headquarters

850 Lincoln Centre Drive Foster City, CA 94404 USA Phone: +1 650.638.5800 Toll Free: +1 800.345.5224 Fax: +1 650.638.5884

Technical Support

For technical support:

Toll Free: +1 800.831.6844 ext 23

Fax: +1 650.638.5891

www.appliedbiosystems.com



PE Corporation is committed to providing the world's leading technology and information for life scientists. PE Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Printed in the USA, 09/2000 Part Number 4308318B