

# VIASURE MULTIPLEX


## *BK + JC Virus Real Time PCR Detection Kit*

### Pathogen and product description

**T**he two most commonly known human polyomaviruses, which are ubiquitous viruses belonging to the *Papovaviridae* family, are BK and JC viruses. Most humans acquire it during childhood, usually remaining asymptomatic. BKV is frequently associated with renal system, such as ureteral stenosis in kidney transplant recipients, hemorrhagic cystitis in hematopoietic stem cell transplant recipients and polyomavirus-associated nephropathy. Actually, early diagnosis of BKV has been shown to positively influence the organ survival. JC virus causes central nervous system diseases in immunocompromised patients, including progressive multifocal leukoencephalopathy (PML), JCV granule cell neuronopathy, and JCV encephalopathy.

Real Time PCR assays have been shown to be a sensitive and specific diagnosis tool for the detection of BK Virus and JC Virus. BKV can be mainly detected in urine and plasma (blood and serum specimens as well). JC virus can be mostly identified in cerebrospinal fluid (CSF) and plasma (blood), as well as urine (less usual).

VIASURE *BK+JC Virus* Real Time PCR Detection Kit is designed for the diagnosis of BK Virus and/or JC Virus in clinical samples. After DNA isolation, the identification of BK Virus and JC Virus is performed by the amplification of a conserved region of the *VP1* gene for BK Virus, *T* gene for JC Virus, using specific primers and a fluorescent-labelled probe.



**Analytical sensitivity**

**VIASURE BK + JC Virus Real Time PCR Detection Kit** has a detection limit of  $\geq 10$  DNA copies per reaction for BK and JC Virus. (figures 1 and 2).




Figure 1. Dilution series of JC Virus ( $10^7$ – $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).





Figure 2. Dilution series of BK Virus ( $10^7$ – $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5).

**Components**

Reagent/Material	Description	Colour	Quantity
<b>BK + JC Virus 8-well strips</b>	A mix of enzymes, primers-probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 X 8-well strip
<b>Rehydration Buffer</b>	Solution to reconstitute the stabilized product	Blue	1 vial x 1,8 mL
<b>BK + JC Virus Positive Control</b>	Non-infectious synthetic lyophilized cDNA	Red	1 vial
<b>Negative Control</b>	Non template control	Violet	1 vial x 1 mL
<b>Water RNase/DNase free</b>	Water RNase/DNase free	White	1 vial x 1 mL
<b>Tear-off 8-cap strips</b>	Optical caps for sealing Wells during thermal cycling	Transparent	6/12 x 8-cap strip

**Work Flow**

One-step rehydration of wells and add your extracted DNA




**STEP 1**  
Add 15 µl of rehydration buffer into each well



**STEP 2**  
Add 5 µl of DNA sample / positive control / negative control



**STEP 3**  
Load the strips into the thermocycler and run the specified protocol



**STEP 4**  
Interpretate results

**Kit References**

Reference	Description
VS-BJV106L	Viasure BK + JC Virus Real Time PCR Detection Kit 6 x 8-well strips, low profile
VS-BJV106H	Viasure BK + JC Virus Real Time PCR Detection Kit 6 x 8-well strips, high profile
VS-BJV112L	Viasure BK + JC Virus Real Time PCR Detection Kit 12 x 8-well strips, low profile
VS-BJV112H	Viasure BK + JC Virus Real Time PCR Detection Kit 12 x 8-well strips, high profile
VS-BJV113L	Viasure BK + JC Virus Real Time PCR Detection Kit 96-well plate, low profile
VS-BJV113H	Viasure BK + JC Virus Real Time PCR Detection Kit 96-well plate, high profile



CERTEST BIOTEC, S.L.  
Pol. Industrial Río Gállego II, Calle J, Nº 1,  
50840, San Mateo de Gállego, Zaragoza (SPAIN)  
www.certest.es

