
CerTest
BIOTEC

CerTest Turbilatex

FOB TURBILATEX® Combo

FOB latex turbidimetric assay

Cat.TL-022FB100ED, TL-022FB200ED, TL-022FB400ED, TL-022FB800ED

IU-TL -022FB100ED-200ED-400ED-800ED en-es v.0220

ENGLISH

INTENDED USE

FOB Turbilatex® is a latex turbidimetric assay **for the quantitative detection of haemoglobin (faecal occult blood) in human stool samples.**

This assay is simple and widely applicable. This product is optimized for several automated analyser.

For professional *in vitro* diagnostic use only.

SYNTHESIS

Colorectal cancer is the second leading cause of illness and death in Western world. The screening with faecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps. Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning.

Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases.

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut. FOB Turbilatex® tests detect intact or nearly intact human haemoglobin, being a very specific technique for detecting loss of blood from the lower intestine.

PRINCIPLE

FOB Turbilatex® is a latex turbidimetric assay **for the quantitative detection of haemoglobin in human stool samples.** This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is to make a presumptive diagnosis of faecal occult blood (gastrointestinal bleeding).

FOB latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If

antigen is present in the sample it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at determinate wavelength and is proportional to the quantity of antigen contained in the sample.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- A trained person in Turbidimetric technique and autoanalyzer use is required.
- The analyser should be ready before performing any assay.
- Read and follow the instructions for use provided with the kit.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). Do not freeze.

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

REAGENTS

Materials provided

References	Cat.: <u>TL-022FB100ED</u>	Cat.: <u>TL-022FB200ED</u>	Cat.: <u>TL-022FB400ED</u>	Cat.: <u>TL-022FB800ED</u>
- TL-022FB01 FOB Turbilatex® Reagent 1	1x22mL	2 x 22mL	4 x 22mL	8 x 22mL
- TL-022FB02 FOB Turbilatex® Reagent 2	1 x 8mL	1 x 13mL	2 x 13mL	4 x 13mL

References	Cat.: <u>TL-022FB100ED</u>	Cat.: <u>TL-022FB200ED</u>	Cat.: <u>TL-022FB400ED</u>	Cat.: <u>TL-022FB800ED</u>
- TL-022FB70/ TL-022FB71/ TL-022FB72- TL-022FB73/ TL-022FB74/ TL-022FB75 FOB Turbilatex® Calibrators			6 x 1mL	
- TL-022FB08/ TL-022FB09 FOB Turbilatex® Control 1&2			4 x 1mL	

Materials and instruments that could be required but no provided

- Specimen collection container.
- TL-022UN03E Universal Turbilatex® Sample Diluent (Required for the Procedure A).
- TL-0221213E Microtubes.
- Disposable gloves and laboratory equipment.
- Automated analyser.
- MST-0018MU Universal Turbilatex® Sample Collection Vial (Required for the Procedure B).
- MST-0019U Universal Turbilatex® Sample Collection Vial (Required for the Procedure B).

SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of human stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 3 days prior to testing. If not immediately tested, freeze the stored samples at -20 °C maximum 6 months. In this case, the sample will be totally thawed, and brought to room temperature (15-30°C) before testing. Homogenize stool samples as thoroughly as possible prior to preparation.

To process the collected stool samples:

Procedure A: using catalogue reference TL-022UN03E

Use a microtube for each sample to be tested. Label microtube with name or number of patient.

1. Homogenize the sample. Add 20 mg of sample into microtube.
2. Add 2mL of sample diluent (TL-022UN03E).

3. Shake vigorously the microtube in order to assure good sample dispersion (vortex) until sample is completely dissolved.
4. Centrifuge for 15 minutes at 10000 g or 10 minutes at 15000 g.
5. Take the supernatant to automated analyser vial. Diluted sample can be stored for 7 days in the refrigerator (2-8 °C).

Procedure B: using catalogue reference MST-0018MU/MST-0019U.

Follow the indication of the instruction for use of the product.

ANALYTICAL PROCEDURE

R1 and R2 are ready to use.

Preparation of the calibration curve


For calibration use only the following materials:

TL-022FB70, TL-022FB71, TL-022FB72, TL-022FB73, TL-022FB74 and TL-022FB75, FOB Turbilatex® Calibrator vials: liquid calibrator containing human haemoglobin at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a month is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

FOB Turbilatex® Calibrator vials (calibrator 0 to 5) are ready to use.

Picture 1. Prepared calibrators for performing the calibration curve (5 points and blank)

Quality control


For quality control only use the following materials:

TL-022FB08 FOB Turbilatex® Control 1 and TL-022FB09 FOB Turbilatex® Control 2: liquid controls at two different concentrations of human haemoglobin. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

FOB Turbilatex® Control vials (Control 1 and Control 2) are ready to use.

Picture 2. Prepared controls C1 and C2 (check control concentration on the label)

Analytical procedure, see annex information.

INTERPRETATION OF RESULTS

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 5 µg of hHb/g of stool for diagnostic procedures and 20 µg of hHb/g of stool for screening procedures.

Positive results determine the abnormal presence of human Haemoglobin (hHb) in stool samples.

FOB Turbilatex® cut-off value:

Haemoglobin concentration values lower than 50 ng of hHb/mL of stool are considered normal values and that is not indicative of a bleeding in the gastrointestinal tract.

Haemoglobin concentration values between 50 and 200 ng hHb/mL of stool are considered as an abnormal presence being indicative of slight bleeding, therefore monitoring and follow-up of the patient is recommended.

Haemoglobin concentration values higher than 200 ng of hHb/mL of stool are indicative of a strong bleeding in the gastrointestinal tract.

FOB Turbilatex® cut-off value: 50-200 ng hHb/mL (x100 dilution factor, 5-20 µg hHb/g in stool samples).

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Limit of Detection (LOD): **see annex information.**

Prozone

Studies have been made up to a concentration of 1 mg of hHb/mL of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Clinical sensitivity and specificity

An evaluation was performed comparing a turbidimetric assay (FOB Turbilatex®, CerTest) and another commercial immunoassay (Eiken FOB Latex, EIKEN Chemicals, CO. LTD).

The results were as follows:

	Sensitivity	Specificity
FOB Turbilatex® vs EIKEN FOB Latex	97.5%	100%

The results showed a high sensitivity and specificity to detect human haemoglobin (hHb) using FOB Turbilatex®.

Interferences

An evaluation was performed to determine the interferences of FOB Turbilatex®. No interferences were founded against other substances tested: BSA, Hemin (concentration equal or lower to 0.25 mg/mL), Ascorbic acid.

Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against other faecal markers occasionally present in faeces, such as:

human calprotectin, bovine and pig haemoglobin, human lactoferrin, bovine and human transferrin.

The cross reactions with bovine, pig haemoglobins and bovine transferrin are not present at concentrations below 0.2 µg/mL, 0.9 µg/mL and 250 µg/mL respectively. Take note that these concentrations are higher than the concentrations that could be found on the faecal samples.

LIMITATIONS

1. FOB Turbilatex® should **only be used in human stool samples**. Proper faecal specimens must be obtained.
2. Positive results determine the presence of human haemoglobin in faecal samples. A positive result should be followed up with additional diagnostic invasive procedures to determine the exact cause and source of the blood in the stool.
3. If symptoms or situation still persist, haemoglobin determination should be carried out invasive techniques. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
4. This test may be less sensitive for detecting upper gastrointestinal bleeding because blood degrades as it passes through the gastrointestinal tract.
5. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

EXPECTED VALUES

Colorectal cancer affects both men and women of all racial and ethnic groups, being most often found in people aged 50 years or older. Regarding men, colorectal cancer is the third most common cancer after prostate and lung cancers. Regarding women, colorectal cancer is the third most common cancer after breast and lung cancers.

Some studies established equal or higher 5.1 µg hHb/g faeces as cut-off values to allow detect adult patients with gastrointestinal bleeding that will require to diagnosis

additional diagnostic invasive procedures. In FOB Turbilatex® test, this value is equivalent to 50 ng hHb/mL of diluent.

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1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2)CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO *IN VITRO*

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website: www.certest.es

Note/Nota: Turbilatex® registered brand name/ Turbilatex® marca registrada
Nº: CTM-010179299; community (OAMI)/ Nº: CTM-010179299; comunitaria (OAMI)
Nº: 4,650,387 ; EEUU (UPSTO)

CerTest Turbilatex

One step ahead

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)

Tel. (+34) 976 520 354

Fax (+34) 976 106 268

certest@certest.es | turbilatex@certest.es

www.certest.es