INSTRUCTIONS



ZebaTM Micro Spin Desalting Columns, 7K MWCO

89877 89878

1500.3

Number Description

Zeba Micro Spin Desalting Columns, 75μL, 25 columns
 Zeba Micro Spin Desalting Columns, 75μL, 50 columns

Note: Each column contains 150µL of slurry in 0.03% sodium azide. This product is recommended for

processing compounds > 7,000 Da.

Storage: Upon receipt store at 4°C. Product is shipped at ambient temperature.

Introduction

The Thermo Scientific Zeba Micro Spin Desalting Columns contain a high-performance resin that offers exceptional desalting and protein recovery characteristics. Small sample volumes (2-12 μ L) containing as low as 20 μ g/mL of protein can be processed while providing \geq 95% retention of salts and other small molecules (< 1,000 Da). These columns require no chromatography system or cumbersome column preparation or equilibration. Additionally, the spin-column method eliminates the need to wait for samples to emerge by gravity flow, allowing multiple-sample processing in ~6 minutes.

Desalting Procedure

For maximal protein recovery, add 3μ L of buffer or water above the sample as a stacker. This stacker may be omitted if desired; however, a reduction in protein recovery will result.

A. Additional Materials Required

- Variable-speed bench-top microcentrifuge
- 1.5 microcentrifuge collection tubes

B. Spin Column Preparation

- 1. Remove column's bottom closure and loosen cap (do not remove cap).
- Place column in a 1.5-2.0mL microcentrifuge collection tube and centrifuge at 1,000 × g for 1 minute to remove storage solution.
- 3. Blot bottom of column to remove excess liquid.

C. Sample Loading

- 1. Place column in a new collection tube, remove cap and apply 2-12μL of sample to the top of the compact resin bed.
- 2. Apply a 3µL stacker of ultrapure water or buffer to the top of the resin bed after the sample has absorbed.
- 3. Centrifuge at $1,000 \times g$ for 2 minutes to collect the desalted sample. Discard desalting column after use.

Buffer Exchange Procedure

A. Additional Materials Required

- Variable-speed bench-top microcentrifuge
- 1.5-2.0mL microcentrifuge collection tubes
- Buffer for exchange



B. Spin Column Preparation

- 1. Remove column's bottom closure and loosen cap (do not remove cap).
- 2. Place column in a 1.5-2.0mL microcentrifuge collection tube.
- 3. Centrifuge at $1,000 \times g$ for 1 minute to remove storage solution.
- 4. Add 50μ L of buffer on top of the resin bed. Centrifuge at $1,000 \times g$ for 1 minute to remove buffer.
- 5. Repeat step 4 two to three additional times, discarding buffer from the collection tube.

C. Sample Loading

- 1. Place column in a new collection tube, remove cap and apply 2-12µL of sample to the top of the compact resin bed.
- 2. Apply a 3µL stacker of buffer to the top of the resin bed after the sample has absorbed.
- 3. Centrifuge at $1,000 \times g$ for 2 minutes to collect the sample. Discard desalting column after use.

Troubleshooting

Problem	Possible Cause	Solution
Sample or buffer does not flow through resin	Centrifugation problem	Ensure that centrifuge is in proper working condition
		Ensure bottom closure is removed
Sample contamination	Improper sample loading	Load sample directly to the center of the resin bed; carefully touch pipette tip to resin to expel all sample
		Avoid contact with sides of column
	Centrifugation problem	Do not exceed recommended speed or centrifugation time
Low yield	No stacker used	Apply a 3µL water or buffer stacker after the sample has entered the resin

Related Thermo Scientific Products

89870	PepClean TM C-18 Spin Columns, 25/pkg.
23235	Micro BCA TM Protein Assay Kit, sufficient reagents to perform 480 standard tube assays
23225	BCA TM Protein Assay Kit, sufficient reagents to perform 500 standard tube assays
89879	Zeba Micro Spin Columns (Empty), 50 columns with screw caps

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.