INSTRUCTIONS



Mem-PERTM Plus Membrane Protein Extraction Kit

89842

Number

Description

89842

Mem-PER Plus Membrane Protein Extraction Kit, contains sufficient lysis and extraction reagents for approximately 50 mammalian cell pellet fractions containing 5×10^6 cells each or 25 tissue samples containing 20-40mg of tissue

Kit Contents:

Cell Wash Solution, 225mL; store at 4°C Solubilization Buffer, 25mL; store at 4°C Permeabilization Buffer, 50mL; store at -20°C

Storage: Upon receipt store Cell Wash Solution and Solubilization Buffer at 4°C. Aliquot and store Permeabilization Buffer at -20°C to avoid freeze/thaws. Product is shipped on ice.

Table of Contents

Introduction]
Important Product Information	1
Additional Materials Required	
Procedure for Membrane Protein Extraction from Different Sample Types	
Protocol 1: Adhesion Mammalian Cells	
Protocol 2: Suspension Mammalian Cells	2
Protocol 3: Soft Tissue	
Protocol 4: Hard Tissue	3
Froubleshooting	
Related Thermo Scientific Products	

Introduction

The Thermo ScientificTM Mem-PERTM Plus Membrane Protein Extraction Kit is used for the enrichment of integral membrane proteins and membrane-associated proteins from cultured mammalian cells or tissue. The kit uses a mild detergent-based, selective extraction protocol, which eliminates the hassle of phase separation based on hydrophobicity, allowing better reproducibility and higher throughput.

The cells are permeabilized with a mild detergent to allow the release of soluble cytosolic proteins. A second detergent then solubilizes membrane proteins. Extraction efficiencies will vary depending on the number of times the integral membrane protein(s) of interest spans the lipid bilayer. Membrane proteins with at least 1-2 transmembrane domains are typically extracted with an efficiency of up to 90%. Cross-contamination of cytosolic proteins into the membrane fraction is usually less than 10%.

Important Product Information

- For optimal results, include protease and phosphatase inhibitors (e.g., Product No. 78440) in the Permeabilization and Solubilization Buffers.
- To directly quantify proteins in the cytosolic and/or membrane fraction, use a protein assay such as the Thermo ScientificTM BCA Protein Assay (Product No. 23225), Micro BCATM Protein Assay (Product No. 23235) or PierceTM 660nm Protein Assay (Product No. 22660).



- Both fractions are compatible with SDS-PAGE. Use the Thermo ScientificTM PierceTM SDS-PAGE Sample Prep Kit (Product No. 89888) if the protein of interest is in low abundance and a large volume is required for adequate detection.
- If the detergent in the membrane fraction interferes with downstream applications, use Thermo ScientificTM PierceTM Detergent Removal Spin Columns (Product No. 87776-7) to decrease the amount of detergent in the sample. However, the sample:resin ratio will need to be empirically determined since removing too much detergent may result in protein precipitation.

Additional Material Required

- Protease and phosphatase inhibitors (e.g., Product No. 78440)
- Temperature-controlled mixer
- Cell scraper to remove adherent cells from plates
- 5mL microcentrifuge tubes
- For soft tissues, a 2mL Dounce Tissue Grinder (e.g., Kontes or Wheaton Tenbroeck) is required.
- For hard tissues, a hand-held homogenizer for 0.5-1.5mL samples (e.g., BrinkmannTM PolytronTM PT-1200CL Homogenizer)

Procedure for Membrane Protein Extraction from Different Sample Types

Choose the appropriate protocol for the sample type used.

Protocol 1: Adherent Mammalian Cells

- 1. Resuspend 5×10^6 cells in the growth media by scraping the cells off the surface of the plate with a cell scraper. Centrifuge harvested cell suspension at $300 \times g$ for 5 minutes.
- 2. Wash cell pellet with 3mL of Cell Wash Solution and centrifuge at $300 \times g$ for 5 minutes.
- 3. Carefully remove and discard the supernatant. Resuspend the cells in 1.5mL of Cell Wash Solution and transfer to a 2mL centrifuge tube. Centrifuge at $300 \times g$ for 5 minutes and discard supernatant.
- 4. Add 0.75mL of Permeabilization Buffer to the cell pellet. Vortex briefly to obtain a homogeneous cell suspension. Incubate 10 minutes at 4°C with constant mixing.
- 5. Centrifuge permeabilized cells for 15 minutes at $16,000 \times g$. Carefully remove the supernatant containing cytosolic proteins and transfer to a new tube.
- 6. Add 0.5mL of Solubilization Buffer to the pellet and resuspend by pipetting up and down. Incubate tubes at 4°C for 30 minutes with constant mixing.
- 7. Centrifuge tubes at $16,000 \times g$ for 15 minutes at 4°C. Transfer supernatant containing solubilized membrane and membrane-associated proteins to a new tube.
- 8. Proceed to downstream application. Immediately use cytosolic and membrane fractions stored on ice or store aliquots at -80°C for future use.

Protocol 2: Suspension Mammalian Cells

- 1. Harvest 5×10^6 cells by centrifugation at $300 \times g$ for 5 minutes. Wash harvested cell pellet with 3mL of Cell Wash Solution and centrifuge at $300 \times g$ for 5 minutes.
- 2. Carefully remove and discard the supernatant. Resuspend the cells in 1.5mL of Cell Wash Solution and transfer to a 2mL centrifuge tube. Centrifuge at $300 \times g$ for 5 minutes and discard supernatant.
- 3. Add 0.75mL of Permeabilization Buffer to the cell pellet. Vortex briefly to obtain a homogeneous cell suspension. Incubate 10 minutes at 4°C with constant mixing.
- 4. Centrifuge permeabilized cells for 15 minutes at 16,000 × g. Carefully remove the supernatant containing cytoslic proteins and transfer to a new tube.



- 5. Add 0.5mL of Solubilization Buffer to the pellet and resuspend by pipetting up and down. Incubate tubes at 4°C for 30 minutes with constant mixing.
- 6. Centrifuge tubes at $16,000 \times g$ for 15 minutes at 4°C. Transfer supernatant containing solubilized membrane and membrane-associated proteins to a new tube.
- 7. Proceed to downstream application. Immediately use cytosolic and membrane fractions stored on ice or store aliquots at -80°C for future use.

Protocol 3: Soft Tissue

- 1. Place 20-40mg of soft tissue in a 5mL microcentrifuge tube. Add 4mL of Cell Wash Solution to the tissue, vortex briefly and discard wash.
- 2. Transfer to a 2mL tissue grinder and cut the tissue into small pieces with a pair of scissors. Add 1mL of permeabilization buffer to the tissue and homogenize until an even suspension is obtained (6-10 strokes).
- 3. Add 1mL of Permeabilization Buffer and transfer homogenate to a new tube, incubating for 10 minutes at 4°C with constant mixing.
 - **Optional:** Remove non-homogenized tissue using the protocol provided with the Thermo ScientificTM PierceTM Tissue Strainers (Product No. 87791) before pelleting permeabilized cells.
- 4. Centrifuge at $16,000 \times g$ for 15 minutes at 4°C to pellet permeabilized cells. Carefully remove the supernatant containing cytosolic proteins and transfer to a new tube.
- 5. Resuspend the pellet in 1mL of Solubilization Buffer. Pipette up and down to obtain a homogeneous suspension. Incubate 30 minutes at 4°C with constant mixing.
- 6. Centrifuge tubes at $16,000 \times g$ for 15 minutes at 4°C. Transfer supernatant containing solubilized membrane and membrane-associated proteins to a new tube.
- 7. Proceed to downstream application. Immediately use cytosolic and membrane fractions stored on ice or store aliquots at -80°C for future use.

Protocol 4: Hard Tissue

- 1. Place 20-40mg of hard tissue in a 5mL microcentrifuge tube. Add 4mL of Cell Wash Solution to the tissue, vortex briefly and discard the wash.
- 2. Add 1mL of Permeabilization Buffer to the tissue and cut the tissue into small pieces with a pair of scissors.
- 3. Homogenize minced tissue with a hand-held Polytron Homogenizer. Use a low setting (to prevent foaming) until an even suspension is obtained.
- 4. Add an additional 1mL of Permeabilization Buffer to the homogenate and incubate for 10 minutes on ice with constant mixing.
 - **Optional:** Remove non-homogenized tissue using protocol provided with the Pierce Tissue Strainers before pelleting permeabilized cells.
- 5. Centrifuge at $16,000 \times g$ for 15 minutes at 4°C to pellet permeabilized cells. Carefully remove the supernatant containing cytosolic proteins and transfer to a new tube.
- 6. Resuspend the pellet in 1mL of Solubilization Buffer. Pipette up and down to obtain a homogeneous suspension. Incubate 30 minutes at 4°C with constant mixing.
- 7. Centrifuge tubes at $16,000 \times g$ for 15 minutes at 4°C. Transfer supernatant containing solubilized membrane and membrane-associated proteins to a new tube.
- 8. Proceed to downstream application. Immediately use cytosolic and membrane fractions stored on ice or store aliquots at -80°C for future use.



Troubleshooting

78840 87790

87791

Problem	Possible Cause	Solution
Contamination of membrane fraction with cytosolic proteins	Cells were not permeabilized	Increase incubation time or strokes of homogenization
	Permeabilization Buffer was improperly stored	Store Permeabilization Buffer at -20°C
	Incomplete removal of permeabilized cell supernatant	Ensure complete removal of cytosolic extract
Low membrane protein yield	Membranes were solubilized by Permeabilization Buffer	Decrease incubation time or strokes
	Incomplete membrane protein isolation	Increase incubation time
Low overall protein yield	Not enough cells/tissue	Increase cell number or amount of starting tissue (mg)

Related Thermo Scientific Products

87785	Halt ^{1M} Protease Inhibitor Cocktail, EDTA-Free (100X), 1mL
87786	Halt Protease Inhibitor Cocktail (100X), 1mL
78441	Halt Protease and Phosphatase Inhibitor Cocktail, EDTA-Free (100X), 1mL
78440	Halt Protease and Phosphatase Inhibitor Cocktail (100X), 1mL
78420	Halt Phosphatase Inhibitor Cocktail (100X), 1mL
88665-9	Pierce Protease Inhibitor Tablets
23225	BCA Protein Assay Kit
89888	Pierce SDS-PAGE Sample Prep Kit
87777	Pierce Detergent Removal Spin Columns, 0.5mL

Subcellular Protein Fractionation Kit for Cultured Cells

Subcellular Protein Fractionation Kit for Tissue

Pierce Tissue Strainers, 250µm, 2.5mL, 50 units

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2013 Thermo Fisher Scientific Inc. All rights reserved. Polytron is a trademark of Kinematica AG Company. Unless otherwise indicated, all other trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.