INSTRUCTIONS

EZ-Link Sulfo-NHS-LC-LC-Biotin



Pub. No. MAN0016132 Pub. Part No. 2160236.5 Rev. A.0

<u>21338</u>

Number


21338

Description

EZ-Link Sulfo-NHS-LC-LC-Biotin, 50mg

Molecular Weight: 669.75

Spacer Arm: 30.5Å



Storage: Upon receipt store desiccated at -20°C. Product is shipped at ambient temperature.

Table of Contents

Introduction	1
Important Product Information	2
Additional Materials Required	
Procedure for Biotinylating Proteins in Solution	
Procedure for Biotinylating Cell Surface Proteins	
Related Thermo Scientific Products	4
Cited References	4
Product References	

Introduction

The Thermo ScientificTM EZ-LinkTM Sulfo-NHS-LC-LC-Biotin (sulfosuccinimidyl-6-[biotinamido]-6-hexanamido hexanoate) enables simple and efficient biotin labeling of antibodies, proteins and any other primary amine-containing macromolecules. Specific labeling of cell surface proteins is another common application for these uniquely water-soluble and membrane impermeable reagents.

Biotin is a small, naturally occurring vitamin that binds with high affinity to avidin and streptavidin proteins. Because it is so small (244Da), biotin can be conjugated to many proteins without altering their biological activities. Labeled proteins may be purified from unlabeled proteins using immobilized streptavidin and avidin affinity gels (see Related Thermo Scientific Products), and they may be detected easily in ELISA, dotblot or Western blot applications using streptavidin or avidin-conjugated probes.

N-Hydroxysuccinimide (NHS) esters of biotin are the most popular type of biotinylation reagent. NHS-activated biotins react efficiently with primary amino groups (-NH₂) in pH 7-9 buffers to form stable amide bonds. Proteins, including antibodies, generally have several primary amines in the side chain of lysine (K) residues and the N-terminus of each polypeptide that are available as targets for labeling with NHS-activated biotin reagents. Several different NHS esters of biotin are available, with varying properties and spacer arm lengths. Sulfo-NHS ester reagents are water soluble, enabling reactions to be performed in the absence of organic solvents such as DMSO or DMF.

Cell surface biotinylation is an important tool for studying the expression and regulation of receptors and transporters, differentiation of plasma membrane proteins from those localized to organelle membranes, and distribution of membrane proteins in polarized epithelial cells. The specificity of Sulfo-NHS-LC-LC-Biotin reagents for cell surface labeling has been demonstrated in these applications.^{1,2} Because these molecules dissolve readily in polar solutions and are charged by the sodium sulfonate group on the succinimidyl ring, they cannot permeate the cell membrane. As long as the cell remains intact, only primary amines exposed on the surface will be biotinylated with these Sulfo-NHS-LC-LC-Biotin reagents.



Important Product Information

- Sulfo-NHS-LC-LC-Biotin is moisture-sensitive. Store the vial of biotin reagent at -20°C with desiccant. To avoid moisture condensation onto the product, equilibrate vial to room temperature before opening.
- As directed in the procedure, dissolve the biotin reagent immediately before use. The NHS-ester moiety readily hydrolyzes and becomes non-reactive; therefore, weigh and dissolve only a small amount of the reagent at a time, and do not prepared stock solutions for storage. Discard any unused reconstituted reagent.
- Avoid buffers containing primary amines (e.g., Tris or glycine) as these will compete with the reaction (see Figure 2). If necessary, dialyze or otherwise desalt to exchange the protein sample into an amine-free buffer such as phosphate buffered saline (PBS; see Related Thermo Scientific Products).
- When biotinylating proteins in solution, excess non-reacted biotin and reaction byproducts are easily removed by size exclusion using either desalting columns or dialysis (See Additional Information and Related Thermo Scientific Products). A 10mL desalting column is best suited for processing biotinylation reactions involving 1-10mg of protein in approximately 0.5-2mL. For smaller amounts of protein and/or smaller reaction volumes, both the biotinylation reaction and subsequent buffer exchange may be performed in a single Thermo ScientificTM Slide-A-LyzerTM MINI Dialysis Unit. For larger reaction volumes than can be processed with a desalting column, either split the sample between two columns or use an appropriate Slide-A-Lyzer Dialysis Cassette for buffer exchange steps.

Additional Materials Required

- Phosphate-buffered Saline (PBS) or other amine-free buffer having pH 7-8 for use as reaction buffer (see Important Product Information and Related Thermo Scientific Products)
- Desalting columns or dialysis units for buffer exchange (see Important Product Information and Related Thermo Scientific Products)

Procedure for Biotinylating Proteins in Solution

A. Calculations


The extent of biotin labeling depends on the distribution of amino groups on the protein, and the concentration and amount of protein to be labeled. By using the appropriate molar ratio of biotin to protein, the extent of labeling can be controlled. Compared to reactions involving concentrated protein solutions, labeling reactions with dilute solutions require a greater fold molar excess of biotin reagent to achieve the same incorporation level. Generally, use ≥ 12 -fold molar excess of biotin for proteins at 2-10mg/mL or ≥ 20 -fold molar excess of biotin for proteins at $\leq 2mg/mL$. Adjust the molar ratio of biotin reagent to protein to obtain the level of incorporation desired.

1. Calculate millimoles of biotin reagent to add to the reaction for a 20-fold molar excess:

mL protein
$$\times \frac{\text{mg protein}}{\text{mL protein}} \times \frac{\text{mmol protein}}{\text{mg protein}} \times \frac{20 \text{ mmol Biotin}}{\text{mmol protein}} = \text{mmol Biotin}$$

- 20 = Molar fold excess of biotin for 2 mg/ml protein sample
- 2. Calculate microliters of 10mM biotin reagent solution (prepared in Step B.3) to add to the reaction:

mmol Biotin
$$\times \frac{1,000,000 \,\mu\text{L}}{\text{L}} \times \frac{\text{L}}{10 \,\text{mmol}} = \mu\text{L}$$
 Biotin





B. Biotin Labeling Reaction

- 1. Remove the vial of biotin reagent from freezer and equilibrate it to room temperature before opening in step 3.
- 2. Dissolve 1-10mg protein in 0.5-2.0mL PBS according to the calculation made in section A.

Note: Protein that is already dissolved in amine-free buffer at pH 7.2-8.0 may be used without buffer exchange or dilution with PBS. Proteins in Tris or other amine-containing buffers must be exchanged into a suitable buffer.

- 3. Immediately before use, prepare a 10mM solution of the biotin reagent by dissolving 2.0mg reagent in 300µL of ultrapure water.
- 4. Add the appropriate volume (see Calculations in section A) of 10mM biotin reagent solution to the protein solution.
- 5. Incubate reaction on ice for two hours or at room temperature for 30 minutes.

Note: Other than the possibility of ordinary protein degradation or microbial growth, there is no harm in reacting longer than the specified time.

6. Protein labeling is complete at this point, and although excess non-reacted and hydrolyzed biotin reagent remains in the solution, it is often possible to perform preliminary tests of the labeled protein by ELISA or Western blot. Once proper function and labeling of the protein has been confirmed, the labeled protein may be purified for optimal performance and stability using desalting or dialysis. If the Thermo Scientific[™] Pierce[™] Biotin Quantitation Kit (HABA assay; see Related Thermo Scientific Products) will be performed to determine the level of biotin incorporation, the protein first must be desalted or dialyzed to remove non-reacted biotin.

Procedure for Biotinylating Cell Surface Proteins

Many variations of this procedure exist in the literature (see Product References). Labeling may be performed on cells in suspension or on adherent cells in culture plates. In the latter situation, diffusion of the Sulfo-NHS-LC-LC-Biotin Reagent to all surfaces of the cells will be limited, and labeling will occur predominately on the exposed surface. Culture media must be washed from the cells; otherwise, amine-containing components will compete and quench the reaction to cell surface proteins. Using a more concentrated cell suspension is most effective because less biotin reagent is required in the reaction. Generally, a final concentration of 2-5mM Biotin Reagent is effective. NHS reactions occur more rapidly at higher pH; therefore, pH 8.0 is used in the following example so that labeling can be completed as quickly as possible.

- 1. Wash cells three times with ice-cold PBS (pH 8.0) to remove amine-containing culture media and proteins from the cells.
- 2. Suspend cells at a concentration of $\sim 25 \times 10^6$ cells/mL in PBS (pH 8.0).
- Add 1.0mg of Sulfo-NHS-LC-LC-Biotin reagent per mL of cell suspension (results in ~2mM biotin reagent). Alternatively, add 200µL of the 10mM biotin reagent solution (see step B.3 on previous page) per mL of cell suspension.
- 4. Incubate reaction mixture at room temperature for 30 minutes.

Note: Performing this incubation at 4°C may reduce active internalization of the biotin reagent.

- 5. Wash cells three times with PBS plus 100mM glycine to quench and remove excess biotin reagent and byproducts.
- 6. Lyse and/or analyze biotin-labeled cells as required for the research method.

Additional Information

Please visit the website for additional information including the following:

- Tech Tip #14: Perform labeling and other reactions in Slide-A-Lyzer Dialysis Cassettes
- Tech Tip #43: Protein stability and storage
- HABA Calculator for computing the results associated with the HABA assay measurements



Related Thermo Scientific Products

	28372	BupH TM Phosphate Buffered Saline Packs, 40 pack
	69576	Slide-A-Lyzer MINI Dialysis Unit Kit
	66382	Slide-A-Lyzer Dialysis Cassette Kits
	89890	Zeba TM Spin Desalting Columns, 7K MWCO, 2mL, 25 columns
	89892	Zeba Spin Desalting Columns, 7K MWCO, 5mL, 25 columns
	89894	Zeba Spin Desalting Columns, 7K MWCO, 10mL, 25 columns
	28005	Pierce Biotin Quantitation Kit
	20347	Streptavidin Agarose Resin, 2mL
	20228	Pierce™ Monomeric Avidin Agarose, 5mL
	21126	Streptavidin, Horseradish Peroxidase Conjugated, 1mg
	21445	EZ-Link Sulfo-NHS-SS-Biotinylation Kit
	21945	EZ-Link Micro Sulfo-NHS-SS-Biotinylation Kit
	21925	EZ-Link Micro Sulfo-NHS-Biotinylation Kit
	21935	EZ-Link Micro Sulfo-NHS-LC-Biotinylation Kit
	21326	EZ-Link Sulfo-NHS-Biotin, No-Weigh™ Format, 8 × 1mg microtubes
	21328	EZ-Link Sulfo-NHS-SS-Biotin, No-Weigh Format, 8 × 1mg microtubes
	21329	EZ-Link NHS-PEO₄-Biotin, No-Weigh Format, 8 × 2mg microtubes
	21327	EZ-Link Sulfo-NHS-LC-Biotin, No-Weigh Format, 8 × 1mg microtubes

Cited References

1. Daniels, G.M. and Amara, S.G. (1998). Selective labeling of neurotransmitter transporters at the cell surface. *Methods Enzymol* 296:307-18.

2. Huh, K-H. and Wenthold, R.J. (1999). Turnover analysis of glutamate receptors identifies a rapidly degraded pool of the *N*-methyl-D-aspartate receptor subunit, NR1, in cultured cerebellar granule cells. *J Biol Chem* **274:**151-7.

Product References

- Ali, M.K. and Bergson, C. (2003). Elevated intracellular calcium triggers recruitment of the receptor cross-talk accessory protein calcyon to the plasma membrane. J Biol Chem 278(51):51654-63.
- Borroto, A., *et al.* (2003). Impaired trafficking and activation of tumor necrosis factor-α-converting enzyme in cell mutants defective in protein ectodomain shedding. *J Biol Chem* 278(28):25933-9.
- Chyung, J.H. and Selkoe, D.J. (2003). Inhibition of receptor-mediated endocytosis demonstrates generation of amyloid β-protein at the cell surface. *J Biol Chem* 278(51): 51035-43.
- Frickel, E-M., et al. (2002). TROSY-NMR reveals interaction between Erp57 and the tip of the calreticulin P-domain. PNAS 99(4):1954-9.
- Gimferrer, I., et al. (2003). The accessory molecules CD5 and CD6 associate on the membrane of lymphoid T cells. J Biol Chem 278(10):8564-71.
- Li, H. and Pajor, A.M. (2003). Mutagenesis of the N-glycosylation site of hNaSi-1 reduces transport activity. Am J Physiol Cell Physiol 285:C1188-96.
- Lukashevich, I.S., et al. (2003). Arenavirus-mediated liver pathology: acute lymphocytic choriomeningitis virus infection of rhesus macaques is characterized by high-level interleukin-6 expression and hepatocyte proliferation. J Virol **77**(3):1727-37.
- Ohnishi, T., *et al.* (2003). MD-2 is necessary for the toll-like receptor 4 protein to undergo glycosylation essential for its translocation to the cell surface. *Clin Diagn Lab Immunol* **10(3)**:405-410.



This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2016 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.