

# SuperScript<sup>®</sup> Indirect cDNA Labeling System

# For generating fluorescently labeled cDNA to use in microarray screening

Catalog numbers L1014-01, L1014-02, and L1014-03

Revision date 9 January 2012 Publication Part number 25-0655

MAN0000381

#### **Table of Contents**

Accessory ProductsviOverview.1Detailed Methods.3Isolating RNA.3First-Strand cDNA Synthesis.5Purifying First-Strand cDNA.7Labeling with Fluorescent Dye.9Purifying Labeled cDNA.10Appendix.12Assessing Labeling Efficiency13Technical Support.15Purchaser Notification.16References.17	]	Kit Contents and Storage	iv
Overview.1Detailed Methods.3Isolating RNA.3First-Strand cDNA Synthesis.5Purifying First-Strand cDNA.7Labeling with Fluorescent Dye.9Purifying Labeled cDNA.10Appendix.12Assessing Labeling Efficiency.12Troubleshooting.13Technical Support.15Purchaser Notification.16References.17		Accessory Products	vi
Detailed Methods       .3         Isolating RNA       .3         First-Strand cDNA Synthesis       .5         Purifying First-Strand cDNA       .7         Labeling with Fluorescent Dye       .9         Purifying Labeled cDNA       10         Appendix       .12         Troubleshooting       .13         Technical Support       .15         Purchaser Notification       .16         References       .17	(	Overview	1
Isolating RNA.3First-Strand cDNA Synthesis.5Purifying First-Strand cDNA.7Labeling with Fluorescent Dye.9Purifying Labeled cDNA10Appendix12Assessing Labeling Efficiency12Troubleshooting.13Technical Support.15Purchaser Notification.16References.17	Detailed	Methods	3
First-Strand cDNA Synthesis       .5         Purifying First-Strand cDNA       .7         Labeling with Fluorescent Dye       .9         Purifying Labeled cDNA       10         Appendix       12         Assessing Labeling Efficiency       12         Troubleshooting       13         Technical Support       15         Purchaser Notification       16         References       17	]	Isolating RNA	3
Purifying First-Strand cDNA       .7         Labeling with Fluorescent Dye       .9         Purifying Labeled cDNA       10         Appendix       12         Assessing Labeling Efficiency       12         Troubleshooting       13         Technical Support       15         Purchaser Notification       16         References       17	]	First-Strand cDNA Synthesis	5
Labeling with Fluorescent Dye	]	Purifying First-Strand cDNA	7
Purifying Labeled cDNA       10         Appendix       12         Assessing Labeling Efficiency       12         Troubleshooting       13         Technical Support       15         Purchaser Notification       16         References       17	]	Labeling with Fluorescent Dye	9
Appendix12Assessing Labeling Efficiency12Troubleshooting13Technical Support15Purchaser Notification16References17	]	Purifying Labeled cDNA	10
Assessing Labeling Efficiency	Appendi	х	12
Troubleshooting       13         Technical Support       15         Purchaser Notification       16         References       17		Assessing Labeling Efficiency	12
Technical Support	,	Troubleshooting	13
Purchaser Notification	,	Technical Support	15
References	1	Purchaser Notification	16
	]	References	17

#### Kit Contents and Storage

Kit Sizes	The SuperScript <sup>®</sup> and a Purification the labeling comp	The SuperScript <sup>®</sup> Indirect cDNA Labeling System is supplied with either a Core Module and a Purification Module, or a Core Module only. Note that the Core Module contains the labeling components.			
	<u>Cat no.</u>	Number of Labeling Reactions	Modules		
	L1014-01	10	Core and Purification		
	L1014-02	30	Core and Purification		
	L1014-03	30	Core only		
Shipping and Storage	The Core Module temperature. Upo the components o	The Core Module is shipped on dry ice and the Purification Module is shipped at room temperature. Upon receipt, store the components of the Core Module at $-20^{\circ}$ C and store the components of the Purification Module at room temperature.			

**Core Module** The components of the Core Module should be stored at  $-20^{\circ}$ C.

Item	Components/Concentration	Kit Size		
		10 Rxns	30 Rxns	
SuperScript <sup>®</sup> III Reverse Transcriptase	400 U/µl in: 20 mM Tris-HCl (pH 7.5) 100 mM NaCl 0.1 mM EDTA 1 mM DTT 0.01% (v/v) NP-40 50% (v/v) glycerol	20 µl	60 µl	
5X First-Strand Buffer	250 mM Tris-HCl (pH 8.3, room temp) 375 mM KCl 15 mM MgCl2	1020 µl	1020 µl	
Dithiothreitol (DTT)	0.1 M DTT in water	250 µl	250 µl	
dNTP Mix	dATP, dGTP, dCTP, dTTP, one aminoallyl-modified nucleotide, and one aminohexyl-modified nucleotide at optimal concentrations in DEPC-treated water	20 µl	60 µl	
2X Coupling Buffer	—	100 µl	300 µl	
Anchored Oligo(dT)20 primer	$2.5 \mu g/\mu l$ in DEPC-treated water	20 µl	60 µl	
Random hexamer primers	$0.5 \mu g/\mu l$ in DEPC-treated water	10 µl	30 µl	
Glycogen	20 mg/ml	20 µl	60 µl	
DMSO	—	200 µl	750 µl	
RNaseOUT™	40 U/µl	10 µl	30 µl	
DEPC-treated Water	—	2 ml	6 ml	
3 M Sodium Acetate	pH 5.2	1 ml	2 ml	
Control HeLa RNA	1 μg/μl	20 µl	20 µl	

#### Kit Contents and Storage, continued

#### cDNA Labeling Purification Module

The components of the Purification Module should be stored at room temperature. This module is included with Catalog Numbers L1014-01 and L1014-02.

		Kit	Size		
	Item	10 Rxns	30 Rxns		
	Loading Buffer (you must add 100% isopropanol to create the final buffer; see below)	4.3 ml	13 ml		
	Wash Buffer (you must add 100% ethanol to create the final buffer; see below)	8.33 ml	25 ml		
	S.N.A.P. <sup>™</sup> Columns	22 columns	62 columns		
	Clear Collection Tubes	22 tubes	62 tubes		
	Amber Collection Tubes	11 tubes	31 tubes		
Preparing Loading Buffer with Isopropanol	The Loading Buffer supplied in each Purification Module must be mixed with 100% isopropanol prior to use. The Loading Buffer plus isopropanol is stable for at least six months at room temperature.				
	Add the amount of isopropanol indicated below directly to each bottle of Loading Buffer. Be sure to mark the appropriate checkbox on the bottle to indicate that you have added the isopropanol.				
	10-RxrLoading Buffer4.3 m100% Isopropanol10.0 mTotal Volume14.3 m	Kit         30-Rx           nl (entire bottle)         13           nl         _30           nl         _43	<u>n Kit</u> ml (entire bottle) ml		
Preparing Wash Buffer with Ethanol	The Wash Buffer supplied in each Purification Mc prior to use. The Wash Buffer plus ethanol is stabl temperature.	dule must be mixed e for at least six mon	with 100% ethanol ths at room		
	Add the amount of ethanol indicated below directly to each bottle of Wash Buffer. Be sure to mark the appropriate checkbox on the bottle to indicate that you have added the ethanol.				
	10-Rxn           Wash Buffer         8.33 m           100% Ethanol         25.00 m           Total Volume         33.33 m	Kit         30-Rx           nl (entire bottle)         25           nl         75           nl         100	<u>n Kit</u> ml (entire bottle) <u>ml</u> ml		
Product Use	For research use only. Not intended for any huma uses.	ın or animal diagnos	tic or therapeutic		

#### **Accessory Products**

#### Additional Products

Many of the reagents in the SuperScript<sup>®</sup> Indirect cDNA Labeling System, as well as additional reagents that may be used with this system, are available separately from Invitrogen. Ordering information is provided below.

Product	Quantity	Catalog no.
Alexa Fluor <sup>®</sup> 555 and Alexa Fluor <sup>®</sup> 647 Reactive Dye Decapacks	2 × 10 vials	A-32755
Alexa Fluor <sup>®</sup> 555 Reactive Dye Decapack	10 vials	A-32756
Alexa Fluor <sup>®</sup> 647 Reactive Dye Decapack	10 vials	A-32757
RNase Away <sup>™</sup> Reagent	250 ml	10328-011
Yeast tRNA	25 mg 50 mg	15401-011 15401-029
PureLink <sup>™</sup> Micro-to-Midi Total RNA Purification System	50 reactions	12183-018
PureLink <sup>™</sup> 96 RNA Purification System	384 reactions	12173-011
TRIzol® Reagent	100 ml	15596-026
	200 ml	15596-018
FastTrack <sup>®</sup> 2.0 mRNA Isolation Kit	6 reactions 18 reactions	K1593-02 K1593-03
FastTrack <sup>®</sup> MAG Micro mRNA Isolation Kit	12 reactions	K1580-01
FastTrack® MAG Maxi mRNA Isolation Kit	6 reactions	K1580-02
RNaseOUT <sup>™</sup> Recombinant Ribonuclease Inhibitor	5000 units	10777-019
Human Cot-1 DNA®	500 µg	15279-011
Mouse Cot-1 DNA®	500 µg	18440-016
Random primers	9 A <sub>260</sub> units	48190-011
UltraPure <sup>™</sup> DEPC-treated water	4 × 1.25 ml	10813-012
UltraPure <sup>™</sup> 10% SDS solution	$4 \times 100 \text{ ml}$	15553-027
UltraPure <sup>™</sup> 20X SSC	1 L	15557-044
UltraPure™ 20x SSPE	1 L	15591-043

Introduction	The SuperScript <sup>®</sup> Indirect cDNA Labeling System is a highly efficient system for generating fluorescently labeled cDNA for use on microarrays in gene expression studies (De Risi <i>et al.</i> , 1996; Eisen & Brown, 1999). It uses an aminoallyl-modified nucleotide and an aminohexyl-modified nucleotide together with other dNTPs in a cDNA synthesis reaction with SuperScript <sup>®</sup> III Reverse Transcriptase. After a purification step to remove unincorporated nucleotides, the amino-modified cDNA is coupled with a monoreactive, N-hydroxysuccinimide (NHS)-ester fluorescent dye. A final purification step removes any unreacted dye, and the fluorescently labeled cDNA is ready for hybridization to microarrays. This system uses 5–20 µg of total RNA or 0.4–2 µg of mRNA as starting material, and is compatible with Alexa Fluor <sup>®</sup> 555 and Alexa Fluor <sup>®</sup> 647 fluorescent dyes from Invitrogen, Cy3 <sup>™</sup> and Cy5 <sup>™</sup> dyes from Amersham Biosciences, or other monoreactive NHS-ester dyes from a variety of manufacturers.				
Advantages of the System	<ul> <li>SuperScript<sup>®</sup> III Reverse Transcriptase in the first-strand synthesis reaction ensures high specificity and high yields of cDNA, as well as more full-length cDNA</li> <li>Use of two amino-modified nucleotides in the cDNA synthesis reaction results in a greater incorporation of fluorescent dye and higher signal intensity with small amounts of starting material</li> </ul>				
	• Unbiased incorporation of amino-modified dNTPs and the high efficiency of the coupling reaction result in an even distribution of fluorescent signal and high overall levels of fluorescence, increasing the sensitivity and reproducibility of array hybridizations				
	• System includes all major reagents and materials for preparing fluorescently labeled cDNA, except fluorescent dye				
Experimental Outline	The flow chart below outlines the experimental steps of the system:				
	Perform first-strand cDNA synthesis using SuperScript <sup>™</sup> III RT and amino-modified dNTPs				
	L1014-01 or L1014-02) or your method of choice (catalog no. L1014-03)				
	Perform fluorescent dye coupling reaction				
	or L1014-02) or your method of choice (catalog no. L1014-03)				
	Proceed to hybridization				

## Overview, continued

Advantages of SuperScript <sup>®</sup> III Reverse Transcriptase	SuperScript <sup>®</sup> III Reverse Transcriptase is an engineered version of M-MLV RT with reduced RNase H activity and increased thermal stability. The enzyme can be used to synthesize first-strand cDNA from total RNA or mRNA at temperatures up to 55°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases.		
	The SuperScript <sup>®</sup> III RT in this kit is provided at an optimal concentration and used at an optimal temperature for incorporating amino-modified nucleotides in first-strand cDNA synthesis.		
Anchored Oligo(dT) <sub>20</sub>	Anchored $oligo(dT)_{20}$ primer is a mixture of 12 primers, each consisting of a string of 20 deoxythymidylic acid (dT) residues followed by two additional nucleotides represented by VN, where:		
	• V is dA, dC, or dG		
	• N is dA, dC, dG or dT		
	The VN "anchor" allows the primer to anneal only at the 5' end of the poly(A) tail of mRNA, providing more efficient cDNA synthesis for labeling applications.		
Dye Compatibility	This kit has been developed using Alexa Fluor® fluorescent dyes from Invitrogen and CyDye™ fluorescent dyes from Amersham Biosciences. See page 9 for more information.		
	Other monofunctional, N-hydroxysuccinimide (NHS)-reactive fluorescent dyes are compatible with this system.		
Materials Supplied by the User	In addition to the kit components, you should have the following items on hand before using the SuperScript® Indirect cDNA Labeling System.		
	Monofunctional, NHS-reactive fluorescent dye		
	Vortex mixer		
	Microcentrifuge		
	Aerosol resistant pipette tips		
	Water baths or incubator		
	• 1 N NaOH		
	• 1 N HCl		
	Sterile microcentrifuge tubes		
	• 100% Isopropanol		
	• 100% Ethanol		
	• 75% Ethanol		
Control Reaction	We recommend performing the labeling procedure using the Control HeLa RNA included in the system to determine the efficiency of the labeling reaction. The section on <b>First-Strand cDNA Synthesis</b> (page 5) describes how to set up the control reaction and page 12 has equations for calculating the efficiency of the labeling procedure.		

#### **Detailed Methods**

## **Isolating RNA**

Introduction	High-quality, intact RNA is essential for full-length, high-quality cDNA synthesis. In this step, you isolate total RNA or mRNA using a method of choice.			
Important	The quality of the RNA is critical for successful labeling and hybridization. The presence of contaminants in the RNA may significantly increase background fluorescence in your microarrays. Carefully follow the recommendations below to prevent RNase contamination.			
General Handling of RNA	<ul> <li>When working with RNA:</li> <li>Use disposable, individually wrapped, sterile plasticware.</li> <li>Use aerosol resistant pipette tips for all procedures.</li> <li>Use only sterile, new pipette tips and microcentrifuge tubes.</li> <li>Wear latex gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin.</li> <li>Use proper microbiological aseptic technique when working with RNA.</li> <li>Dedicate a separate set of pipettes, buffers, and enzymes for RNA work.</li> <li>Microcentrifuge tubes can be taken from an unopened box, autoclaved, and used for all RNA work. RNase-free microcentrifuge tubes are available from several suppliers. If it is necessary to decontaminate untreated tubes, soak the tubes overnight in a 0.01% (v/v) aqueous solution of diethylpyrocarbonate (DEPC-treated), rinse the tubes with sterile distilled water, and autoclave the tubes.</li> <li>You can use RNase AWAY<sup>™</sup> Reagent, a non-toxic solution available from Invitrogen (see page vi), to remove RNase contamination from surfaces. For further information on controlling RNase contamination, see Ausubel et al. and Sambrook, et al (Ausubel et al.</li> </ul>			
Isolating RNA	1994; Sambrook <i>et al.</i> , 1989). This system is designed for use with 5–20 µg total RNA or 0.4–2 µg of mRNA. To isolate total RNA, we recommend the PureLink <sup>™</sup> Micro-to-Midi Total RNA Purification System, TRIzol® Reagent, or (for high-throughput applications) the PureLink <sup>™</sup> 96 RNA Purification System. To isolate mRNA, we recommend the FastTrack <sup>®</sup> 2.0 mRNA Isolation Kits or the FastTrack <sup>®</sup> MAG mRNA Isolation Kits. Ordering information is provided on page vi. After you have isolated the RNA, check the quality of your RNA preparation as described on the following page.			

## Isolating RNA, continued

Checking the RNA Quality	To check RNA quality, analyze 500 ng of RNA by agarose/ethidium bromide gel electrophoresis. You can use a regular 1% agarose gel or a denaturing agarose gel (Ausubel <i>et al.</i> , 1994). For total human RNA using a regular agarose gel, mRNA will appear as a smear from 0.5 to 9 kb, and 28S and 18S rRNA will appear as bands at 4.5 kb and 1.9 kb, respectively. The 28S band should be twice the intensity of the 18S band. If you are using a denaturing gel, the rRNA bands should be very clear and sharp.		
	If you do not load enough RNA, the 28S band may appear to be diffuse. A smear of RNA or a lower intensity 28S band with an accumulation of low molecular weight RNA on the gel are indications that the RNA may be degraded, which will decrease the labeling efficiency. If you do not detect any RNA, you will need to repeat RNA isolation. Refer to the <b>Troubleshooting</b> section on page 13.		
Storing RNA	After preparing the RNA, we recommend that you proceed directly to <b>First-Strand cDNA Synthesis</b> on page 5. Otherwise, store the RNA at -80°C.		

# First-Strand cDNA Synthesis

Introduction	After you have isolated RNA and checked the quality of your RNA preparation, you are ready to synthesize cDNA.		
Before Starting	<ul> <li>The following materials are supplied by the user:</li> <li>5-20 μg total RNA or 0.4-2 μg mRNA</li> <li>1 N NaOH</li> <li>1 N HCl</li> <li>Water baths, heating block, or incubator set at 46°C and 70°C</li> <li>Ice</li> <li>0.5-ml or 1.5-ml RNase-free microcentrifuge tubes</li> <li>The following materials are supplied in the kit:</li> <li>Anchored Oligo(dT)<sub>20</sub> primer</li> <li>Random hexamers (for mRNA starting material only)</li> <li>dNTP mix, including amino-modified nucleotides</li> <li>5X First-Strand buffer</li> <li>0.1 M DTT</li> <li>RNaseOUT<sup>™</sup></li> <li>SuperScript<sup>®</sup> III RT</li> <li>DEPC-treated water</li> <li>Control HeLa RNA, optional</li> <li>3 M Sodium Acetate, pH 5.2</li> </ul>		
Control HeLa RNA	Control HeLa RNA is included in the kit to help you determine the efficiency of the labeling procedure. We strongly recommend that you perform the control reaction if you are a first-time user of the SuperScript <sup>®</sup> Indirect cDNA Labeling System. Instructions are provided on the next page to set up separate cDNA synthesis reactions for your sample and the Control HeLa RNA. Equations for calculating the amount of coupled dye in the control reaction are provided on page 12.		
Note	RNaseOUT <sup>™</sup> Recombinant RNase Inhibitor has been included in the system to safeguard against degradation of target RNA due to ribonuclease contamination of the RNA preparation.		

# First-Strand cDNA Synthesis, continued

First-Strand cDNA Synthesis	The mR	following procedure is designed to convert NA into first-strand cDNA.	5–20 µg of tot	tal RNA or 0.4–2 μg of
Reaction	<b>Not</b> 10 µ	<b>Note:</b> If you are setting up a control reaction (recommended for first-time users), use 10 $\mu$ l of the Control HeLa RNA (1 $\mu$ g/ $\mu$ l) supplied in the kit.		
	1.	Mix and briefly centrifuge each component	before use.	
	2.	Prepare reaction(s) as follows, using 0.5-ml	or 1.5-ml RNa	ase-free tubes:
		Component 5–20 μg total RNA or 0.4–2 μg mRNA Control HeLa RNA (1 μg/μl) Anchored Oligo(dT) <sub>20</sub> Primer (2.5 μg/μl) Random hexamers (only if using mRNA)	<u>Sample</u> X μl 2 μl 1 μl *	<u>Control</u> — 10 μl 2 μl
		DEPC-treated water	to 18 µl	to 18 µl
		*For mRNA, use <b>both</b> anchored $oligo(dT)_{20}$ use <b>only</b> 2 µl of anchored $oligo(dT)_{20}$ .	and random ł	hexamers. For total RNA,
	3.	Incubate tubes at 70°C for 5 minutes, and the	en place on ic	ce for at least 1 minute.
	4.	Add the following to each tube (sample and control) on ice:		
		Component 5X First-Strand buffer 0.1 M DTT dNTP mix (including amino-modified nucleotic RNaseOUT <sup>™</sup> (40 U/μl) <u>SuperScript<sup>®</sup> III RT (400 U/μl)</u> Total Volume	des)	<u>Volume</u> 6 μl 1.5 μl 1.5 μl 1 μl <u>2 μl</u> 30 μl
	5.	Mix gently and collect the contents of each tube at 46°C for 2–3 hours. <b>Note:</b> A 3-hour i cDNA yield than a 2-hour incubation.	tube by brief c incubation res	centrifugation. Incubate sults in 20–30% higher
	Afte	er incubation, proceed directly to <b>Alkaline F</b>	Iydrolysis and	d Neutralization, below.
Alkaline Hydrolysis and	Aft deg	er cDNA synthesis, above, immediately perf rade the original RNA:	orm the follow	wing hydrolysis reaction to
Neutralization	1.	Add 15 µl of 1 N NaOH to each reaction tul	be from Step 5	5, above. Mix thoroughly.
	2.	Incubate tube at $70^{\circ}$ C for 10 minutes.		
	3.	Add 15 µl of 1 N HCl to neutralize the pH a	and mix gently	у.
	4.	Add 20 µl 3 M Sodium Acetate, pH 5.2, and	l mix gently.	

# Purifying First-Strand cDNA

Introduction	After you have generated cDNA with amino-modified nucleotides, you need to purify the cDNA to remove unincorporated dNTPs and hydrolyzed RNA.			
	Catalog Numbers L1014-01 and L1014-02 include a Purification Module developed for use with the system. Use the <b>S.N.A.P.™ Column Purification</b> procedure on the next page to purify your amino-modified cDNA using this Purification Module.			
	Catalog Number L1014-03 does not included a Purification Module. Use your preferred method of purification instead of the S.N.A.P.™ Column Purification procedure, and proceed to the Ethanol Precipitation procedure on the next page .			
Important	You must perform the <b>Ethanol Precipitation</b> step on page 8 even if you are using your own purification procedure.			
Before Starting	The following materials are supplied by the user:			
_	Microcentrifuge			
	• 1.5-ml microcentrifuge tube			
	• 100% Ethanol			
	• 75% Ethanol			
	The following materials are supplied in the Core Module:			
	• 2X Coupling Buffer			
	• 3 M Sodium Acetate, pH 5.2			
	• Glycogen (20 mg/ml)			
	The following materials are supplied in the Purification Module (Catalog Numbers L1014-01 and L1014-02):			
	• DEPC-treated water			
	• S.N.A.P. <sup>™</sup> column(s) and clear collection tube(s)			
	• Loading Buffer <b>plus isopropanol</b> (see page v for preparation)			
	• Wash Buffer <b>plus ethanol</b> (see page v for preparation)			
Important	The pellet should be completely dry at the end of the purification procedure to ensure complete removal of the ethanol. The presence of ethanol can inhibit the labeling reaction.			

## Purifying First-Strand cDNA, continued

S.N.A.P. <sup>™</sup> Column Purification	Use Mo	e the following procedure to purify the cDNA using the components of the Purification odule (Cat. nos. L1014-01 and L1014-02).
	If y to <b>l</b>	rou are using Cat. No. L1014-03, purify using your method of choice and then proceed Ethanol Precipitation below.
	<b>No</b> eth	<b>te:</b> Before starting the procedure, be sure to add isopropanol to the Loading Buffer and anol to the Wash Buffer as described on page v.
	1.	Add 500 µl of Loading Buffer prepared as directed on page v to the neutralized cDNA (from <b>Alkaline Hydrolysis and Neutralization,</b> Step 4, previous page). Mix well by vortexing.
	2.	Place a S.N.A.P. <sup>™</sup> Column on a collection tube and load your sample on the S.N.A.P. <sup>™</sup> Column.
	3.	Centrifuge at 14,000 $\times$ g at room temperature for 60 seconds. Remove the collection tube and discard the flow-through.
	4.	Place the S.N.A.P. <sup><math>TM</math></sup> Column onto the same collection tube and add 700 µl of Wash Buffer prepared as directed on page v,
	5.	Centrifuge at 14,000 × $g$ at room temperature for 60 seconds. Remove the collection tube and discard the flow-through.
	6.	Repeat Steps 4–5 above.
	7.	Centrifuge one more time at $14,000 \times g$ at room temperature for 60 seconds. Remove the collection tube and discard the flow-through.
	8.	Place the S.N.A.P. <sup>™</sup> Column onto a new 1.5-ml microcentrifuge tube.
	9.	Add 50 µl of DEPC-treated water to the S.N.A.P. <sup>TM</sup> Column and incubate at room temperature for 1 minute. Centrifuge at $14,000 \times g$ at room temperature for 1 minute.
	10.	Repeat Step 9, using the same microcentrifuge tube. Proceed directly to <b>Ethanol Precipitation</b> on the next page.
Ethanol Precipitation	In t cDl	the second part of the purification procedure, perform an ethanol precipitation of the NA:
•	1.	Add 10 µl of 3 M Sodium Acetate, pH 5.2, to the tube from Step 10, previous page.
	2.	Add 2 $\mu$ l of 20 mg/ml glycogen to the tube and mix.
	3.	Add 300 $\mu l$ of ice-cold 100% ethanol, and incubate the tube at –20°C for at least 30 minutes.
		Note: You can incubate the tube overnight if desired.
	4.	Spin the tube at 14,000 × $g$ at 4°C for 10–20 minutes. Carefully remove and discard the supernatant.
		<b>Note:</b> You can spin at room temperature if a refrigerated centrifuge is unavailable; however, the yield may be slightly less.
	5.	Add 250 µl of ice-cold 75% ethanol and spin the tube at 14,000 × $g$ for 2 minutes. Carefully remove and discard the supernatant.
	6.	Air dry the sample for 10 minutes.
	7.	Resuspend the sample in 5 µl of 2X Coupling Buffer.

# Labeling with Fluorescent Dye

Introduction	After cDNA synthesis and purification, you are ready to label the amino-modified cDNA with fluorescent dye.			
Dye Information	This kit has been validated with the following dyes and dye packs:			
	Alexa Fluor <sup>®</sup> 555 Reactive Dye Decapack (10 vials) (A-32756) Alexa Fluor <sup>®</sup> 647 Reactive Dye Decapack (10 vials) (A-32757) Alexa Fluor <sup>®</sup> 555 and Alexa Fluor <sup>®</sup> 647 Reactive Dye Decapacks (10 vials each dye) (Cat. no. A-32755)			
	CyDye Post-Labeling Reactive Dye Pack (12 vials each Cy3 <sup>™</sup> and Cy5 <sup>™</sup> ) (Amersham Biosciences, #RPN 5661) Cy3 <sup>™</sup> Mono-Reactive Dye Pack (Amersham Biosciences, #PA 23001) Cy5 <sup>™</sup> Mono-Reactive Dye Pack (Amersham Biosciences, #PA 25001)			
	This kit is also compatible with other monofunctional, NHS-reactive fluorescent dyes.			
Before Starting	The following items will be used in the following procedure:			
_	• DMSO (supplied in the kit)			
	• DEPC-treated water (if using Alexa Fluor <sup>®</sup> dves)			
	• Fluorescent dye(s)			
Important	Fluorescent dyes are sensitive to photobleaching. When preparing the reaction, be careful to minimize exposure of the dye solution to light. The dye coupling reaction must be incubated in the dark.			
CAUTION	DMSO is hygroscopic and will absorb moisture from the air. Water absorbed from the air will react with the NHS ester of the dye and significantly reduce the coupling reaction efficiency. Keep the DMSO supplied in the kit in an amber screw-capped vial at $-20^{\circ}$ C, and let the vial warm to room temperature before opening to prevent condensation.			
Labeling Reaction	Follow the steps below to couple fluorescent dye to your amino-modified first-strand cDNA. Use only the DMSO provided with this kit.			
	1. <b>Alexa Fluor<sup>®</sup> dye</b> vials — Add 2 μl of DMSO directly to each dye vial.			
	<b>Cy3<sup>™</sup> or Cy5<sup>™</sup> dye</b> vials — Individual reaction size (RPN5661): Add 5 µl DMSO directly to each dye vial. Large size (PA23001 and PA25001): Add 45 µl DMSO directly to each dye vial. Use 5 µl of this DMSO/dye solution in the next step.			
	Dye from another manufacturer: Prepare as directed.			
	2. Add the DMSO/dye solution to the tube from <b>Ethanol Precipitation</b> , Step 7, previous page. If you are using Alexa Fluor <sup>®</sup> dyes, add 3 μl of DEPC-treated water to the tube to bring the total volume to 10 μl.			
	3. Mix well and incubate the tube at room temperature in the dark for 1 hour. Reaction can be stored overnight if necessary. Store any unused dye solution according to manufacturer's directions.			

# Purifying Labeled cDNA

Introduction	In this step, you purify the labeled cDNA to remove any unreacted dye. Cat nos. L1014-01 and L1014-02 include a Purification Module developed for use with the system. Follow the procedure below to purify your labeled cDNA.			
	Cat lab	no. L1014-03 does not include a Purification Module. Use your preferred method of eled cDNA purification, and then continue to hybridization.		
Before Starting	The following items are supplied by the user:			
	•	Microcentrifuge		
	The L1(	e following items are supplied in the Purification Module (Cat nos. L1014-01 and )14-02):		
	•	3 M Sodium Acetate, pH 5.2		
	٠	DEPC-treated water (supplied in the kit)		
	•	S.N.A.P. <sup><math>m</math></sup> column(s) and collection tube		
	٠	Amber collection tubes		
	•	Loading Buffer <b>plus isopropanol</b> (see page v for preparation)		
	•	Wash Buffer <b>plus ethanol</b> (see page v for preparation)		
S.N.A.P. <sup>™</sup> Column Purification	Use	e the following procedure to purify the cDNA:		
	1.	Add 20 µl of 3 M Sodium Acetate, pH 5.2, to the dye-coupled cDNA solution from Step 3, page 9.		
	2.	Add 500 µl of Loading Buffer plus isopropanol to the cDNA solution. Mix well by vortexing.		
	3.	Place a S.N.A.P. <sup>™</sup> Column onto a clear collection tube and load the cDNA/buffer solution onto the S.N.A.P. <sup>™</sup> Column.		
	4.	Centrifuge at 14,000 × $g$ at room temperature for 60 seconds. Remove the collection tube and discard the flow-through.		
	5.	Place the S.N.A.P. <sup>™</sup> Column on the same collection tube and add 700 µl of Wash Buffer plus ethanol to the column.		
	6.	Centrifuge at 14,000 × $g$ at room temperature for 60 seconds. Remove the collection tube and discard the flow-through.		
	7.	Repeat Steps 5–6 of this procedure, using the same collection tube.		
	8.	Centrifuge one more time at $14,000 \times g$ at room temperature for 60 seconds. Remove the collection tube and discard the flow-through.		
	9.	Place the S.N.A.P. <sup>™</sup> Column onto a new <b>amber</b> collection tube.		
	10.	Add 50 µl of DEPC-treated water to the S.N.A.P. <sup><math>TM</math></sup> Column and incubate at room temperature for 1 minute.		
	11.	Centrifuge at 14,000 × g at room temperature for 1 minute and collect the flow-through. The flow-through contains your purified dye-coupled cDNA.		
	The free Lat	e sample can stored at –20° C for up to one week prior to hybridization. Avoid eze/thawing. To determine the efficiency of the labeling reaction, proceed to <b>Assessing</b> peling Efficiency on page 12.		

**Hybridization** After purification, you are ready to use the labeled cDNA in any application of choice, including glass microarray hybridization. Follow the preparation and hybridization instructions for your specific application.

# Appendix

#### Assessing Labeling Efficiency

Introduction	You Flue exp	ı can use the following proc or®-labeled or CyDye™-label ected amount of labeled cD	edure and led cDNA NA using	l formulas to measure and determine the ei the Control HeLa RN	e the amount of Alexa ficiency of the reaction. The JA is noted below.	
Absorption Wavelengths and Baselines	The following table shows the absorbance and baseline wavelengths for CyDyes $^{\mbox{\tiny TM}}$ and Alexa Fluor $^{\mbox{\tiny B}}$ dyes:					
	Label		Absorbance Wavelength		<b>Baseline Wavelength</b>	
	Ale Ale	xa Fluor® 555 or Cy3™ xa Fluor® 647 or Cy5™		550 nm 650 nm	650 nm 750 nm	
Calculating the Amount of Coupled Dye	To calculate the amount of coupled dye:					
	1.	Transfer the <b>undiluted</b> sample from Step 11, page 10, into a clean cuvette, and scan at 240–800 nm using a UV/visible spectrophotometer. If you are using a 100-µl cuvette, transfer the entire sample; for smaller cuvettes, transfer an appropriate amount of sample.				
		<b>Note:</b> The labeled cDNA must be purified as described on page 10 before scanning, as any unreacted dye will interfere with the detection of labeled cDNA.				
	2.	Calculate the amount of lal	beled cDN	JA using the formula	below:	
		$cDNA (ng) = (A_{260} - A_{320}) \times 37 ng/\mu l \times Volume in \mu l$				
	The amount of cDNA generated from 10 $\mu$ g of Control HeLa RNA should be >250 ng. If it is <250 ng, see <b>Troubleshooting</b> on page 13.					
	3.	Calculate the amount of flu	ıorescentl	y labeled dye using a	formula below:	
		Alexa Fluor <sup>®</sup> 555 (pmole) = $(A_{550}-A_{650})/0.15 \times 50$ (elution volume)				
		Alexa Fluor <sup>®</sup> 647 (pmole) = $(A_{650}-A_{750})/0.24 \times 50$ (elution volume)				
		$Cy3^{TM}$ (pmole) = (A <sub>550</sub> -A <sub>650</sub> )	/0.15 × 50	) (elution volume)		
		$Cy5^{TM}$ (pmole) = (A <sub>650</sub> -A <sub>750</sub> )	/0.25 × 50	) (elution volume)		

## Troubleshooting

Observation	Cause	Solution
28S and 18S bands are not observed after isolation of total RNA and agarose gel electrophoresis	Too little RNA loaded on the gel	Be sure to load at least 250 ng of RNA for analysis.
	RNA is degraded due to RNase activity	Follow the guidelines on page 3 to avoid RNase contamination. Use a fresh sample for RNA isolation.
Yield of cDNA from the first-strand synthesis reaction is low	Temperature too high during cDNA synthesis	Perform the cDNA synthesis at 46°C.
	Incorrect reaction conditions used	Verify that all reaction components are included in the reaction and use reagents provided in the system. Verify the reaction conditions using the control RNA provided in the kit.
	Concentration of template RNA is too low	Increase the concentration of template RNA. Use at least 5 $\mu$ g of total RNA or 0.4 $\mu$ g of mRNA.
	Poor quality RNA used or RNA is degraded	Check the quality of your RNA preparation (see page 4). If RNA is degraded, use fresh RNA.
	RNase contamination	Use the RNaseOUT <sup>™</sup> included in the kit to prevent RNA degradation.
	RT inhibitors are present in your RNA sample	Inhibitors of RT include SDS, EDTA, guanidinium chloride, formamide, sodium phosphate and spermidine (Gerard, 1994). Remove inhibitors from your RNA sample by performing an additional 70% ethanol wash after ethanol precipitation during RNA isolation and purification.
		Test for the presence of inhibitors by mixing 1 $\mu$ g of control RNA with 25 $\mu$ g total RNA or 1 $\mu$ g mRNA and compare the yields of first-strand synthesis.
	Improper storage of SuperScript <sup>®</sup> III RT	Store the enzyme at $-20^{\circ}$ C.
Yield of labeled cDNA from the control reaction is low	cDNA has been lost in the purification step following cDNA synthesis	Follow the S.N.A.P. <sup>™</sup> Column purification procedure and ethanol precipitation procedure without modifications. Overnight ethanol precipitation may result in higher yields.
	cDNA has been lost in the purification step after labeling	Measure the amount of labeled cDNA in the control reaction before and after purification. Follow the purification procedure without modifications.

## Troubleshooting, continued

Observation	Cause	Solution
Amount of coupled dye in the control reaction is low (<40 pmoles) and/or fluorescence of labeled cDNA is low	Reaction tubes have been exposed to light	Avoid direct exposure of the labeling reaction to light. Use an amber tube for collection of the final product.
	Dye solution has been exposed to light	Repeat labeling reaction with fresh mixture of dye, being careful to avoid direct exposure to light.
	DMSO used to prepare dye mixture was contaminated with water.	Prepare a new mixture of dye using fresh DMSO. Carefully follow the instructions for storing and handling DMSO in the <b>Caution</b> on page 9.
	Inefficient labeling due to improper purification	Follow all purification steps carefully and without modification.
	2X Coupling Buffer was not stored properly	Store 2X Coupling Buffer at –20°C.

## **Technical Support**

Obtaining Support	For the latest services and support information for all locations, go to <b>www.lifetechnologies.com</b>			
	At the website, you can:			
	• Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities			
	• Search through frequently asked questions (FAQs)			
	Submit a question directly to Technical Support (techsupport@lifetech.com)			
	• Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents			
	Obtain information about customer training			
	Download software updates and patches			
Safety Data Sheets (SDS)	Safety Data Sheets (SDSs) are available at <b>www.lifetechnologies.com/sds</b>			
Certificate of Analysis	The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to <b>www.lifetechnologies.com/support</b> and search for the Certificate of Analysis by product lot number, which is printed on the box.			
Limited Warranty	Life Technologies and/or its affiliate(s) warrant their products as set forth in the Life Technologies General Terms and Conditions of Sale found on the Life Technologies web site at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies.			
	Life Technologies and/or its affiliate(s) disclaim all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Life Technologies and/or its affiliate(s) be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.			

Limited Use Label License No. 358: Research Use Only	The purchase of this product conveys to the purchaser the limited, non- transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining
	activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing
	Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

\_\_\_\_\_

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994) *Current Protocols in Molecular Biology*, Greene Publishing Associates and Wiley-Interscience, New York
- De Risi, J., Penland, L., Brown, P. O., Bittner, M. L., Meltzer, P. S. R., M., Chen, Y., Su, Y. A., and Trent, J. M. (1996) Use of a cDNA microarray to analyse gene expression patterns in human cancer. *Nature Genet.*, 14, 457-460
- Eisen, M. B., and Brown, P. O. (1999) DNA arrays for analysis of gene expression. *Methods Enzymol.*, 303, 179-205
- Gerard, G. F. (1994) Inhibition of SuperScript II Reverse Transcriptase by Common Laboratory Chemicals. *FOCUS*, *16*, 102-103
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Second Ed., Cold Spring Harbor Laboratory Press, Plainview, New York

©2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. CyDye, Cy3, and Cy5 are trademarks of Amersham Biosciences.

#### Notes

#### Notes

#### Notes

