

Revised: 29-July-2004

CandyCane™ Glycoprotein Molecular Weight Standards (C21852)

Quick Facts

Storage upon receipt:

- -20°C (long term)
- 4°C (short term)

Concentration: Each protein ~0.5 µg/µL

Introduction

CandyCaneTM glycoprotein molecular weight standards contain a mixture of glycosylated and non-glycosylated proteins ranging from 14,000 to 180,000 molecular weight. When separated by polyacrylamide gel electrophoresis, the standards appear as alternating bands corresponding to glycosylated and nonglycosylated proteins (Figure 1). Thus, these standards serve both as molecular weight markers and as positive and negative controls for methods that detect glycoslyated proteins, such as the methods provided in our Pro-Q[®] Glycoprotein Stain Kits (see *Product List*, below).




Figure 1. Glycosylated and nonglycosylated proteins in the Candy-Cane glycoprotein molecular weight standards. The standards were electrophoresed through two identical 13% polyacrylamide gels. Both lanes contain ~0.5 μ g of protein in each band. The left lane is stained with SYPRO[®] Ruby protein gel stain to detect all eight marker proteins. The right lane is stained using Pro-Q Emerald 300 glycoprotein detection reagent. The protein molecular weights are indicated, and the identities of the proteins can be found in Table 1.

Protein	Molecular Weight	Glycoprotein?	Con A	WGA	GS-II
α_2 -Macroglobulin ¹	180,000	yes *	+	+	_
Phosphorylase b ²	97,000	no	-	-	-
Glucose oxidase ³	82,000	yes	+	+	+
Bovine serum albumin ⁴	66,000	no †	(+)	(+)	-
α_1 -Acid glycoprotein ⁴	42,000	yes	(+)	+	-
Carbonic anhydrase ⁵	29,000	no	-	-	(+)
Avidin ⁶	18,000	yes *	+	+	-
Lysozyme ⁶	14,000	no	_	_	_

Table 1. Characteristics of proteins included in the CandyCane glyco-

Con A = concanavalin A; WGA = wheat germ agglutinin; GS-II = *Griffonia simplicifolia* lectin II. Staining with the indicated lectins is indicated by either + strong signal, (+) weak signal, or – no signal. 1. from human plasma;
2. from rabbit muscle; 3. from *Aspergillis niger*; 4. from bovine serum;
5. from bovine erythrocyte; 6. from chicken egg. * These glycoproteins may not transfer efficiently to nitrocellulose or poly(vinylidene fluoride) (PVDF) membranes. † A comigrating glycoprotein contaminant may be visible with wheat germ agglutinin– and concanavalin A–detection strategies.

Materials

Contents

The CandyCane glycoprotein molecular weight standards are sold in a unit size of 400 μ L, sufficient volume for approximately 200 gel lanes. Each protein is present at ~0.5 μ g/ μ L. The proteins are provided in a storage solution that contains 10 mM Tris, pH 7.5, 200 mM NaCl, 15 mM dithiothreitol, 2 mM EDTA, 3 mM sodium azide and 35% glycerol. Characteristics of the proteins and staining patterns with various lectins are listed in Table 1.

Storage

For long-term storage, the standards should be kept at -20°C to prevent microbial contamination.

Application

For 8 cm \times 10 cm gels, mix 0.5 µL of the protein standard with 7.5 µL of SDS gel-loading buffer, heat at 95°C for 4 minutes and pipet into the well reserved for the standard. This will result in ~250 ng of each protein per lane. For larger gels, increase the amount of standards. Figure 1 shows the standards and their molecular weights.

Product List Current prices may be obtained from our Web site or from our Customer Service Department.

Cat #	Product Name	Unit Size
C21852 P21855 P21857 P21875 P21872	CandyCane [™] glycoprotein molecular weight standards *200 gel lanes* Pro-Q [®] Emerald 300 Glycoprotein Gel Stain Kit *with SYPRO [®] Ruby protein gel stain* *10 minigels* Pro-Q [®] Emerald 300 Glycoprotein Gel and Blot Stain Kit *10 minigels or mingel blots* Pro-Q [®] Emerald 488 Glycoprotein Gel and Blot Stain Kit *10 minigels or minigel blots* Pro-Q [®] Glycoprotein Blot Stain Kit #5 *with <i>Griffonia simplicifolia</i> lectin II (GS-II) and DDAO phosphate* *>20 minigel blots*	1 kit 1 kit

Contact Information

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Leiden, the Netherlands. All others should contact our Technical Assistance Department in Eugene, Oregon.

Please visit our Web site - www.probes.com - for the most up-to-date information

Molecular Probes, Inc.	Molecular Probes Europe BV
29851 Willow Creek Road, Eugene, OR 97402	Poortgebouw, Rijnsburgerweg 10
Phone: (541) 465-8300 • Fax: (541) 335-0504	2333 AA Leiden, The Netherlands
	Phone: +31-71-5233378 • Fax: +31-71-5233419
Customer Service: 6:00 am to 4:30 pm (Pacific Time)	
Phone: (541) 335-0338 • Fax: (541) 335-0305 • order@probes.com	Customer Service: 9:00 to 16:30 (Central European Time)
	Phone: +31-71-5236850 • Fax: +31-71-5233419
Toll-Free Ordering for USA and Canada:	eurorder@probes.nl
Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228	
	Technical Assistance: 9:00 to 16:30 (Central European Time)
Technical Assistance: 8:00 am to 4:00 pm (Pacific Time)	Phone: +31-71-5233431 • Fax: +31-71-5241883
Phone: (541) 335-0353 • Toll-Free: (800) 438-2209	eurotech@probes.nl
Fax: (541) 335-0238 ● tech@probes.com	

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. Our products are not available for resale or other commercial uses without a specific agreement from Molecular Probes, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to busdev@probes.com. All names containing the designation [®] are registered with the U.S. Patent and Trademark Office.

Copyright 2004, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.