

Leo

Precision in your hands



okolab



Features

- User-friendly touch screen interface
- Compatible with any incubator on the market
- Diffusion or Aspiration gas sampling modes
- Pump suction rate 100 ml/min
- Sampling tube and de-humidification kit included
- Long term logging
- Micro USB for data download
- Calibration:
 - a) Removable sensors can be shipped to certified laboratories to be calibrated. Calibration is stored in sensor on-board memory. Leo reads calibration upon sensor installation
 - b) In situ sensor calibration via LEO calibration

CO₂-O₂ Module - Always included



CO₂-O₂ module can be easily extracted from LEO and sent out for calibration.

	Type	Range	Accuracy
CO ₂	Non Dispersive InfraRed (NDIR) dual wave length detector with pressure and temperature compensation	0-20%	±(1.0% of full scale +2% scaling) at 6% = ±0.32%
O ₂	Fluorescence based optical sensor	0-25%	±(1.0% of full scale) at 6% = ±0.25%

The module stores calibration data in its own memory, that Leo reads when the module is put in place.



LEO is a battery operated, handheld analyzer designed to measure CO₂, O₂, Temperature, VOC, Relative Humidity and pH.

Includes CO₂ and O₂ sensors.

Optional modules: Temperature, Humidity, VOC and pH



CO₂ - O₂ Sampling

Diffusion Mode

For incubator featuring gas sampling port with spontaneous outlet



Gas flows into LEO through a dedicated diffusion inlet. Suction pump is not activated.

Aspiration Mode

For incubator featuring gas sampling port without a spontaneous outlet





LEO activates its **suction pump** to sample gas from the incubator. LEO releases the sampled gas through its output port. Gas can be VOC filtered and reinjected into the incubator.



T Modules


Leo features two external temperature modules, which can operate at the same time:

- T1 is a small, flexible immersible thermocouple, ideal for measurements in liquids, such as the culture media in the dish.
- T2 is a thermistor, ideal for measurements of incubator temperature, accessed through the sampling port.




	Type	Range	Accuracy
T1	K-type Thermocouple. PFA-insulated; Length: 2 m; OD: 0.13 mm	20-45 °C	± 0.1 °C
T2	PT1000 Probe length: 150 mm; OD: 1.5 mm; Cable length: 1 m	20-45 °C	± 0.1 °C

Each sensor module stores calibration data in its own local memory, that LEO reads upon connection.



T1



T2



	Type
LEO	Handheld, touch screen, battery operated analyzer. Includes CO ₂ -O ₂ module, sampling tube and de-humidification kit.
CO ₂ -O ₂ module (included in Leo)	CO ₂ : Non Dispersive InfraRed (NDIR) dual wave length detector O ₂ : fluorescence based optical sensor
T1 module	K-type Thermocouple; OD: 0.13 mm
T2 module	PT1000 Thermoresistance - Probe length: 150 mm; OD: 1.5 mm
H module	Linearized and temperature compensated sensor. Range: 0-100% Accuracy: ±1.5% RH
VOC module	Photoionization detector (PID). Range:0-2ppm Accuracy: ±3% of reading Minimum detection limit: 0.5 ppb
pH module	Automatic temperature compensation probe Accuracy @ 20°C: pH: ±0.01 pH Range: 0-14
CO ₂ /O ₂ calibration kit	6 mm O tygon tube with a calibrated orifice deliveing a flow of 100 ml/min to Leo when connected to a gas source at 1 barg (14.7 psig)
Hard Case	Hard travel case for Leo and its accessories