Keep this rotor instruction manual and the centrifuge manual in the file.

Swing Rotors for Refrigerated Centrifuge

R4SS

Before using these rotors, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

Contents

1. Specifications	_____	1
2. Tubes/Bottles	_____	3
3. How to use	_____	5
4. Maintenance	_____	7
5. Packing list	_____	9
6. Spare parts	_____	9
7. Decontamination	_____	10
8. Rotor Retirement	_____	10
(Appendix) K Factors	_____	11

Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.

These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

WARNING


- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor body).
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density, kinds of microplates, etc.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

CAUTION

- Balance the tubes/bottles and the samples including the caps or adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use the tubes/bottles within their actual capacities.
- Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge. The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature. Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/bottles if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it. Contact an authorized sales or service representative.

1. Specifications

1.1 Construction

Over-speed adapter(Aluminum alloy)
Red : for CR-NX series Blue : for other than CR-NX series
Figure1-1 Construction of R4SS rotor

1.2 Specifications

Applicable centrifuge*1

(Some centrifuges are on sale in only Japan.
Visit to our website.)

Max. speed*2: 4,000 rpm

Max. RCF: 2,900 × g

Allowable mean density of sample: 1.2 g/mL

Nominal rotor capacity: 15 mL × 40 tubes = 600 mL

Rotor Size: ϕ 348 × 154 mm

Weight: 4.9 kg(including buckets)


Material of rotor: Stainless steel

Material of bucket: Aluminum alloy

*1: For applicable centrifuges, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”.

*2: The maximum allowable speed varies with centrifuge models and so on. For further details, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”.

1.3 Cross-sectional view of rotor

*: Rmin. is 5.1 cm when you use 15 mL culture tube.

Figure1-2 Cross-sectional view of R4SS rotor

1.4 Characteristic of rotor

Table1-1 Characteristics

Rotor speed (rpm)	RCF (× g)		K factor ※
	R = Rmin. (5.1 cm)	R = Rmax. (16.2 cm)	
1,000	57	181	292,760
1,500	128	408	130,116
2,000	228	724	73,190
2,500	356	1,130	46,842
3,000	513	1,630	32,529
3,500	698	2,220	23,899
4,000	912	2,900	18,298

※For details about K factor, refer to “(Appendix) K Factors”

Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)

N: speed (rpm)

1.5 Allowable speed of the rotor

WARNING :

Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic, kinds of tubes/bottles, and centrifuge model.
Do not exceed the allowable speed of the rotor.

The speed must be limited depending on the total weight of contents (adapter, tubes/bottles, and sample) in the bucket. When the total weight is over 300 g, calculate the allowable speed according to the following equation.

(1) Allowable speed depending on the total weight of contents in the bucket.

WARNING :

the total weight of contents : adapter, tubes/bottles, and sample

$$\text{Allowable speed (rpm)} = 4,000 \text{ (rpm)} \times \sqrt{\frac{300 \text{ (g)}}{\text{Total weight of contents in the bucket (g)}}}$$

1.6 Rotor registration data for the CR-G II series centrifuges

If this rotor is not registered in the CR-G II series centrifuges, register this rotor by entering the following data according to “3-5-5 Entry of New Rotor” on the instruction manual of the CR-G II series centrifuges.

Item/No.	DATA	Item/No.	DATA	Item/No.	DATA	Item/No.	DATA
ROTOR No.	79	6	2088	12	83083	18	0
1	49000	7	4176	13	32032	19	0
2	15164	8	100	14	32032	20	0
3	30148	9	18000	15	15160	21	8024
4	6084	10	11000	16	1244	22	1255
5	4086	11	82052	17	0	23	111241

1.7 Registration of Rotor data for the CR-NX series centrifuges

If this rotor is not registered in the CR-NX series centrifuges, Contact an authorized sales or service representative.

1.8 Rotor Compatibility

The CR-NX series are compatible with the CR-NX rotors only. For further details about applicable centrifuges, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”.

2. Tubes/Bottles

2.1 Applicable tubes/bottles

- Adapters, tubes and bottles are not included in this rotor. Please order applicable ones selecting from the following list (designate the part numbers) or purchase applicable ones on the market.

CAUTION:

15TCX6S and 50TCX2S adapters are specifically designed for conical culture tubes. Do not use round-bottom or self-standing culture tubes.

- When using tubes or bottles on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer and do not use any tube or bottle whose allowable speed or RCF is not specified. Otherwise the tubes or bottles may be broken during operation. We will check the applicability of tubes or bottles upon request. Contact an authorized sales or service representative.
- Clean and sterilize tubes or bottles on the market according to the instructions of the manufacturer.
- Check the size of tubes before using the tubes (see the following list). When the tube length is longer than the tube length shown in the following list, you cannot use the tube in your adapter.
- Maximum speeds listed are guidelines only. Because of variances in user methodologies, no guarantee of performance is expressed or implied.


Table 2-1 Applicable tubes/bottles

Adapter			Tube				Max, speed (Max. RCF) (Max. radius)
PART NO.	Name	Nominal Capacity x Q'ty	PART NO.	Name	Actual Capacity	Size (mm)	
S308546A	15 x 10S adapter (4 pcs.) (Material: *1)	15mL x 40	(Marketed)	Blood collection tube	*4	(D) 15.4 – 16.0×(L) 100	4,000rpm(2,860 × g) (Rmax.:16.0cm)
			464976A	15 glass tube (50 pcs.)	13mL	(D) 16.5×(L) 105	3,000rpm(1,610 × g) (Rmax.:16.0cm)
S307334A	15 x 8S adapter (4 pcs.) (Material: *2)	15mL x 32	(Marketed)	Blood collection tube	*4	(D) 15.4 – 16.0 ×(L) 100	4,000rpm(2,860 × g) (Rmax.:16.0cm)
			464976A	15 glass tube (50 pcs.)	13mL	(D) 16.5×(L) 105	3,600rpm(2,320 × g) (Rmax.:16.0cm)
S308695A	10 x 10SB adapter (4 pcs.) (Material: *2) *3	10mL x 40	(Marketed)	Blood collection tube	*4	(D) 12.4 – 15.6 ×(L) 75 – 100	4,000rpm(2,800 × g) (Rmax.:15.7cm)
S307333A	10 x 8S adapter (4 pcs.) (Material: *2)	10mL x 32	(Marketed)	Blood collection tube	*4	(D)12.4 – 13.2 ×(L) 100	4,000rpm(2,860 × g) (Rmax.:16.0cm)
			464977A	10 glass tube (50 pcs.)	7.5mL	(D) 12.0×(L) 105	3,600rpm(2,320 × g) (Rmax.:16.0cm)
S307332A	5 x 8S adapter (4 pcs.) (Material: *2)	5mL x 32	(Marketed)	Blood collection tube	*4	(D) 12.4 – 13.2 ×(L) 75 – 78	4,000rpm(2,380 × g) (Rmax.:13.3cm)
S308844A	50 x 2S adapter (4 pcs.) (Material: *2)	50 mL x 8	473866A	50 glass tube (10 pcs.)	45mL	(D) 35.0×(L) 100	3,500rpm(2,120 × g) (Rmax.:15.5cm)
S307335A	15TC x 6S adapter (4 pcs.) (Material: *2)	15 mL x 24	(Marketed)	15 mL culture tube (conical bottom)	15mL	(D) 17.0×(L) 120	4,000rpm(2,900 × g) (Rmax.:16.2cm)
S307336A	50TC x 2S adapter (4 pcs.) (Material: *2)	50 mL x 8	(Marketed)	50 mL culture tube (conical bottom)	50mL	(D) 30.0×(L) 115	4,000rpm(2,860 × g) (Rmax.:16.0cm)
		Unit (mm)					Unit (mm)

*1: Polypropylene, stainless steel *2: Polypropylene, nitrile rubber

*3: Check the size of the applicable blood collection tubes before using the blood collection tubes (see the following).

*4: Actual capacities of blood collection tubes vary depending on the applications. Refer to the manufactures catalog.

3. How to use

3.1 Preparation of the rotor

- (1) Check that there is no crack, corrosion or deformation on the rotor and the buckets.
 - Check that the drive pin of the drive hole (crown hole) is not deformed.
 - Check that there is no corrosion, foreign matters, dirt, and scratches in the inside bottom of the bucket and in the drive hole (crown hole).

CAUTION:

If there is any abnormality such as corrosion or cracks, stop use of the rotor and the buckets and contact an authorized sales or service representative.

- (2) Check that there is no foreign matter in the buckets.

3.2 Preparation of the sample and the tubes/bottles

Inject the sample into the tubes/bottles.

WARNING :

- These rotors and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

CAUTION:

- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor (including buckets) and adapters. Using such a sample could corrode or deteriorate them.
- Chemical resistance and the strength of tubes/bottles varies with speed, temperature, and so on. Before using sample, fill tubes/bottles with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the actual capacity specified in the manufacturer. Otherwise the sample leakage may occur during operation.
- Use the tubes/bottles that are same products. Otherwise, the tubes/bottles may be damaged due to their different form or the imbalance operation may occur. And it may result in the damage to the centrifuge and the rotor.
- Use the tubes/bottles at a temperature between 4°C and 25°C. Otherwise, tubes/bottles may be broken or deformed during operation.
- Before using tubes/bottles on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.
- It takes much time for this rotor to reach the desired temperature. We recommend that the rotor, buckets, adapters, tubes/bottles, and samples are warmed or cooled to the desired temperature before operation.

3.3 Balancing the microplates

The allowable imbalance of this rotor is 2g (Allowable imbalance: Difference in weight of contents (adapters, tubes/bottles, and samples) in the symmetrically arranged buckets).




CAUTION:

- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise, the precipitation levels may be different by centrifugation and there may be variations in position of center of gravity. And it may result in damage to the centrifuge and the rotor due to the imbalance operation.
- Balance between symmetrically arranged microplates filled with the same-density sample. Variations in sample density cause imbalance operation. And it may result in the damage to the centrifuge and the rotor.

Use of tubes with caps is recommended to prevent contamination of the sample with dust or water.

Check that there is no foreign matter in the adapters then load the tubes. When loading tubes less than the number of tube cavities in the adapters, load the tubes symmetrically with respect to the center of the adapters. Examples of incorrect tube setting are shown below.

Center of adapter

●: Tubes are loaded.

○: Tubes are not loaded.



When it is impossible to load the tubes symmetrically for lack of a tube, prepare a dummy tube filled with the same-density solution and load the tubes symmetrically.

3.4 Mounting the rotor and the buckets

Mount the balanced buckets with respect to the center of the rotor. And Mount the rotor to the drive shaft.

WARNING:

- Be sure to set all the buckets regardless of the presence or absence of sample. Operation with only two buckets is very dangerous because the rotor may be deformed and the buckets may come off. Never use any bucket or adapter that is not designated for this rotor by us.
- Be sure to set all the buckets of the same serial number (The serial number is shown on the upper side of the bucket). Otherwise, it may cause the imbalance operation and it may result in the damage to the centrifuge and the rotor.
- Securely hook the buckets on the rotor pins. Otherwise the buckets may come off while the rotor is running.

CAUTION:

- Do not load adapters and tubes/bottles asymmetrically. Asymmetrically loading may cause imbalance operation and the centrifuge may be damaged.

Check that the buckets swing smoothly. If the swinging is not smooth, clean both the rotor pins and the pin holes of the buckets and apply lubricant for screw (standard accessory) lightly.

3.5 Operation

Run the rotor according to the procedures specified in the centrifuge instruction manual.

3.6 Taking out samples

- (1) After the operation, carefully remove the rotor from the centrifuge.
- (2) Remove the adapters from the buckets. Then take out the separated sample from the tubes/bottles.

4. Maintenance

4.1 Corrosion resistance of the rotor, buckets, and adapters

⚠ WARNING:
Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor, buckets, and adapters.

The material of the adapters is polypropylene and nitrile rubber. The rotor is made of stainless steel. The buckets are made of aluminum alloy. The buckets made of aluminum alloy have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode them and decrease their strength. Use samples that will not affect them referring to the attached chemical resistance chart.

4.2 Maintenance of the rotor and the buckets

(1) Maintenance of the rotor and the buckets

After the operation, properly maintain the rotor and the buckets to prevent corrosion that can cause them breakage.

- Normal maintenance

⚠ CAUTION:
Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor and buckets with tap water or a dilute solution of neutral detergent and rinse them out with distilled water. Wipe them with a soft cloth. Dry the rotor which the buckets and adapters are removed from. Check that the rotor is completely dried and then put a light coat of silicone grease (vacuum grease)(standard accessory of the centrifuge) on the rotor. Store the rotor in a dry place that is not exposed to the direct sunlight.

- Maintenance after use of a corrosive sample

Wash the rotor and buckets with tap water immediately after the operation. Then perform the normal maintenance.

- Maintenance when foreign substances are adhered to the rotor and buckets

⚠ CAUTION:
Use a neutral detergent having a pH between 5 and 9, otherwise the rotor and buckets can be discolored or corroded.

Soak the rotor and the buckets in warm water for one or two hours and wash the inside of the buckets and the crown hole at the rotor bottom with a soft brush to remove foreign substances. Do not forget to remove foreign substances adhered in the crown hole at the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

(2) Maintenance of the adapters

- In the event sample leaked from the tube to the inside of the adapter, tilt the adapter to drain sample from the weep hole, and rinse it out with tap water (Fig.1). Turn the adapter upside down to dry it well.
- In the event a tube is broken while the rotor is running and the fragments or foreign matter remains in the adapter, disassemble the adapter according to the following procedure to remove the fragments or foreign matter (Fig.2).




① Insert a flat-blade screwdriver into the groove between the adapter(U) and the adapter(L).

② Turn the flat-blade screwdriver about 90° to raise the adapter(U).

Perform the above steps ① and ② for the opposite side.

③ Remove the adapter(U) from the adapter(L).

④ Clean the inside of the adapter completely and reinstall the adapter(U) and the adapter(L).

4.3 Sterilizing

WARNING:

Never sterilize rotor, buckets, and adapters by boiling, otherwise they will be brittle.

(1) Sterilizing the rotor

Sterilize the rotor in any of the following methods; gas sterilization (ethylene oxide or formaldehyde), chemicals sterilization (ethanol (70%), hydrogen peroxide (3%) or formalin (3%)), or ultraviolet rays sterilization(200-300nm).

CAUTION:

- Never sterilize the rotor by autoclaving. Otherwise the strength of the rotor may be significantly decreased.
- Do not dip the rotor in the formalin (3%) solution more than 2 hours.

(2) Sterilizing the buckets

Sterilize the buckets in any of the following methods; gas sterilization (ethylene oxide or formaldehyde), chemicals sterilization (ethanol (70%), hydrogen peroxide (3%) or formalin (3%)), or ultraviolet rays sterilization(200-300nm). Turn the buckets upside down to dry well.

CAUTION:

- Never sterilize the buckets by autoclaving. Otherwise the strength of the buckets may be significantly decreased.
- Do not dip the buckets in the formalin (3%) solution more than 2 hours.

(3) Sterilizing the adapters

Sterilize the adapters in any of the following method; chemicals sterilization (ethanol (70%) or hydrogen peroxide (3%)). The adapters can not be sterilized by autoclaving due to their specifications.

4.4 Inspecting rotor

Periodically check every 100 hours or once a month whether the rotor and buckets are not corroded. Especially be careful of the inside bottom of the bucket and the crown hole of the rotor bottom, because the rotor and the buckets can be brittle if these portions are corroded. If their surfaces are discolored, dented or cracked, they are corroded. Do not use such a corroded rotor or a bucket. Contact an authorized sales or service representative immediately for further inspection.


CAUTION:

If there is any abnormality such as corrosion or cracks, stop use of the rotor and the buckets and contact an authorized sales or service representative.

5. Packing list

Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	PART No. of spare parts (Qty)
①	Rotor	1	—
②	Bucket	4	—
③	Lubricant for screw	1	84810601
④	Rotor instruction manual	1	S999925

6. Spare parts

Figure	PART No.	Name	Qty
	84810601	Lubricant for screw	1
	S999925	Rotor Instruction manual	1

7. Decontamination

 **WARNING:**

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

8. Rotor Retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

(Appendix) K Factors

• Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit(S).

Namely, the time taken for particles of 1 S to settle from Rmin. to Rmax. is K factor.
The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R_{\max}) - \ln(R_{\min})}{\omega^2} \cdot \frac{10^{13}}{3,600}$$

ω : Angular velocity $\omega = 2\pi/60 \times N$

N: Speed (rpm)

Rmin.: Distance from axis of rotation to top of solution (cm)

Rmax.: Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t(in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, when using the 15 mL culture tube, the K factor at the maximum speed of 4,000 rpm of Model R4SS rotor is 18,298. The time taken for particles of 1,000 S to be settled from R min. to R max. using this rotor is as follows:

$$t = \frac{18,298}{1,000} = 18.3 \text{ hours}$$

The time taken for settling of particles of 1,000 S by the same rotor at a speed of 3,000 rpm is:

$$K = \left[\frac{4,000}{3,000} \right]^2 \times 18,298 = 32,530 \quad t = \frac{32,530}{1,000} = 32.5 \text{ hours}$$

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Copyright© 2021 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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