# Growth and Maintenance of $Flp-In^{TT}$ cell lines USER GUIDE

Catalog Numbers R750-07, R752-07, R758-07, R760-07, R761-07, and R762-07 Document Part Number 250323 Publication Number MAN0000143 Revision 4.0





Life Technologies Corporation | 3175 Staley Road | Grand Island, NY 14072 For descriptions of symbols on product labels or product documents, go to **thermofisher.com/symbols-definition**.

The information in this guide is subject to change without notice.

**DISCLAIMER**: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

#### Revision history: Pub. No. MAN0000143

Revision	Date	Description
4.0	-	Converted to CCMS. Rebranded. Updated contents to state liquid nitrogen vapor phase. Added the Important guidelines for thawing and storing cells topic.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

TRADEMARKS: All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2020 Thermo Fisher Scientific Inc. All rights reserved.

## Contents

Growth and maintenance of Flp-In <sup><math>^{\text{M}}</math></sup> cell lines
Product information
Contents and storage
Description of the system
Methods 10
Culturing Flp-In <sup>™</sup> cell lines
General cell handling
Before starting
Thawing adherent cells
Thawing suspension cells
Passaging adherent cells 12
Passaging suspension cells 12
Freezing cells
Introduction
Preparing freezing medium
Freezing the cells
Transfecting cells
Introduction
Transfection methods 14
Important
Generating stable expression cell lines
Note
Polyclonal selection of isogenic cell lines
Select individual cell lines
Note
APPENDIX A Ordering information
Accessory products
Introduction
Cell culture reagents 17
Additional reagents 17
Flp-In <sup>™</sup> products
References

APPENDIX B Safety	20
Chemical safetyBiological hazard safety	
Documentation and support	23
Customer and technical support	



storage

# Growth and maintenance of Flp-In<sup>™</sup> cell lines

## **Product information**

#### Contents and Shipping and storage

All cell lines are shipped on dry ice. Store in liquid nitrogen vapor-phase upon receipt.

#### Important guidelines for thawing and storing cells

- Upon receipt, immediately thaw cells or place into vapor-phase liquid nitrogen storage until ready to use. **Do not store the cells at -80°C**.
- Avoid short-term extreme temperature changes. When storing cells in liquid nitrogen after shipping on dry ice, allow the cells to remain in liquid nitrogen for 3-4 days before thawing.
- Prior to starting experiments, ensure you have established cells and have frozen stocks on hand. Upon receipt, grow and freeze multiple vials of cells to ensure that you have an adequate supply of early-passage cells.

#### Contents

This manual is supplied with the following cell lines:

Cell Line	Catalog no.
Flp-In <sup>™</sup> -293	R750-07
Flp-In <sup>™</sup> -CV-1	R752-07
Flp-In <sup>™</sup> -CHO <sup>™</sup>	R758-07
Flp-In <sup>™</sup> -BHK	R760-07
Flp-In <sup>™</sup> -3T3	R761-07
Flp-In <sup>™</sup> -Jurkat	R762-07

All cell lines are supplied as one vial containing  $1 \times 10^7$  frozen cells in 1 mL of Freezing Medium.



# Description of the<br/>systemIntroductionThe Flp-In™ cell lines stably express the lacZ-Zeocin™ fusion gene and are designed<br/>for use with the Flp-In™ System (Cat. nos. K6010-01 and K6010-02). Each cell line<br/>contains a single integrated Flp Recombination Target (FRT) site from pFRT/lacZeo or<br/>pFRT/lacZeo2 as confirmed by Southern blot analysis. See "Parental cell lines" on

pFRT/*lac*Zeo2 as confirmed by Southern blot analysis. See "**Parental cell lines**" on **page 6** and "Description of Flp-In<sup>™</sup> cell lines" on page 7 for information about the generation of the Flp-In<sup>™</sup> cell lines. For more information about the Flp-In<sup>™</sup> System and its components, refer to the Flp-In<sup>™</sup> System manual, visit **www.lifetechnologies.com**, or contact Technical Support (see ). The Flp-In<sup>™</sup> System manual is also available from our website.

Generation of Flp-In<sup>™</sup> expression cell lines requires cotransfection of the Flp-In<sup>™</sup> cell line with a Flp-In<sup>™</sup> expression vector containing your gene of interest and the Flp recombinase expression plasmid, pOG44 (O'Gorman et al., 1991). Flp recombinase mediates insertion of your Flp-In<sup>™</sup> expression construct into the genome at the integrated FRT site through site-specific DNA recombination (O'Gorman et al., 1991; Sauer, 1994). Stable cell lines expressing your gene of interest from the Flp-In<sup>™</sup> expression vector can be generated by selection using hygromycin B. For more information about FRT sites and Flp recombinase-mediated DNA recombination, refer to the Flp-In<sup>™</sup> System manual.

#### Parental cell lines

The table below provides a brief description of the source of the parental cell line used to generate each Flp-In<sup>™</sup> cell line. The parental cell lines were obtained from the American Type Culture Collection (ATCC<sup>®</sup>). The ATCC<sup>®</sup> number for each cell line is included. For further information about the parental cell lines, refer to the ATCC<sup>®</sup> website (www.atcc.org).


Cell Line	Characteristic	Source	ATCC <sup>®</sup> Number
293	Adherent	Human embryonic kidney (Graham et al., 1977)	CRL-1573
CV-1	Adherent	African Green Monkey kidney (Kit et al., 1965)	CCL-70
СНО <sup>™</sup> -К1	Adherent	Chinese Hamster ovary (Kao and Puck, 1968)	CCL-61
внк	Adherent	Baby hamster kidney (Talavera and Basilico, 1977)	CCL-10
NIH/3T3	Adherent	Mouse (NIH Swiss) embryonic fibroblast (Jainchill et al., 1969)	CRL-1658
Jurkat	Suspension	Human T-cell leukemia (Weiss et al., 1984)	TIB-152



## Description of Flp-In<sup>™</sup> cell lines


All of the Flp-In<sup>TM</sup> cell lines (except Flp-In<sup>TM</sup>-CHO<sup>TM</sup>; see the following section) contain a single integrated FRT site and stably express the *lacZ*-Zeocin<sup>TM</sup> fusion gene from the pFRT/*lacZ*eo plasmid under the control of the SV40 early promoter (see the following diagram). The location of the FRT site in each Flp-In<sup>TM</sup> cell line has not been mapped, but is presumed to have integrated into a transcriptionally active genomic locus as determined by generation of a Flp-In<sup>TM</sup> expression cell line containing the pcDNA<sup>TM</sup>5/FRT/CAT or pEF5/FRT/GW-CAT control plasmid. The Flp-In<sup>TM</sup> cell lines should be maintained in medium containing Zeocin<sup>TM</sup> Selection Antibiotic (see **"Media for cell lines" on page 8**).

For more information about pFRT/*lac*Zeo, refer to the Flp-In<sup>™</sup> System manual. For information about pcDNA<sup>™</sup>5/FRT/CAT or pEF5/FRT/GW-CAT, refer to the pcDNA<sup>™</sup>5/FRT or pEF5/FRT/V5-DEST<sup>™</sup> manuals, respectively.



## Flp-In<sup>™</sup>-CHO<sup>™</sup> cell line

The Flp-In<sup>TM</sup>-CHO<sup>TM</sup> cell line contains a single integrated FRT site and stably expresses the *lacZ*-Zeocin<sup>TM</sup> fusion gene from the pFRT/*lacZeo2* plasmid. Note that pFRT/*lacZeo2* contains a mutated SV40 early promoter ( $P_{sv40\Delta}$ ) which is severely abrogated in its activity. The SV40 $\Delta$  early promoter in pFRT/*lacZeo2* exhibits approximately 60-fold less activity than the wild-type SV40 early promoter in pFRT/*lacZeo.* Because of the minimal activity of the SV40 $\Delta$  promoter, we expect that stable transfectants expressing the *lacZ*-Zeocin<sup>TM</sup> gene from pFRT/*lacZeo2* should contain FRT sites which have integrated into the most transcriptionally active genomic loci. The location of the FRT site in the Flp-In<sup>TM</sup>-CHO<sup>TM</sup> cell line has not been mapped, but has been demonstrated to have integrated into a highly transcriptionally active genomic locus as determined by generation of a Flp-In<sup>TM</sup> expression cell line containing the pcDNA<sup>TM</sup>5/FRT/luc (luciferase-expressing) control plasmid. The Flp-In<sup>TM</sup>-CHO<sup>TM</sup> cell line should be maintained in medium containing Zeocin<sup>TM</sup> Selection Antibiotic (see "**Media for cell lines**" on page 8). For more information about pFRT/*lacZeo2* and the SV40 $\Delta$  early promoter, refer to the pFRT/*lacZeo2* manual.





#### Media for cell lines

The following table lists the recommended complete medium, freezing medium, and antibiotic concentration required to maintain and culture each  $Flp-In^{TM}$  cell line.

Cell Line	Complete Medium	[Antibiotic]	Freezing Medium
Flp-In <sup>™</sup> -293	D-MEM <sup>™</sup> (high glucose)	100 µg/mL Zeocin <sup>™</sup>	90% complete medium
	10% FBS <sup>[1]</sup>	Selection Antibiotic	10% DMS0
	2 mM L-glutamine		
	1% Pen-Strep (optional)		
Flp-In <sup>™</sup> -CV-1	D-MEM <sup>™</sup> (high glucose)	100 µg/mL Zeocin <sup>™</sup>	90% complete medium
	10% FBS <sup>[2]</sup>	Selection Antibiotic	10% DMS0
	2 mM L-glutamine		
	1% Pen-Strep (optional)		
Flp-In <sup>™</sup> -CHO <sup>™</sup>	Ham's F12	100 µg/mL Zeocin <sup>™</sup>	90% complete medium
	10% FBS <sup>[3]</sup>	Selection Antibiotic	10% DMS0
	2 mM L-glutamine		
	1% Pen-Strep (optional)		
Flp-In <sup>™</sup> -BHK	D-MEM <sup>™</sup> (high glucose)	100 µg/mL Zeocin <sup>™</sup>	90% complete medium
	10% FBS <sup>[4]</sup>	Selection Antibiotic	10% DMS0
	2 mM L-glutamine		
	1% Pen-Strep (optional)		
Flp-In <sup>™</sup> -3T3	D-MEM <sup>™</sup> (high glucose)	100 µg/mL Zeocin <sup>™</sup>	90% complete medium
	10% donor calf serum	Selection Antibiotic	10% DMS0
	2 mM L-glutamine		
	1% Pen-Strep (optional)		
Flp-In <sup>™</sup> -Jurkat	RPMI 1640	100 µg/mL Zeocin <sup>™</sup>	90% complete medium
	10% FBS <sup>[5]</sup>	Selection Antibiotic	10% DMS0
	2 mM L-glutamine		
	1% Pen-Strep (optional)		

<sup>[1]</sup> FBS = fetal bovine serum

<sup>[2]</sup> FBS = fetal bovine serum

[3] FBS = fetal bovine serum

[4] FBS = fetal bovine serum

<sup>[5]</sup> FBS = fetal bovine serum



#### Important guidelines

- FBS does not need to be heat inactivated for use with these cell lines.
- Cell lines should be maintained in medium containing Zeocin<sup>™</sup> Selection Antibiotic at the concentrations listed in the previous section.
- If adherent cells (e.g., Flp-In<sup>™</sup>-293, Flp-In<sup>™</sup>-CV-1, Flp-In<sup>™</sup>-CHO<sup>™</sup>, Flp-In<sup>™</sup>-3T3, Flp-In<sup>™</sup>-BHK) are split at a 1:5 to 1:10 dilution, they will generally reach 80–90% confluence in 3–4 days.
- Suspension Flp-In<sup>™</sup>-Jurkat cells will demonstrate optimal growth characteristics if maintained at a cell density between 1 × 10<sup>5</sup> cells/mL and 1 × 10<sup>6</sup> cells/mL.
- When maintaining Flp-In<sup>™</sup>-Jurkat cells in suspension culture, do not allow the medium to turn yellow; this indicates that cells have reached too high a density or that the medium is depleted of nutrients. If this occurs, either add fresh complete media to the cells or passage them.

# Methods



# Culturing Flp-In<sup>™</sup> cell lines

General cell	Follow the guidelines below to successfully grow and maintain your cells.
handling	• All solutions and equipment that come in contact with the cells must be sterile. Always use proper sterile technique and work in a laminar flow hood.
	<ul> <li>Before starting experiments, be sure to have cells established and also have some frozen stocks on hand. We recommend that you always use early-passage cells for your experiments. Upon receipt of the cells from Life Technologies<sup>™</sup>, grow and freeze multiple vials of the particular cell line to ensure that you have any adequate supply of early-passage cells.</li> </ul>
	• Cells should be at the appropriate confluence (approximately 60%) and >90% viability prior to transfection (see "Transfecting cells" on page 14).
	<ul> <li>For general maintenance of cells, pass all cell lines when they are 80–90% confluent (for adherent cells) or when they reach a density of 2–4 × 10<sup>6</sup> cells/mL (for suspension cells).</li> </ul>
	<ul> <li>Use trypan blue exclusion to determine cell viability. Log phase cultures should be &gt;90% viable.</li> </ul>
Before starting	Be sure to have the following solutions and supplies available:
5	• 15-mL sterile, conical tubes
	• 5-, 10-, and 25-mL sterile pipettes
	Cryovials
	<ul> <li>Phosphate-Buffered Saline (see "Accessory products" on page 17 for ordering information)</li> </ul>
	• 0.4% Trypan blue in PBS and hemacytometer (for counting cells)
	<ul> <li>Reagents to prepare the appropriate complete medium (see "Media for cell lines" on page 8)</li> </ul>
	• Freezing Medium (see "Media for cell lines" on page 8 and "Preparing freezing medium" on page 13)
	Table-top centrifuge
	• 75-cm <sup>2</sup> flasks, 175-cm <sup>2</sup> flasks and other appropriately-sized tissue culture flasks or plates
	Trypsin/versene (EDTA) solution or other trypsin solution

Thawing adherent cells	The following protocol is designed to help you thaw adherent cells to initiate cell culture. All cell lines are supplied in vials containing $1 \times 10^7$ cells in 1 mL of Freezing Medium.
	<ol> <li>Remove the vial of cells from the liquid nitrogen and thaw quickly in a 37℃ water bath.</li> </ol>
	<ol> <li>Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol, and transfer the cells to a T-75 flask containing 12 mL of complete medium without Zeocin<sup>™</sup> Selection Antibiotic.</li> </ol>
	<b>3.</b> Incubate the flask at $37^{\circ}$ C for 2–4 hours to allow the cells to attach to the bottom of the flask.
	<ol> <li>Aspirate off the medium and replace with 12 mL of fresh, complete medium without Zeocin<sup>™</sup> Selection Antibiotic.</li> </ol>
	<b>5.</b> Incubate cells overnight at 37℃.
	<ol> <li>The next day, aspirate off the medium and replace with fresh, complete medium containing Zeocin<sup>™</sup> Selection Antibiotic (at the recommended concentration listed in "Media for cell lines" on page 8).</li> </ol>
	<b>7.</b> Incubate the cells and check them daily until the cells are 80–90% confluent (2– 7 days).
	<b>8.</b> Proceed to "Passaging adherent cells" on page 12.
Thawing suspension cells	The following protocol is designed to help you thaw suspension cells to initiate cell culture. All cell lines are supplied in vials containing $1 \times 10^7$ cells in 1 mL of Freezing Medium.
	1. Remove the vial of cells from the liquid nitrogen and thaw quickly in a 37℃ water bath.
	<ol> <li>Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol, and transfer the cells to a sterile, conical tube containing 5 mL of complete medium without Zeocin<sup>™</sup> Selection Antibiotic.</li> </ol>
	<b>3.</b> Centrifuge for 3 minutes at $750 \times g$ at room temperature.
	<ol> <li>Aspirate off the medium and resuspend the cells in 12 mL of fresh, complete medium without Zeocin<sup>™</sup> Selection Antibiotic.</li> </ol>
	<b>5.</b> Transfer <sup><math>^{\text{TM}}</math></sup> the cells to a T-75 flask and incubate cells overnight at 37°C.
	6. The next day, add Zeocin <sup>™</sup> Selection Antibiotic to the cells (at the recommended concentration listed in "Media for cell lines" on page 8).

**7.** Incubate the cells and count them daily until the cells reach a density ranging from  $2-4 \times 10^6$  cells/mL (2–7 days).

**Note:** You may add fresh, complete medium containing  $\text{Zeocin}^{\text{TM}}$  Selection Antibiotic to the cells every few days.

8. Proceed to "Passaging suspension cells" on page 12.

#### 1. When cells are $\sim$ 80–90% confluent, remove all medium from the flask.

- **2.** Wash cells once with 10 mL PBS to remove excess medium and serum. Serum contains inhibitors of trypsin.
- **3.** Add 5 mL of trypsin/versene (EDTA) solution to the monolayer and incubate 1– 5 minutes at room temperature until cells detach. Check the cells under a microscope and confirm that most of the cells have detached. If cells are still attached, incubate a little longer until most of the cells have detached.
- 4. Add 5 mL of complete medium to stop trypsinization.
- 5. Briefly pipet the solution up and down to break up clumps of cells.
- 6. To maintain cells in 75-cm<sup>2</sup> flasks, transfer 1 mL of the 10 mL cell suspension from Step 5 on page 12 to a new 75-cm<sup>2</sup> flask and add 15 mL fresh, complete containing Zeocin<sup>™</sup> Selection Antibiotic. If you want the cells to reach confluency sooner, split the cells at a lower dilution (e.g., 1:4).

**Note:** To expand cells into  $175\text{-cm}^2$  flasks, add 28 mL of fresh, complete medium containing Zeocin<sup>TM</sup> Selection Antibiotic to each of three  $175\text{-cm}^2$  flasks, then transfer 2 mL of the cell suspension to each flask to obtain a total volume of 30 mL.

- **7.** Incubate flasks in a humidified, 37°C, 5% CO<sub>2</sub> incubator.
- **8.** Repeat "Passaging adherent cells" on page 12–"Passaging adherent cells" on page 12 as necessary to maintain or expand cells.
- 1. Passage suspension cells when they reach a density of  $2-4 \times 10^6$  cells/mL.
- To maintain cells in 75-cm<sup>2</sup> flasks, transfer 1–1.5 mL of cell suspension from Step 1 on page 12 to a new 75-cm<sup>2</sup> flask containing 13–14 mL of fresh, complete medium with Zeocin<sup>™</sup> Selection Antibiotic.

Note: You may split the cells at a lower dilution (e.g., 1:4), if desired.

**3.** To expand cells into 175-cm<sup>2</sup> flasks, add 28 mL of fresh, complete medium containing Zeocin<sup>™</sup> Selection Antibiotic to each of three 175-cm<sup>2</sup> flasks, then transfer 2 mL of the cell suspension to each flask to obtain a total volume of 30 mL.

Note: You may also expand cells into a spinner flask, if desired.

- **4.** Incubate flasks in a humidified, 37°**C**, CO<sub>2</sub> incubator.
- 5. Repeat Steps 1 on page 12–4 on page 12 as necessary to maintain or expand cells.

## Passaging adherent cells

Passaging

suspension cells

## Freezing cells

Introduction Preparing freezing medium	<ul> <li>When freezing the Flp-In<sup>™</sup> cell lines, we recommend the following:</li> <li>Freeze cells at a density of at least 3 × 10<sup>6</sup> cells/mL.</li> <li>Use a freezing medium composed of 90% complete medium and 10% DMSO. Complete medium is medium containing serum.</li> <li>Guidelines to prepare freezing medium and freeze cells are provided in this section.</li> <li>Freezing medium should be prepared fresh immediately before use.</li> <li>In a sterile, conical centrifuge tube, mix together the following reagents for every 1 mL of freezing medium needed:</li> </ul>		
	Fresh complete medium	0.9 mL	
	DMSO	0.1 mL	
	<b>2.</b> Place the tube on ice. Discard any remain	ning freezing medium after use.	
Freezing the cells	Before starting, label cryovials and prepare from <b>medium</b> " on page 13). Keep the freezing med		
	1. To collect cells, perform the following:	o collect cells, perform the following:	
	<ul> <li>For adherent cells, follow Steps 1 on adherent cells" on page 12.</li> </ul>	page 13–5 on page 13 of <b>"Passaging</b>	
	• For suspension cells, transfer cells to	a sterile, conical centrifuge tube.	
	<b>2.</b> Count the cells.		
	<b>3.</b> Pellet cells at $250 \times g$ for 5 minutes in a ta and carefully aspirate off the medium.	ble top centrifuge at room temperature	
	<b>4.</b> Resuspend the cells at a density of <b>at lea</b> medium.	st $3 \times 10^6$ cells/mL in chilled freezing	
	<b>5.</b> Place vials in a microcentrifuge rack and each cryovial.	aliquot 1 mL of the cell suspension into	
	<ol> <li>Freeze cells in an automated or manual, following standard procedures. For idea should be a decrease of 1℃ per minute.</li> </ol>	0 11	
	<b>7.</b> Transfer <sup><math>^{\text{TM}} vials to liquid nitrogen for long</math></sup>	g-term storage.	
	<b>Note:</b> You may check the viability and restoring cryovials in liquid nitrogen by fo <b>"Thawing adherent cells" on page 11</b> or <b>page 11</b> , as appropriate.	llowing the procedure outlined in	



## Transfecting cells

Introduction	To generate stable Flp-In <sup>™</sup> expression cell lines, you will cotransfect your Flp-In <sup>™</sup> expression construct and the pOG44 plasmid into the Flp-In <sup>™</sup> cell line and select for stable transfectants using hygromycin B. General guidelines and recommendations for transfection are provided in this section. We recommend that you read through this section before beginning.
Transfection methods	<ul> <li>The Flp-In<sup>™</sup> cell lines are generally amenable to transfection using standard methods including calcium phosphate precipitation (Chen and Okayama, 1987; Wigler et al., 1977), lipid-mediated transfection (Felgner et al., 1989; Felgner and Ringold, 1989), and electroporation (Chu et al., 1987; Shigekawa and Dower, 1988). Although other transfection reagents may be suitable, we recommended using Lipofectamine<sup>™</sup> 2000 Reagent (see "Additional reagents" on page 17 for ordering information) to introduce Flp-In<sup>™</sup> expression constructs into the following Flp-In<sup>™</sup> cell lines:</li> <li>Flp-In<sup>™</sup>-293</li> <li>Flp-In<sup>™</sup>-CV-1</li> <li>Flp-In<sup>™</sup>-CHO<sup>™</sup></li> <li>Flp-In<sup>™</sup>-Jurkat</li> <li>Flp-In<sup>™</sup>-BHK</li> </ul>
Important	We have observed down-regulation of the viral CMV promoter and subsequent loss of gene expression when pcDNA <sup>™</sup> 5/FRT-based expression constructs are introduced into Flp-In <sup>™</sup> -3T3 or Flp-In <sup>™</sup> -BHK cells. This behavior is not observed with pEF5/FRT- based expression constructs. If you are generating Flp-In <sup>™</sup> expression cell lines using the Flp-In <sup>™</sup> -3T3 or Flp-In <sup>™</sup> -BHK cell line, we recommend that you express your gene

of interest from a pEF5/FRT-based plasmid (e.g., pEF5/FRT/V5-DEST<sup>™</sup>).

### Generating stable expression cell lines

To generate Flp-In<sup>™</sup> expression cell lines, cotransfect your Flp-In<sup>™</sup> expression construct and the pOG44 plasmid into the Flp-In<sup>™</sup> cell line of choice, and select for stable transfectants using hygromycin B. Before transfection, you may want to test the sensitivity of the Flp-In<sup>™</sup> cell line to hygromycin B to more accurately determine the hygromycin B concentration to use for selection. A suggested range of hygromycin B concentrations to use for selection of your Flp-In<sup>™</sup> expression vector is listed below. For more information, refer to the Flp-In<sup>™</sup> System manual. Hygromycin B may be obtained from Life Technologies<sup>™</sup> (see "Additional reagents" on page 17 for ordering information).

**IMPORTANT!** Following cotransfection, your Flp-In<sup>TM</sup> expression clones should become sensitive to Zeocin<sup>TM</sup> Selection Antibiotic; therefore, your selection medium should NOT contain Zeocin<sup>TM</sup> Selection Antibiotic.

Cell Line (After Transfection with Flp-In <sup>™</sup> Expression Vector)	Estimated Hygromycin B Concentration (µg/mL)
Flp-In <sup>™</sup> -293	100-200
Flp-In <sup>™</sup> -CV-1	100-200
Flp-In <sup>™</sup> -CHO <sup>™</sup>	500-600
Flp-In <sup>™</sup> -BHK	100–200
Flp-In <sup>™</sup> -3T3	100–200
Flp-In <sup>™</sup> -Jurkat	200-400

Note

When transfecting Flp-In<sup>™</sup>-CHO<sup>™</sup> cells, we recommend following these guidelines:

- 48 hours after transfection, split the cells directly into medium containing the appropriate concentration of hygromycin B.
- Split the cells such that they are no more than 25% confluent. If the cells are too dense, the antibiotic will not kill the cells. Antibiotics work best on actively dividing cells.

Polyclonal selection of isogenic cell lines

Because every Flp-In<sup>™</sup> cell line contains a single integrated FRT site, all of the hygromycin-resistant foci that you obtain after cotransfection with the Flp-In<sup>™</sup> expression vector and pOG44 should be isogenic (i.e., the Flp-In<sup>™</sup> expression vector should integrate into the same genomic locus in every clone; therefore, all clones should be identical). To obtain stable expression cell lines, you may perform polyclonal selection and screening of your hygromycin-resistant cells. After hygromycin B selection, simply pool the hygromycin-resistant foci and screen the entire population of cells for the following phenotypes:

- Zeocin<sup>™</sup> Selection Antibiotic sensitivity
- Lack of β-galactosidase activity
- Expression of the gene of interest



Select individual cell lines	If desired, single hygromycin-resistant foci can be isolated and expanded to generate individual clonal cell lines. To isolate individual clones, simply pick 5–20 hygromycin- resistant foci and expand the cells. You may verify that your Flp-In <sup>TM</sup> expression construct has integrated into the FRT site by testing each clone for Zeocin <sup>TM</sup> Selection Antibiotic sensitivity and lack of $\beta$ -galactosidase activity. Select those clones that are hygromycin-resistant, Zeocin <sup>TM</sup> Selection Antibiotic-sensitive, and lack $\beta$ -galactosidase activity, and assay for expression of your gene of interest.
Note	Note that in rare instances, it is possible to generate a Flp-In <sup>TM</sup> expression cell line in which the Flp-In <sup>TM</sup> expression plasmid has undergone both Flp recombinase-mediated integration into the FRT site and random integration into a second genomic site. In this case, clones will still exhibit hygromycin B resistance. To test for these second site integrants, transfect the cells with the pOG44 plasmid and select for Zeocin <sup>TM</sup> Selection Antibiotic resistance. The Flp recombinase should mediate excision of the Flp-In <sup>TM</sup> expression plasmid at the FRT site and restore the <i>lacZ</i> -Zeocin <sup>TM</sup> fusion gene. The resulting cells should exhibit $\beta$ -galactosidase activity, Zeocin <sup>TM</sup> Selection Antibiotic resistance, and continued expression of the gene of interest. Alternatively, you may perform Southern blot analysis to identify second site integrants if suitable restriction enzymes are selected.



# **Ordering information**

## Accessory products

# **Introduction** The products listed in this section are intended for use with the Flp-In<sup>™</sup> Cell Lines and the Flp-In<sup>™</sup> System. For more information, refer to **www.lifetechnologies.com** or contact Technical Support (see ).

Cell cultureA large variety of Gibco<sup>™</sup> cell culture products are available from Life Technologies<sup>™</sup><br/>to facilitate growth and maintenance of the Flp-In<sup>™</sup> cell lines. For more information<br/>about the products listed below, refer to www.lifetechnologies.com or contact<br/>Technical Support (see ). Note: Reagents are available in other sizes.

Item	Amount	Catalog no.
Dulbecco's Modified Eagle Medium (D-MEM <sup>™</sup> )	500 mL	11960-044
Ham's F-12	500 mL	11765-054
RPMI Medium 1640	500 mL	11875-093
Fetal Bovine Serum	500 mL	16000-044
Donor Calf Serum	500 mL	16030-074
L-Glutamine, 200 mM (100X)	100 mL	25030-081
Penicillin-Streptomycin	100 mL	15070-063
0.05% Trypsin-EDTA (1X), Phenol Red	100 mL	25300-054

#### Additional reagents

The products listed below may be used with the Flp-In<sup>TM</sup> Cell Lines. Zeocin<sup>TM</sup> Selection Antibiotic is available for maintenance and growth of the Flp-In<sup>TM</sup> cell lines. For more information, refer to **www.lifetechnologies.com** or contact Technical Support (see ).

Item	Amount	Catalog no.
Zeocin <sup>™</sup> Selection Antibiotic	8 × 1.25 mL	R250-01
	50 mL	R250-05
Hygromycin B, 50 mg/mL	20 mL	10687-010
Lipofectamine <sup>™</sup> 2000 Reagent	0.75 mL	11668-027
	1.5 mL	11668-019
Lipofectamine <sup>™</sup> Reagent	1 mL	18324-012



ltem	Amount	Catalog no.
PLUS <sup>™</sup> Reagent	0.85 mL	11514-015
Phosphate-Buffered Saline (PBS), pH 7.4	500 mL	10010-023

## Flp-In<sup>™</sup> products

The plasmids required to generate Flp-In<sup>™</sup> host cell lines and expression cell lines are available separately from Life Technologies<sup>™</sup>. For more information about the features of each vector, refer to **www.lifetechnologies.com** or contact Technical Support (see ). Ordering information is provided below.

Product	Amount	Catalog no.
pFRT/ <i>lac</i> Zeo	20 µg	V6015-20
pFRT/ <i>lac</i> Zeo2	20 µg	V6022-20
pOG44	20 µg	V6005-20
pcDNA <sup>™</sup> 5/FRT	20 µg	V6010-20
pcDNA <sup>™</sup> 5/FRT/V5-His TOPO <sup>™</sup> TA Expression Kit	20 reactions	K6020-01
pSecTag/FRT/V5-His TOPO <sup>™</sup> TA Expression Kit	20 reactions	K6025-01
pEF5/FRT/V5-DEST <sup>™</sup> Gateway <sup>™</sup> Vector Pack	6 µg	V6020-20

## References

Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. Molec. Cell. Biol. 7, 2745-2752.

Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. Nucleic Acids Res. *15*, 1311-1326.

Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. Proc. West. Pharmacol. Soc. *32*, 115-121.

Felgner, P. L. a., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. Nature 337, 387-388.

Graham, F. L., Smiley, J., Russell, W. C., and Nairn, R. (1977). Characteristics of a Human Cell Line Transformed by DNA from Human Adenovirus Type 5. J. Gen. Virol. *36*, 59-74.

Jainchill, J. L., Aaronson, S. A., and Todaro, G. J. (1969). Murine Sarcoma and Leukemia Viruses: Assay Using Clonal Lines of Contact-Inhibited Mouse Cells. J. Virol. 4, 549-553.

Kao, F. T., and Puck, T. T. (1968). Genetics of Somatic Mammalian Cells, VII. Induction and Isolation of Nutritional Mutants in Chinese Hamster Cells. Proc. Natl. Acad. Sci. USA *60*, 1275-1281.



Kit, S., Dubbs, D. R., DeTorres, R. A., and Melnick, J. L. (1965). Enhanced Thymidine Kinase Activity Following Infection of Green Monkey Kidney Cells by Simian Adenoviruses, Simian Papovavirus SV40, and an Adenovirus-SV40 Hybrid. Virology 27, 453-457.

O'Gorman, S., Fox, D. T., and Wahl, G. M. (1991). Recombinase-Mediated Gene Activation and Site-Specific Integration in Mammalian Cells. Science 251, 1351-1355.

Sauer, B. (1994). Site-Specific Recombination: Developments and Applications. Curr. Opin. Biotechnol. *5*, 521-527.

Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. BioTechniques *6*, 742-751.

Talavera, A., and Basilico, C. (1977). Temperature Sensitive Mutants of BHK Cells Affected in Cell Cycle Progression. J. Cell. Physiol. *92*, 425-436.

Weiss, A., Wiskocil, R. L., and Stobo, J. D. (1984). The Role of T3 Surface Molecules in the Activation of Human T Cells: A Two-stimulus Requirement for IL2 Production Reflects Events Occurring at a Pre-translational Level. J. Immunol. *133*, 123-128.

Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer<sup>™</sup> of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. Cell *11*, 223-232.

# Safety





**WARNING!** GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, and so on). To obtain SDSs, see the "Documentation and Support" section in this document.

## **Chemical safety**



**WARNING!** GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below. Consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the "Documentation and Support" section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with sufficient ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if needed) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

**WARNING!** HAZARDOUS WASTE (from instruments). Waste produced by the instrument is potentially hazardous. Follow the guidelines noted in the preceding General Chemical Handling warning.



**WARNING!** 4L Reagent and Waste Bottle Safety. Four-liter reagent and waste bottles can crack and leak. Each 4-liter bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position.



## **Biological hazard safety**

**WARNING!** Potential Biohazard. Depending on the samples used on this instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.



**WARNING! BIOHAZARD.** Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Conduct all work in properly equipped facilities with the appropriate safety equipment (for example, physical containment devices). Safety equipment can also include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

• U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 5th Edition, HHS Publication No. (CDC) 21-1112, Revised December 2009; found at:

https://www.cdc.gov/labs/pdf/

CDC-BiosafetymicrobiologicalBiomedicalLaboratories-2009-P.pdf

World Health Organization, *Laboratory Biosafety Manual*, 3rd Edition, WHO/CDS/CSR/LYO/2004.11; found at:

www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf

# **Documentation and support**

## **Customer and technical support**

Visit **thermofisher.com/support** for the latest service and support information.

- Worldwide contact telephone numbers
- Product support information
  - Product FAQs
  - Software, patches, and updates
  - Training for many applications and instruments
- Order and web support
- Product documentation
  - User guides, manuals, and protocols
  - Certificates of Analysis
  - Safety Data Sheets (SDSs; also known as MSDSs)

**Note:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at **www.thermofisher.com/us/en/home/global/terms-and-conditions.html**. If you have any questions, please contact Life Technologies at **www.thermofisher.com/support**.

