

Iodoacetamide, Single-Use

MAN0011661

Rev. B.0

Pub. Part No. 2162099

A39271

Number	Description
A39271	Iodoacetamide, Single-Use , 30 × 9.3 mg Synonyms: 2-iodoacetamide, IAM Molecular Weight: 184.96 Formula: C ₂ H ₄ INO CAS: 144-48-9 Melting Point: 94°C Appearance: White solid in amber, screw capped vials

Storage: Upon receipt store product at room temperature.

Note: Product labels have been provided for your convenience. Please label the vials using one of the labels provided in the Al foil pouch to avoid any confusion as you work with this No-Weigh reagent.

Introduction

Iodoacetamide is a sulfhydryl-reactive alkylating reagent used to block reduced cysteine residues for protein characterization and peptide mapping. Alkylation with iodoacetamide after cystine reduction results in the covalent addition of a carbamidomethyl group (57.07 Da) and prevents the formation of disulfide bonds. Reducing agents added after alkylation will react with excess iodoacetamide.

Alkylation Procedure

Note: Iodoacetamide is unstable and light-sensitive. Prepare solutions immediately before use and perform alkylation in the dark. If iodoacetamide is present in limiting quantities and a slightly alkaline pH, cysteine modification will be the exclusive reaction. Excess iodoacetamide or non-buffered iodoacetamide reagent can also alkylate amines (lysine, N-termini), thioethers (methionine), imidazoles (histidine) and carboxylates (aspartate, glutamate).

1. Add 5 µl of 2% SDS and 45 µl of 200 mM ammonium bicarbonate (pH 8.0) to 20-100 µg of protein sample. Adjust volume to 100 µl with ultrapure water.
2. Add 5 µl of 200 mM Tris(2-carboxyethyl) phosphine hydrochloride (TCEP•HCl, Product No. 20490) and incubate sample at 55°C for 1 hour.
3. Immediately before use, dissolve one tube of iodoacetamide (9.3 mg) with 132 µl of 200 mM ammonium bicarbonate (pH 8.0) to make 375 mM iodoacetamide. Protect solution from light. **The maximum useable volume of the vial is 500µL.**
4. Add 5 µl of the 375 mM iodoacetamide to the sample and incubate for 30 minutes protected from light.
5. Proceed to proteolytic digestion before MS analysis or other processing.

Troubleshooting

Problem	Possible Cause	Solution
Sulfhydryls not blocked	Iodoacetamide hydrolysis	Make iodoacetamide solutions immediately before each use and dispose of excess reconstituted reagent
Sulfhydryls partially blocked	Insufficient iodoacetamide used	Use at least a 10-fold excess of iodoacetamide to sulfhydryls
	Incorrect reaction buffer	Avoid buffers that contain sulfhydryls or that are not at a slightly alkaline pH
	Insufficient reaction time	Allow reaction to proceed for 30 minutes at room temperature
Amines or other functional groups labeled	Incorrect pH of the reaction buffer	Maintain the reaction buffer pH at 7.5-8.0
	Excess reagent or time	Reduce the amount of reagent or incubation time
		Acetone precipitate the alkylated sample to limit exposure during digestion

Related Thermo Scientific Products

A39255	DTT, No-Weigh™ Format, 48 × 7.7mg
20290	DTT, 5g
20490	TCEP•HCl [Tris(2-carboxyethyl) phosphine hydrochloride], 1g
23011	Methyl Methanethiosulfonate (MMTS), 200mg
23030	N-Ethylmaleimide (NEM), 25g
77720	Bond-Breaker™ TCEP Solution, Neutral pH, 5mL
77712	Immobilized TCEP Disulfide Reducing Gel, 5 mL
77701	Immobilized Reductant Column, 1 × 2 mL
90051	Lys-C Endoproteinase, MS Grade, 20 µg
90053	Asp-N Endoproteinase, MS Grade, 2 µg
90054	Glu-C Endoproteinase, MS Grade, 5 × 10 µg
90055	Trypsin Endoproteinase, modified, TPCK treated, MS Grade, 5 × 20 µg
90056	Chymotrypsin Endoproteinase, TLCK treated, MS Grade, 4 × 25 µg

General References

- Gurd, F.R.N. (1967). Carboxymethylation. *In* Methods in Enzymology, (C.H.W. Hirs, ed.), Vol. 11, p 532. Academic Press, New York.
- Vithayathil, P.J. and Richards, F.M. (1960). Modification of the methionine residue in the peptide component of ribonuclease S. *J Biol Chem* **235**:2343-51.
- Cole, R.D, *et al.* (1958). On the cysteine content of hemoglobin. *J Biol Chem* **233**:1359-63.

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2018 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.