

TMT Mass Tagging Kits and Reagents

Pub. No. MAN0011639
Pub. Part No. 2162073.14
Rev B.0

90060-90068

Number	Description
90060	TMTduplex Isotopic Label Reagent Set , sufficient reagents for 5 duplex isotopic experiments Contents: TMT⁰ Label Reagent , 5 × 0.8mg TMT⁶-127 Label Reagent , 5 × 0.8mg
90061	TMTsixplex Isobaric Label Reagent Set , sufficient reagents for 1 sixplex isobaric experiment Contents: TMT⁶-126 Label Reagent , 1 × 0.8mg TMT⁶-127 Label Reagent , 1 × 0.8mg TMT⁶-128 Label Reagent , 1 × 0.8mg TMT⁶-129 Label Reagent , 1 × 0.8mg TMT⁶-130 Label Reagent , 1 × 0.8mg TMT⁶-131 Label Reagent , 1 × 0.8mg
90062	TMTsixplex Isobaric Label Reagent Set , sufficient reagents for 2 sixplex isobaric experiments Contents: TMT⁶-126 Label Reagent , 2 × 0.8mg TMT⁶-127 Label Reagent , 2 × 0.8mg TMT⁶-128 Label Reagent , 2 × 0.8mg TMT⁶-129 Label Reagent , 2 × 0.8mg TMT⁶-130 Label Reagent , 2 × 0.8mg TMT⁶-131 Label Reagent , 2 × 0.8mg
90063	TMTduplex Isobaric Mass Tagging Kit , sufficient reagents for 5 duplex isobaric experiments Contents: TMT⁰ Label Reagent , 5 × 0.8mg TMT²-126 Label Reagent , 5 × 0.8mg TMT²-127 Label Reagent , 5 × 0.8mg Dissolution Buffer (1 M triethyl ammonium bicarbonate), 5mL Denaturing Reagent (10% SDS), 1mL Reducing Reagent (0.5M TCEP), 1mL Iodoacetamide , 12 × 9mg Quenching Reagent (50% hydroxylamine), 1mL Pierce™ Trypsin Protease, MS Grade , 2 × 20μg Trypsin Storage Solution , 250μL Albumin, Bovine , 2.5mg

90064 TMTsixplex Isobaric Mass Tagging Kit, sufficient reagents for 5 sixplex isobaric experiments

Contents:

TMT⁰ Label Reagent, 5 × 0.8mg
TMT⁶-126 Label Reagent, 5 × 0.8mg
TMT⁶-127 Label Reagent, 5 × 0.8mg
TMT⁶-128 Label Reagent, 5 × 0.8mg
TMT⁶-129 Label Reagent, 5 × 0.8mg
TMT⁶-130 Label Reagent, 5 × 0.8mg
TMT⁶-131 Label Reagent, 5 × 0.8mg
Dissolution Buffer (1M triethyl ammonium bicarbonate), 5mL
Denaturing Reagent (10% SDS), 1mL
Reducing Reagent (0.5 M TCEP), 1mL
Iodoacetamide, 12 × 9mg
Quenching Reagent (50% hydroxylamine), 1mL
Pierce Trypsin Protease, MS Grade, 5 × 20µg
Trypsin Storage Solution, 250µL
Albumin, Bovine, 2.5mg

90065 TMTduplex Isobaric Label Reagent Set, sufficient reagents for 5 duplex isobaric experiments

Contents:

TMT²-126 Label Reagent, 5 × 0.8mg
TMT²-127 Label Reagent, 5 × 0.8mg

90066 TMTsixplex Label Reagent Set, sufficient reagents for 5 sixplex isobaric experiments

Contents:

TMT⁶-126 Label Reagent, 5 × 0.8mg
TMT⁶-127 Label Reagent, 5 × 0.8mg
TMT⁶-128 Label Reagent, 5 × 0.8mg
TMT⁶-129 Label Reagent, 5 × 0.8mg
TMT⁶-130 Label Reagent, 5 × 0.8mg
TMT⁶-131 Label Reagent, 5 × 0.8mg

90067 TMTzero Label Reagent, 5 × 0.8mg, sufficient reagents for 5 samples

90068 TMTsixplex Label Reagent Set, sufficient reagents for 12 sixplex isobaric experiments

Contents:

TMT⁶-126 Label Reagent, 2 × 5mg
TMT⁶-127 Label Reagent, 2 × 5mg
TMT⁶-128 Label Reagent, 2 × 5mg
TMT⁶-129 Label Reagent, 2 × 5mg
TMT⁶-130 Label Reagent, 2 × 5mg
TMT⁶-131 Label Reagent, 2 × 5mg

Storage: Upon receipt store at -20°C. Reagents are shipped with dry ice.

Note: These products are for research use only – do not use for diagnostic procedures.

Contents

Introduction	3
Procedure Summary.....	4
Important Product Information	4
Additional Materials Required.....	4
Material Preparation	5
Preparing and Labeling Peptides with the TMT Isobaric Mass Tags	5
Troubleshooting.....	6
Additional Information	6
A. Data Acquisition Methods	6
B. Data Analysis and Quantitation	7
C. Information Available from our Website.....	8
Related Thermo Scientific Products	8
General References.....	8

Introduction

The Thermo Scientific™ TMT™ Isobaric Mass Tagging Kits and Reagents enable multiplex relative quantitation by mass spectrometry (MS). Each mass-tagging reagent within a set has the same nominal mass (i.e., isobaric) and chemical structure composed of an amine-reactive NHS-ester group, a spacer arm and an MS/MS reporter (Figure 1). The reagent sets can be used to label two or six peptide samples prepared from cells or tissues. For each sample, a unique reporter in the low mass region of the MS/MS spectrum (i.e., 126-127Da for TMT² and 126-131Da for TMT⁶ Isobaric Label Reagents) is used to measure relative protein expression levels during peptide fragmentation.

The TMTduplex™ Isotopic Label Reagent Set contains TMTzero™ and one of the TMTsixplex™ Reagents (TMT⁶-127) to be used as “light” and “heavy” tags for MS-level peptide quantitation similar to duplex isotopic metabolic labeling (e.g., SILAC) or isotopic dimethylation labeling. These isotopic pairs can also be used in targeted quantitation strategies, including selective reaction monitoring (SRM, see the Additional Information Section). Advantages of the TMTduplex and TMTsixplex Isobaric Label Reagents include increased sample multiplexing for relative quantitation, increased sample throughput and fewer missing quantitative channels among samples.




Figure 1. Chemical structure of the TMT Label Reagents. **A.** Functional regions of the reagent structure, including MS/MS fragmentation sites by higher energy collision dissociation (HCD) and electron transfer dissociation (ETD). **B.** TMTduplex Reagent structures and isotope positions (*); only HCD differentiates between these two reporters. **C.** TMTsixplex Reagent structures and isotope positions (*).

Procedure Summary

Protein extracts isolated from cells or tissues are reduced, alkylated and digested overnight. Samples are labeled with the TMT Reagents and then mixed before sample fractionation and clean-up. Labeled samples are analyzed by high resolution Orbitrap LC-MS/MS before data analysis to identify peptides and quantify reporter ion relative abundance (Figure 2).




Figure 2. Schematic for using the Thermo Scientific TMTsixplex Isobaric Mass Tagging Reagents.

Important Product Information

- The TMT Reagents are moisture-sensitive. To avoid moisture condensation onto the product, vial must be equilibrated to room temperature before opening.
- Anhydrous acetonitrile is the recommended solvent to dissolve reagents. Stock solutions are stable for one week when stored at -20°C . For long term storage of unused reagent, remove all solvent by drying and store with desiccant at -20°C . Anhydrous ethanol can be used as an alternative solvent to dissolve reagents but is not recommended for stock solution storage.
- The TMT Reagents are amine-reactive and modify lysine residues and the peptide N-termini. All amine-containing buffers and additives must be removed before digestion and labeling.
- All samples must be digested, labeled and then mixed equally before desalting, fractionation and LC-MS/MS. For optimal results, use 25-100 μg of peptide for each labeling reaction.
- To avoid contamination of MS samples, always wear gloves when handling samples and gels. Use ultrapure MS-grade reagents. Perform sample preparation in a cleaned work area.
- The TMTzero Label Reagent can be used to optimize methods before multiplexed analysis of samples with the TMTduplex or TMTsixplex Reagent Set.

Additional Materials Required

- Microcentrifuge tubes
- Anhydrous acetonitrile (Thermo Scientific™ Acetonitrile HPLC grade, Product No. 51101)
- Water, LC-MS Grade (Product No. 51140)
- Chilled (-20°C) acetone
- Protein assay (e.g., Thermo Scientific™ BCA Protein Assay Kit, Product No. 22235)
- 75-300 μm capillary C_{18} reversed-phase column
- High-resolution Orbitrap Mass Spectrometer, ion trap or time-of-flight (TOF) mass spectrometer with online or offline liquid chromatography (LC) system
- Data analysis software such as Thermo Scientific™ Proteome Discoverer™ or Mascot™ Software (Matrix Science, Ltd.)
- Optional: C_{18} spin tips or columns (e.g., Thermo Scientific™ Pierce™ C_{18} Spin Columns, Product No. 89870 or Pierce™ C_{18} Tips, Product No. 87784)

Material Preparation

Note: The 50% hydroxylamine and 10% SDS stock solutions provided with the kit may precipitate during storage. Warm both solutions to room temperature and vortex before use. The amounts listed below are sufficient for preparing and labeling 6 samples.

100mM TEAB (triethyl ammonium bicarbonate)	Add 500µL of the Dissolution Buffer (1M TEAB) to 4.5mL of ultrapure water.
Lysis Buffer	Add 200µL of the Denaturing Reagent (10% SDS) to 1.8mL of 100mM TEAB.
200mM TCEP	Add 70µL of the Reducing Reagent (0.5M TCEP) to 70µL of ultrapure water. Then add 35µL of the Dissolution Buffer (1M TEAB).
5% Hydroxylamine	Add 50µL of the Quenching Reagent (50% hydroxylamine) to 450µL of 100mM TEAB.

Preparing and Labeling Peptides with the TMT Isobaric Mass Tags

Note: BSA can be used as a control sample for method optimization. Dissolve BSA to 1mg/mL using 100mM TEAB. Use 25-100µg of protein per labeling reaction. The Thermo Scientific™ Pierce™ Mass Spec Sample Prep Kit for Cultured Cells can also be used to prepare peptide digests for TMT reagent labeling.

A. Preparing Whole Cell Protein Extracts

- Culture cells to harvest at least 100µg of protein per condition. For best results, culture a minimum of 2×10^6 cells.
Note: Rinse cells 2-3 times with 1X PBS to remove cell culture media. Pellet cells using low-speed centrifugation (i.e., $< 1000 \times g$) to prevent premature cell lysis.
- Lyse the cells by adding five cell-pellet volumes of Lysis Buffer (i.e., 100µL of Lysis Buffer for a 20µL cell pellet).
Note: Lysis buffers such as 8M urea (Product No. 29700) in 50mM TEAB or HEPES buffer, pH 8 may be used as alternative denaturing cell lysis buffers. For urea-based lysis buffer, protein samples must be diluted to $< 1M$ urea before digestion, and the final C18 desalting step (C.6) is not optional. Addition of protease and/or phosphatase inhibitors during lysis is optional and may interfere with MS analysis.
Note: Depending on the Lysis Buffer used it may be necessary to reduce sample viscosity by shearing DNA using a microtip sonicator or addition of a nuclease (e.g., Thermo Scientific™ Pierce™ Universal Nuclease for Cell Lysis, Product No. 88700).
- Centrifuge lysate at $16,000 \times g$ for 10 minutes at 4°C.
- Carefully separate the supernatant and transfer into a new tube.
- Determine the protein concentration of the supernatant using established methods such as the BCA Protein Assay Kit (Product No. 23227).
Note: Use samples at $\geq 2\text{mg/mL}$. Less concentrated samples may be used; however, it might be necessary to use larger volumes of reducing/alkylating reagents.
- Transfer 100µg per condition (two for the TMTduplex or six for the TMTsixplex Label Reagents) into a new microcentrifuge tube and adjust to a final volume of 100µL with 100mM TEAB.
- Add 5µL of the 200mM TCEP and incubate sample at 55°C for 1 hour.
- Immediately before use, dissolve one tube of iodoacetamide (9mg) with 132µL of 100mM TEAB to make 375mM iodoacetamide. Protect solution from light.
- Add 5µL of the 375mM iodoacetamide to the sample and incubate for 30 minutes protected from light at room temperature.
- Add six volumes (~600µL) of pre-chilled (-20°C) acetone and freeze at -20°C. Allow the precipitation to proceed for at least 4 hours up to overnight.
Note: Methanol/chloroform is the recommended solvent for precipitation of proteins derived from tissue extracts.
- Centrifuge the samples at $8000 \times g$ for 10 minutes at 4°C. Carefully invert the tubes to decant the acetone without disturbing the white pellet. Allow the pellet to dry for 2-3 minutes.

B. Protein Digestion

- Resuspend 100µg of acetone-precipitated (or lyophilized) protein pellets with 100µL of 50mM TEAB.
Note: An acetone-precipitated pellet might not completely dissolve; however, after proteolysis at 37°C, all the protein (peptides) will be solubilized.
- Immediately before use, add 20µL of the Trypsin Storage Solution to the bottom of the trypsin glass vial and incubate for 5 minutes. Store any remaining reagent in single-use volumes at -80°C (e.g., 2.5µg of trypsin per 100µg of protein).
- Add 2.5µL of trypsin (i.e., 2.5µg) per 100µg of protein. Digest the sample overnight at 37°C.

C. Peptide Labeling

- Immediately before use, equilibrate the TMT Label Reagents to room temperature. For the 0.8mg vials, add 41µL of anhydrous acetonitrile to each tube. For the 5mg vials, add 256µL of solvent to each tube. Allow the reagent to dissolve for 5 minutes with occasional vortexing. Briefly centrifuge the tube to gather the solution.
Note: Reagents dissolved in anhydrous acetonitrile are stable for one week when stored at -20°C. Anhydrous ethanol can be used as an alternative solvent to dissolve reagents but is not recommended for stock solution storage.
- Optional: Measure protein digest concentration using Thermo Scientific™ Pierce™ Quantitative Fluorescent Peptide Assay (Product No. 23290) or Thermo Scientific™ Pierce™ Quantitative Colorimetric Peptide Assay (Product No. 23275).
- Carefully add 41µL of the TMT Label Reagent to each 100µL sample (25-100µg protein digest). Alternatively, transfer the reduced and alkylated protein digest to the TMT Reagent vial.
- Note:** Labeling more than 100µg of protein digest per reaction requires additional TMT Label Reagent.
- Incubate the reaction for 1 hour at room temperature.
- Add 8µL of 5% hydroxylamine to the sample and incubate for 15 minutes to quench the reaction.
- Combine samples at equal amounts in new microcentrifuge tube and store at -80°C.
Note: TMT-labeled peptide concentration can be measured using Thermo Scientific™ Pierce™ Quantitative Colorimetric Peptide Assay. The Thermo Scientific™ Pierce™ Quantitative Fluorescent Peptide Assay cannot be used to measure TMT-labeled peptide concentrations.
- Optional: Clean-up samples with C18 spin tips (Product No. 87784) or columns (Product No. 89870) before LC-MS analysis. Peptide clean up is recommended before LC-MS analysis but is not required. Fractionation of labeled peptides using Thermo Scientific™ Pierce™ High pH Reversed-Phase Peptide Fractionation Kit (Product No. 84868) is recommended before LC-MS analysis to increase the number of peptide identifications.

Troubleshooting

Problem	Possible Cause	Solution
Poor labeling	An amine-based buffer was used	Use a non-amine-based buffer
	Incorrect buffer pH	Make sure the buffer pH is ~8.0
	Too much sample was used	Label 25-100µg per sample
Protein precipitation	Lack of detergent present	Add detergent, such as 0.05% SDS to the preparation
	pH decreased	Make sure the pH is > 7.5

Additional Information

A. Data Acquisition Methods

Quantitation of peptides labeled with Thermo Scientific™ Tandem Mass Tag™ Reagents requires a mass spectrometer capable of MS/MS fragmentation, such as an ion trap, quadrupole time of flight, time of flight-time of flight (TOF-TOF) or triple quadrupole instrument. Higher energy collision dissociation (HCD) is recommended for TMT reporter ion fragmentation. Optimal HCD fragmentation energy is instrument-dependent and can be optimized using TMTzero Reagents.

Electron transfer dissociation (ETD) may be used as an alternative fragmentation method for peptide identification and quantitation. The choice of MS/MS fragmentation method(s) depends on the instrument capabilities such as collisionally induced dissociation (CID), pulsed-Q dissociation (PQD), higher energy collisional dissociation (HCD), or electron transfer dissociation (ETD). TMT Reagent reporter ions are not visible in ion traps following traditional CID fragmentation.

Table 1. Instruments and MS/MS fragmentation options for peptide identification and quantitation with Thermo Scientific TMT Reagents.

<u>Instrument</u>	<u>Fragmentation Method</u>	<u>Reference(s)</u>
Thermo Scientific Orbitrap™ Fusion™ Tribrid™ Mass Spectrometer	HCD/SPS-MS3	McAllister, G.C., <i>et al.</i> (2014), Viner, <i>et al.</i> (2013)
Thermo Scientific Orbitrap Elite™ Mass Spectrometer	HCD/MS3	McAllister, G.C., <i>et al.</i> (2012), Viner, <i>et al.</i> (2012)
Thermo Scientific Q Exactive™ Mass Spectrometer	HCD/MS2	Wühr, <i>et al.</i> (2012)
Thermo Scientific Orbitrap Velos Pro™, LTQ-Orbitrap™ XL, or MALDI-Orbitrap™ XL Mass Spectrometer	HCD/MS2	Ting, <i>et al.</i> (2011), Wenger, <i>et al.</i> (2011), Schirle, <i>et al.</i> (2012), Lee, <i>et al.</i> (2011), Xiong, <i>et al.</i> (2011), Strupat, <i>et al.</i> (2008)
Thermo Scientific™ Velos Pro™ ion trap	Trap HCD/MS2	Biringer, <i>et al.</i> (2011)
Thermo Scientific Orbitrap Elite ETD, Velos Pro ETD, LTQ-OrbitrapXL ETD	HCD/MS2 or ETD/MS2	Viner, <i>et al.</i> (2009)
Q-TOF	CID	Van Ulsen, <i>et al.</i> (2009)
TOF-TOF	CID	Dayon, <i>et al.</i> (2008)
Triple Quadrupole	CID/SRM	Stella, <i>et al.</i> (2011), Byers, <i>et al.</i> (2009)

B. Data Analysis and Quantitation

The masses for peptide modification by the TMT zero, duplex, and sixplex reagents are present in the UNIMOD database (www.unimod.org) and are listed below. Several software packages directly support the modifications by TMT Reagents and the relative quantitation of reporter ions released from labeled peptides, including Thermo Scientific™ Proteome Discoverer™ 1.1 and above, Matrix Science Mascot™ 2.1 and above, and Proteome Software Scaffold™ Q+. For data acquired using a combination of fragmentation methods (i.e., HCD/MS3 or HCD/ETD), Proteome Discoverer may be necessary to merge spectra for identification and quantitation.

Table 2. Modification masses of the Thermo Scientific TMT Label Reagents.

<u>Label Reagent</u>	<u>Reagent Reporter Ion</u>	<u>Modification Mass (monoisotopic)</u>	<u>Modification Mass (average)</u>	<u>HCD Monoisotopic Reporter Mass*</u>	<u>ETD Monoisotopic Reporter Mass**</u>
TMT ⁰ -126	126	224.152478	224.2994	126.127726	114.127725
TMT ² -126	126	225.155833	225.2921	126.127726	114.127725
TMT ² -127	127C	225.155833	225.2921	127.131081	114.127725
TMT ⁶ -126	126	229.162932	229.2634	126.127726	114.127725
TMT ⁶ -127	127N	229.162932	229.2634	127.124761	115.124760
TMT ⁶ -128	128C	229.162932	229.2634	128.134436	116.134433
TMT ⁶ -129	129N	229.162932	229.2634	129.131471	117.131468
TMT ⁶ -130	130C	229.162932	229.2634	130.141145	118.141141
TMT ⁶ -131	131	229.162932	229.2634	131.138180	119.138176

* HCD is a collisional fragmentation method that generates six unique reporter ions from 126 to 131Da.

**ETD is a non-ergodic fragmentation method that generates six unique reporter ions from 114 to 119Da.

C. Information Available from our Website

- Tech Tip Protocol #49: Acetone precipitation of proteins
- Tech Tip Protocol #19: Remove detergent from protein samples

Related Thermo Scientific Products

90110	TMT10plex™ Isobaric Label Reagent Set, 10 × 0.8mg
90113	TMT10plex Isobaric Mass Tag Labeling Kit
90406	TMT10plex Isobaric Label Reagent Set, 10 × 5mg
90114	1M Triethylammonium bicarbonate (TEAB), 50mL
90115	50% Hydroxylamine, 5mL
90100	iodoTMTzero™ Label Reagent, 5 × 0.2mg
90101	iodoTMTsixplex™ Label Reagent Set, 1 × 0.2mg
90103	iodoTMTsixplex Isobaric Mass Tag Labeling Kit
90076	Immobilized Anti-TMT Antibody Resin
90075	Anti-TMT Antibody, 0.1mL
90104	TMT Elution Buffer, 20mL
84840	Pierce™ Mass Spec Sample Prep Kit for Cultured Cells
23227	BCA Protein Assay Kit
23275	Pierce Quantitative Colorimetric Peptide Assay
23290	Pierce Quantitative Fluorescent Peptide Assay
90057	Pierce Trypsin Protease, MS Grade
90051	Lys-C Protease, MS Grade
88300	Fe-NTA Phosphopeptide Enrichment Kit
88301	Pierce TiO ₂ Phosphopeptide Enrichment and Clean-up Kit
84868	Pierce High pH Reversed-Phase Peptide Fractionation Kit
88321	Pierce Peptide Retention Time Calibration Mixture, 200µL
87784	Pierce C18 Tips, 100µL bed, 96 tips
89870	Pierce C18 Spin Columns, 25 columns
28904	Trifluoroacetic Acid, Sequanal Grade

General References

- Altelaar A.F., et al. (2012). Benchmarking stable isotope labeling based quantitative proteomics. *J Proteomics* Oct 22. pii: S1874-3919(12)00704-X. doi: 10.1016/j.jprot.2012.10.009.
- Bantscheff, M., et al. (2008). Robust and sensitive iTRAQ quantification on an LTQ Orbitrap Mass Spectrometer. *Mol Cell Proteomics* 7:1702-13.
- Biringer, R.G., et al. (2011). Quantitation of TMT-Labeled Peptides Using Higher-Energy Collisional Dissociation on the Velos Pro Ion Trap Mass Spectrometer. Application note # 520. www.lifetechologies.com
- Byers, H.L. (2009). Candidate verification of iron-regulated *Neisseria meningitidis* proteins using isotopic versions of tandem mass tags (TMT) and single reaction monitoring. *J Prot* 73(2):231-9.
- Dayon, L., et al. (2008). Relative quantification of proteins in human cerebrospinal fluids by MS/MS using 6-plex isobaric tags. *Anal Chem* 80(8):2921-31.
- Dillon, R., et al. (2011). Discovery of a Novel B-Raf Fusion Protein Related to c-Met Drug Resistance. *J Proteome Res* 10(11):5084-94.
- Erikson, B.K., et al. (2015). Evaluating multiplexed quantitative phosphopeptide analysis on a hybrid quadrupole mass filter/linear ion trap/orbitrap mass spectrometer. *Anal Chem* 87(2):1241-9.
- Keshishian, H., et al. (2015). Multiplexed, quantitative workflow for sensitive biomarker discovery in plasma yields novel candidates for early myocardial injury. *Mol Cell Proteomics*. 2015 Feb 27. pii: mcp.M114.046813
- Lee, M.V., et al. (2011). A dynamic model of proteome changes reveals new roles for transcript alteration in yeast. *Mol Syst Biol* 7:514.
- McAllister, G.C., et al. (2014). MultiNotch MS3 enables accurate, sensitive, and multiplexed detection of differential expression across cancer cell line proteomes. *Anal Chem* 86(14):7150-8.

- McAllister, G.C., *et al.* (2012). Increasing the multiplexing capacity of TMTs using reporter ion isotopologues with isobaric masses. *Anal Chem* **84(17)**:7469-78.
- Murphy, J.P., *et al.* (2014). Combining amine metabolomics and quantitative proteomics of cancer cells using derivatization with isobaric tags. *Proteomics* **86(7)**:3585-93.
- Paulo, J.A., *et al.* (2014). A comprehensive proteomic and phosphoproteomic analysis of yeast deletion mutants of 14-3-3 orthologs and associated effects of rapamycin. *Nature* **(2-3)**:474-86.
- Ross, P.L., *et al.* (2004). Multiplexed protein quantitation in *Saccharomyces cerevisiae* using amine-reactive isobaric tagging reagents. *Mol Cell Proteomics* **3(12)**:1154-69.
- Savitski, M.M., *et al.* (2014). Tracking cancer drugs in living cells by thermal profiling of the proteome. *Science* **346(6205)**:1255784
- Schirle, M., *et al.* (2012). Kinase inhibitor profiling using chemoproteomics. *Methods Mol Biol* **795**:161-77.
- Schwartz, J. *et al.* (2008). Relative quantitation of protein digests using tandem mass tags and pulsed-Q dissociation (PQD). Application note # 452. www.lifetechnologies.com
- Stella, R., *et al.* (2011). Relative Quantification of Membrane Proteins in Wild-type and PrP-knockout Cerebellar Granule Neurons. *J Proteome Res* doi: 10.1021/pr200759m. <http://dx.doi.org>
- Strupat K., *et al.* (2008). Accurate MS and MSn Analysis with the Thermo Scientific MALDI LTQ Orbitrap. Application note # 30150. www.lifetechnologies.com
- Ting, L., *et al.* (2011). MS3 eliminates ratio distortion in isobaric multiplexed quantitative proteomics. *Nature Methods* **8**: 937–940.
- Van Ulsen, P., *et al.* (2009). Identification of proteins of *Neisseria meningitidis* induced under iron-limiting conditions using the isobaric tandem mass tag (TMT) labeling approach. *Proteomics* **9(7)**:1771-81.
- Viner, R.I., *et al.* (2013). Increasing the multiplexing of protein quantitation from 6- to 10-Plex with reporter ion isotopologues. PN_ASMS_W617_RViner_R1.
- Viner, R.I., *et al.* (2012). Relative quantitation of TMT-labeled proteomes – Focus on sensitivity and precision. Application note #566.
- Viner, R.I., *et al.* (2009). Quantification of post-translationally modified peptides of bovine α -crystallin using tandem mass tags and electron transfer dissociation. *J Proteomics* **72(5)**:874-85.
- Wenger, C.D., *et al.* (2011). Gas-phase purification enables accurate, multiplexed proteome quantification with isobaric tagging. *Nat Methods* **8(11)**:933-5.
- Xiong, L., *et al.* (2011). Mass spectrometric studies on epigenetic interaction networks in cell differentiation. *J Biol Chem* **286(15)**:13657-68.
- Zhang, T., *et al.* (2010). Improving quantitation of TMT-labeled peptides using stepped higher-energy collisional dissociation. Application note # 483 www.lifetechnologies.com

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts (“Documentation”). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER’S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER’S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. Tandem Mass Tag and TMT are trademarks of Proteome Sciences plc. iTRAQ is a trademark of AB Sciex Pte. Ltd. Mascot is a trademark of Matrix Science. Scaffold is a trademark of Proteome Software. Unless otherwise indicated, all other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.