

Paclitaxel and Fluorescent Paclitaxel Conjugates

Quick Facts

Storage upon receipt:

- ≤–20°C
- Desiccate
- Protect from light

Ex/Em: See Table 1

Solvent for stock: Methanol, ethanol, DMSO

Introduction

Molecular Probes offers research-grade paclitaxel (P3456), BODIPY® FL paclitaxel (P7500), Oregon Green® 488 paclitaxel (P22310, also known as Flutax-2) and a red-orange fluorescent BODIPY® 564/570 paclitaxel (P7501) (Table 1). Paclitaxel, formerly referred to as taxol in some scientific literature, is the approved generic name for TAXOL®, a registered trademark of Bristol-Myers Squibb Co. for an anticancer pharmaceutical.

The diterpenoid paclitaxel is a potent anti-neoplastic agent^{1,2} originally isolated from the bark, and more recently from the needles, of the western yew tree, *Taxus brevifolia*.^{3,4} The anti-mitotic and cytotoxic action of paclitaxel is related to its ability to promote tubulin assembly into stable aggregated structures that cannot be depolymerized by dilution, calcium ions, cold or a number of microtubule disrupting drugs.^{5,6} Paclitaxel also decreases the critical concentration of tubulin required for microtubule assembly. Cultured cells treated with paclitaxel are blocked in the G2 (the “gap” between DNA synthesis and mitosis) and M (mitosis) phases of the cell cycle.⁷

Table 1. Molecular Probes' paclitaxel and fluorescent paclitaxel conjugates.

Cat #	Label	Ex *	Em *	MW
P3456	Unlabeled	NA	NA	854
P7500	BODIPY FL	505	515	1024
P7501	BODIPY 564/570	565	571	1099
P22310	Oregon Green 488	496	524	1319

* Approximate fluorescence excitation (Ex) and emission (Em) maxima, in nm. NA, not applicable.

The green-fluorescent BODIPY FL and Oregon Green 488 paclitaxels and the red-orange-fluorescent BODIPY 564/570 paclitaxel are demonstrated tools for studying microtubules in solution.⁸⁻¹⁰

For the BODIPY® FL (Figure 1A) and BODIPY® 564/570 (Figure 1B) paclitaxels the *N*-benzoyl substituent of the 3-phenyl-isoserine portion of native paclitaxel is replaced by a BODIPY® propionyl substituent.

The fluorescent label on Oregon Green® 488 paclitaxel is attached by derivatization of the 7-β-hydroxy group of native paclitaxel (Figure 1C), a strategy that permits selective binding of the probe to microtubules.¹¹

Though Oregon Green® 488 paclitaxel has been successfully used in the past for live-cell tubulin labeling, Molecular Probes now offers TubulinTracker™ Green, an acetylated version of Oregon Green® 488 paclitaxel which is more soluble and cell-permeant, allowing more uniform and consistent cell labeling, superior selectivity for polymerized tubulin in live cells, and increased ease of use.




Figure 1. Structures of A) BODIPY FL, B) BODIPY 564/570 paclitaxel analogs, and C) Oregon Green 488.

Storage and Handling

Unlabeled paclitaxel is supplied as a powder in a unit size of 5 mg. BODIPY® FL and BODIPY® 564/570 paclitaxel are supplied lyophilized in a unit size of 10 µg, and the Oregon Green® 488 paclitaxel, in a unit size of 100 µg. These products should be stored frozen at ≤-20°C, desiccated, and protected from light. Concentrated stock solutions can be prepared in methanol, ethanol, or dimethylsulfoxide (DMSO), aliquoted and stored

for up to one month frozen at ≤-20°C, protected from light. Aqueous working solutions should be prepared immediately before use. A typical working concentration for experimentation with unlabeled paclitaxel is 0.1 µM and with labeled paclitaxel is 1 µM.

Warning: Paclitaxel is potentially harmful if inhaled, ingested or allowed to contact skin or eyes. Avoid prolonged or repeated exposure. Do not get on skin, in eyes, or on clothing. Wear gloves when using these products and wash thoroughly after handling.

References

1. Invest New Drugs 9, 59 (1991); 2. J Am Chem Soc 93, 2325 (1971); 3. J Am Chem Soc 110, 5917 (1988); 4. Tetrahedron 42, 4451 (1986); 5. J Cell Biol 112, 1177 (1991); 6. Pharmacol Ther 25, 83 (1984); 7. Cancer Treat Rep 62, 1219 (1978); 8. BioTechniques 25, 188 (1998); 9. Biochemistry 40, 11975 (2001); 10. Biochemistry 41, 12436 (2002); 11. Int J Mol Med 2, 161 (1998).

Product List *Current prices may be obtained from our Web site or from our Customer Service Department.*

Cat #	Product Name	Unit Size
P3456	paclitaxel (Taxol equivalent) *for use in research only*	5 mg
P7501	paclitaxel, BODIPY® 564/570 conjugate (BODIPY® 564/570 Taxol)	10 µg
P7500	paclitaxel, BODIPY® FL conjugate (BODIPY® FL Taxol)	10 µg
P22310	paclitaxel, Oregon Green® 488 conjugate (Oregon Green® 488 Taxol)	100 µg
T34075	TubulinTracker™ Green (Oregon Green® 488 Taxol, bis-acetate) *for live-cell imaging*	1 set

Contact Information

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Please visit our website — probes.invitrogen.com — for the most up-to-date information.

Molecular Probes, Inc.

29851 Willow Creek Road, Eugene, OR 97402
Phone: (541) 465-8300 • Fax: (541) 335-0504

Customer Service: 6:00 am to 4:30 pm (Pacific Time)

Phone: (541) 335-0338 • Fax: (541) 335-0305 • probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

Technical Service: 8:00 am to 4:00 pm (Pacific Time)

Phone: (541) 335-0353 • Toll-Free (800) 438-2209
Fax: (541) 335-0238 • probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100 • Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License

For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic,

vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402. Tel: (541) 465-8300. Fax: (541) 335-0504.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2005, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.