## INSTRUCTIONS



# TMB Substrate Kit

34021 0204.3

Number Description
34021 TMB Substrate Kit

**Kit Contents:** 

TMB Solution (0.4g/L), 200mL

Peroxide Solution, 200mL, contains 0.02% hydrogen peroxide in citric acid buffer

**Storage:** Upon receipt store kit at 4°C. Kit is shipped at ambient temperature.

## Introduction

The Thermo Scientific TMB Substrate Kit contains soluble TMB (3,3',5,5'-tetramethylbenzidine) that is ideal for sensitive ELISA-based applications. TMB is a chromogenic substrate for horseradish peroxidase. Once oxidized by the enzyme, this substrate yields a blue product that absorbs at 370nm and 652nm. Adding sulfuric or phosphoric acid to stop the reaction changes the product to yellow that absorbs at 450nm.

## **Procedure for using TMB for ELISAs**

Note: Color variations of the TMB Solution are normal and will not interfere with ELISA results.

- 1. Immediately before use, mix equal volumes of the TMB Solution and the Peroxide Solution.
- 2. Add 100μL of the TMB substrate solution to each microplate well.
- 3. Incubate plate at room temperature for 15-30 minutes or until the desired color develops.

**Note:** High concentrations of HRP yield a greenish solution. Stop the reaction before any wells display a green product. A precipitated product indicates the presence of too much HRP and the need to optimize experimental conditions.

- 4. Stop reaction with 100μL of 2M sulfuric acid or 8M acetic acid with 1M sulfuric acid.
- 5. Measure the absorbance of each well at 450nm.

### Related Thermo Scientific Products

| 15075 | Reagent Reservoirs, 200/pkg.                          |
|-------|---|
| 15082 | Microtube Racked System, 10 racks                     |
| 15036 | Sealing Tape for 96-Well Plates, 100/pkg.             |
| 21130 | Pierce® High Sensitivity Streptavidin-HRP, 0.5mL      |
| 31030 | Pierce High Sensitivity NeutrAvidin®-HRP, 0.5mL       |
| 21124 | Streptavidin, Horseradish Peroxidase Conjugated, 2mg  |
| 31002 | NeutrAvidin, Horseradish Peroxidase Conjugated, 2mg   |
| 37538 | StartingBlock <sup>TM</sup> (PBS) Blocking Buffer, 1L |
| 37542 | StartingBlock (TBS) Blocking Buffer, 1L               |

#### **General References**

Bos, E.E., et al. (1981). 3,3'5,5'-Tetramethylbenzidine as an Ames test negative chromogen for horseradish peroxidase in enzyme immunoassay. *J Immunoassay* 2:187-204.

Josephy, P., et al. (1982). The horseradish peroxidase-catalyzed oxidation of 3, 5,3',5' Tetramethylbenzidine. J Biol Chem 257(7):3669-75.



This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at <a href="www.thermoscientific.com/pierce">www.thermoscientific.com/pierce</a>. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.