SeeBlue® Plus2 Pre-Stained Protein Standard Protocol

## Protocol

novex

by *life* technologies

sure the solution is	s homogeneous.
Applicaton	
Electrophoresis	Blotting
10 µL	3 µL
20 µL	5 µL
20 µL	7 µL
30 µL	10 µL
	Electrophoresis   10 μL   20 μL   20 μL

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support

## Important licensing information

This products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

## Disclaimer

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT, TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT. INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Buffer

kDa

198

98

62 49

38

28

17

14

6

3

10 µL of SeeBlue®

Plus2 Pre-Stained

Protein Standard

Bis-Tris Gel with

MES SDS Running

on a NuPAGE® 10%

For support, visit www.lifetechnologies.com/support.

Publication Part No. LC5925.pps Publication No. MAN0007923 Rev B.0 **Catalog Numbers** Package Size: LC5925 500 uL contents Storage Store at 2°C to 8°C. conditions

Required materials

Selection

Product

auide

 Polyacrylamide gel(s) Electrophoresis apparatus for protein analysis

Pre-Stained Protein Standards Go online to view related products.

SeeBlue<sup>®</sup> Pre-Stained Protein Standard allows you to accurately determine molecular weight ranges and evaluate western transfer efficiency during electrophoresis.



This standard consists of 10 pre-stained protein bands lescription (8 blue and 2 contrasting colors) ranging in molecular weight from ~3-200 kDa.

> Storage buffer: Tris-HCl, Formamide, SDS, and Phenol Red.



Use with NuPAGE<sup>®</sup>, Tris-Glycine, and Tricine gels. Click here to view our online Protein Gel Selection Guide.

Important quidelines

DO NOT apply heat or add a reducing agent. SeeBlue<sup>®</sup> Plus2 Pre-Stained Standard is supplied ready-to-use.



Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support.

For Research Use Only. Not for use in diagnostic procedures. 18 March 2014

