CERTEST

Calprotectin+

Lactoferrin


ONE STEP human calprotectin
and human lactoferrin COMBO
CARD TEST

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ENGLISH

INTENDED USE

CerTest Calprotectin+Lactoferrin one step combo card test is a coloured chromatographic immunoassay (cut-off qualitative test) for the simultaneous qualitative determination of human calprotectin (hCp) and human lactoferrin (hLf) in stool samples that may reflect gastrointestinal inflammation caused by several pathologies (inflammatory bowel disease, colorectal cancer and some enteropathies).

CerTest Calprotectin+Lactoferrin combo card test offers a simple, a highly sensitive and non-invasive screening assay to determine intestinal inflammatory activity, to monitor treatment response and to predict risk of relapse.

INTRODUCTION

Calprotectin is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool during bowel inflammation. The stability of the protein to degradation keeps it stable in faeces for up to 7 days at room temperature, making it an ideal analyte. Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s). Calprotectin inhibits zinc-dependent enzyme systems, as a result kills microbes and induces apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is a remarkably resistant to proteolytic degradation and so is stable in stools kept at room temperature for 7 days.

Lactoferrin is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response and is released from faecal leukocytes. This protein is resistant to proteolysis in the faeces and may serve as a marker of inflammation in the intestine. The major cause of faecal neutrophils in patients with chronic diarrhoea is chronic inflammatory intestine disease of the colon (i.e., Crohn's Disease and Ulcerative Colitis). Lactoferrin has been also studied as a predictor of infection with invasive enteropathogens in children with diarrhoea. Bacterial inflammatory diarrhoea may be caused by *Shigella*, *Salmonella*, *Campylobacter* and *Clostridium difficile*.

The hCp and hLf are non-invasive markers of intestinal inflammation (for example in Ulcerative Colitis (UC) and Crohn's Disease (CD)).

Suspicion of gastrointestinal inflammatory disorders is raised in patients with persistent (≥ 4 weeks) or recurrent (≥ 2 episodes in six months) abdominal pain and diarrhoea. Additionally, rectal bleeding, weight loss or anemia increase the probability of the condition. Endoscopic evaluation with histopathological sampling is generally considered indispensable to diagnosis patients with gastrointestinal inflammation disorder suspicion.

TEST PRINCIPLE

CerTest Calprotectin+Lactoferrin is based on the principle of a cut-off qualitative immunochromatographic assay for the determination of human calprotectin and human lactoferrin in stool samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human calprotectin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human calprotectin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human lactoferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human lactoferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is human calprotectin positive, the antigens of the diluted sample reacts with the red-coloured conjugate complex (anti-human calprotectin monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is human lactoferrin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human lactoferrin monoclonal antibodies-red polystyrene microspheres) in the strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human calprotectin antibodies present on the membrane of strip A (test line) and the anti-human lactoferrin antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no human calprotectin and human lactoferrin presence and yet, the antigens may be present in a concentration lower than the cut-off value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-human calprotectin and anti-human lactoferrin antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate

complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest Calprotectin+Lactoferrin test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1); making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

MATERIALS

MATERIALS PROVIDED

- CerTest Calprotectin+Lactoferrin combo card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

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

TEST PROCEDURE

Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.


1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest Calprotectin+Lactoferrin combo card test from its sealed bag just before using it.
3. Take the stool collection tube, open the cap 2 (4) and dispense 3 drops in the circular window marked with the letter A (5) and 3 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
4. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.


Calprotectin Strip A-procedure

Add 3 drops in A-circular window

Lactoferrin Strip B-procedure


Add 3 drops in B-circular window

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

<p>1.</p> <p>A: Green → Negative hCp B: Green → Negative hLf</p>	<p>2.</p> <p>A: Green/Red → Positive hCp B: Green/Red → Positive hLf</p>	<p>3.</p> <p>A: Green/Red → Positive hCp B: Green → Negative hLf</p>	<p>4.</p> <p>A: Green → Negative hCp B: Green/Red → Positive hLf</p>
--	--	--	--

INVALID RESULTS

	A (Calprotectin (hCp))	B (Lactoferrin (hLf))	Interpretation of the results
1.	- GREEN	- GREEN	There is no human calprotectin and/or human lactoferrin presence which might mean neither active gastrointestinal inflammation, nor risk of relapse (CD or UC relapse).
2.	+ GREEN-RED	+ GREEN-RED	There is human calprotectin and human lactoferrin presence which might mean an inflammatory intestinal disease (caused by ulcerative colitis (UC), Crohn's disease (CD), some type of carcinomas or an enteropathy caused by <i>Shigella</i> , <i>Salmonella</i> , <i>Campylobacter</i> or <i>Clostridium difficile</i>) or risk of relapse in clinical remission.
3.	+ GREEN-RED	- GREEN	There is human calprotectin presence which might mean active gastrointestinal inflammation or risk of relapse in clinical remission.
4.	- GREEN	+ GREEN-RED	There is human lactoferrin presence which might mean an inflammatory intestinal disease.
5.	Any other result	Any other result	Invalid result, either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the result lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay using a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of human calprotectin and human lactoferrin present in the specimen.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- The intensity of test line may vary from very strong at high human calprotectin/lactoferrin concentration to faint when the human calprotectin/lactoferrin concentration is close to cut-off value of the test.
- CerTest Calprotectin+Lactoferrin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of human calprotectin and/or human lactoferrin in faecal samples; nevertheless, it can be due to several causes (Inflammatory bowel disease, colorectal cancer and some other enteropathies). A positive result should be followed up with additional diagnostic invasive procedures, a colonoscopy and a biopsy in order to confirm the diagnosis and to establish the inflammation extent.

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- A negative result is not meaningful because of it is possible the human calprotectin and/or human lactoferrin concentration in the stool sample is lower than cut-off values. If symptoms or situation still persist, calprotectin and lactoferrin determination should be carried out invasive techniques Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that involve mainly monocuclear inflammation.
- Stool from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show positive calprotectin result.
- Lactoferrin is a component of breast milk; the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.
- Neonatal faecal calprotectin levels have been reported higher than normal children with a median of 167µg/g (range 22-860µg/g).
- Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's Disease and Ulcerative Colitis, would be positive for faecal calprotectin (the level of calprotectin in Crohn's Disease is higher than in Ulcerative Colitis) and would be positive for faecal lactoferrin. CerTest Calprotectin+Lactoferrin could be sensitive for this diagnosis in patients with chronic diarrhea.
- Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

EXPECTED VALUES

Some studies established equal or higher 50µg hCp/g faeces and 10µg hLf/g faeces as cut-off values to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Cut-off values of CerTest Calprotectin+Lactoferrin are 500ng/mL (50µg hCp/g faeces) for human calprotectin and 100ng/mL (10µg hLf /g faeces) for human lactoferrin.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest Calprotectin+Lactoferrin, CerTest), other commercial immunoassays (Calprest®, Eurospital) for the Strip A and other commercial ELISA tests (Human Lactoferrin ELISA kit, ICL) for the Strip B. The results were as follows:

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Data for cut-off value cut-off 10 µg hLf/g

IC test: CerTest Calprotectin+Lactoferrin (Strip B: Lactoferrin)	ELISA test: Human Lactoferrin ELISA kit			Total
		+	-	
+		27	13	40
-		3	47	50
Total		30	60	90

CerTest Calprotectin+Lactoferrin (Strip B: Lactoferrin) vs Human Lactoferrin ELISA kit		
	Mean Value	95% CI (confidence interval)
Sensitivity	90.0%	73.5 – 97.9%
Specificity	78.3%	65.8 – 87.9%
PPV	67.5%	50.9 - 81.4%
NPV	94.0%	83.5 –98.7%

Specificity for a cut-off value of 2.5 µg hLf/g: 100.0% (mean value); 83.9-100.0% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human calprotectin and human lactoferrin using CerTest Calprotectin+Lactoferrin.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest Calprotectin+Lactoferrin; no cross reactivity against other gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.
5. Gisbert, J.P. et al. Role of Biological Markers in Inflammatory Bowel Disease. *Gastroenterol. Hepatol.* May 2007, 30(3), p. 117-29.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/

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