





CA(PEG)_n Reagents

Carboxy-amine PEGylation reagents

Pub. No. MAN0011634

Rev. B.0

Pub. Part No. 2162066

Number	Description	
26120 26121	CA(PEG) ₄ , 100mg CA(PEG) ₄ , 1 g Form: White to pale yellow solid Molecular Weight: 265.30 Spacer Arm: 18.1 Å	 <p>Carboxy-PEG₄-Amine</p>
26122 26123	CA(PEG) ₈ , 100mg CA(PEG) ₈ , 1 g Form: White to pale yellow solid Molecular Weight: 441.71 Spacer Arm: 33.6 Å	 <p>Carboxy-PEG₈-Amine</p>
26124 26125	CA(PEG) ₁₂ , 100mg CA(PEG) ₁₂ , 1 g Form: White to pale yellow solid Molecular Weight: 617.72 Spacer Arm: 46.8 Å	 <p>Carboxy-PEG₁₂-Amine</p>
26126 26127	CA(PEG) ₂₄ , 100mg CA(PEG) ₂₄ , 1 g Form: White to pale yellow waxy solid Molecular Weight: 1146.35 Spacer Arm: 89.8 Å	 <p>Carboxy-PEG₂₄-Amine</p>

Storage: Upon receipt store desiccated at 4°C. Product is shipped at ambient temperature.

Introduction

The carboxy-PEG_n-amine (CA[PEG]_n) PEGylation reagents are zwitterionic, amino acid derivatives that are used for modifying proteins or surfaces such as beads, nanoparticles and self-assembled monolayers. Modification of proteins adds polyethylene glycol (PEG) spacers, which impart increased water solubility, reduced immunogenicity of the labeled molecule and enhanced *in vivo* stability in solution.¹ Functionalization of solid surfaces, such as quantum dots, self-assembled monolayers and nanoparticles, with polyethylene glycol spacers significantly reduces nonspecific protein binding.²⁻⁷ CA(PEG)_n reagents used with MA(PEG)_n reagents in surface modification can form a hydrophilic “lawn” of methyl ether-terminated PEGs with periodic exposed carboxy-terminated PEGs (Figure 1). The exposed carboxy groups can be coupled to affinity ligands using the carbodiimide coupling reaction with EDC and sulfo-NHS.

Typical PEGylation reagents contain heterogeneous mixtures of different PEG chain lengths; however, our PEGylation reagents are homogenous compounds of defined molecular weight and spacer length, providing precision in optimizing modification applications.




Figure 1. Two-step modification of surfaces with CA(PEG)_n-MA(PEG)_n reagents.

Important Product Information

- The CA(PEG)_n reagents are low-melting solids that are difficult to weigh and dispense. To facilitate handling, make a stock solution by dissolving the reagent with water or methylene chloride. Ethyl acetate and tetrahydrofuran (THF) are also effective and compatible solvents.
- Use the MA(PEG)_n reagents (see Related Thermo Scientific Products) in combination with CA(PEG)_n reagents to modify surfaces and minimize nonspecific binding.
- Use non-amine-containing buffers at pH 7-9 such as PBS (20mM sodium phosphate, 150mM NaCl; pH 7.4) (Product No. 28372); 20mM HEPES; 100mM carbonate/bicarbonate; or 50mM borate. Do not use buffers that contain primary amines, such as Tris or glycine, which compete with acylation.
- The CA(PEG)_n-to-MA(PEG)_n ratio and the reagent mixture-to-surface carboxylic acid molar ratio in the reaction affects the number of carboxy groups modified on the surface and the number of new carboxylic acid residues available for further modification. Optimize these ratios to obtain the modification level needed for the specific application.

Procedure for Coupling CA(PEG)_n–MA(PEG)_n Mixtures to Carboxylated Surfaces

The following protocol, adapted from a procedure described by Grabarek and Gergely⁸ is a two-step coupling reaction using EDC and NHS or Sulfo-NHS. The CA(PEG)_n–MA(PEG)_n mixture is coupled to a carboxylated surface without exposing the carboxylic acids on CA(PEG)_n to EDC. The activation reaction requires quenching with a thiol-containing compound.

The activation reaction with EDC and Sulfo-NHS is most efficient at pH 4.5-7.2; however, the reaction of Sulfo-NHS-activated molecules with primary amines is most efficient at pH 7-8. For best results, perform the first reaction in MES buffer (or other non-amine, non-carboxy buffer) at pH 5-6, then raise the pH to 7.2-7.5 with phosphate buffer (or other non-amine buffer) immediately before reacting with the CA(PEG)_n–MA(PEG)_n mixture. Use DTT to quench the activation reaction. The conjugation reaction is quenched using hydroxylamine, Tris or glycine.

Materials Required

- Water, methylene chloride, ethyl acetate or tetrahydrofuran (THF) for preparing the reagent stock solution
- Small-volume, non-coring syringes for dispensing the reagent stock solution while minimizing exposure to air
- EDC (Product No. 77149)
- Activation Buffer: MES-buffered saline (0.1M MES, 0.5M NaCl; pH 6.0 or 0.1M MES, 0.9% NaCl; pH 4.7; Product No. 28390)
- Conjugation Buffer: Phosphate-buffered saline, PBS (20mM sodium phosphate, 0.15M NaCl; pH 7.2, Product No. 28372)
- NHS or Sulfo-NHS (Product No. 24500 and 24510, respectively)
- Dithiothreitol (DTT; Product No. 20290 or 20291)
- Hydroxylamine•HCl (Product No. 26103)

Procedure

1. Equilibrate EDC, NHS or sulfo-NHS, CA(PEG)_n, and MA(PEG)_n to room temperature before opening bottles.
2. Prepare CA(PEG)_n and MA(PEG)_n stock solutions by dissolving 100mg of each reagent (~100μL) in the desired amount of solvent (e.g., water, methylene chloride).
3. Add appropriate amounts of EDC and NHS or sulfo-NHS to the appropriate amount of carboxylated surface in Activation Buffer and react for 15 minutes at room temperature.
4. Add DTT to quench the EDC.

Note: For surfaces that can be easily washed, the quenching step can be skipped and the surface washed with Coupling Buffer to remove any remaining EDC and NHS.

5. Add the CA(PEG)_n–MA(PEG)_n mixture prepared in Conjugation Buffer to the activated surface and react for 2 hours at room temperature.
6. To quench the reaction, add hydroxylamine or another amine-containing buffer. Hydroxylamine hydrolyzes non-reacted NHS on the solid surface and results in hydroxamate formation. Other quenching methods involve adding Tris, lysine, glycine or ethanolamine; however, these primary amine-containing compounds modify carboxyls.

Note: The newly introduced carboxy groups can be further modified by repeating steps 4 and 5.

7. Add the desired amine-containing substrate, prepared in Coupling Buffer, to the activated surface and react for 2 hours at room temperature.
8. Quench the reaction as described in step 7.

Related Thermo Scientific Products

26110	MA(PEG) ₄ , 100mg
26111	MA(PEG) ₄ , 1g
26112	MA(PEG) ₈ , 100mg
26113	MA(PEG) ₈ , 1g
26114	MA(PEG) ₁₂ , 100mg
26115	MA(PEG) ₁₂ , 1g
26116	MA(PEG) ₂₄ , 100mg
26117	MA(PEG) ₂₄ , 1g
28390	BupH™ MES Buffered Saline, 10 packs, makes 5 L
28372	BupH Phosphate Buffered Saline, 40 packs, makes 20 L
77149	EDC, 10mg
24500	NHS (<i>N</i> -hydroxysuccinimide), 25g
24510	Sulfo-NHS (sulfo <i>N</i> -hydroxy succinimide), 500mg
20290	DTT, 5g
20291	DTT, No-Weigh™ Format, 48 tubes × 7.7mg
26103	Hydroxylamine, 25g

Cited References

1. Morar, A.S., *et al.* (2006). PEGylation of proteins: A structural approach. *BioPharm. Int.* April 34-46.
2. Prime, K.L. and Whitesides, G.M. (1991). Self-assembled organic monolayers: model systems for studying absorption of proteins at surfaces. *Science* **252**:1164.
3. Bentzen, E.L., *et al.* (2005). Surface modification to reduce non-specific binding of quantum dots in live cell assays. *Bioconjugate Chem.* **16**:1488-94.
4. Lin, P-C., *et al.* (2006). Ethylene glycol-protected magnetic nanoparticles for a multiplexed immunoassay in human plasma. *Small* **2(4)**:485-9.
5. Zheng, M., *et al.* (2003). Ethylene glycol monolayer protected nanoparticles for eliminating nonspecific binding with biological molecules. *J. Am. Chem. Soc.* **125**:7790-1.
6. Verma, A. and Rotello, V.M. (2005). Surface recognition of biomacromolecules using nanoparticle receptors. *Chem. Commun.* **3**:303-12.
7. Kidambi, S., *et al.* (2004). Selective depositions on polyelectrolyte multilayers: self-assembled monolayers of m-dPEG acid as molecular template. *J. Am. Chem. Soc.* **126**:4697-03.
8. Grabarek, Z. and Gergely, J. (1990). Zero-length crosslinking procedure with the use of active esters. *Anal. Biochem.* **185**:131-5.

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2017 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all other trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.