

96 well Magnetic-Ring Stand

Catalog # (P/N): AM10050

One Unit Includes: One 96 well Magnetic-Ring Stand

Precautions: The magnets in this stand are extremely strong. Do not allow the unit to come in contact with metal objects or other

magnets. Damage will occur to magnetized media, such as diskettes or credit cards, near the stand. Damage may

also occur to computers and CRT-based monitors near the stand.

Cleaning: If solutions are spilled onto the 96 well Magnetic-Ring Stand, rinse the unit with tap water and allow it to air dry

completely.

USER INFORMATION

Product Description: The Ambion 96 well Magnetic-Ring Stand is designed for paramagnetic bead precipitation from standard 96-well

U-bottom microplates. The stand has a plastic base in standard SBS 96-well microtiter plate footprint. The use of 96 extremely strong permanent magnetic rings enables fast and easy magnetic separation. With this magnetic stand, paramagnetic beads in a U-bottom plate will pellet in a donut shape ~6 mm in diameter, which enables easy removal of supernatant by robot or with multichannel pipettor. It sits flat on a lab bench, as well as in a microtiter plate adapter on robotic workstations. The 96 well Magnetic-Ring Stand is ideal for use with Ambion MagMAX™-96 magnetic bead-based RNA isolation kits (P/N AM1830, AM1835, AM1836, AM1837, AM1839). For more information, visit

www.ambion.com.

Applications: If you are using this magnetic stand with a multichannel pipettor, especially an electric pipettor, we strongly

recommend refining your pipetting technique to ensure you do not disrupt the pellet. A normal ring-shaped pellet will form when liquid covers the bead pellet. The pellet may take on an irregular (non-ring) shape after the liquid is

removed; this is normal and will not affect performance.

If any magnets come out of place, carefully re-insert them so that adjacent magnets have opposite poles in contact with the 96-well plate. Loose magnets may accelerate towards one another or other metal objects nearby causing damage. Replace loose magnets immediately to avoid damage. For additional information, contact the Ambion Technical Services Department (1-800-888-8804, option 2 or TechServ@ambion.com).

QUALITY CONTROL

Each lot is functionally tested for magnetic separation of beads in a 96-well U-bottom microplate. **Functional Testing:**

OTHER INFORMATION

Material Safety Data Sheets (MSDSs) can be printed or downloaded from product-specific links on our website at the following **Material Safety Data Sheets:**

address: www.ambion.com/techlib/msds. Alternatively, e-mail your request to MSDS Inquiry CCRM@appliedbiosystems.com. Specify the catalog or part number(s) of the product(s), and we will e-mail the associated MSDSs unless you specify a preference for fax delivery. For customers without access to the internet or fax, our technical service department can fulfill MSDS requests placed by

telephone or postal mail. (Requests for postal delivery require 1–2 weeks for processing.)

Warranty and Liability: For research use only. Not for use in diagnostic procedures.

Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service

and technical support teams. Warranty information for the accompanying consumable product is available at www.ambion.com/info/warranty in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions,

and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about

post-warranty support.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document.

Applied Biosystems disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Applied Biosystems be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection

with or arising from this document, including but not limited to the use thereof.

Trademarks, Patents, and Licensing:

Applied Biosystems, AB (Design), and Ambion are registered trademarks, and MagMAX is a trademark of Applera Corporation or its subsidiaries in the US and/or certain other countries. All other trademarks are the sole property of their respective owners.

© 2008 Ambion, Inc. All rights reserved. 4386456A