Protocol Pub. No. MAN0007855 Rev. 1.0

BLOCK-iT™ Alexa Fluor® Red Fluorescent Control

Packag Conten	e Catalog Number Size 14750-100 20 μM stock in solution • 1.75 mL Nuclease-free Water
Storage Conditi	 Store at or below -20°C. Do not store in a frost-free freezer. (Dried oligonucleotides are shipped at room temperature.)
Require Materia	 • RNase-free reagents • Transfection reagent e.g. Lipofectamine[®] RNAiMAX
O Timing	Transfection preparation: 15 minutes Final incubation: 1–3 days
Selecti Guide	on siRNAs Go online to view related products.
Produc Descrip	 t BLOCK-iTTM Alexa Fluor[®] Red Fluorescent Control is an Alexa Fluor[®] 555-labeled, double-stranded, RNA (dsRNA) duplex for assessing lipid-mediated transfection for RNAi experiments. Alexa Fluor[®] 555 dye can be detected using standard filter sets designed for Cy[®]3, DsRed, Texas Red[®], or rhodamine fluorophores.
Importa Guideli	 Handling instructions: RNA oligonucleotides are susceptible to degradation by exogenous ribonucleases introduced during handling. Wear gloves when handling this product. Use RNase-free reagents, tubes, and barrier pipette tips. Transfection efficiency varies according to the cell type and transfection agent used. If you are transfecting your cell line for the first time, we recommend starting with a 10 nM final concentration of the BLOCK-iTTM Alexa Fluor[®] Red Fluorescent Control to determine the optimal amount for a strong fluorescence signal. Transfect BLOCK-iTTM Alexa Fluor[®] Red Fluorescent Control using the same methodology as for your
Online	experimental siRNA duplexes. Visit our product page for additional information and protocols. For support,

visit www.lifetechnologies.com/support.

For Research Use Only. Not for use in diagnostic procedures.



This control is designed for use in RNAi analysis to facilitate assessment and optimization of cationic lipid-mediated delivery or electroporation of dsRNA oligonucleotides into mammalian cells. It has the following characteristics:

invitrogen

by *life* technologies"

- Chemical modifications that enhance the stability and allow assessment of fluorescence signal for a significantly longer time than is obtained with other unmodified, fluorescently labeled RNA. It has the same length, charge, and configuration as standard siRNA.
- Proven correlation of transfection efficiency with siRNA molecules.
- Localization primarily to the nucleus upon uptake, designed strictly for use as a tool for siRNA uptake assessment.
- No significant sequence similarity to mouse, rat, or human transcript sequences and has been tested in multiple cell lines and shown to have no significant impact on cell proliferation, apoptosis, or cell morphology.

RNAi Transfection Protocol

See page 2 to view guidelines for transfecting siRNAs using Lipofectamine[®] RNAiMAX Reagent.

Transfection Amounts per Well

Use 10 nM siRNA duplex as a starting point.

	96-well	24-well	6-well
Final siRNA	1 pmol	5 pmol	25 pmol
Final Lipofectamine® RNAiMAX	0.3 µL	1.5 µL	7.5 µL

Limited Product Warranty and Disclaimer Details Limited Use Label License

RNAi Transfection Protocol

This procedure is designed for one RNA amount combined with one amount of Lipofectamine[®] RNAiMAX. The prepared mix is enough to have triplicates (96-well), duplicates (24-well), and single well (6-well) transfections, and account for pipetting variations.

	Timeline		Steps		Procedure Details				
Day 0			Seed cells to be 60-80% confluent at transfection		Component	96-well	24-well	6-well	
	1				Adherent cells	$1-4 \times 10^{4}$	$0.5-2 \times 10^{5}$	$0.25 - 1 \times 10^{6}$	
Day 1			Dilute Lipofectamine®		Opti-MEM [®] Medium	25 µL	50 µL	150 µL	
	2		RNAiMAX Reagent in Opti-MEM [®] Medium		Lipofectamine [®] RNAiMAX Reagent	1.5 µL	3 µL	9 µL	
		>	Dilute ciDNA in		Opti-MEM [®] Medium	25 µL	50 µL	150 µL	
	3		Opti-MEM [®] Medium	siRNA (10 µM)	0.5 μL (5 pmol)	1 μL (10 pmol)	3 μL (30 pmol)		
			Add diluted siRNA to		Diluted siRNA	25 μL	50 µL	150 µL	
	4		RNAiMAX Reagent (1:1 ratio)	RNAiMAX Reagent (1:1 ratio)		Diluted Lipofectamine® RNAiMAX Reagent	25 μL	50 μL	150 µL
	5	5	Incubate		Incubate for 5 minutes at room temperature.				
	6		Add siRNA-lipid complex to cells		Component	96-well	24-well	6-well	
		1		siRNA-lipid complex per well	10 µL	50 µL	250 µL		
				Final siRNA used per well	1 pmol	5 pmol	25 pmol		
					Final Lipofectamine [®] RNAiMAX used per well	0.3 µL	1.5 µL	7.5 µL	
Day 2-4	7		Visualize/analyze transfected cells		Incubate cells for 1–3 days at 37°C. Then, analyze transfected cells.				