# **INSTRUCTIONS**



Pub. No. MAN0011588

Pub Part No. 2161935.1

Rev. B.0

# Low-Fluorescence PVDF Transfer Membrane

22860

Number

22860

#### Description

Low-Fluorescence PVDF Transfer Membrane, 0.2 µm, 7 cm × 8.4 cm, 10 sheets

**Storage:** Store membranes flat at ambient temperature and away from chemical vapors. Some solvent vapors may partially dissolve the membranes, which will disrupt the pore structure.

#### Introduction

Polyvinylidene difluoride (PVDF) membranes are hydrophobic and have high binding affinity for proteins and nucleic acids. The Low-Fluorescence PVDF Membrane has been optimized for use with fluorescent probes. This membrane has lower background levels and increased sensitivity for fluorescent probing than regular PVDF or nitrocellulose membranes.

## **Important Product Information**

- Bromophenol blue can emit a strong fluorescent signal that interferes with protein detection. Either do not use bromophenol blue in the gel loading buffer or make sure that the dye front has migrated from the gel before transfer.
- Do not use a ballpoint pen to write on membranes because ink emits a fluorescent signal and can smear across the membrane causing high background levels. If a membrane needs to be marked, use a pencil.
- Once the membrane has been hydrated, only use forceps to move the membrane. Grip only the edges and avoid scratching or bending the membrane, which increases background levels.
- Stripping and reprobing fluorescent Western blots is not recommended because results are typically inconsistent.

## **Example Procedure for Transferring Proteins to a PVDF Membrane**

**Note:** Always wear gloves when handling PVDF membranes because oils from fingers may prevent proper wetting. Proteins from hands may also bind to the membranes, causing background.

#### A. Materials required

- Completed gel separation of proteins: Use any suitable protocol to separate proteins by electrophoresis
- 100% Methanol
- Cold transfer buffer: Use any suitable buffer for transferring proteins from acrylamide gels to membranes (see the Related Products Section)
- Western blotting filter paper (see Related Products Section)
- Transfer unit (including fiber pads) and power supply

#### Method


1. Remove gel from cassette and equilibrate in cold transfer buffer for 10-15 minutes.

Note: Incubation time is based on 1.0 mm thick gels. Increase incubation time for thicker gels.

- 2. For orientation with the gel, cut a notch in a corner of the membrane.
- 3. Wet the membrane in 100% methanol for 15 seconds and rinse with ultrapure water for 2 minutes with agitation.
- 4. Equilibrate the membrane in cold transfer buffer for 10-15 minutes.



- 5. Wet the filter paper and fiber pads in cold transfer buffer for 5 minutes.
- 6. Assemble the components in the following manner for transfer:



7. Attach the leads and perform the transfer at 200 mA for 90 minutes.

**Note:** Transfer time and efficiency will vary depending on polyacrylamide concentration, gel thickness, the presence of SDS or methanol, pH and ionic strength of the transfer buffer, and the molecular weight of the protein. Determine the optimal transfer conditions empirically.

- 8. When the transfer is complete, disconnect the leads from the transfer unit and disassemble the transfer components to remove the membrane.
- 9. Keep membrane moist until ready to use.

#### **Related Products**

Please see the catalog or website for a complete listing of products for Western blotting.

XP04200BOX	Novex <sup>TM</sup> Tris-glycine protein gels (see <u>thermofisher.com/proteingels</u> for a complete listing)
NW04120BOX	Bolt <sup>TM</sup> Bis-Tris Plus protein gels (see <u>thermofisher.com/proteingels</u> for a complete listing)
LC2676	Novex™ Tris-Glycine SDS Sample Buffer
LC2675	Novex Tris-Glycine SDS Running Buffer
88600	Western Blotting Filter Paper, 100 sheets
LC5615	iBright™ Prestained Protein Ladder
A32727	Goat Anti-Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor Plus 555
A32732	Goat Anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor Plus 555
A32728	Goat Anti-Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor Plus 647
A32733	Goat Anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor Plus 647

For research use only. Not for use in diagnostic procedures.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current versions of product instructions are available at thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2017 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.