invitrogen detection technologies

Revised: 20-June-2005

Labeled aha-dUTP and aha-dCTP

Quick Facts

Storage upon receipt:

• ≤-20°C

Concentration: 1 mM nucleotide solution in TE buffer (pH 8.0)

Introduction

Labeled aha-dUTP and aha-dCTP (5-aminohexylacrylamido-dUTP and -dCTP) nucleotides are modified with a unique hexylacrylamide linker at the C-5 position of uridine and cytosine, respectively. This modification serves as a spacer between the nucleotide and the dye and reduces interactions between these two portions of the molecule, resulting in brighter conjugates and increased hapten accessibility for secondary detection reagents. Labeled nucleotides can be used to generate labeled nucleic acid hybridization probes for many molecular biology and molecular cytogenetics applications, including multicolor techniques.¹⁻⁵

The Alexa Fluor® 555 and Alexa Fluor® 647 dyes used to label these nucleotides are compatible with commonly used microarray scanners, and provide greater signal correlation (R²) values than the spectrally similar Cy™3 and Cy™5 dye pair, improving the resolution of two-color microarray gene expression assays to 1.3-fold changes in expression. The exceptionally bright and photostable Alexa Fluor® dyes are also essentially insensitive to pH and are highly water soluble.

Biotin- or fluorescein-labeled nucleic acid probes are easily detected with streptavidin conjugates or labeled anti-fluorescein or anti-biotin antibodies, respectively, and have been used in two-color microarray assays, Southern and Northern blots, colony and plaque hybridizations, DNA sequencing, primer extension, DNA and RNA amplification, and bead-based separation techniques. In these applications, the labeled samples are generally detected with streptavidin- or antibody-conjugated enzymes in conjunction with chemiluminescent, fluorescent, or colorimetric substrates, such as those employed in our Tyramide Signal Amplification (TSATM) Kits. Please consult our *Handbook of Fluorescent Probes and Research Products* or visit our Web site (www.probes.com) for more information.

Table 1. Labeled 5-aminohexylacrylamido-dUTP and -dCTP nucleotides.

Label	Ex/Em*	Color	Catalog Number
Biotin	NA	Nonfluorescent	B32766, B32772
Fluorescein	495/525	Green fluorescent	F32767, F32773
Alexa Fluor® 555	555/570	Orange fluorescent	A32762, A32770
Alexa Fluor® 647	650/670	Far-red fluorescent	A32763, A32771

^{*} Excitation (Ex) and fluorescence emission (Em) maxima, in nm.

Materials

Each labeled 5-aminohexylacrylamido-dUTP nucleotide is supplied as a 1 mM nucleotide solution in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0). These nucleotides are >95% pure as determined by HPLC and spectrophotometric analysis. Nucleotides should be stored frozen at ≤−20°C and protected from light. When stored properly, these products will be stable for up to two years. AVOID REPEATED FREEZE/THAW CYCLES.

Spectral Properties and Filter Selection

For optimal performance of the nucleotides in your application, use a high-quality optical filter set that closely matches the spectral characteristics of the conjugates. Please note that these nucleotides may undergo small spectral shifts upon incorporation into polynucleotides. Table 1 lists the approximate excitation and emission maxima determined for each of the fluorescently labeled nucleotides in 50 mM phosphate buffer (pH 7.0). Full spectra for the labeled 5-aminohexylacrylamido-dUTP and -dCTP nucleotides are available at our Web site (www.probes.com).

Enzymatic Incorporation Protocol

We have optimized the following reverse transcription protocol using the SuperScriptTM Direct cDNA Labeling Kit and SuperscriptTM III reverse transcriptase (Invitrogen Corp.). Modification of this protocol may be required for use with other

systems; the most useful parameter to modify is the ratio of labeled aha-dNTP to dNTP (i.e., the ratio of aha-dUTP to dTTP or aha-dCTP to dCTP, depending on the aha-dNTP used in the reaction). For cDNA synthesis using labeled aha-dUTP, optimal results are obtained using 500 μ M final concentration of dGTP, dATP, and dCTP, and a 3:1 mixture of labeled aha-dUTP to unlabeled dTTP to give a final combined nucleotide concentration of 300 μ M (i.e. 200 μ M labeled aha-dUTP and 100 μ M dTTP). For cDNA synthesis using labeled aha-dCTP, optimal results are obtained using 500 μ M final concentration of dGTP, dATP, and dTTP, and a 3:1 mixture of labeled aha-dCTP to unlabeled dCTP to give a final combined nucleotide concentration of 300 μ M (i.e. 200 μ M labeled aha-dCTP and 100 μ M dCTP).

Materials Required

- Labeled aha-dNTP, 1 mM solution (provided)
- 10-40 µg of total RNA
- SuperScript[™] III Reverse transcriptase (400 U/μL)
- d(GACT) mixture containing:
 7.5 mM each of dGTP, dATP, dCTP, and 1.5 mM of dTTP (for reactions incorporating labeled aha-dUTP)

OR

7.5 mM each of dGTP, dATP, dTTP, and 1.5 mM of dCTP (for reactions incorporating labeled aha-dCTP)

- Anchored Oligo (dT)20 Primer, 2.5 mg/mL
- 5X first-strand buffer for reverse transcriptase (250 mM Tris, 375 mM KCl, 15 mM MgCl, pH 8.3)
- 100 mM DTT
- RNaseOUTTM (40 U/μL)

Labeling Reaction

1. Add dH_2O to the sample of RNA to bring the final volume to $8~\mu L$.

- 2. Add 2 µL of primer and mix thoroughly.
- **3.** Heat the template and primer mixture to 70°C for 10 minutes and then transfer immediately to ice.
- 4. To the tube containing the denatured template and primers, add:
- $6 \mu L 5X$ first-strand buffer
- 3 µL 100 mM DTT
- 2 μL d(GACT)
- 1 μL RNaseOUTTM
- 6 µL Labeled aha-dNTP (1 mM solution)
- 2 μL SuperScriptTM III
- **5.** Mix gently and incubate the reaction at 25°C for 10 minutes, followed by a three-hour incubation at 46°C.
- **6.** After cDNA synthesis, above, immediately perform the following hydrolysis reaction to degrade the RNA. Place the reverse transcription reaction at 95°C for 5 minutes to inactivate the reverse transcriptase and denature the RNA:cDNA hybrids. Snap cool by immediately placing the reaction vial into an ice bath. Add 15 μL of 0.1M NaOH to the reaction, mix thoroughly, and incubate at 65°C for 15 minutes. Add 15 μL of 0.1M HCl to neutralize the pH and mix gently.

Purification of RT-Labeled DNA

Purify the labeled DNA using a Qiagen QIAQuick™ PCR Purification Kit using the manufacturer's protocol, followed by ethanol precipitation.

References

1. Proc Natl Acad Sci U S A 64, 600 (1969); 2. Nature 223, 582 (1969); 3. Nature 265, 472 (1977); 4. Exp Cell Res 138, 485 (1980); 5. Proc Natl Acad Sci U S A 83, 2934 (1986); 6. Anal Biochem 331, 243 (2004); 7. J Histochem Cytochem 47, 1179 (1999); 8. Nonisotopic DNA Probe Techniques, Larry J. Kricka, Ed., Academic Press, Inc., San Diego, CA (1992).

Product List Current prices may be obtained from our Web site or from our Customer Service Department.

Cat #	Product Name	Unit Size
A32770	Alexa Fluor® 555-aha-dCTP *1 mM in TE buffer*	50 μl
A32762	Alexa Fluor® 555-aha-dUTP *1 mM in TE buffer*	50 μI
A32771	Alexa Fluor® 647-aha-dCTP *1 mM in TE buffer*	50 μl
A32763	Alexa Fluor® 647-aha-dUTP *1 mM in TE buffer*	50 μl
B32772	biotin-aha-dCTP *1 mM in TE buffer*	25 µl
B32766	biotin-aha-dUTP *1 mM in TE buffer*	25 µl
F32773	fluorescein-aha-dCTP *1 mM in TE buffer*	50 μl
F32767	fluorescein-aha-dUTP *1 mM in TE buffer*	50 μl

Contact Information

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Please visit our website — probes.invitrogen.com — for the most up-to-date information.

Molecular Probes, Inc.

29851 Willow Creek Road, Eugene, OR 97402 Phone: (541) 465-8300 • Fax: (541) 335-0504

Customer Service: 6:00 am to 4:30 pm (Pacific Time)

Phone: (541) 335-0338 • Fax: (541) 335-0305 • probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

Technical Service: 8:00 am to 4:00 pm (Pacific Time) Phone: (541) 335-0353 • Toll-Free (800) 438-2209 Fax: (541) 335-0238 • probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd. 3 Fountain Drive Inchinnan Business Park Paisley PA4 9RF, UK

Phone: +44 (0) 141 814 6100 • Fax: +44 (0) 141 814 6260

Email: euroinfo@invitrogen.com

Technical Services: eurotech@invitrogen.com

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License

For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full r

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2005, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.