Elite pH and pH Spear Tester Operating Instructions

68X001326 Rev.0

June, 2017

This instruction manual is applicable for Elite pH tester and Elite pH Spear tester.

Replacing the Batteries

The Elite pH tester uses four AAA 1.5 V batteries.

1. To remove the battery cover, See Figure 1. Clear the front catch and then the back catch, before sliding the cover off.
2. To remove the battery plate, push the center tab towards the front of the tester as show in Figure 2. Once unlocked, remove the plate to access the batteries.
3. Invert the tester upside down to remove the batteries. Each side uses two AAA batteries. Orient each battery with positive terminal facing downward.
4. To lock the battery plate, align the small tabs into the guide ribs on the housing and then press down. See Figure 4.


Figure 1: Removing battery cover

Figure 5: Battery plate unlock position

Figure 2: Push to unlock

Figure 3: Align tabs

Figure 4: Push down to lock


Figure 6: Battery plate lock position


Getting Started

The Elite pH tester has been factory calibrated and usually works well out of the box. However, after extended periods of non-use, it is best to ensure that the sensor is hydrated. Fill the cap with electrode storage solution or pH buffer and soak the sensor as needed to get a stable pH reading. A brief rinse with deionized (DI) water is OK, but avoid soaking or storing in deionized water as this will shorten pH electrode life. Prior to taking the measurements, periodic calibration with certified standards is recommended for best accuracy.

Buffer Set Selection


Select either USA (pH 4.01, pH 7.00 and pH 10.01) or NIST (pH 4.01, pH 6.86, and pH 9.18) to match your calibration standard values.

1. Press to power on.
2. Press to enter setup window. Press to select Settings. The display shows Buffer and Backlight.
3. Press to select Buffer. The display shows USA and NIST.
4. Press to select USA or scroll down by pressing to toggle between the two buffer standards.
5. The display shows the selected buffer with a .


Backlight Settings

1. Press to enter setup window. Press to select Settings. The display shows Buffer and Backlight.
2. Scroll down by pressing to toggle between Buffer and Backlight. Press to select Backlight. The display shows ON and OFF.
3. Scroll down by pressing to toggle between ON and OFF. Backlight ON increases readability in low light conditions.
4. Press to select the desired backlight option. The display shows the selected backlight option with a .

Temperature Settings


1. Press to enter setup window. Scroll down by pressing to select Temp Set. Press to select Temp set. The display shows Set °C/°F and Temp Cal.
2. Press to select Set °C/°F. The display automatically shows °C and °F.
3. Scroll down by pressing to toggle between °C and °F.
4. Press to select the desired temperature unit.

Temperature Calibration

1. Press to enter setup window. Scroll down by pressing to highlight Temp Set. Press to select Temp set. The display shows Set °C/°F and Temp Cal.
2. Scroll down by pressing to toggle between Set °C/°F and Temp Cal. Press to select Temp Cal.
3. The lower display shows the current measured temperature reading based on the last set offset and the upper display shows the current measured temperature reading based on factory default calibration.
4. Dip the tester into a solution of known temperature and allow time for the temperature sensor to stabilize.
5. Press to adjust the temperature value, then press to confirm the new temperature value of the solution.

Note: To exit this program without confirming the calibration, press .


pH Calibration

Calibration should be done regularly, preferably once a week or more if using the tester daily. You can calibrate up to three points using standards from either the USA or the NIST buffer groups. For best results, use a minimum of two pH buffers such as pH 4 and pH 7.

1. Press  to power the tester on if needed.
2. Dip sensor in about 2 cm to 3 cm into the pH standard buffer solution. If using the supplied cap, fill to the "MAX" line. For best results, allow the reading to stabilize if needed before proceeding.
3. Press the  to enter calibration mode. The CAL indicator will be displayed. The upper display will show the unadjusted factory default reading while the lower display will search for the corresponding value from the selected buffer group.

Note: All testers have dual display during calibration mode. To abort calibration, press .

4. The timer icon will blink while the reading stabilizes (allow 1 to 2 minutes). Once the reading is stabilized, the timer stops blinking. The tester will automatically confirm the calibration value and return to the measurement window when successful.
5. Rinse electrode and repeat with other buffers if necessary.


Measurement

1. Press  to power the tester on if needed.
2. Dip the sensor in about 2 cm to 3 cm into the sample. If using the supplied cap, fill to the "MAX" line. Stir and let the reading stabilize. The timer icon will blink during this time. Once the reading is stable,  will appear.

CAUTION: Testing dry samples is not recommended and can lead to sensor damage. Soils must be wet and free from particulates that may scratch the glass sensor. Excessive force into dry samples can cause glass breakage.

3. Note the pH value or press  to freeze the reading. To release the reading, press  again.
4. Press  for 3 seconds to turn off. If you do not press a button for 8.5 minutes, the tester will automatically shut off to conserve batteries.

User Reset

You can reset the pH calibration to the user's default settings by using the user reset function. Buffer selection and temperature user calibration are not affected by the user reset function.


1. Press  to enter setup window. Scroll down by pressing  to select Reset. Press  to select Reset. The display shows User Reset and Fact Reset.
2. Press  to select User Reset. The display automatically shows No and Yes.
3. Scroll down by pressing  to toggle between No and Yes.
4. Press  to confirm either No or Yes. The display shows the user reset option with a .

Factory Reset

You can reset the pH calibration to the factory reset default by using the factory reset function.

1. Press to enter setup window. Scroll down by pressing the to select Reset. Press to select Reset. The display shows User Reset and Fact Reset.
2. Scroll down by pressing the to toggle between the resets. Press to select Fact Reset. The display automatically shows No and Yes.
3. Scroll down by pressing to toggle between No and Yes.
4. Press to confirm either No or Yes. The display shows the factory reset option with a ✓.


MTS - Manual Temperature Setting (applicable to pH Spear tester only)

The MTS range is 0°C to 60.0°C (32.0°F to 140.0°F). User reset will set temperature to default value 25°C or 77°F.

1. Press to select setup window. Scroll down by pressing to select Temp Set. Press to select Temp set. The display shows Set °C/°F and Temp Cal.
2. Scroll down by pressing to toggle between Set °C/°F and Temp Cal. Press to select Temp Cal.
3. The lower display shows the current measured temperature reading based on the last set offset and the upper display shows the current measured temperature reading based on factory default calibration.
4. Dip the tester into a solution of known temperature and allow time for the in built temperature sensor to stabilize.
5. Press to adjust the temperature value or press the to confirm the calibrated value as new temperature value of the solution.



Note: To exit this program without confirming the calibration, press before the automatic confirmation takes place.

HOLD

This feature help to freeze the display reading for a delayed observation.

1. Press button to freeze the measurement.
2. Press again to release the measurement.

Sensor Maintenance


1. Before measuring soil pH with the pH Spear, wet your soil sample with distilled water and ensure that the soil is free of particulates. Even though the pH Spear is ruggedized, hard surfaces such as stones and pebbles can still cause breakage.
2. Rinse the sensor with sensor storage solution after each measurement. Care has to be taken not to damage the sensor's glass electrode especially while rinsing the pH Spear penetrating tip electrode.
3. In aggressive chemicals, dirty or viscous solutions, and solutions with heavy metals or proteins, take readings quickly and rinse electrode immediately afterward. For the pH Spear, the remnants of the semi solid samples on the penetrating electrode can be removed by rubbing it with some table salt and then rinsing. Mild detergent and warm water can be used to wash the penetrating electrode clean.
4. For best results, keep pH electrode stored upright in protective storage cap with electrode storage solution to keep the sensor hydrated. For pH Spear, keep tip moistened in sleeve with electrode storage solution. Never store pH and pH Spear in deionized water for extended periods. Use deionized water for rinsing only.

Sensor Replacement

You can replace the sensor module at the fraction of the cost of a new tester. When the tester fails to calibrate or gives fluctuating readings in calibration standards, you need to change the sensor.

1. With dry hands, grip the ring with sensor facing you. Twist the ring clockwise. Save the ring for later use.
2. Pull the old sensor module away from the tester.
3. Align the four tabs on the new module so that they match the four slots on the tester.
4. Gently push the module onto the slots to sit it in position. Push the smaller O-ring fully onto the new sensor module. Push the other O-ring over the module and thread it into place by firmly twisting counter clockwise.

Note: It is necessary that you recalibrate your tester prior to measurement after a sensor replacement.

Applications

The Elite pH tester can be used for the following:

- Aquariums
- Aquaculture
- Boilers
- Car washes
- Ecology
- Food sectors
- Hydroponics
- Labs
- Sanitation plants
- Spas
- Steam generators
- Studies
- Water and wastewater treatment
- Water quality testing in pools and more

Warranty

This instrument is supplied with a warranty against manufacturing defects for a period of one year from the date of purchase.

Return of Items

Authorization must be obtained from your distributor before returning items for any reason. When applying for authorization, please include information regarding the reason the item(s) are to be returned.

We reserve the right to make improvements in design, construction and appearance of products without notice. Prices are subject to change without notice.

Self Diagnostic Messages

| | |
|--|---|
|  | Batteries are weak and need replacement soon. |
| stable error | Appears when calibration is attempted but the reading is not yet stable. Wait for the reading to stabilize or manually confirm the calibration by pressing enter. |
| buffer error | The buffer is outside of the calibration range. |
| slope error | The 2 nd or 3 rd calibration point is not within 80% to 120% slope range. |
| over range | The reading is above the measuring range of tester. |
| under range | The reading is below the measuring range of tester. |

Accessories

| Ordering Code | Product Description |
|--|---|
| pH Tester | |
| ELITEPH | Elite pH tester with batteries |
| ELITEPHSPEAR | Elite pH Spear tester with batteries |
| pH Buffer Solutions and Sachets | |
| ECBU4BS | pH 4.01 (NIST traceable), box of 20 x 20 mL sachets |
| ECBU7BS | pH 7.00 (NIST traceable), box of 20 x 20 mL sachets |
| ECBU10BS | pH 10.01 (NIST traceable), box of 20 x 20 mL sachets |
| ECRINWT | pH deionized water (NIST traceable), box of 20 x 20 mL sachets |
| ECBU4710R | pH buffer sachet set, pH 4.01, 7.00, 10.01, deionized water, 5 x 20 mL each |
| Storage and Cleaning Solutions | |
| ECRE005 | Storage solution for pH sensor, 480 mL |
| ECDPCBT | Protein remover solution, 480 mL |
| Other Accessories | |
| ELITECAP | Replacement cap for Elite pH |
| TSLANYARD | Tester lanyard |

| Specification | Elite pH | Elite pH Spear |
|--------------------------------|---|--|
| pH | | |
| pH range | -1.00 to 15.0 pH | -1.00 to 15.0 pH |
| Resolution | 0.01 pH | 0.01 pH |
| Relative accuracy | ± 0.01 pH | ± 0.01 pH |
| Calibration points | Up to 3 points | Up to 3 points |
| Buffer set standard selection | USA 4.01/7.00/10.01 NIST 4.01/6.86/9.18 | USA 4.01/7.00/10.01 NIST 4.01/6.86/9.18 |
| Calibration window | ±1.00 pH | ±1.00 pH |
| Calibration type | Point to Point | Point to Point |
| Temperature | | |
| Temperature range | 0°C to 60°C 32.0°F to 140.0°F | 0°C to 60°C 32.0°F to 140.0°F |
| Temperature resolution | 0.1°C / 0.1°F | 0.1°C / 0.1°F |
| Temperature accuracy | from 0°C to 50°C (±0.5°C / ±0.9°F + 1 LSD); from 50°C to 60°C (±1.0°C / ±1.8 °F + 1 LSD) | NA |
| Temperature compensation | Yes (Automatic Temperature Compensation) | Yes (Manual Temperature Compensation) |
| Temperature calibration window | ± 5.0°C / ± 9.0°F | 0°C to 60°C 32.0°F to 140.0°F |
| General | | |
| Display | Graphics, Dot Matrix 80X100 Pixel | Graphis, Dot Matrix 80X100 Pixel |
| Auto off | 8.5 minutes (From Last Key press) | 8.5 minutes (From Last Key press) |
| Reset | User / Factory | User / Factory |
| Power requirement | Four AAA 1.5 V batteries | Four AAA 1.5 V batteries |
| Batthey life | >250 hours | >250 hours |
| Water proofing | IP67 | IP67 |
| Ambient operating temperature | 5°C to 45°C | 5°C to 45°C |
| Relative humidity | 5% to 85% Non-Condensing | 5% to 85 % Non-Condensing |
| Storage temperature | -20°C to 60°C | -20°C to 60°C |
| Storage humidity | 5% to 85% Non-Condensing | 5% to 85% Non-Condensing |

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