

**MOLECULAR PROBES°** 

# PRODUCT INSERT

## STREPTAVIDIN CONJUGATES

Product code	Form	Volume	Excitation (nm)	Peak Emission (nm)
SA1001	FITC	1.0 ml	488	525
SA1010	Cy3	0.5 ml	488	565
SA1004-1	R-PE	0.5 ml	488	575
SA1004-4	R-PE	2.0 ml		
SA1017	$PE-TR^{\dagger}$	0.5 ml	488	615
SA1006	$TC^{\ddagger}$	0.5 ml	488	670
SA1018	PE-Cy5.5	0.5 ml	488	694
SA1012	PE-Cy7	0.5 ml	488	767
SA1005	APC	0.5 ml	600-650	660
SA1011	Cy5	0.5 ml	600-650	670
SA1019	APC-Cy5.5	0.5 ml	600-650	694
SA1014	APC-Cy7	0.5 ml	600-650	767
SA1027	APC-Alexa Fluor <sup>®</sup> 750	1.0 ml	600-650	775

### PRODUCT DESCRIPTION

Streptavidin

Clone: N/A

Isotype: N/A

Lot No.: See label Expiration: See label

Concentration: See label

Buffer: Phosphate buffered saline (PBS)

**Preservative:** 0.1% *sodium azide*. Sodium azide is an extremely toxic and dangerous compound, particularly when combined with acids or metals. Solutions containing sodium azide should be disposed of properly.

**Stabilizer:** A highly purified grade of BSA has been added as a stabilizing agent.

### **STORAGE & HANDLING**

Store reagents at 2-8°C. Light exposure should be avoided. Use dim light during handling, incubation with cells and prior to analysis. It is recommended that cells be analyzed within 18 hours of staining. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

## PRODUCT CHARACTERIZATION

Streptavidin binds to biotin with high affinity  $(10^{-15} \text{ M})$ . When conjugated to fluorochromes, streptavidin is a potent second step reagent used to detect primary, biotinylated antibodies.

### PRODUCT QUALITY CONTROL

Every lot of fluorochrome-conjugated streptavidin is tested by flow cytometry using biotinylated primary antibodies. From this testing it is recommended that between 0.05 and 0.2  $\mu$ g of streptavidin be used per 1 x 10<sup>6</sup> cells in a 100  $\mu$ l staining volume. Because conditions may vary, each investigator should determine the optimal amount of streptavidin to be used for each application.

- <sup>†</sup> TR, Texas Red<sup>®</sup>
- <sup>‡</sup> TC, TRI-COLOR<sup>®</sup>, PE-Cy<sup>®</sup>5

The efficiency of energy transfer in tandem dyes can be significantly decreased by exposure to visible light. We recommend that longer wavelength fluorochrome conjugates, e.g. PE-Cy<sup>®</sup>7, PE-Alexa Fluor<sup>®</sup> 700, be protected from light during staining and while awaiting analysis, e.g. cover with aluminum foil.

The Texas Red<sup>®</sup>, Alexa Fluor<sup>®</sup> and Pacific Blue<sup>TM</sup> dye conjugates in this product are sold under license from Molecular Probes, Inc., for research use only or as analyte specific reagents, except for use in combination with microarrays or high content screening, and are covered by pending and issued patents.

Explanation of symbols					
Symbol	Description	Symbol	Description		
REF	Catalogue Number	LOT	Batch code		
RUO	Research Use Only	IVD	In vitro diagnostic medical device		
X	Use by	ł	Temperature limitation		
***	Manufacturer	EC REP	European Community authorised representative		
[-]	Without, does not contain	[+]	With, contains		
even Light	Protect from light	Â	Consult accompanying documents		
[]i	Directs the user to consult instructions for use (IFU), accompanying the product.				

(Rev 06/10) DCC-10-1766

Cy<sup>®</sup> is a trademark of GE/Amersham Biosciences.

For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

PI L11017 (Page 1 of 1)

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, <u>www.invitrogen.com</u>). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.