

Pierce Chromatin Prep Module

Pub. No. MAN0016151
Pub. Part No. 2162217.6
Rev. A.0**26158****Number****Description**

26158

Pierce Chromatin Prep Module**Kit Contents:****Membrane Extraction Buffer**, 10mL; store at 4°C**Nuclear Extraction Buffer**, 10mL; store at 4°C**MNase Digestion Buffer**, 5mL; store at 4°C**MNase Stop Solution**, 0.5mL; store at 4°C**Halt™ Protease and Phosphatase Inhibitor Cocktail EDTA-free (100X)**, 4 × 90µL microtubes; store at 4°C**Glycine Solution (10X)**, 15mL; store at 4°C**PBS (20X)**, 15mL; store at 4°C**Micrococcal Nuclease (ChIP Grade) (10U/µL)**, 25µL; store at -20°C**Proteinase K**, (20mg/mL), 0.25mL; store at -20°C**DTT, Lyophilized**, 1 vial; store at room temperature or 4°C**Storage:** Upon receipt store contents of yellow pouch at -20°C. Store all other components at 4°C. The kit is shipped at ambient temperature.**Introduction**

The Thermo Scientific™ Pierce™ Chromatin Prep Module provides a simple, fast and reproducible method to isolate chromatin for use in assays that monitor transcription regulation through histone modification¹ (epigenetics) or transcription factor-DNA binding interactions, such as chromatin immunoprecipitation (ChIP). The module contains all reagents necessary to isolate crosslinked chromatin from mammalian cells. The specially titrated and tested micrococcal nuclease digests the DNA, eliminating variable results caused by the traditional method of sonication. The advantages of enzymatic digestion include digestion reproducibility, reaction control, and easy titration of the enzyme for specific cell types. The included lysis reagents extract and solubilize the crosslinked DNA:protein complexes, providing simple and convenient isolation of chromatin-bound DNA without a dounce homogenizer. Because there is less than 15% contamination from other cellular compartments, the protein of interest can be efficiently enriched in downstream applications.

Additional Materials Required

- Mammalian cell culture reagents and equipment
- Chemical fume hood
- Cell scrapers
- Thermomixer or heat block
- 50 and 15mL conical tubes
- 16% Formaldehyde (Product No. 28906)
- DMSO (optional)
- Sodium chloride (5M)
- Nuclease-free water

Protocol for Chromatin Prep Module

Note: Read the entire protocol before beginning the assay.

A. Experimental Design

- When designing a ChIP experiment, consider the number of immunoprecipitations desired from the chromatin sample and the number of cell culture conditions (e.g., drug-treated vs. non-treated). For accurate comparison, immunoprecipitate the same target proteins and controls from each cell culture condition.
- Performing a ChIP using normal rabbit IgG is an effective negative control. An abundant and ubiquitous DNA-binding protein, such as histone H3, can be used as a positive control.
- If performing multiple ChIPs from a single cell culture condition, crosslinked chromatin may be prepared in bulk. Scale the reagent amounts according to the number of ChIPs desired (Table 1) and sub-divide the chromatin sample after diluting with the IP Wash Buffer. The total input control obtained from the cell culture condition is common to this set of ChIPs.

Table 1. Reagent volumes to use based on the number of ChIPs.

Example number of HeLa cells	2×10^6	6×10^6	2×10^7	6×10^7
Number of ChIPs	1	3	10	30
1X PBS	2.1mL	6.3mL	21mL	63mL
Lysis Buffer 1	0.1mL	0.3mL	1mL	3mL
Lysis Buffer 2	0.05mL	0.15mL	0.5mL	1.5mL
MNase Digestion Buffer Working Solution	0.1mL	0.3mL	1mL	3mL
1X IP Dilution Buffer	0.45mL	1.2mL	4.5mL	12mL

B. Material Preparation

Note: The amounts listed below are for one ChIP. If you are preparing chromatin in bulk or performing multiple ChIPs, multiply the reagent amounts by the number of ChIPs being performed.

Glycine Solution (10X)	If precipitate formed, warm solution to 37°C in a water bath for 30 minutes and vortex before use.
PBS	If precipitate formed, warm solution to 37°C in a water bath and vortex before use. Dilute PBS (20X) to 1X with nuclease-free water.
1M DTT	Add 50µL of ultrapure water to the DTT. Store this solution at -20°C for up to 6 months.
Lysis Buffer 1	Add 1µL of the Halt Cocktail to 100µL of Membrane Extraction Buffer in a microcentrifuge tube and place on ice.
Lysis Buffer 2	Add 0.5µL of the Halt Cocktail to 50µL of Nuclear Extraction Buffer in a microcentrifuge tube and place on ice.

MNase Digestion Buffer Add 0.1 μ L of 1M DTT to 100 μ L of 1X MNase Digestion Buffer and place tube at room temperature.
Working Solution

C. Crosslinking and Cell Pellet Isolation

Note 1: Perform all centrifugations at room temperature.

Note 2: This procedure is for one CHIP. Multiply reagent amounts by the number of CHIPs being performed (see Table 1).

1. Culture adherent mammalian cells and treat as desired.

Optional: If you are unfamiliar with cell type being used, culture an extra dish of cells for determining cell number. Before crosslinking, trypsinize and determine the cell number of the extra dish of cells.

2. To each dish containing cell culture media, add sufficient quantity of 16% formaldehyde to obtain a final concentration of 1% formaldehyde.
Caution: Formaldehyde is a skin irritant and the fumes are toxic. Use proper personal protective, laboratory safety and disposal equipment.
3. Mix well by gently swirling the dish. Incubate at room temperature for 10 minutes in a chemical fume hood.
4. To each dish containing cell culture media and formaldehyde, add Glycine Solution (10X) to a final concentration of 1X. Mix well by gently swirling the dish. Incubate at room temperature for 5 minutes in the chemical fume hood.
5. Aspirate formaldehyde/glycine-containing media in the fume hood. Properly dispose the formaldehyde-containing waste.
6. Wash the cells twice with one media volume of ice-cold PBS, removing each wash by aspiration.
7. Add 10 μ L of the Halt Cocktail to 1mL of ice-cold PBS. Add the solution to the cells, and detach cells by scraping. Transfer the cell suspension to a 1.5mL microcentrifuge tube using a pipette.
8. Centrifuge tubes at 3000 \times g for 5 minutes.
9. Remove the PBS. Store the cell pellet(s) at -80°C, or proceed directly to Section C: Lysis and MNase Digestion.

D. Lysis and MNase Digestion

Note: For best results, empirically determine the optimal crosslinking time and Micrococcal Nuclease digestion conditions for each cell type (see Appendix A).

1. Use the crosslinked HeLa cells prepared above. Thaw stored cells on ice.
2. Add 100 μ L Lysis Buffer 1 containing protease inhibitors to the cell pellet and pipette up and down to break up the pellet. Vortex the tube for 15 seconds and incubate on ice for 10 minutes.
3. Centrifuge at 9000 \times g for 3 minutes and remove the supernatant.
4. Resuspend nuclei in 100 μ L of MNase Digestion Buffer Working Solution.
5. Add 0.25 μ L of Micrococcal Nuclease (CHIP Grade) (10U/ μ L), vortex the tube and incubate in a 37°C water bath for 15 minutes, mixing by inversion every 5 minutes.
6. Add 10 μ L MNase Stop Solution to stop the reaction, vortex briefly and incubate on ice for 5 minutes.
7. Centrifuge at 9000 \times g for 5 minutes to recover the nuclei. Remove the supernatant.
8. Resuspend nuclei in 50 μ L of Lysis Buffer 2 containing protease/phosphatase inhibitors and incubate on ice for 15 minutes, vortexing for 15 seconds every 5 minutes.
9. Centrifuge at 9000 \times g for 5 minutes and transfer the supernatant, containing the digested chromatin, to a new 1.5mL tube. Proceed to the immunoprecipitation or store samples at -80°C.

Troubleshooting

Problem	Possible Cause	Solution
Chromatin fragments are too large (> 1000 base pairs)	Crosslinking time is too long	Shorten crosslinking time
	Cell to Micrococcal Nuclease (MNase) ratio is too high	Increase amount of MNase or decrease cell number (refer to the MNase digestion optimization protocol in Appendix A)
Chromatin fragments are too small (< 200 base pairs)	Cell to MNase ratio is too low	Decrease amount of MNase or increase cell number (refer to the MNase digestion optimization protocol in Appendix A)

Appendix A: Optimization of micrococcal nuclease digestion

1. Prepare crosslinked cell pellet from 10^7 cells as described in Section C.
2. Follow the lysis procedure in Section D steps 1-4; however, use 0.5mL of the Lysis Buffer 1, Membrane Extraction Buffer and MNase Digestion Buffer Working Solution.
3. Transfer 100 μ L of the resuspended nuclei into five tubes (2×10^6 nuclei/tube), labeling the tubes 0, 1, 2, 4 and 6.
4. Add 2 μ L of the Micrococcal Nuclease (ChIP Grade) to 18 μ L of MNase Digestion Buffer (1:10 dilution).
5. Add 0, 1, 2, 4 and 6 μ L of the diluted Micrococcal Nuclease to the corresponding tubes of nuclei and incubate for 15 minutes in a 37°C water bath mixing by inversion every 5 minutes.
6. Add 20 μ L of MNase Stop Solution to each tube and incubate on ice for 5 minutes.
7. Centrifuge at $9000 \times g$ for 5 minutes to recover the nuclei and remove the supernatant.
8. Resuspend nuclei in 50 μ L of Lysis Buffer 2 and incubate on ice for 15 minutes, vortexing for 15 seconds every 5 minutes.
9. Centrifuge at $9000 \times g$ for 5 minutes.
10. Prepare five 1.5mL microcentrifuge tubes labeled 0, 2, 4 and 6 each containing 6.6 μ L Nuclease Free Water, 2.4 μ L of 5M NaCl, and 1 μ L Proteinase K (20mg/mL).
11. Transfer 50 μ L of each nuclear supernatant, containing digested chromatin, to the corresponding tube.
12. Vortex each tube for 10 seconds and incubate the tubes at 65°C in a heat block for 1.5 hours.
13. Analyze 15-25 μ L of each sample and DNA size marker by agarose gel electrophoresis. Optimal digestion should yield fragments from 200 to 1000 base pairs with a more intense ladder of bands occurring at approximately 160, 320 and 480 base pairs, which corresponds to the 1, 2, and 3 nucleosome units.

Related Thermo Scientific Products

See our website for a listing of ChIP-validated antibodies.

28908	16% Formaldehyde (w/v), Methanol-free, 10 \times 1mL
20688	Dimethylsulfoxide (DMSO), Sequanal grade, 950mL
26158	Pierce Chromatin Prep Module
26159	ChIP Grade Protein A/G Plus, 0.65mL
26160	Proteinase K Solution (20mg/mL), 0.25mL
28372	BupH™ Phosphate Buffered Saline Pack, 40 packs
78425	Halt Protease Inhibitor Cocktail, EDTA-Free (100X), 24 vials
78443	Halt Protease and Phosphatase Inhibitor Cocktail, EDTA-Free (100X), 24 vials
78833	NE-PER Nuclear and Cytoplasmic Extraction Kit
20148	LightShift™ Chemiluminescent EMSA Kit
23225	Pierce BCA Protein Assay Kit
22662	Pierce 660nm Protein Assay Kit

AB1219	Absolute™ QPCR SYBR™ Green Fluorescein Mix
AB-0900/w	Thermo-Fast™ 96 Semi-Skirted PCR Plate, opaque white
AB-1170	ABsolute QPCR Seal

Cited References

1. Lee, D.Y., *et al.* (2009). Role of protein methylation in regulation of transcription. *Endocrine Rev* **26(2)**:147-70.

General References

- Aparicio, O., *et al.* (2004). Chromatin immunoprecipitation for determining the association of proteins with specific genomic sequences *in vivo*. *Curr Protoc Cell Biol* Chapter 17:Unit 17.7.
- Das, P.M., *et al.* (2004). Chromatin immunoprecipitation assay. *Biotechniques*. **6**:961-9.
- Wells, J. and Farnham, P.J. (2002). Characterizing transcription factor binding sites using formaldehyde crosslinking and immunoprecipitation. *Methods* **1**:48-56.
- Yuichi Y. and Hughes D.E., (2006). C/EBP α and HNF6 protein complex formation stimulates HNF6-dependent transcription by CBP coactivator recruitment in HepG2 cells. *Hepatology* **43(2)**:276-86.
- Zeng, P.Y., *et al.* (2006). *In vivo* dual cross-linking for identification of indirect DNA-associated proteins by chromatin immunoprecipitation. *Biotechniques* **41**:694-8.

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, *in vitro* diagnostic uses, *ex vivo* or *in vivo* therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2013 Thermo Fisher Scientific Inc. All rights reserved. SYBR is a trademark of Life Technologies Corporation. All (other) trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.