

# VIASURE MULTIPLEX

## *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit

### Pathogen and product description

**C***ampylobacter* species are gram-negative, nonspore forming, spiral, or curved-shaped bacteria. Among the more than 26 species currently classified in the *Campylobacter* genus, most human diseases are attributed to three major food-borne species: *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari*. These three species can be isolated from poultry and are of greatest concern to the poultry industry.

*Campylobacter* is considered one of the most common causes of diarrheal illness worldwide. Poultry is a major reservoir and source of transmission of *Campylobacter* to humans. In particular *C. jejuni* is the common species found in beef, *C. coli* is often isolated from pork, and *C. lari* is predominant in shorebirds. Other risk factors include consumption of animal products and water, contact with animals, and even person-to-person transmission (fecal-oral or via fomites).

Infection with *Campylobacter* causes gastroenteritis characterised by fever, vomiting, headaches, and abdominal pain with watery or bloody diarrhea, for a median duration of 6 days. Besides gastroenteritis, these three species can cause periodontitis, septicemia, and second trimester intrauterine growth restriction. Furthermore, *C. jejuni* infection may lead to autoimmune conditions such as Guillain-Barré syndrome (GBS) and Miller Fisher syndrome (MFS).

*Campylobacter* enteritis is usually self-limiting and typically does not require antimicrobial therapy. In

these cases, maintenance of proper hydration and electrolyte balance is the most important tenets of treatment. However, in severe and prolonged cases of enteritis, bacteremia, or other extraintestinal infection, prompt antimicrobial treatment is indicated.

Traditional microbiological methods for *Campylobacter* identification include enrichment, culturing, isolation, and phenotypic characterization. The procedures are laborintensive, time consuming, and with a relatively narrow differentiation spectrum among target species. These factors present challenges for the identification of *Campylobacter* from patient samples or contaminated food. Fortunately, methodologies based on molecular biology have been developed to improve laboratory approaches, such as Real Time PCR. Multiplex qPCR can detect several targets in one reaction, saving time, effort, and sample. This qPCR assay can be used for identifying *C. coli*, *C. lari* and *C. jejuni* or secondary screening for confirmation of *Campylobacter* bacteria to the species level.

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit is designed for the diagnosis of *Campylobacter coli*, *C. lari* and/or *C. jejuni* in clinical samples. After DNA isolation, the identification of *C. coli*, *C. lari* and *C. jejuni* is performed by the amplification of a conserved region of the *hip0* gene for *Campylobacter jejuni*, *Gyrasa A* gene for *Campylobacter lari* and *CeuE* gene for *Campylobacter coli*, using specific primers and a fluorescent-labelled probe.



## Analytical sensitivity

VIASURE *Campylobacter coli*, *C. lari* & *C. jejuni* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *Campylobacter coli*, *Campylobacter lari* and *Campylobacter jejuni* (Figures 1, 2 and 3).




Figure 1. Dilution series of *Campylobacter jejuni* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).




Figure 2. Dilution series of *Campylobacter lari* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).




Figure 3. Dilution series of *Campylobacter coli* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5).

## Components


Reagent/Material	Description	Colour	Quantity
<b>Campylobacter coli, C. lari &amp; C. jejuni 8-well strips</b>	A mix of enzymes, primers-probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
<b>Rehydration Buffer</b>	Solution to reconstitute the stabilized product	Blue	1 vial x 1,8 mL
<b>Campylobacter coli, C. lari &amp; C. jejuni Positive Control</b>	Non-infectious synthetic lyophilized cDNA	Red	1 vial
<b>Extraction Control</b>	Non-infectious nucleic acid lyophilized	Green	1 vial
<b>Negative Control</b>	Non template control	Violet	1 vial x 1 mL
<b>Water RNase/DNase free</b>	RNase/DNase free water	White	1 vial x 1 mL
<b>Tear-off 8-cap strips</b>	Optical caps for sealing wells during thermal cycling	Transparent	6/12 x 8-cap strip

## Kit References

Reference	Description
VS-CLJ106L	VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile
VS-CLJ106H	VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile
VS-CLJ112L	VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile
VS-CLJ112H	VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile
VS-CLJ113L	VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit 96-well plate, low profile
VS-CLJ113H	VIASURE <i>Campylobacter coli</i> , <i>C. lari</i> & <i>C. jejuni</i> Real Time PCR Detection Kit 96-well plate, high profile

## Work Flow

One-step rehydration of wells and add your extracted DNA



CERTEST BIOTEC, S.L.  
Pol. Industrial Río Gállego II, Calle J, Nº 1,  
50840, San Mateo de Gállego, Zaragoza (SPAIN)  
www.certest.es

