eBioscience™ Super Bright Complete Staining Buffer

Catalog Numbers SB-4401-42 and SB-4401-75

Pub. No. MAN0018634 Rev. B.0



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **thermofisher.com/support**.

Product description

eBioscience Super Bright Complete Staining Buffer is designed for use as a supplement to Flow Cytometry Staining Buffer in immunofluorescent staining protocols of cells in suspension. eBioscience Super Bright Complete Staining Buffer is only necessary when using more than one polymer dye-conjugated antibody in the same sample to prevent nonspecific polymer interactions, which can result in data appearing under-compensated. eBioscience Super Bright Complete Staining Buffer is provided in a convenient 5 μ L/test format and is compatible with traditional fluorochromes, Brilliant Violet dyes, and standard flow cytometry protocols.

eBioscience [™] Super Bright Complete Staining Buffer is a fully compatible replacement for Super Bright Staining Buffer (Cat. Nos. SB-4400-42, SB-4400-75). Additionally, it shows decreased background in the 450/40 channel on the violet laser.

Product specifications

Concentration ^[1]	5 μL/test
Storage	Store at 2–8°C.
Application	Flow cytometry.
Testing	This product is tested by flow cytometric analysis of normal human peripheral blood cells or mouse splenocytes.
Batch code	See product label.
Use by	See product label.
Related product	Flow Cytometry Staining Buffer (Product No. 00-4222).

 $^{^{[1]}\,}$ A test is defined as the amount of buffer to be used in a final volume of 100 $\mu L.$

Important product information

- eBioscience Super Bright Complete Staining Buffer is not compatible with UltraComp eBeads Compensation Beads (Cat. No. 01-2222). If using UltraComp eBeads Compensation Beads as a compensation tool, solely use Flow Cytometry Stain Buffer (Cat. No. 00-4222) for any antibody dilutions.
- eBioscience [™] Super Bright Complete Staining Buffer is provided in a convenient 5 μL/test format.
- eBioscience Super Bright Complete Staining Buffer is compatible with traditional fluorochromes and Live/Dead and Fixable Viability eFluor dyes.
- Staining in the presence of eBioscience Super Bright Complete Staining Buffer will result in less background signal in the 450/50 bandpass detector on the violet laser (used to collect signal for eFluor 450 and Super Bright 436) than Super Bright Staining Buffer (Cat. No. SB-4400).
- eBioscience Super Bright Complete Staining Buffer is compatible with RBC lysis protocols, such as 1-step Fix/Lyse (Cat. No. 00-5333) and 10X RBC Lysis Buffer (multi-species) (Cat. No. 00-4300).
- eBioscience Super Bright Complete Staining Buffer can also be used at the appropriate test concentration when preparing bulk (multi-test) antibody cocktails.



Workflow

Materials required

- eBioscience [™] Super Bright Complete Staining Buffer (Cat. No. SB-4401-42)
- 12 × 75 mm round-bottom test tubes
- Flow Cytometry Staining Buffer (Cat. No. 00-4222)
- Primary antibodies (directly conjugated to fluorochromes)

Procedure

- 1. Add 5 μL of Super Bright Complete Staining Buffer to each tube. Staining buffer can be added directly to tubes or to previously aliquoted cells in tubes. If adding to cells, mix well by pipetting up and down or gently vortexing the sample.
- 2. Add appropriate amounts of each fluorochrome-conjugated antibody, including Super Bright and traditional fluorochrome-conjugated antibodies, to the tubes containing Super Bright Complete Staining Buffer.
- 3. Mix well after addition of each antibody by pipetting up and down or gently vortexing the sample.
 - Note: If a cocktail of antibodies is prepared in bulk, it should be used fresh to minimize nonspecific polymer dye interactions.
- 4. If cells were not previously added to the tubes, aliquot $100 \mu L$ of cells to the buffer-antibody cocktail promptly.
- 5. Mix samples well by pipetting up and down or gently vortexing.
- 6. Incubate for 30 minutes in the dark at 2-8°C.
- 7. Wash the cells by adding 2 mL/tube of Flow Cytometry Staining Buffer. Centrifuge at 400-600 × g for 5 minutes. Discard supernatant.
- 8. Repeat step 7.
- 9. Resuspend cells in an appropriate volume of Flow Cytometry Staining Buffer.
- 10. Analyze samples by flow cytometry or, if staining for intracellular targets, proceed with *Best Protocols: Staining Intracellular Antigens for Flow Cytometry* (available on our website).

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Life Technologies Corporation | 5781 Van Allen Way | Carlsbad, CA 92008 For descriptions of symbols on product labels or product documents, go to **thermofisher.com/symbols-definition**.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses. Super Bright Polymer Dyes are sold under license from Becton, Dickinson and Company.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.



23 April 2019