

# pTrcHis and pTrcHis2 TOPO<sup>®</sup> TA Expression Kits

# Five-minute cloning of *Taq* polymerase-amplified PCR products for expression in *E. coli*

Catalog no. K4410-01 (pTrcHis TOPO®)

Catalog nos. K4400-01 and K4400-40 (pTrcHis2 TOPO<sup>®</sup>)

**Version J** 23 February 2006 25-0209

A Limited Label License covers this product (see Purchaser Notification). By use of this product, you accept the terms and conditions of the Limited Label License.

**User Manual** 

### **Table of Contents**

Table of Contents	iii
Kit Contents and Storage	iv
Methods	1
Overview	1
Designing PCR Primers for pTrcHis-TOPO <sup>®</sup>	6
Designing PCR Primers for pTrcHis2-TOPO <sup>®</sup>	8
Producing PCR Products	
TOPO <sup>®</sup> Cloning and Transformation	11
Expression of the PCR Product	
Troubleshooting	
Appendix	19
Recipes	
Purifying PCR Products	21
Addition of 3' A-Overhangs Post-Amplification	
TOPO TA Cloning <sup>®</sup> Control Reactions	
Map and Features of pTrcHis-TOPO <sup>®</sup>	
Map and Features of pTrcHis2-TOPO <sup>®</sup>	
pTrcHis-TOPO <sup>®</sup> /lacZ	
pTrcHis2-TOPO <sup>®</sup> /lacZ	
Technical Service	
Purchaser Notification	
Qualifying the Product	
References	

#### **Kit Contents and Storage**

The pTrcHis and pTrcHis2 TOPO<sup>®</sup> TA Expression Kits are shipped on dry ice. Each kit Shipping and contains a box with TOPO TA Cloning® reagents (Box 1) and a box with TOP10 One Storage Shot<sup>®</sup> competent cells (Box 2). Store Box 1 at -20°C and Box 2 at -80°C. The pTrcHis TOPO<sup>®</sup> TA Expression Kit contains the pTrcHis-TOPO<sup>®</sup> vector, which pTrcHis and allows you to clone in frame with an N-terminal tag. The pTrcHis2 TOPO® TA pTrcHis2 TOPO<sup>®</sup> Expression Kit contains the pTrcHis2-TOPO® expression vector, which allows you to **TA Expression** clone in frame with a C-terminal tag. See the table below for ordering information. Kits Product **Pack Size** Catalog no. pTrcHis TOPO<sup>®</sup> TA Expression Kit K4410-01 20 (containing pTrcHis-TOPO<sup>®</sup> vector) pTrcHis2 TOPO® TA Expression Kit 20 K4400-01 (containing the pTrcHis2-TOPO<sup>®</sup> vector) K4400-40 40

#### TOPO TA Cloning<sup>®</sup> Reagents

The TOPO TA Cloning<sup>®</sup> reagents (Box 1) for both kits are listed below. **Please note that the user must supply** *Taq* **polymerase.** Store Box 1 at -20°C.

Item	Concentration	Amount
pTrcHis-TOPO <sup>®</sup> OR	10 ng/µl plasmid DNA in:	1 tube of 25
pTrcHis2-TOPO <sup>®</sup> vector	50% glycerol	μl
	50 mM Tris-HCl, pH 7.4 (at 25°C)	
	1 mM EDTA	
	2 mM DTT	
	0.1% Triton X-100	
	100 µg/ml BSA	
	phenol red	
10X PCR Buffer	100 mM Tris-HCl, pH 8.3 (at 42°C)	100 µl
	500 mM KCl	
	25 mM MgCl <sub>2</sub>	
	0.01% gelatin	
dNTP Mix	12.5 mM dATP; 12.5 mM dCTP; 12.5 mM dGTP; 12.5 mM dTTP	10 µl
	neutralized at pH 8.0 in water	

#### Kit Contents and Storage, continued

Item	Concentration	Amount
1 M IPTG	1 M IPTG in sterile water	1 ml
Xpress <sup>™</sup> Forward Sequencing Primer	0.1 μg/µl in TE Buffer	20 µ1
(for pTrcHis-TOPO <sup>®</sup> )		
pTrcHis Forward Sequencing Primer	0.1 µg/µl in TE Buffer	20 µ1
(for pTrcHis2-TOPO <sup>®</sup> )		
pTrcHis Reverse Sequencing Primer	0.1 µg/µl in TE Buffer	20 µ1
(for both vectors)		
Control PCR Primers	0.1 µg/µl each in TE Buffer	10 µ1
Control PCR Template	0.1 µg/µl in TE Buffer	10 µ1
Sterile Water		1 ml
Expression Control Plasmid	0.5 µg/µl in TE Buffer	10 µ1
(pTrcHis-TOPO <sup>®</sup> / <i>lacZ</i> or pTrcHis2- TOPO <sup>®</sup> / <i>lacZ</i> )		

#### pTrcHis TOPO TA Cloning<sup>®</sup> Reagents, continued

#### Sequences of Primers

The table below provides the sequences and pmoles supplied of the Forward and Reverse sequencing primers. Two micrograms of each primer are supplied.

Primer	Sequence	pMoles Supplied
Xpress <sup>™</sup> Forward	5'-TATGGCTAGCATGACTGGT-3'	342
pTrcHis Forward	5'-GAGGTATATATTAATGTATCG-3'	309
pTrcHis Reverse	5'-GATTTAATCTGTATCAGG-3'	363

#### Kit Contents and Storage, continued

#### One Shot<sup>®</sup> Reagents

The table below describes the items included in the One Shot® competent cell kit. Store at -80°C.

Item	Composition	Amount
SOC Medium	2% Tryptone	6 ml
(may be stored at room	0.5% Yeast Extract	
temperature or $+4^{\circ}C$ )	10 mM NaCl	
	2.5 mM KCl	
	10 mM MgCl <sub>2</sub>	
	10 mM MgSO <sub>4</sub>	
	20 mM glucose	
TOP10 cells		21 x 50 µl
pUC19 Control DNA	10 pg/µl in 5 mM Tris-HCl, 0.5 mM EDTA, pH 8.0	50 µl

#### Genotype

**TOP10**: Use this strain for general cloning and expression of PCR products in pTrcHis-TOPO<sup>®</sup> or pTrcHis2-TOPO<sup>®</sup>. Please note that this strain cannot be used for single-strand rescue of DNA.

 $\label{eq:star} F^{-}mcrA\ \Delta(mrr-hsdRMS-mcrBC)\ \Phi 80lacZ\Delta M15\ \Delta lacX74\ recA1\ araD139\ \Delta(araleu)7697\ galU\ galK\ rpsL\ (Str^R)\ endA1\ nupG$ 

### Methods

Overview	
Introduction	The pTrcHis and pTrcHis2 TOPO <sup>®</sup> TA Expression Kits provide a highly efficient, rapid cloning strategy ("TOPO <sup>®</sup> Cloning") for direct insertion of <i>Taq</i> polymerase-amplified PCR products into a plasmid vector for expression in <i>E. coli</i> . No ligase, post-PCR procedures, or PCR primers containing specific sequences are required. Once the PCR product of interest is TOPO <sup>®</sup> Cloned into either pTrcHis-TOPO <sup>®</sup> or pTrcHis2-TOPO <sup>®</sup> , the construct is transformed into <i>E. coli</i> and expression induced with isopropyl $\beta$ -thiogalactoside (IPTG).
Description of the Vectors	pTrcHis-TOPO <sup>®</sup> and pTrcHis2-TOPO <sup>®</sup> are designed to facilitate expression of eukaryotic proteins in <i>E. coli</i> . Both vectors contain the following:
	• The <i>trc</i> promoter, a hybrid promoter containing the -35 region from the <i>trpB</i> promoter and the -10 region from the <i>lacUV5</i> promoter for high-level expression in <i>E. coli</i> (Brosius <i>et al.</i> , 1985; Egon <i>et al.</i> , 1983; Mulligan <i>et al.</i> , 1985).
	• The <i>lacO</i> sequence for binding the Lac repressor encoded by the <i>lacI<sup>q</sup></i> gene. In the absence of IPTG, Lac repressor binds to the <i>lacO</i> sequence, repressing transcription. Upon addition of IPTG, expression is induced (Jacob and Monod, 1961; Müller-Hill <i>et al.</i> , 1968).
	• <i>rrnB</i> antitermination sequence that reduces premature transcription termination (Li <i>et al.</i> , 1984).
	• T7 gene 10 translational enhancer sequence for more efficient translational initiation (Olins <i>et al.</i> , 1988).
	• A minicistron containing nucleotides that are efficiently translated in prokaryotic cells for enhanced translational efficiency (Schoner <i>et al.</i> , 1986).
pTrcHis-TOPO <sup>®</sup>	In addition to the features above, pTrcHis-TOPO <sup>®</sup> contains the following additional elements:
	• An N-terminal peptide containing the HisG epitope, the Xpress <sup>™</sup> epitope and a 6xHis tag for detection and purification of the recombinant protein.
	• An enterokinase recognition site for removal of the N-terminal peptide
pTrcHis2-TOPO <sup>®</sup>	In addition to the features listed in <b>Description of Vectors</b> , above, pTrcHis2-TOPO <sup>®</sup> encodes a C-terminal peptide containing the <i>c-myc</i> epitope and a 6xHis tag for detection and purification of the recombinant protein.

**TOPO<sup>®</sup> Cloning** Both pTrcHis-TOPO<sup>®</sup> and pTrcHis2-TOPO<sup>®</sup> are supplied linearized with:

- Single 3'-thymidine (T) overhangs for TA Cloning®
- Topoisomerase I covalently bound to the vector (this is referred to as "activated" vector)

*Taq* polymerase has a nontemplate-dependent terminal transferase activity that adds a single deoxyadenosine (A) to the 3' ends of PCR products. The linearized vector supplied in this kit has single, overhanging 3' deoxythymidine (T) residues. This allows PCR inserts to ligate efficiently with the vector.

TOPO<sup>®</sup> Cloning exploits the ligation activity of topoisomerase by providing an "activated", linearized TA vector using proprietary technology (Shuman, 1994). Ligation of the vector with a PCR product containing 3´A-overhangs is very efficient and occurs spontaneously within 5 minutes at room temperature. The TOPO<sup>®</sup> Cloning reaction can be transformed into chemically competent cells (provided) or electroporated directly into electrocompetent cells.



# Induction of Expression

The strong *trc* promoter regulates expression in *E. coli*. The product of the *lacI*<sup>q</sup> gene encoded in both vectors represses this promoter. To induce expression, IPTG is added to a final concentration of 1 mM, and the culture monitored for expression of the protein of interest.



#### Detection of Recombinant Proteins

Expression of your recombinant protein can be detected using an antibody to the protein itself or to the appropriate epitope. The table below describes the antibodies available for use with pTrcHis-TOPO<sup>®</sup> or pTrcHis2-TOPO<sup>®</sup>. Horseradish peroxidase (HRP)-conjugated antibodies allow one-step detection using colorimetric or chemiluminescent detection methods.

Vector	Epitope	Antibody	Catalog No.
pTrcHis-TOPO <sup>®</sup>	Xpress <sup>™</sup>	Anti-Xpress <sup>™</sup>	R910-25
		Anti-Xpress <sup>™</sup> -HRP	R911-25
	HisG	Anti-HisG	R940-25
		Anti-HisG-HRP	R941-25
pTrcHis2-TOPO®	с-тус	Anti-myc	R950-25
		Anti-myc-HRP	R951-25
	C-terminal polyhistidine tag	Anti-His(C-term)	R930-25
		Anti-His(C-term)-HRP	R931-25

#### Purification of Recombinant Protein

The metal binding domain encoded by the 6xHis tag allows simple, easy purification of your recombinant protein by Immobilized Metal Affinity Chromatography (IMAC) using Invitrogen's ProBond<sup>™</sup> Resin (see below). To purify proteins expressed using pTrcHis-TOPO<sup>®</sup> or pTrcHis2-TOPO<sup>®</sup>, the ProBond<sup>™</sup> Purification System is available separately. Additional ProBond<sup>™</sup> resin is available in bulk. See the table below for ordering information.

Product	Quantity	Catalog no.
ProBond <sup>™</sup> Metal-Binding Resin	50 ml	R801-01
(precharged resin provided as a 50% slurry in 20% ethanol)		
	150 ml	R801-15
Purification Columns	50	R640-50
(10 ml polypropylene columns)		
ProBond <sup>™</sup> Purification System	6 purifications	K850-01
(includes six 2 ml precharged, prepacked ProBond <sup>™</sup> resin columns and buffers for native and denaturing purification)		

Reagents Available Separately Some of the reagents in this kit are available separately. Use the table below for ordering information.

Product	Amount	Catalog no.
10 mM dNTPs	1 ml, 2.5 mM of each dNTP	R725-01
Xpress <sup>™</sup> Forward Primer	2 µg, lyophilized	N576-02

# Designing PCR Primers for pTrcHis-TOPO<sup>®</sup>

Introduction	It is very important to design your PCR primers to ensure you obtain the recombinant protein you need for your studies. Please use the information below and the diagram on the next page to design your PCR primers.	
Special Considerations	pTrcHis-TOPO <sup>®</sup> is designed with the initiati ribosome binding site to ensure maximum tra and is contained in the unique <i>Nco</i> I site.	on ATG correctly spaced from an optimized anslation. This ATG is located at bp 413-415
	The N-terminal peptide can be cleaved off from partially purified or purified recombinant protein using enterokinase. Please note that you will have at least four extra amino acids at the N-terminus of your protein (Asp-Pro-Thr-Leu-). The exact number of additional amino acids will depend on your PCR product.	
Primer Design	Suggestions for primer design are provided i	n the table below.
	If you wish to	Then
	clone in frame with the DNA encoding the N-terminal peptide	the forward PCR primer must be designed to ensure that your ORF is cloned in frame with the DNA encoding the N-terminal peptide.
	remove the N-terminal leader (for expression of truly native protein) <b>Note:</b> Proteins with N-terminal leaders	the forward PCR primer can be designed to include a unique <i>Nco</i> I site which contains the first ATG of your protein.
	tend to express better in <i>E. coli</i> . You may	Ex. 1. 5´-ACC <u>ATG</u> G
	wish to prepare constructs with and without the leader and compare expression.	After TOPO <sup>®</sup> Cloning your PCR product, the vector can be digested with <i>Nco</i> I and religated, assuming there are no internal <i>Nco</i> I sites in your PCR product.
	include the native stop codon for your protein	include the native sequence containing the stop codon in the reverse primer or make sure the stop codon is upstream from the reverse PCR primer binding site.



Do not add 5' phosphates to your primers for PCR. This will prevent ligation into  $pTrcHis-TOPO^{\$}$ .

### **Designing PCR Primers for pTrcHis-TOPO<sup>®</sup>**, continued

pTrcHis-TOPO<sup>®</sup> The diagram below is supplied to help you design appropriate PCR primers to correctly clone and express your PCR product. Restriction sites are labeled to indicate the actual **Cloning Site** cleavage site. The complete sequence is available by downloading from our Web site (www.invitrogen.com) or by calling Technical Service (see page 33). -10 -35 lac operator (lacO) TGTTGACAAT TAATCATCCG GCTCGTATAA TGTGTGGAAT TGTGAGCGGA TAACAATTTC ACACAGGAAA 191 rrnB antitermination sequence 261 CAGCGCCGCT GAGAAAAAGC GAAGCGGCAC TGCTCTTTAA CAATTTATCA GACAATCTGT GTGGGCACTC pTrcHis Forward priming site gene 10 translational enhancer RBS Minicistron 331 GACCGGAATT ATCGATTAAC TTTATTATTA AAAATTAAAG AGGTATATAT TA ATG TAT CGA TTA Met Tyr Arg Leu HisG epitope RBS 6xHis tag Nco I 395 AAT AAG GAG GAA TAA ACC **ATG** GGG GGT TCT CAT CAT CAT CAT CAT CAT' GGT ATG Met Gly Gly Ser His His His His His Gly Met Asn Lys Glu Glu \* \* \* Nhe I Xpress<sup>™</sup> epitope Xpress<sup>™</sup> Forward priming site 449 GCT AGC ATG ACT GGT GGA CAG CAA ATG GGT CGG GAT CTG TAC GAC GAT GAC GAT Ala Ser Met Thr Gly Gly Gln Gln Met Gly Arg Asp Leu Tyr Asp Asp Asp Asp EK recognition sequence Bam HI Eco RI Bst BI HindIII 503 AAG GAT CCA ACC CTT AAG GGCGAATTCA ATTCGAAGCT TGGCTGTTTT TTC CTA GGT TGG GAA Product TTC CCGCTTAAGT TAAGCTTCGA ACCGACAAAA Lys Asp Pro Thr Leu >>> ⊸ EK cleavage site pTrcHis Reverse priming site 551 GGCGGA<u>TGA</u>G AGAAGATTTT CAGCC<u>TGA</u>TA CAGAT<u>TAA</u>AT CAGAACGCAG AAGCGGTCTG ATAAAACAGA rrnB T<sub>1</sub> and T<sub>2</sub> transcription termination sequence ATTTGCCTGG CGGCAGTAGC GCGGTGGTCC CACCTGACCC CATGCCGAAC TCAGAAGTGA AACGCCGTAG 621 691 CGCCGATGGT AGTGTGGGGT CTCCCCATGC GAGAGTAGGG AACTGCCAGG CATCAAATAA AACGAAAGGC 761 TCAGTCGAAA GACTGGGCCT TTCGTTTTAT CTGTTGTTTG TCGGTGAACG CTCTCCTGAG TAGGACAAAT

#### **Designing PCR Primers for pTrcHis2-TOPO®**

initiation codon

native protein)

expression.

remove the small N-terminal

leader (for expression of truly

Note: Proteins with N-terminal

leaders tend to express better in

without the leader and compare

E. coli. You may wish to

prepare constructs with and

Introduction	It is very important to design your PCR primers to ensure you obtain the recombinant protein you need for your studies. Please use the information below and the diagram on the next page to design your PCR primers.		
Special Considerations	pTrcHis2-TOPO <sup>®</sup> is designed wit ribosome binding site to ensure op and is contained in the unique <i>Ncc</i> encoded in the DNA between the	h the initiation ATG correctly spaced from the optimized otimum translation. This ATG is located at bp 413-415 o I site. Please note that there are two more amino acids initiation codon and the TOPO <sup>®</sup> Cloning site.	
	Please note that the C-terminal tag protein without the C-terminal tag	g cannot be cleaved off. If you wish to express your s, see the table below.	
Primer Design	Suggestions for primer design are	provided in the table below.	
	If you wish to	Then	
	include the <i>c-myc</i> epitope and polyhistidine region	the reverse PCR primer must be designed to remove the native stop codon in the gene of interest and preserve the reading frame through the C-terminal tag.	
	<b>NOT</b> include the <i>c-myc</i> epitope and polyhistidine region	include the native sequence containing the stop codon in the reverse primer or make sure the stop codon is upstream from the reverse PCR primer binding site.	
	clone in frame with the	the forward PCR primer must be designed to ensure	

	~
$\approx$	Note

Do not add 5' phosphates to your primers for PCR. This will prevent ligation into pTrcHis2-TOPO<sup>®</sup>.

product.

the protein.

Ex. 1. 5'-ACC ATG G ....

continued on next page

that your ORF is in frame with the initiation codon.

the forward PCR primer can be designed to include

a unique Nco I site which contains the first ATG of

After TOPO<sup>®</sup> Cloning your PCR product, the vector

can be digested with Nco I and religated, assuming

there are no internal Nco I sites in your PCR

### **Designing PCR Primers for pTrcHis2-TOPO<sup>®</sup>**, continued

# **pTrcHis2 TOPO**<sup>®</sup> The diagram below is supplied to help you design appropriate PCR primers to correctly clone and express your PCR product. Restriction sites are labeled to indicate the actual cleavage site. The complete sequence is available by downloading from our Web site (www.invitrogen.com) or by calling Technical Service (see page 33).

	35	-10		lac O	
181	<b>31</b> TGAAATGAGC TGTTGACAAT TAATCATCCG GCTCG	TATAA TG	STGTGGAAT TG	TGAGCGGA TAACAA	ATTTC ACACAGGAAA
	rmB antite	ermination sequ	uence		
261	51 CAGCGCCGCT GAGAAAAAGC GAAGCGGCAC TGCTC	CTTTAA CA	ATTTATCA GA	CAATCTGT GTGGGG	CACTC GACCGGGAATT
	۲q	TrcHis Forward	priming site		
	T7 gene 10 translational enhancer region RBS			Minicistron RBS	S Nco I
341	1 ATCGATTAAC TTTATTATTA AAAATTAA'AG AGGTA	ATATAT TA	ATG TAT CG Met Tyr Ard	A TTA AAT AAG ( g Leu Asn Lys (	GAG GAA TAA ACC Glu Glu ***
Ini	Initiation codon EcoR I BstB   Hin	nd III	SnaB I	myc epitor	De
413	13 ATG GCC CTT PCR AAG GGC GAA TTC GA	AA GCT TA	C'GTA GAA CA	AA AAA CTC ATC	TCA GAA GAG GAT
	Met Ala Leu Lys Gly Glu Phe Gl	lu Ala Ty	r Val Glu G	ln Lys Leu Ile	Ser Glu Glu Asp
	6xHis tag	J	_		
473	73 CTG AAT AGC GCC GTC GAC CAT CAT CAT CA Leu Asn Ser Ala Val Asp His His His Hi	AT CAT CA Ls His Hi	AT TGA GTTTAX .s ***	AACG GTCTCCAGC	T TGGCTGTTTT
	pTrcHis Reverse pr	rimina site			
541	1 ggcggatgag agaagatttt cagcctgata cagat	TAAAT CA	GAACGCAG AA	GCGGTCTG ATAAAA	ACAGA ATTTGCCTGG
	rmB T <sub>1</sub> an	d T <sub>2</sub> transcriptio	on terminators		
621	21 CGGCAGTAGC GCGGTGGTCC CACCTGACCC CATGC	CGAAC TC	CAGAAGTGA AA	CGCCGTAG CGCCGA	ATGGT AGTGTGGGGT
701	1 CTCCCCATGC GAGAGTAGGG AACTGCCAGG CATCA	AATAA AA	CGAAAGGC TC	AGTCGAAA GACTGO	GCCT TTCGTTTTAT

### **Producing PCR Products**

Introduction	This section describes a procedure for PCR using the primers you designed. Please note that other procedures are suitable.				
Materials Supplied	Yo	You will need the following reagents and equipment.			
by the User	<ul><li><i>Taq</i> polymerase</li><li>Thermocycler</li></ul>				
	•	DNA template			
	•	Primers for PCR product			
Producing PCR Products	1.	Set up the following 50 µl PCR r your primers and template. Be su ensure that all PCR products are for plasmids and 100 ng template	eaction. Use the cycling parameters suitable for re to include a 7 to 30 minute extension at 72°C to full length and 3' adenylated. Use 10 ng template of or genomic DNA.		
		DNA Template	10-100 ng		
		10X PCR Buffer	5 µl		
		50 mM dNTPs	0.5 µl		
		Primers (100-200 ng/µl each)	1 µl each		
		Sterile water	add to a final volume of 49 µl		
		<i>Taq</i> Polymerase (1 unit/µl)	1 µl		
		Total Volume	50 μl		
	2.	Check the PCR product by agaro discrete band. If not, see the <b>Not</b>	se gel electrophoresis. You should see a single, e, below.		
	3.	Use the PCR product immediatel store the product at -20°C until re for about 1 week. Long-term stor from your PCR product. This will Cloning reaction.	y in a TOPO <sup>®</sup> Cloning reaction (next page) or eady for use. PCR products may be stored at -20°C age may result in removal of the 3´A-overhangs l decrease cloning efficiency in the TOPO <sup>®</sup>		
Note	If y ens PC	ou do not see a single, discrete band from your PCR, use one of the options below to are a single PCR species in your TOPO <sup>®</sup> Cloning reaction. Please note that small R products will clone preferentially over larger ones.			
	•	Gel-purify your fragment before u Cloning <sup>®</sup> Kit (see page 21). Take contamination.	sing either the pTrcHis or pTrcHis2 TOPO TA special care to avoid sources of nuclease		
	•	Optimize your PCR to eliminate n	nultiple bands and smearing (Innis <i>et al</i> 1990)		

• Optimize your PCR to eliminate multiple bands and smearing (Innis *et al.*, 1990). The PCR Optimizer<sup>™</sup> Kit (Catalog no. K1220-01) from Invitrogen can help you optimize your PCR. Please call Technical Service for more information (page 33).

# **TOPO<sup>®</sup> Cloning and Transformation**

Introduction	TOPO <sup>®</sup> Cloning technology allows you to produce your PCR products, ligate them into the appropriate TOPO <sup>®</sup> vector, and transform the recombinant vector into <i>E. coli</i> all in one day. It is important to have everything you need set up and ready to use to ensure the best possible results. If this is the first time you have TOPO <sup>®</sup> Cloned, you may wish to perform the control reactions on pages 24-26 in parallel with your samples.
Materials Supplied by the User	In addition to microbiological supplies ( <i>i.e.</i> plates and spreaders), you will need the following reagents and equipment.
-	• 42°C water bath
	• 37°C shaking and non-shaking incubator
	• Two LB plates containing 50 µg/ml ampicillin and 0.5% glucose per transformation (see page 19 for recipe)
	<b>Tip</b> : If you already have LB plates containing 50 $\mu$ g/ml ampicillin only, you may spread 20 $\mu$ l of a 2 M (or a 50%) glucose solution onto the plate. Please note that the concentration of glucose does not have to be exact.
	We recommend that you include glucose (25 mM, 0.5%) in the selection medium to ensure stability of your insert. Promoters based on the <i>lac</i> promoter ( <i>i.e. trc</i> ) tend to have higher basal levels of transcription. If your insert is toxic to <i>E. coli</i> , DNA rearrangement may occur. Glucose represses basal level transcription to stabilize your construct.
Mechanism of Glucose Repression	A transcriptional activator protein called CAP (catabolite activator protein) normally binds upstream of the <i>trc</i> promoter and activates transcription. This protein requires cAMP to bind to the DNA. Adding glucose to the medium can reduce intracellular cAMP levels. Supplementing LB medium and agar plates with glucose will repress basal level transcription from the <i>trc</i> promoter.
Preparation	For each transformation, you will need one vial of competent cells and one or two selective plates.
	• Equilibrate a water bath to 42°C.
	• Thaw the vial of SOC medium from Box 2 and bring to room temperature.
	<ul> <li>Warm LB plates containing 50 µg/ml ampicillin and 0.5% glucose at 37°C for 30 minutes.</li> </ul>
	• Thaw <u>on ice</u> 1 vial of One Shot <sup>®</sup> cells for each transformation.

# **TOPO<sup>®</sup> Cloning and Transformation, continued**

TOPO <sup>®</sup> Cloning Reaction	In general, 0.5 to 4 $\mu$ l of a typical PCR sample (10-20 ng/ $\mu$ l) with an average insert length of 400 to 1000 bp will give the proper insert:vector ratio for TOPO <sup>®</sup> Cloning.			
	1.	Set up the following 5 µl TOPO	<sup>9</sup> Cloning reaction.	
		Fresh PCR product	0.5 to 4 µl	
		Sterile Water	add to a final volume of 4 µl	
		TOPO <sup>®</sup> vector	<u>1 μl</u>	
		Final Volume	5 µl	
	2. Mix gently and incubate for <b>5 minutes</b> at room temperature (~25°C). Fo possible results, do not leave for more than <b>5 minutes or the transfor efficiencies may decrease.</b>			
	3.	If needed, the TOPO <sup>®</sup> Cloning reaction may be stored on ice or frozen at -20°C for up to 24 hours. You may see a decrease in the transformation efficiency, but the cloning efficiency should remain high. We recommend that you proceed immediately to <b>Transformation</b> , below.		
One Shot <sup>®</sup> Transformation	1.	Add 2 µl of the TOPO <sup>®</sup> Cloning gently. <b>Do not mix by pipetting</b>	reaction into a vial of One Shot <sup>®</sup> cells and mix <b>up and down.</b>	
Reaction	2.	Incubate on ice for 30 minutes.		
	3.	Heat shock the cells for 30 second	ds at 42°C without shaking.	
	4.	Immediately transfer the tubes to	ice.	
	5.	Add 250 µl of room temperature	SOC medium.	
	6.	Cap the tube tightly and shake th	e tube horizontally at 37°C for 30 minutes.	
	7.	Spread 10-50 µl from each transfincubate overnight at 37°C. We resure well-spaced colonies. For ensure even spreading.	Formation on a prewarmed selection plate and recommend that you plate two different volumes to plating smaller volumes, add 20 $\mu$ l of SOC to	
	8.	An efficient TOPO <sup>®</sup> Cloning rea ~10 colonies for analysis.	ction will produce hundreds of colonies. Pick	
Transformation by Electroporation	Use the '	ONLY electrocompetent cells fo FOP10 One Shot <sup>®</sup> chemically co	r electroporation to avoid arcing. Do not use mpetent cells for electroporation.	

# **TOPO<sup>®</sup> Cloning and Transformation, continued**

Analysis of Positive Clones	1.	Take the 10 colonies and culture them overnight in LB medium containing $50 \mu g/ml$ ampicillin and 0.5% glucose. <b>Note</b> : If you use a rich broth like SOC or Terrific Broth, grow the cells for 4 hours before performing a miniprep; do not grow overnight. We obtain less DNA with overnight growth.	
	2.	Isolate plasmid DNA using your method of choice. If you need ultra-pure plasmid DNA for automated or manual sequencing, we recommend the S.N.A.P. <sup>™</sup> MiniPrep Kit (Catalog no. K1900-01) or the S.N.A.P. <sup>™</sup> MidiPrep Kit (K1910-01).	
	3.	Analyze the plasmids by restriction analysis or by sequencing. For pTrcHis-TOPO <sup>®</sup> use the Xpress <sup>™</sup> Forward and the pTrcHis Reverse sequencing primers for sequencing. For pTrcHis2-TOPO <sup>®</sup> use the pTrcHis Forward and Reverse sequencing primers to sequence your insert. For the sequence surrounding the TOPO <sup>®</sup> Cloning site, please refer to the diagram on page 7 for pTrcHis-TOPO <sup>®</sup> or page 9 for pTrcHis2-TOPO <sup>®</sup> .	
	If you need help with setting up restriction enzyme digests or DNA sequencing, please refer to general molecular biology texts (Ausubel <i>et al.</i> , 1994; Sambrook <i>et al.</i> , 1989).		
Alternative Method of Analysis	You Forv prim perfo Falso cont prote	may wish to use PCR to directly analyze positive transformants. Use either the vard or Reverse sequencing primer and a primer that hybridizes to your insert as PCR ters. If this is the first time you have used this technique, we recommend that you form restriction analysis in parallel to confirm that PCR gives you the correct result. The positive and false negative results can be obtained because of mispriming or aminating template. The following protocol is provided for your convenience. Other pools are suitable.	
	<b>Note</b> : Additional primers and nucleotides are available separately. See page 5 for ordering information.		
	1.	Prepare a PCR cocktail consisting of PCR buffer, dNTPs, primers, and <i>Taq</i> polymerase. Use a 20 $\mu$ l reaction volume. Multiply by the number of colonies to be analyzed ( <i>e.g.</i> 10).	
	2.	Pick 10 colonies and resuspend them individually in 20 $\mu$ l of the PCR cocktail. Remember to patch colonies to a separate plate to preserve the colonies.	
	3.	Incubate the reaction for 10 minutes at 94°C to lyse the cells and inactivate nucleases.	
	4.	Amplify for 20 to 30 cycles (94°C for 1 minute, 55°C for 1 minute, and 72°C for 1 minute).	
	5.	For the final extension, incubate at 72°C for 10 minutes. Hold at +4°C.	
	6.	Visualize by agarose gel electrophoresis.	
		continued on next page	

### **TOPO<sup>®</sup> Cloning and Transformation, continued**



Long-Term Storage If you have problems obtaining transformants or the correct insert, please see pages 24-26. Control reactions are described using reagents supplied in the kit. This will help you troubleshoot your experiment. Please perform the control reactions before calling Technical Service.

Once you have identified the correct clone, be sure to purify the colony and make a glycerol stock for long term storage.

- 1. Streak the original colony out for single colonies on LB plates containing 50  $\mu g/ml$  ampicillin and 0.5% glucose.
- 2. Isolate a single colony and inoculate into 1-2 ml of LB containing  $50 \ \mu g/ml$  ampicillin and 0.5% glucose. Grow until culture reaches mid-log (OD<sub>600</sub> = 0.5-0.7)
- 3. Mix 0.85 ml of culture with 0.15 ml of sterile glycerol, transfer to a cryovial, and store at -80°C.

### **Expression of the PCR Product**

Introduction	Since each recombinant protein has different characteristics that may affect optimum expression, it is helpful to run a time course of expression to determine the best conditions for maximum expression of your particular protein. Use the positive control vector included in each kit as an expression control (see pages 31-32). TOP10 cells may be used as a general host for expression.		
Note	Remember that inclusion of the N- or C-terminal tag will increase the size of your recombinant protein by 3 to 4 kDa.		
Before Starting	Be sure to have the following reagents, solutions, and equipment on hand before starting the experiment:		
	<ul> <li>Positive control (TOP10 cells containing pTrcHis-TOPO<sup>®</sup>/lacZ or pTrcHis2- TOPO<sup>®</sup>/lacZ, see below)</li> </ul>		
	• Negative control (TOP10 cells only, see below)		
	• SOB or LB containing 50 µg/ml ampicillin (see <b>Recipes</b> , page 19) <b>Note</b> : For expression, you generally do not need to include glucose (see the next page for more information)		
	<ul> <li>37°C shaking incubator</li> <li>Thaw 1 M IPTG stock</li> </ul>		
	• 1X SDS-PAGE sample buffer		
	• Reagents and apparatus for SDS-PAGE gel		
Positive and Negative Controls	Details of each positive control vector are provided on pages 31-32. Transform the plasmid into TOP10 One Shot <sup>®</sup> cells as you did for your construct. TOP10 cells that do not contain any vector are used as a negative control.		
Pilot Expression	<ol> <li>For each strain, inoculate 2 ml of SOB or LB containing 50 µg/ml ampicillin with a single recombinant <i>E. coli</i> colony.</li> </ol>		
	2. Grow overnight at 37°C with shaking (225-250 rpm).		
	<ol> <li>The next day, inoculate 10 ml of SOB or LB containing 50 μg/ml ampicillin with 0.2 ml of the overnight culture.</li> </ol>		
	4. Grow the culture at 37°C with vigorous shaking to an $OD_{600} = 0.6$ (the cells should be in mid-log phase).		
	5. Remove a 1 ml aliquot of cells, centrifuge at maximum speed in a microcentrifuge for 30 seconds, and aspirate the supernatant. Freeze the cell pellet at -20°C. This is the zero time point sample.		
	6. Add IPTG to a final concentration of 1 mM (9 μl of a 1 M IPTG stock to 9 ml) and grow at 37°C with shaking.		
	<ol> <li>Take 1 ml samples every hour for 5 hours (or more) and treat as described in Step 5. Label each tube to correspond to the number of hours postinduction.</li> </ol>		

### Expression of the PCR Product, continued

Preparation of	Before starting, prepare SDS-PAGE gels to analyze all the time points you collected.			
Time Point Samples	<ol> <li>When all the time points have been collected, resuspend each pellet in 100 µl of 1X SDS-PAGE sample buffer.</li> </ol>			
	2. Boil 5 minutes and centrifuge briefly. If solution is viscous, sonicate briefly and centrifuge again.			
	3. Analyze 5 µl of each sample on an SDS-PAGE gel.			
Analysis of Time Point Samples	1. Stain the gel with Coomassie blue and look for a band of increasing intensity in the expected size range for the recombinant protein.			
	<b>Note</b> : The tags contribute ~3 to 4 kDa to your protein.			
	2. Use the negative control to distinguish recombinant proteins from background proteins.			
	3. Use the positive control to confirm that growth and induction was done properly. pTrcHis-TOPO <sup>®</sup> / <i>lacZ</i> should yield a 40 kDa protein and pTrc-His2-TOPO <sup>®</sup> / <i>lacZ</i> should yield a 120 kDa protein with maximum expression occurring between 3-4 hours.			
	4. You should be able to determine the optimal time point for maximum expression. If you do not see your protein of interest, please see the <b>Troubleshooting</b> section, page 18.			
Levels of Recombinant Protein	high levels of basal transcription. Use glucose to further repress transcription (see page 11). Supplementing LB medium with 25 mM glucose (0.5% w/v) prior to induction will repress basal level transcription from the <i>trc</i> promoter. To induce expression, pellet the cells, and resuspend them in LB without glucose and add IPTG to induce.			
Scale-Up of Expression	Use the conditions determined previously to grow and induce 50 ml of cells. This is the largest culture volume to use with the 2 ml prepacked columns included in the ProBond <sup>™</sup> Purification System. If you need to purify larger amounts of recombinant protein, you may need more ProBond <sup>™</sup> resin. See page 4 for ordering information.			
	1. Inoculate 2 ml of SOB or LB containing 50 μg/ml ampicillin with a single recombinant <i>E. coli</i> colony.			
	2. Grow overnight at 37°C with shaking (225-250 rpm).			
	3. The next day, inoculate 50 ml of SOB or LB containing 50 $\mu$ g/ml ampicillin with 1 ml of the overnight culture.			
	4. Grow the culture at 37°C with vigorous shaking to an $OD_{600} = 0.6$ (the cells should be in mid-log phase).			
	5. Add IPTG to a final concentration of 1 mM (50 µl of 1 M IPTG stock to 50 ml).			
	6. Grow at 37°C with shaking until the optimal time point is reached. Harvest the cells by centrifugation (3000 x g for 10 minutes at $+4^{\circ}$ C).			
	<ol> <li>At this point, you may proceed directly to purification (please refer to the ProBond<sup>™</sup> Purification System manual) or store the cell pellet at -80°C for future use.</li> </ol>			

### Expression of the PCR Product, continued

Enterokinase Cleavage	For recombinant proteins expressed from pTrcHis-TOPO <sup>®</sup> , you may wish to remove the N-terminal tag from your partially pure or purified protein. We recommend using EnterokinaseMax <sup>™</sup> (Catalog no. E180-01), a recombinant form of bovine enterokinase. For more information, please see our Web site (www.invitrogen.com) or contact
	Technical Service (page 33).

#### Troubleshooting

# Troubleshooting Table

If you have trouble expressing your protein, try some of the suggestions listed below. Please be sure to include the positive and negative controls when testing for expression of your protein to ensure that the cells were grown and induced correctly. If you find that the positive control did not express, it may be that the IPTG solution is too old. Prepare fresh IPTG solution (see page 20).

Problem	Possible Cause	Solution
Recombinant protein is not detected on a Coomassie stained gel.	Low expression of recombinant protein.	Use western blot analysis to detect recombinant protein expression. You may use antibody to your own protein or the appropriate antibody listed on page 4.
Low expression of recombinant protein.	Recombinant plasmid is unstable or protein is slightly toxic	Include glucose in the growth medium to reduce basal levels of transcription. Be sure to check plasmid to ensure that no rearrange-ments have occurred. See page 16 for details.
	<i>Taq</i> polymerase may introduce mutations.	Sequence your construct. If you find mutations, redo your PCR using a proofreading polymerase and add 3' A overhangs using the method on page 23.
No expression/detection of protein	PCR product is out of frame with the N-terminal peptide or the initiation codon and/or C-terminal peptide.	Sequence your construct to confirm the protein is in frame with the N-terminal peptide or the initiation codon and the C-terminal tag (if the tag is desired).
	<i>Taq</i> polymerase may introduce mutations.	Sequence your construct. If you find mutations, redo your PCR using a proofreading polymerase and add 3' A overhangs using the method on page 23.

# Appendix

### Recipes

LB (Luria-Bertani)	Composition:			
Medium and Plates	1.0% Tryptone 0.5% Yeast Extract 1.0% NaCl +/- 0.5% glucose (dextrose)			
	1.	For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 ml deionized water.		
	2.	Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.		
	3.	Prepare a 50% solution of glucose (dextrose). Filter-sterilize or autoclave as described below. <b>Note</b> : Solution may turn yellowish. This is normal.		
	4.	Autoclave on liquid cycle for 20 minutes at 15 psi. Allow solutions to cool to 55°C. Add antibiotic to the medium if needed. Add glucose to a final concentration of 0.5%.		
	5.	Store at room temperature or at +4°C. Shelf life with ampicillin is 1-2 weeks.		
	LB agar plates			
	1.	Prepare LB medium as above, but add 15 g/L agar before autoclaving.		
	2.	Autoclave on liquid cycle for 20 minutes at 15 psi.		
	3.	After autoclaving, cool to ~55°C, add antibiotic (50 $\mu$ g/ml of ampicillin) and glucose to 0.5%, if desired. Pour into 10 cm plates.		
	4.	Let harden, then invert and store at $+4^{\circ}$ C, in the dark. Shelf life with ampicillin is 1-2 weeks.		
SOB Medium (with Ampicillin)	<b>SOI</b> 2% 2 0.5% 0.05 2.5 n 10 m	B (per liter) Tryptone 6 Yeast Extract % NaCl mM KCl nM MgCl <sub>2</sub>		
	1.	Dissolve 20 g tryptone, 5 g yeast extract, and 0.5 g NaCl in 950 ml deionized water.		
	2.	Make a 250 mM KCl solution by dissolving 1.86 g of KCl in 100 ml of deionized water. Add 10 ml of this stock KCl solution to the solution in Step 1.		
	3.	Adjust pH to 7.5 with 5 M NaOH and add deionized water to 1 liter.		
	4.	Autoclave this solution, cool to ~55°C, and add 10 ml of sterile 1 M MgCl <sub>2</sub> . You may also add ampicillin to 50 $\mu$ g/ml.		
	5.	Store at +4°C. Medium is stable for only 1-2 weeks.		

### **Recipes**, continued

#### 1 M IPTG

- 1. To prepare a 1 M stock solution, dissolve 2.38 g of IPTG in 10 ml of deionized water.
- 2. Filter-sterilize and store in 1 ml aliquots at -20°C.

### **Purifying PCR Products**

Introduction	Smearing, multiple banding, primer-dimer artifacts, or large PCR products (>1 kb) may necessitate gel purification. If you intend to purify your PCR product, be extremely careful to remove all sources of nuclease contamination. There are many protocols to isolate DNA fragments or remove oligonucleotides. Please refer to Current Protocols in Molecular Biology, Unit 2.6 for the most common protocols (Ausubel <i>et al.</i> , 1994). Two simple protocols are provided below that work for most people. Please note that cloning efficiency may decrease with purification of the PCR product. You may wish to optimize your PCR to produce a single band (see <b>Producing PCR</b> <b>Products</b> , page 10).		
Note			
Using the S.N.A.P. <sup>™</sup> MiniPrep Kit	The S.N.A.P. <sup>™</sup> MiniPrep Kit (Catalog no. K1900-01) allows you to rapidly purify PCR products from regular agarose gels. You will need to prepare a 6 M sodium iodide, 10 mM sodium sulfite solution in sterile water before starting.		
	1. Electrophorese amplification reaction on a 1 to 5% regular TAE agarose gel.		
	<b>Note</b> : Do not use TBE to prepare agarose gels. Borate interferes with the sodium iodide step, below.		
	2. Cut out the gel slice containing the PCR product and melt it at 65°C in 2 volumes of 6 M NaI.		
	3. Add 1.5 volumes Binding Buffer (provided in the S.N.A.P. <sup>™</sup> MiniPrep Kit).		
	4. Load solution (no more than 1 ml at a time) from Step 3 onto a S.N.A.P. <sup>™</sup> column. Centrifuge 1 minute at 3000 x g in a microcentrifuge and discard the supernatant.		
	5. If you have solution remaining from Step 3, repeat Step 4.		
	6. Add 900 μl of the Final Wash Buffer (provided in the S.N.A.P. <sup>™</sup> MiniPrep Kit).		
	7. Centrifuge 1 minute at 3000 x g in a microcentrifuge and discard the supernatant.		
	8. Centrifuge again at maximum speed for 1 minute to fully dry the resin.		
	<ol> <li>Elute the purified PCR product in 40 μl of TE or sterile water. Use 4 μl for the TOPO<sup>®</sup> Cloning reaction and proceed as described on page 8.</li> </ol>		
Quick S.N.A.P. <sup>™</sup> Method	An even easier method is to simply cut out the gel slice containing your PCR product, place it on top of the S.N.A.P. <sup><math>T</math></sup> column bed, and centrifuge at full speed for 10 seconds. Use 1-2 µl of the flow-through in the TOPO <sup>®</sup> Cloning reaction (page 12). Be sure to make the gel slice as small as possible for best results.		

### Purifying PCR Products, continued

Low-Melt Agarose Method	If y that effi	ou gel-purify your PCR product in low-melt agarose, use the procedure below. Please note gel purification will result in a dilution of your PCR product and decreased cloning ciencies.
	1.	Electrophorese as much as possible of your PCR reaction on a low-melt agarose gel (0.8 to 1.2%) in TAE buffer.
	2.	Visualize the band of interest and excise the band.
	3.	Place the gel slice in a microcentrifuge tube and incubate the tube at 65°C until the gel slice melts.
	4.	Place the tube at 37°C to keep the agarose melted.
	5.	In a fresh tube, mix together 4 $\mu$ l of the melted agarose containing your PCR product and 1 $\mu$ l of TOPO <sup>®</sup> vector.
	6.	Incubate at 37°C for 5 to 10 minutes.
	7.	Transform 2 to 4 µl directly into TOP10 One Shot <sup>®</sup> cells using the method on page 12.

# Addition of 3<sup>´</sup> A-Overhangs Post-Amplification

Introduction	Direct cloning of DNA amplified by <i>Vent</i> <sup>®</sup> or <i>Pfu</i> polymerases into TOPO TA Cloning <sup>®</sup> vectors is often difficult because of very low cloning efficiencies. Proofreading polymerases lack the terminal transferase activity that adds the 3' A-overhangs necessary for TA Cloning <sup>®</sup> . Invitrogen has developed a simple method to clone these blunt-ended fragments.		
Before Starting	You will need the following items:		
-	• <i>Taq</i> polymerase		
	• A heat block equilibrated to 72°C		
	• Phenol-chloroform		
	• 3 M sodium acetate		
	• 100% ethanol		
	• 80% ethanol		
	• TE buffer		
Procedure	This is just one method for adding 3´ adenines. Other protocols may be suitable.		
	1. After amplification with <i>Vent</i> <sup>®</sup> or <i>Pfu</i> polymerase, place vials on ice and add 0.7-1 unit of <i>Taq</i> polymerase per tube. Mix well. It is not necessary to change the buffer.		
	2. Incubate at 72°C for 8-10 minutes (do not cycle).		
	3. Place the vials on ice. Use immediately in a TOPO <sup>®</sup> Cloning reaction.		
	<b>Note:</b> If you plan to store your sample(s) overnight before proceeding with TOPO <sup>®</sup> Cloning, you may want to extract your sample(s) with phenol-chloroform to remove the polymerases. After phenol-chloroform extraction, precipitate the DNA with ethanol and resuspend the DNA in TE buffer to the starting volume of the amplification reaction.		
Note	You may also gel-purify your PCR product after amplification with $Vent^{\text{®}}$ or $Pfu$ (see previous page). After purification, add <i>Taq</i> polymerase buffer, dATP, and 0.5 unit of <i>Taq</i> polymerase and incubate 10-15 minutes at 72°C. Use 4 µl in the TOPO <sup>®</sup> Cloning reaction.		

*Vent*<sup>®</sup> is a registered trademark of New England Biolabs.

### **TOPO TA Cloning® Control Reactions**

Introduction	If you have trouble obtaining transformants or vector containing insert, please perform the following control reactions to help troubleshoot your experiment. Performing the control reactions involves producing a 750 bp control PCR product and TOPO <sup>®</sup> Cloning it using the reagents included in the kit. Successful TOPO <sup>®</sup> Cloning of the control PCR product will yield >85% recombinants.				
Before Starting	Be su perfo	are to prepare LB plates prming the control react	s containing 50 μg/ml ai ion. Use two plates per	npicillin and 0.5% gluc transformation.	ose before
	<b>Tip</b> : If you already have LB plates containing 50 $\mu$ g/ml ampicillin only, 20 $\mu$ l of a 2 M (or a 50%) glucose solution onto the plate. Incubate the p 30 minutes to allow the glucose to diffuse into the plate. Please note that tration of glucose does not have to be exact.				u may spread e at 37°C for e concen-
Producing Control	1.	To produce the 750 bp	control PCR product, s	et up the following 50 μ	1 PCR:
PCR Product		Control DNA Templat	e (100 ng)	1 µl	
		10X PCR Buffer		5 µl	
		50 mM dNTPs		0.5 µl	
		Control PCR Primers (0.1 µg/µl each)		1 µl	
		Sterile Water		41.5 µl	
		Taq Polymerase (1 uni	t/μl)	1 μl	
		Total Volume		50 µl	
	<ol> <li>Overlay with 70 μl (1 drop) of mineral oil.</li> <li>Amplify using the following cycling parameters:</li> </ol>				
			rs:		
		Step	Time	Temperature	Cycles
		Initial Denaturation	2 minutes	94°C	1X
		Denaturation	1 minute	94°C	
		Annealing	1 minute	55°C	25X
		Extension	1 minute	72°C	
		Final Extension	7 minutes	72°C	1X
	4.	Remove 10 $\mu$ l from the discrete 750 bp band s <b>Reactions</b> , next page.	e reaction and analyze b hould be visible. Procee	y agarose gel electropho d to the <b>Control TOPC</b>	oresis. A <sup>®</sup> Cloning

#### **TOPO TA Cloning® Control Reactions, continued**

#### Control TOPO<sup>®</sup> Cloning Reactions

Using the control PCR product produced on the previous page and the pTrcHis-TOPO<sup>®</sup> or pTrcHis2-TOPO<sup>®</sup> vectors, set up two 5  $\mu$ l TOPO<sup>®</sup> Cloning reactions as described below.

1. Set up control TOPO<sup>®</sup> Cloning reactions:

	Reagent	"Vector Only"	"Vector + PCR Insert"
	Control PCR Product		1 µl
	Sterile Water	4 μ1	3 µ1
	TOPO <sup>®</sup> vector	1 µl	1 µl
2.	Incubate at 25°C (room temperature) for <b>5 minutes and place on ice. Do not incu- bate longer than 5 minutes at room temperature.</b>		
3.	Transform 2 $\mu$ l of each reac (page 12).	tion into separate vials of TC	DP10 One Shot <sup>®</sup> cells

- 4. Spread 10-50 μl of each transformation mix onto LB plates containing 50 μg/ml ampicillin and 0.5% glucose. We recommend that you plate two different volumes of the transformation reaction to ensure well-spaced colonies. For plating small volumes, add 20 μl SOC to ensure even spreading.
- 5. Incubate overnight at 37°C.

Analysis of<br/>ResultsHundreds of colonies from the "Vector + PCR Insert" reaction should be produced.<br/>Select 10 colonies and culture in LB containing 50 μg/ml ampicillin and 0.5% glucose.<br/>Isolate plasmid DNA and analyze by restriction enzyme digestion.<br/>For pTrcHis-TOPO<sup>®</sup>, use *Eco*R I and *Bam*H I.<br/>For pTrcHis2-TOPO<sup>®</sup>, use *Nco* I and *Eco*R I.<br/>Digestion of each recombinant vector will yield 3 fragments: 100 bp, 650 bp, and 4.4 kb<br/>(vector backbone). Greater than 85% of the transformants will yield this digestion<br/>pattern.

The "Vector Only" plate should contain very few colonies (<10% of the number of colonies on the "Vector + PCR Insert" plate).

Transformation<br/>ControlpUC19 plasmid is included to check the transformation efficiency of the One Shot®<br/>competent cells. Transform with 10 pg per 50  $\mu$ l of cells using the protocol on page 12.<br/>Plate 10  $\mu$ l of the transformation reaction plus 20  $\mu$ l SOC on an LB plate containing<br/>50  $\mu$ g/ml ampicillin. The transformation efficiency should be ~1 x 10° cfu/µg DNA.

#### **TOPO TA Cloning® Control Reactions, continued**

#### Factors Affecting Cloning Efficiency

Please note that lower cloning efficiencies will result from the following variables. Most of these are easily correctable, but if you are cloning large inserts, you may not obtain the expected 90% (or more) cloning efficiency.

Variable	Solution
TOPO <sup>®</sup> Cloning reactions longer than 5 minutes at room temperature	Be sure to incubate for only 5 minutes. Incubations longer than 5 minutes will decrease transformation efficiency.
pH>9 in PCR amplification reaction	Check the pH of the PCR amplification reaction and adjust with 1 M Tris-HCl, pH 8.
Incomplete extension during PCR	Be sure to include a final extension step of 7 to 30 minutes during PCR. Longer PCR products will need a longer extension time.
Cloning large inserts (>1 kb)	Increase amount of insert. Or gel-purify as described on page 21. For cloning large inserts (3 to 10 kb), try the TOPO XL PCR Cloning <sup>™</sup> Kit (Catalog no. K4700-01).
Excess (or overly dilute) PCR product	Reduce (or concentrate) the amount of PCR product.
Cloning blunt-end PCR products	Add 3' A-overhangs by incubating with <i>Taq</i> polymerase (page 23).
PCR cloning artifacts ("false positives")	TOPO <sup>®</sup> Cloning is very efficient for small fragments (< 100 bp) present in certain PCR reactions. Gel-purify your PCR product (page 21) or optimize your PCR.
PCR product does not contain sufficient 3' A-overhangs even though you used <i>Taq</i> polymerase	<i>Taq</i> polymerase is less efficient at adding a nontemplate 3' A next to another A. <i>Taq</i> is most efficient at adding a nontemplate 3' A next to a C. You may have to redesign your primers so that they contain a 5' G instead of a 5' T (Brownstein <i>et al.</i> , 1996).

#### Map and Features of pTrcHis-TOPO<sup>®</sup>

#### pTrcHis-TOPO<sup>®</sup> Map

The map below shows the features of pTrcHis-TOPO<sup>®</sup>. For the full sequence of the vector, you may download it from our Web site (www.invitrogen.com) or call Technical Service (page 33).



# Map and Features of pTrcHis-TOPO<sup>®</sup>, continued

# Features of pTrcHis-TOPO<sup>®</sup>

The important elements of pTrcHis-TOPO<sup>®</sup> (4390 bp) are described in the following table. All features have been functionally tested.

Feature	Benefit
trc promoter region	Provides high-level, inducible expression of recombinant proteins in <i>E. coli</i> . It is a hybrid promoter consisting of the - 35 region from <i>trpB</i> and the -10 region from the <i>lacUV5</i> promoter (Egon <i>et al.</i> , 1983).
<i>lac</i> operator ( <i>lacO</i> )	Binding site of the <i>lac</i> repressor to provide regulated expression of the <i>trc</i> promoter (Jacob and Monod, 1961).
<i>rrnB</i> antitermination signal	Sequence from the <i>rrnB</i> gene that reduces the level of premature transcription termination (Li <i>et al.</i> , 1984).
T7 gene 10 translational enhancer	Sequence from bacteriophage T7 gene 10 that optimizes translation initiation (Olins <i>et al.</i> , 1988).
Minicistron, RBS, and Initiation ATG	A short open reading frame containing nucleotide sequences that is efficiently translated in prokaryotic cells. A ribosome binding site (RBS) is present within the coding sequence 5' to the translation termination codon. This RBS and termination codon are positioned in frame and three nucleotides upstream from the translation initiation codon used to express the fusion protein of interest. Following translation of the open reading frame of the minicistron, ribosomes efficiently reinitiate translation at the second initiation site (Schoner <i>et al.</i> , 1986).
6xHis tag	Forms metal-binding site for affinity purification of recombinant fusion protein on metal-chelating resin ( <i>i.e.</i> $ProBond^{TM}$ ).
HisG epitope	Allows detection of the fusion protein by the Anti-HisG Antibodies (see page 4).
Xpress <sup>™</sup> Forward priming site	Permits sequencing of your insert from the 5' end.
Xpress <sup>™</sup> epitope	Allows detection of the fusion protein by the Anti-Xpress <sup><math>TM</math></sup> Antibodies (see page 4).
Enterokinase recognition site	Encodes the binding site for bovine enterokinase to permit removal of the N-terminal peptide from your protein.
TOPO <sup>®</sup> Cloning site	Allows fast insertion of your PCR product for expression.
pTrcHis Reverse priming site	Permits sequencing of your insert from the 3 <sup>-</sup> end.
<i>rrn</i> B transcription termination region	Strong transcription termination region (Orosz et al., 1991).
Ampicillin resistance gene (β-lactamase)	Allows selection of the plasmid in <i>E. coli</i> .
pBR322-derived origin	Medium copy replication and growth in E. coli.
<i>lacI</i> <sup>q</sup> gene	Encodes the <i>lac</i> repressor for regulation of the <i>trc</i> promoter (Müller-Hill <i>et al.</i> , 1968).

#### Map and Features of pTrcHis2-TOPO®

#### pTrcHis2-TOPO<sup>®</sup> Map

The map below shows the features of pTrcHis2-TOPO<sup>®</sup>. For the full sequence of the vector, you may download it from our Web site (www.invitrogen.com) or call Technical Service (page 33).



# Map and Features of pTrcHis2-TOPO<sup>®</sup>, continued

# Features of pTrcHis2-TOPO<sup>®</sup>

The important elements of pTrcHis2-TOPO<sup>®</sup> (4381 bp) are described in the following table. All features have been functionally tested.

Feature	Benefit
<i>trc</i> promoter region	Provides high-level, inducible expression of recombinant proteins in <i>E. coli</i> . It is a hybrid promoter consisting of the - 35 region from <i>trpB</i> and the -10 region from the <i>lacUV5</i> promoter (Egon <i>et al.</i> , 1983).
<i>lac</i> operator ( <i>lacO</i> )	Binding site of the <i>lac</i> repressor to provide regulated expression of the <i>trc</i> promoter (Jacob and Monod, 1961).
<i>rrnB</i> antitermination signal	Sequence from the <i>rrnB</i> gene that reduces the level of premature transcription termination (Li <i>et al.</i> , 1984).
T7 gene 10 translational enhancer	Sequence from bacteriophage T7 gene 10 that optimizes translation initiation (Olins <i>et al.</i> , 1988).
pTrcHis Forward priming site	Permits sequencing of your insert from the 5' end.
Minicistron, RBS, and Initiation ATG	A short open reading frame containing nucleotide sequences that is efficiently translated in prokaryotic cells. A ribosome binding site (RBS) is present within the coding sequence 5' to the translation termination codon. This RBS and termination codon are positioned in frame and three nucleotides upstream from the translation initiation codon used to express the fusion protein of interest. Following translation of the open reading frame of the minicistron, ribosomes efficiently reinitiate translation at the second initiation site (Schoner <i>et al.</i> , 1986).
TOPO <sup>®</sup> Cloning site	Allows fast insertion of your PCR product for expression.
C-terminal <i>myc</i> epitope (optional)	Allows detection of the fusion protein by the Anti- <i>Myc</i> Antibody (Catalog no. R930-25) (Evan <i>et al.</i> , 1985).
C-terminal polyhistidine region (optional)	Forms metal-binding site for affinity purification of recombinant fusion protein on metal-chelating resin ( <i>i.e.</i> $ProBond^{TM}$ ).
	In addition, it allows detection of the recombinant protein with Anti-His (C-term) Antibody (see page 4).
pTrcHis2 Reverse priming site	Permits sequencing of your insert from the 3' end.
<i>rrn</i> B transcription termination region	Strong transcription termination region (Orosz et al., 1991).
Ampicillin resistance gene (β-lactamase)	Allows selection of the plasmid in <i>E. coli</i> .
pBR322-derived origin	Medium copy replication and growth in E. coli.
<i>lacI</i> <sup>q</sup> gene	Encodes the <i>lac</i> repressor for regulation of the <i>trc</i> promoter (Müller-Hill <i>et al.</i> , 1968).

### pTrcHis-TOPO<sup>®</sup>/lacZ



### pTrcHis2-TOPO<sup>®</sup>/lacZ

#### Description

pTrcHis2-TOPO<sup>®</sup>/*lac*Z is a 7450 bp control vector containing the gene for  $\beta$ -galactosidase fused to the C-terminal peptide. It was constructed by amplifying the *lac*Z gene from pTrcHis2/*lac*Z and TOPO<sup>®</sup> Cloned into pTrcHis2-TOPO<sup>®</sup>. It yields a 120 kDa expression product.

Map of Control Vector The figure below summarizes the features of the pTrcHis2-TOPO<sup>®</sup>/*lacZ* vector. For the full sequence of the vector, you may download it from our Web site (www.invitrogen.com) or call Technical Service (page 33).



### **Technical Service**

Web Resources	Visit the Invitrogen Web site at <u>www.invitrogen.com</u> for:			
	• Technical resources, including manuals, vector maps and sequences, application notes, MSDSs, FAOs, formulations, citations, handbooks, etc.			
	Complet	te technical service contact inform	nation	
	Access to	o the Invitrogen Online Catalog		
	Addition	nal product information and spec	cial offers	
Contact Us	For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page ( <u>www.invitrogen.com</u> ).			
Corporate Headquarters: Invitrogen Corporation 1600 Faraday Avenue Carlsbad, CA 92008 USA Tel: 1 760 603 7200 Tel (Toll Free): 1 800 955 6288 Fax: 1 760 602 6500 E-mail: tech_service@invitrogen.com		Japanese Headquarters: Invitrogen Japan LOOP-X Bldg. 6F 3-9-15, Kaigan Minato-ku, Tokyo 108-0022 Tel: 81 3 5730 6509 Fax: 81 3 5730 6519 E-mail: jpinfo@invitrogen.com	European Headquarters: Invitrogen Ltd Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF, UK Tel: +44 (0) 141 814 6100 Tech Fax: +44 (0) 141 814 6117 E-mail: <u>eurotech@invitrogen.com</u>	
Material Data Safety Sheets (MSDSs)	MSDSs are available on our Web site at <u>www.invitrogen.com</u> . On the home page, click on <b>Technical Resources</b> and follow instructions on the page to download the MSDS for your product.			
Limited Warranty	download the MSDS for your product. Invitrogen is committed to providing our customers with high-quality goods a services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about Invitrogen product or service, contact our Technical Service Representatives. Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, fr of charge, any product that does not meet those specifications. <u>This warranty</u> <u>limits Invitrogen Corporation's liability only to the cost of the product</u> . No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing pr to acceptance of the order. Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefored Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives. Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expresses or implied, including any warranty of merchantability or fitness for a particular purpose.			

#### **Purchaser Notification**

#### Limited Use Label License No: 5 Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: outlicensing@invitrogen.com

#### Limited Use Label License No: 22 Vectors and Clones Encoding Histidine Hexamer

This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

#### **Qualifying the Product**

#### **Restriction Digest**

Supercoiled pTrcHis, pTrcHis2, pTrcHis-TOPO<sup>®</sup>/*lacZ*, and pTrcHis2-TOPO<sup>®</sup>/*lacZ* are qualified by restriction digest. The table below lists the restriction enzymes and the expected fragments.

<b>Restriction Enzyme</b>	pTrcHis	pTrcHis2	pTrcHis-TOPO <sup>®</sup> /lacZ	pTrcHis2-TOPO <sup>®</sup> /lacZ
BamH I	4400 bp (linearizes)	N/A	N/A	N/A
Hind III	4400 bp (linearizes)	N/A	N/A	4371, 3079 bp
Xmn I/Nco I	3900, 700, 160 bp	N/A	N/A	N/A
Apa I	N/A	4405 bp (linearizes)	N/A	N/A
SnaB I	N/A	4405 bp (linearizes)	N/A	N/A
Cla I	N/A	N/A	4558, 707, 46 bp	N/A
BamH I/EcoR I	N/A	N/A	4371, 940 bp	N/A
Nco I	N/A	N/A	N/A	7450 bp (linearizes)

TOPO <sup>®</sup> Cloning Efficiency	Once pTrcHis and pTrcHis2 have been adapted with topoisomerase I, they are lot- qualified using the control reagents included in the kit. Under conditions described on pages 24-26, a 750 bp control PCR product was TOPO <sup>®</sup> -Cloned into each vector and subsequently transformed into the One Shot <sup>®</sup> competent <i>E. coli</i> included with the kit. Each lot of vector should yield greater than 85% cloning efficiency.	
Primers	All primers have been lot-qualified by DNA sequencing experiments using the dideoxy chain termination technique.	
One Shot <sup>®</sup> TOP10 Competent <i>E. coli</i>	All competent cells are tested for transformation efficiency using the control plasmid. Transformed cultures are plated on LB plates containing 100 $\mu$ g/ml ampicillin and the transformation efficiency is calculated. Test transformations are performed in duplicate. Transformation efficiency should be ~1 x 10 <sup>9</sup> cfu/ $\mu$ g DNA for chemically competent cells. In addition, untransformed cells are tested for appropriate antibiotic sensitivity and lack of phage contamination.	

#### References

Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). Current Protocols in Molecular Biology (New York: Greene Publishing Associates and Wiley-Interscience).

Brosius, J., Erfle, M., and Storella, J. (1985). Spacing of the -10 and -35 Regions in the *tac* Promoter. J. Biol. Chem 260, 3539-3541.

Brownstein, M. J., Carpten, J. D., and Smith, J. R. (1996). Modulation of Non-Templated Nucleotide Addition by *Taq* DNA Polymerase: Primer Modifications that Facilitate Genotyping. BioTechniques 20, 1004-1010.

Egon, A., Brosius, J., and Ptashne, M. (1983). Vectors Bearing a Hybrid *trp-lac* Promoter Useful for Regulated Expression of Cloned Genes in *Escherichia coli*. Gene 25, 167-178.

Evan, G. I., Lewis, G. K., Ramsay, G., and Bishop, V. M. (1985). Isolation of Monoclonal Antibodies Specific for *c-myc* Proto-oncogene Product. Mol. Cell. Biol. *5*, 3610-3616.

Innis, M. A., Gelfand, D. H., Sninsky, J. J., and White, T. S. (1990) PCR Protocols: A Guide to Methods and Applications. Academic Press, San Diego, CA.

Jacob, F., and Monod, J. (1961). Genetic Regulatory Mechanisms in the Synthesis of Proteins. J. Mol. Biol. *3*, 318-328.

Li, S. C., Squires, C. L., and Squires, C. (1984). Antitermination of *E. coli* rRNA Transcription is Caused by a Control Region Segment Containing Lambda *nut*-like Sequences. Cell *38*, 851-860.

Müller-Hill, B., Crapo, L., and Gilbert, W. (1968). Mutants That Make More *lac* Repressor. Proc. Natl. Acad. Sci. USA *59*, 1259-1262.

Mulligan, M. E., Brosius, J., and Clure, W. R. (1985). Characterization *in vitro* of the Effect of Spacer Length on the Activity of *Escherichia coli* RNA Polymerase at the *tac* Promoter. J. Biol. Chem. 260, 3539-3538.

Olins, P. O., Devine, C. S., Rangwala, S. H., and Kavka, K. S. (1988). T7 Phage Gene 10 Leader RNA, a Ribosome-binding Site the Dramatically Enhances the Expression of Foreign Genes in *Escherichia coli*. Gene 73, 227-235.

Orosz, A., Boros, I., and Venetianer, P. (1991). Analysis of the Complex Transcription Termination Region of the *Escherichia coli rrnB* Gene. Eur. J. Biochem. 201, 653-659.

Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). Molecular Cloning: A Laboratory Manual, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).

Schoner, B. E., Belagaje, R. M., and Schoner, R. G. (1986). Translation of a Synthetic Two-cistron mRNA in *Escherichia coli*. Proc. Natl. Acad. Sci. USA *83*, 8506-8510.

Shuman, S. (1994). Novel Approach to Molecular Cloning and Polynucleotide Synthesis Using Vaccinia DNA Topoisomerase. J. Biol. Chem. 269, 32678-32684.

©1999-2006 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

# invitrogen®

Corporate Headquarters Invitrogen Corporation 1600 Faraday Avenue Carlsbad, CA 92008 T: 1 760 603 7200 F: 1 760 602 6500 E: tech.service@invitrogen.com

For country-specific contact information visit our web site at **www.invitrogen.com**