APEX™ Antibody Labeling Kits

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage*	Stability
Reactive fluorescent label (Component A)	5 vials	Notapplicable		
APEX™ antibody labeling tips (Component B)	5 each	ποταρρικασιε		
wash buffer (Component C)	1.8 mL	0.1 M phosphate buffered saline (PBS), pH 7.5, 2 mM azide · 2–6°C · Desiccate		When stored as directed this
dimethylsulfoxide (DMSO) (Component D)		Not applicable	Protect from light Do not freeze	kit is stable for 6 months.
labeling buffer (Component E)	100 μL	50 mM borate buffer, pH 8.3		
neutralization buffer (Component F)		1 M Tris, pH 9.0		
elution buffer (Component G)	1 mL	0.2 M acetic acid, pH 3.3		
elution syringe (Component H)	1 each	Not applicable		

*These storage conditions are appropriate when storing the entire kit upon receipt. For optimal storage conditions for each component, see individual component labels.

Number of assays: Sufficient material is supplied for 5 labelings of 10–20 µg of IgG antibody based upon the protocol below.

Approximate fluorescence excitation/emission maxima: See Table 2.

Introduction

The APEX[™] Antibody Labeling Kits provide a convenient method for covalently attaching a fluorophore to small amounts of IgG antibody (approximately, 10–20 µg). A primary antibody directly labeled with a fluorophore often produces lower background fluorescence and less nonspecific binding. Further, multiple primary antibodies of the same isotype or derived from the same species can easily be used in the same experiment if they are directly labeled with compatible fluorophores. Many IgG antibodies are often available only in small quantities and packaged with stabilizing proteins, such as BSA, or other contaminants which can interfere with the amine-reactive labeling reagents commonly used to covalently attach the fluorophore to the antibody. Removal of these contaminants often result in significant loss of the IgG antibody.

The APEX[™] Antibody Labeling Kits utilize a solid-phase labeling technique that captures the IgG antibody on the resin inside the APEX[™] antibody labeling tip (Figure 1). Any contaminants, including stabilizing proteins or amine-containing buffers are eluted through the tip. After applying the amine-reactive, fluorescent label to the IgG antibody on the resin, a fluorescent IgG conjugate is formed and is eluted from the resin using elution buffer. The fluorescent IgG conjugate is ready for use in an imaging or flow cytometry assay in as little as 2.5 hours with minimal hands on time. The typical yield of labeled antibody using this method ranges from 40–80%.

The APEX[™] Antibody Labeling Kits includes all reagents required to perform 5 separate labeling reactions of 10–20 µg of IgG antibody with one of Molecular Probes' superior fluorophores or biotin-XX. The biotin-XX provides a 14-atom spacer between the antibody and the biotin moiety.



Table 2. Spectral characteristics and applications of the labels available in the APEX[™] Antibody Labeling Kits.

Fluorophore or biotin label	Cat. no.	Excitation (nm)	Emission (nm)	Application
Alexa Fluor® 488	A10468	496	519	Fluorescent label for use in imaging or flow cytometry; Hapten for signal amplification with anti-Alexa Fluor® 488 antibodies.
Alexa Fluor® 555	A10470	555	565	
Alexa Fluor® 568	A10494	578	603	Fluorescent label for use in imaging.
Alexa Fluor® 594	A10474	590	617	
Alexa Fluor® 647	A10475	650	665	Fluorescent label for use in imaging or flow cytometry.
Oregon Green® 488	A10476	496	524	Fluorescent label for use in imaging or flow cytometry; hapten for signal amplification with anti-fluorescein/Oregon Green® antibodies.
Pacific Blue™ dye	A10478	416	451	Fluorescent label for use in imaging or flow cytometry.
Biotin-XX	A10495	NA*	NA*	Streptavidin, avidin, and neutravidin binding partner for use in imaging or flow cytometry.
*NA = not applicable.	·	·	•	

Materials Required but Not Provided	• 100–200 μL pipette			
	• 10–20 μ g IgG antibody in \leq 10 μ L in neutral pH buffer such as phosphate buffered saline (PBS), Tris-buffered saline (TBS), Tris-HCl, HEPES, borate, or equivalent (the sample can contain serum or other stabilizing proteins)			
	Microcentrifuge tubes (2 per reaction)			
Caution	DMSO (Component D), provided as a solvent in this kit, is known to facilitate the entry of organic molecules into tissues. Handle reagents containing DMSO using equipment and practices appropriate for the hazards posed by such materials. Dispose of the reagents in compliance with all pertaining local regulations.			
Spectral Characteristics and Applications	For the best results in experiments, it is important to match the light source, excitation filters, and emission filters to the spectral characteristics of the dye. Refer to Table 2 for details.			

Experimental Protocols

Antibody Labeling Procedure	We recommend using a 200 μ L pipette to dispense the buffers and antibody solution into the APEX TM antibody labeling tip. Use the elution syringe (Component H) to push the buffers into or through the labeling tip. Unlike the pipette, upon insertion into the labeling tip and depression of the plunger, the elution syringe supplies constant pressure to the APEX TM antibody labeling tip and reduces the risk of disturbing the resin bead. Because the elution syringe does not come into contact with the buffers or the resin, you can use it for all five labeling reactions.
	Add buffers to the upper surface of the resin in the APEX [™] antibody labeling tip using a gel- loading tip to minimize gel disturbance and introduction of air into the column (Figure 2A). After applying a buffer to the top of the resin, attach the APEX [™] antibody labeling tip to a 200 µL-pipette (Figure 2B) or the elution syringe (Component H) and push the buffer into the resin until the top level of the buffer just enters the top of the resin bed. When pushing a buffer into or through the APEX [™] antibody labeling tip, do not introduce air into the resin bed.
1.1	Gently tap the APEX [™] antibody labeling tip (Component B) on a hard surface to settle all resin at the bottom of the tip.
1.2	Remove both caps from the APEX [™] antibody labeling tip and place the labeling tip into a clean microcentrifuge tube.
1.3	B Hydrate the APEX [™] antibody labeling tip (Component B) by applying 100 µL of wash buffer (Component C) to the resin in the APEX [™] antibody labeling tip (Figure 2A) with a gel-loading tip.
	Note: If you are using a pipette instead of the elution syringe (Component H), remove the pipette before releasing the plunger to ensure that the resin bed is not disturbed. Apply the APEX ^{max} antibody labeling tip directly to the pipette and gently push the wash buffer through the tip into the microcentrifuge tube (Figure 2B). The hydrated resin bed volume is 10–15 µL.



Figure 2. Using the APEX[™] antibody labeling tip. Panel A: Applying solutions to the resin in the tip. Panel B: Pushing solutions onto the resin in the tip by attaching the APEX[™] antibody labeling tip to a pipette.

- 1.4 Apply 10–20 µg of IgG antibody solution to the top of the resin in the APEX[™] antibody labeling tip, then gently push the antibody solution onto the resin using the elution syringe (Component H). The antibody volume should not exceed 10 µL. After the antibody solution is pushed onto the column a drop may elute from the tip. Discard this eluent as waste.
- 1.5 To the vial of reactive dye (Component A), add the following:
 - 2 µL DMSO (Component D); pipet up and down to dissolve
 - 18 µL Labeling buffer (Component E); pipet up and down to dissolve
- $1.6 \ \ \ \ Apply 10 \ \mu L \ of the reactive dye from step 1.5 to the top of the resin, then gently push the solution onto the tip. A small amount of dye may elute from the tip. Discard this eluent as waste.$
- 1.7 Incubate the tip for 2 hours at room temperature or overnight at 4°C.
- **1.8** Wash the APEX[™] antibody labeling tip twice with 50 μL each with wash buffer (Component C) by applying 50 μL to the top of the resin, then pushing through the tip into the microcentrifuge tube.
- 1.9 To a clean microcentrifuge tube, add 10 µL neutralization buffer (Component F).
- 1.10 Position the APEX[™] antibody labeling tip on the microcentrifuge tube containing the neutralization buffer and apply 40 µL elution buffer (Component G) to the top of the resin. Push through the tip to elute the labeled antibody into the microcentrifuge tube containing neutralization buffer.

Note: The ratio of neutralization buffer to elution buffer must remain 1:4 to ensure the correct pH. The elution can be performed with $20-40 \ \mu\text{L}$ elution buffer. Reducing the elution volume may increase the antibody concentration but result in a reduced total antibody yield.

1.11	Mix the labeled antibody solution to ensure neutralization. The final eluate volume is ${\sim}50\mu\text{L}.$
	Cap the microcentrifuge tube containing the labeled antibody solution and place the tube on ice until use. The labeled antibody is ready for use in your imaging or flow cytometry application or store the antibody (see below).
	Discard the APEX ^{\cong} antibody labeling tip as biohazardous waste. Do not reuse the APEX ^{\cong} antibody labeling tip.
Labeled Antibody Storage	The labeled antibody solution can be stored in the elution/neutralization buffer at 4°C for short-term storage (up to 2 weeks). For long-term storage, exchange the storage buffer with PBS or equivalent buffer by dialysis or gel filtration and store at -20 °C. You may add other stabilization agents such as BSA to the labeled antibody solution, if desired.
Labeling Kits	Invitrogen offers several other antibody and protein labeling kits optimized for labeling of smaller amounts of IgG antibody, or larger amounts of IgG antibody or proteins >30 kDa (Table 3).

Table 3. Antibody and protein labeling kits from Invitrogen.

lgG amount	Product	Features
<1–20 µg	Zenon [®] IgG Labeling Kit	 Labeled antibodies ready to use in 10 minutes Isotype-specific labeling Fast, noncovalent attachment of label Labeling compatible with stabilizing proteins such as BSA
10–20 µg	APEX™ Antibody Labeling Kit	 Labeled antibodies ready to use in 2.5 hours (~15 minutes hands on time) Covalent attachment of label Labeling compatible with stabilizing proteins such as BSA
20–100 µg	Microscale Protein Labeling Kit	 Labeled antibodies ready to use in 2 hours (~30 minutes hands on time) Covalent attachment of label Optimized for proteins between 10–150 kDa, including IgG antibodies (~150 kDa) Stabilizing proteins must be removed from sample before labeling
100 µg	Monoclonal Antibody Labeling Kit	 Labeled antibodies ready to use in 90 minutes (~15 minutes hands on time) Covalent attachment of label Optimized for IgG antibodies (~150 kDa) Stabilizing proteins must be removed from sample before labeling Designed to label polyclonal and monoclonal IgG antibodies
1 mg	Protein Labeling Kit	 Labeled antibodies ready-to-use in 2 hours (~30 minutes hands on time) Covalent attachment of label Optimized for IgG antibodies (~150 kDa) Stabilizing proteins must be removed from sample before labeling

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
A10468	APEX [™] Alexa Fluor [®] 488 Antibody Labeling Kit	1 kit
A10470	APEX [™] Alexa Fluor [®] 555 Antibody Labeling Kit	1 kit
A10474	APEX [™] Alexa Fluor [®] 594 Antibody Labeling Kit	1 kit
A10475	APEX [™] Alexa Fluor [®] 647 Antibody Labeling Kit	1 kit
A10476	APEX™ Oregon Green® 488 Antibody Labeling Kit	1 kit
A10478	APEX [™] Pacific Blue [™] Antibody Labeling Kit	1 kit
A10494	APEX [™] Alexa Fluor [®] 568 Antibody Labeling Kit	1 kit
A10495	APEX [™] Biotin-XX, Antibody Labeling Kit.	1 kit
Related Produ	icts	
A10196	Qdot [®] 625 Streptavidin conjugate *1µM solution*	200 µL
Q10101MP	Qdot [®] 605 streptavidin conjugate *1 µM solution*	200 µL
Q10111MP	Qdot [®] 585 streptavidin conjugate *1 µM solution*	200 µL
Q10121MP	Qdot [®] 655 streptavidin conjugate *1 µM solution*	200 µL
Q10131MP	Qdot [®] 565 streptavidin conjugate *1 µM solution*	200 µL
Q10141MP	Qdot [®] 525 streptavidin conjugate *1 µM solution*	200 µL
Q10151MP	Qdot [®] Streptavidin Sampler Kit *1 μM solutions*	1 kit
Q10161MP	Qdot [®] 705 streptavidin conjugate *1 μM solution*	200 µL
Q10171MP	Qdot [®] 800 streptavidin conjugate *1 µM solution*	200 µL
S866	streptavidin, R-phycoerythrin conjugate (SAPE) *1 mg/mL*	1 mL
S11223	streptavidin, Alexa Fluor® 488 conjugate	1 mg
S11227	streptavidin, Alexa Fluor [®] 594 conjugate	1 mg
S11237	streptavidin, Alexa Fluor® 430 conjugate	1 mg
S21374	streptavidin, Alexa Fluor® 647 conjugate	1 mg
S32354	streptavidin, Alexa Fluor [®] 488 conjugate *2 mg/mL*	0.5 mL
S32357	streptavidin, Alexa Fluor [®] 647 conjugate *2 mg/mL*	0.5 mL
S32362	streptavidin, allophycocyanin conjugate *premium grade* *1 mg/mL*	250 μL
SNN2004	streptavidin-HRP conjugate *ELISA Grade*	1 mg

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road Eugene, OR 97402 Phone: (541) 465-8300 Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time) Phone: (541) 335-0338 Fax: (541) 335-0305 probesorder@invitrogen.com

Toll-Free Ordering for USA: Order Phone: (800) 438-2209

Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time) Phone: (541) 335-0353 Toll-Free (800) 438-2209 Fax: (541) 335-0238 probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd. 3 Fountain Drive Inchinnan Business Park Paisley PA4 9RF, UK Phone: +44 (0) 141 814 6100 Fax: +44 (0) 141 814 6200 Email: euroinfo@invitrogen.com Technical Services: eurotech@invitrogen.com

For country-specific contact information, visit www.invitrogen.com.

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation [®] are registered with the U.S. Patent and Trademark Office.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Copyright 2010, Life Technologies Corporation. All rights reserved. This information is subject to change without notice.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.