

## Click-iT® EdU Flow Cytometry Assay Kits

\*50 tests\*

Catalog numbers C10418, C10424, C10425

**Table 1** Contents and storage

Material	Amount	Concentration	Storage*	Stability
EdU (Component A)	10 mg	NA	<ul style="list-style-type: none"> <li>• 2–6°C</li> <li>• Desiccate</li> <li>• Protect from light</li> <li>• Do not freeze</li> </ul>	When stored as directed, this kit is stable for up to 1 year after receipt.
Alexa Fluor® 488 azide (Cat. no. C10425), Alexa Fluor® 647 azide (Cat. no. C10424), or Pacific Blue™ azide (Cat. no. C10418) (Component B)	1 vial	NA		
Dimethylsulfoxide (DMSO) (Component C)	4.25 mL	NA		
Click-iT® fixative (Component D)	5 mL	4% paraformaldehyde in PBS		
Click-iT® saponin-based permeabilization and wash reagent (Component E)	50 mL	10X solution		
CuSO <sub>4</sub> (Component F)	0.5 mL	100 mM aqueous solution		
Click-iT® EdU buffer additive (Component G)	400 mg	NA		

\*These storage conditions are appropriate when storing the entire kit upon receipt. For optimal storage conditions for each component, see vial labels. NA = Not applicable.

**Number of assays:** Sufficient material is supplied for 50 reactions, based on the protocol below.


**Approximate fluorescence excitation/emission maxima:** Pacific Blue™ azide: 410/455 nm; Alexa Fluor® 647 azide: 650/670 nm; Alexa Fluor® 488 azide: 495/519 nm.

## Introduction


Measuring a cell's ability to proliferate is a fundamental method for assessing cell health, determining genotoxicity, and evaluating anti-cancer drugs. The most accurate method is by directly measuring DNA synthesis. Initially, this was performed by incorporation of radioactive nucleosides, i.e., <sup>3</sup>H-thymidine. This method was replaced by antibody-based detection of the nucleoside analog bromo-deoxyuridine (BrdU). The Click-iT® EdU Flow Cytometry Assay Kits are novel alternatives to the BrdU assay. EdU (5-ethynyl-2'-deoxyuridine) is a nucleoside analog to thymidine and is incorporated into DNA during active DNA synthesis. Detection is based on a click reaction,<sup>1-4</sup> a copper catalyzed covalent reaction between an azide and an alkyne. In this application, the alkyne is found in the ethynyl moiety of EdU, while the azide is coupled to Pacific Blue™ dye, Alexa Fluor® 647 dye, or Alexa Fluor® 488 dye. Standard flow cytometry methods are used for determining the percentage of S-phase cells in the population (Figure 1, page 2).

The advantage of Click-iT® EdU labeling is that the small size of the dye azide allows for efficient detection of the incorporated EdU using mild conditions. Standard aldehyde-based fixation and detergent permeabilization are sufficient for the Click-iT® detection reagent to gain access to the DNA. This is in contrast to BrdU assays that require DNA denaturation (using acid, heat, or digestion with DNase) to expose the BrdU so that it may be detected with an anti-BrdU antibody. Sample processing for the BrdU assay can result in signal alteration of the cell cycle distribution as well as the destruction of antigen recognition sites when using the acid denaturation method. In contrast, the EdU cell proliferation kit is compatible with cell cycle dyes (Figure 2). The EdU assay can also be multiplexed with antibodies against surface and intracellular markers. However, some reagents or antibody conjugates may not be compatible with the Click-iT® EdU detection reaction and may need some additional steps to ensure compatibility (see Table 2, page 3 for details).

**Figure 1** Fluorescence signal from Pacific Blue™, Alexa Fluor® 488, and Alexa Fluor® 647 Click-iT® EdU Flow Cytometry Assay Kits. Jurkat (human T-cell leukemia) cells were treated with 10 µM EdU for 2 hours and detected according to the recommended staining protocol. The figures show a clear separation of proliferating cells which have incorporated EdU and nonproliferating cells which have not. Panel A shows data from cells labeled with Pacific Blue™ azide analyzed on an Attune® Acoustic Focusing Cytometer using 405 nm excitation with a 450/40 nm bandpass emission filter; Panel B shows data from cells labeled with Alexa Fluor® 488 azide analyzed on an Attune® Acoustic Focusing Cytometer using 488 nm excitation and a 530/30 nm bandpass emission filter; Panel C shows data from cells labeled with Alexa Fluor® 647 azide analyzed on a flow cytometer using 633 nm excitation and a 660/20 nm bandpass emission filter.



**Figure 2** Dual parameter plot of Click-iT® EdU Alexa Fluor® 488 and FxCycle™ Violet. Jurkat (human T-cell leukemia) cells were treated with 10 µM EdU for 2 hours and detected according to the recommended staining protocol. Data was collected and analyzed using an Attune® Acoustic Cytometer using 488 nm excitation and a 530/30 bandpass for detection of the EdU Alexa Fluor® 488 azide and 405 nm excitation and a 450/40 bandpass for detection of the FxCycle™ Violet fluorescence. This figure combines DNA content with EdU; cells that are positive for both labels are in S-phase of the cell cycle.



**Table 2** Click-iT® EdU detection reagent compatibility

Fluorescent molecule	Compatibility*
Qdot® nanocrystals	Use Qdot® nanocrystals <b>after</b> the Click-iT® detection reaction.
Fluorescent proteins (GFP)	Use anti-GFP antibodies <b>before</b> the Click-iT® detection reaction or use organic dye-based reagents for protein expression detection.
Organic dyes such as Alexa Fluor® dyes, fluorescein (FITC)	Compatible
PerCP, allophycocyanin (APC) and APC-based tandems (i.e., Alexa Fluor® 680-APC)	Compatible
R-phycoerythrin (R-PE) and R-PE based tandems (i.e., Alexa Fluor® 610-RPE)	Use R-PE and R-PE-based tandems <b>after</b> the Click-iT® detection reaction.
TC-FIAsh™/TC-ReAsH™ reagents	Detect the tetracysteine (TC) tag with FIAsh™ or ReAsH™ reagents <b>before</b> the Click-iT® detection reaction.

\*Compatibility indicates whether fluorescent molecule itself or the detection methods involve components that are unstable in the presence of copper catalyst used for the Click-iT® EdU detection reaction. Not all anti-GFP antibodies recognize the same antigen site. Rabbit and chicken anti-GFP antibodies perform well, but the mouse monoclonal antibodies tested do not generate an acceptable amount of fluorescence and are not recommended for this application.

## Before Starting

### Materials Required but Not Provided

- 1% bovine serum albumin (BSA) in phosphate buffered saline (PBS), pH 7.1–pH 7.4
- Buffered saline solution, such as PBS, D-PBS, or TBS
- Deionized water or 18 MΩ purified water
- 12 × 75-mm tubes, or other flow cytometry tubes

### Cautions

- DMSO (Component C), provided as a solvent in this kit, is known to facilitate the entry of organic molecules into tissues. Handle reagents containing DMSO using equipment and practices appropriate for the hazards posed by such materials. Dispose of the reagents in compliance with all pertaining local regulations.
- Click-iT® fixative (Component D) contains paraformaldehyde, which is harmful. Use with appropriate precautions.
- Click-iT® saponin-based permeabilization and wash reagent (Component E) contains sodium azide, which yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

### Preparing Reagents

- 1.1 Allow vials to warm to room temperature before opening.
- 1.2 To prepare a 10 mM solution of EdU, add 4 mL of DMSO (Component C) or aqueous solution (PBS) to Component A and mix well. After use, store any remaining stock solution at ≤−20°C. When stored as directed, the stock solution is stable for up to 1 year.
- 1.3 To prepare a working solution of Pacific Blue™ azide (Cat. no. C10418), Alexa Fluor® 647 azide (Cat. no. C10424), or Alexa Fluor® 488 azide (Cat. no. C10425), add 130 μL of DMSO to Component B and mix well. After use, store any remaining working solution at ≤−20°C. When stored as directed, this working solution is stable for up to 1 year.

- 1.4** To prepare 500 mL of 1X Click-iT<sup>®</sup> saponin-based permeabilization and wash reagent, add 50 mL of Component E to 450 mL of 1% BSA in PBS. Smaller amounts can be prepared by diluting a volume of Component E 1:10 with 1% BSA in PBS. After use, store any remaining solutions at 2–6°C. When stored as directed, the 1X solution is stable for 6 months and the 10X solution is stable for 12 months after receipt.

**Note:** Component E contains sodium azide (see Cautions, page 3).


- 1.5** To make a 10X stock solution of the Click-iT<sup>®</sup> EdU buffer additive (Component G), add 2 mL of deionized water to the vial and mix until the Click-iT<sup>®</sup> EdU buffer additive is fully dissolved. After use, store any remaining stock solution at ≤–20°C. When stored as directed, the stock solution is stable for up to 1 year.

## Experimental Protocols

---

The following protocol was developed with Jurkat cells, a human T cell line, and using an EdU concentration of 10 μM, and can be adapted for any cell type. Growth medium, cell density, cell type variations, and other factors may influence labeling. In initial experiments, we recommend testing a range of EdU concentrations to determine the optimal concentration for your cell type and experimental conditions. If currently using a BrdU based assay for cell proliferation, a similar concentration to BrdU is a good starting concentration for EdU. If using whole blood as the sample, we recommend heparin as the anticoagulant for collection.

**Figure 3** Workflow diagram for the Click-iT<sup>®</sup> EdU Flow Cytometry Assay Kits.



## Labeling Cells with EdU

- 2.1 Suspend the cells in an appropriate tissue culture medium to obtain optimal conditions for cell growth. Disturbing the cells by temperature changes or washing prior to incubation with EdU slows the growth of the cells during incorporation.
- 2.2 Add EdU to the culture medium at the desired final concentration and mix well. We recommend a starting concentration of 10  $\mu\text{M}$  for 1–2 hours. For longer incubations, use lower concentrations. For shorter incubations, higher concentrations may be required. For a negative staining control, include cells from the same population that have not been treated with EdU.
- 2.3 Incubate under conditions optimal for cell type for the desired length of time. Altering the amount of time the cells are exposed to EdU or subjecting the cells to pulse labeling with EdU allows the evaluation of various DNA synthesis and proliferation parameters. Effective time intervals for pulse labeling and the length of each pulse depend on the cell growth rate.
- 2.4 Harvest cells and proceed immediately to step 3.1 if performing antibody surface labeling; otherwise continue to step 4.1.

## Staining Cell-Surface Antigens with Antibodies (Optional)

- 3.1 Wash cells once with 3 mL of 1% BSA in PBS, pellet cells by centrifugation, and remove supernatant.
- 3.2 Dislodge the pellet and resuspend cells at  $1 \times 10^7$  cells/mL in 1% BSA in PBS.
- 3.3 Add 100  $\mu\text{L}$  of cell suspension or whole blood sample to flow tubes.
- 3.4 Add surface antibodies and mix well (Table 2, page 3).

**Note:** Do not use PE, PE-tandem, or Qdot<sup>®</sup> antibody conjugates before performing the click reaction; wait until step 6.1 for labeling with these fluorophores.

- 3.5 Incubate for the recommended time and temperature, protected from light.
- 3.6 Proceed to step 4.1 for cell fixation.

## Fixation and Permeabilization

The Click-iT<sup>®</sup> saponin-based permeabilization and wash reagent can be used with whole blood or cell suspensions containing red blood cells, as well as with cell suspensions containing more than one cell type. This permeabilization and wash reagent maintains the morphological light scatter characteristics of leukocytes while lysing red blood cells.

- 4.1 Wash the cells once with 3 mL of 1% BSA in PBS, pellet the cells, and remove the supernatant.
- 4.2 Dislodge the pellet, add 100  $\mu\text{L}$  of Click-iT<sup>®</sup> fixative (Component D), and mix well.
- 4.3 Incubate the cells for 15 minutes at room temperature, protected from light.
- 4.4 Wash the cells with 3 mL of 1% BSA in PBS, pellet the cells, and remove the supernatant. Repeat the wash step if red blood cells or hemoglobin are present in the sample. Remove all residual red blood cell debris and hemoglobin before proceeding.
- 4.5 Dislodge the cell pellet and resuspend the cells in 100  $\mu\text{L}$  of 1X Click-iT<sup>®</sup> saponin-based permeabilization and wash reagent (prepared in step 1.4), and mix well. Incubate the cells for 15 minutes or proceed directly to step 5.1 for click labeling.

## Click-iT® Reaction

- 5.1 Prepare 1X Click-iT® EdU buffer additive by diluting the 10X stock solution (prepared in step 1.5) 1:10 in deionized water.
- 5.2 Prepare the Click-iT® reaction cocktail according to Table 3.

**Table 3** Click-iT® reaction cocktails

Reaction components	Number of reactions						
	1	2	5	10	15	30	50
PBS, D-PBS, or TBS	438 µL	875 µL	2.19 mL	4.38 mL	6.57 mL	13.2 mL	21.9 mL
CuSO <sub>4</sub> (Component F)	10 µL	20 µL	50 µL	100 µL	150 µL	300 µL	500 µL
Fluorescent dye azide (prepared in step 1.3)	2.5 µL	5 µL	12.5 µL	25 µL	37.5 µL	75 µL	125 µL
Reaction Buffer Additive (prepared in step 5.1)	50 µL	100 µL	250 µL	500 µL	750 µL	1.5 mL	2.5 mL
Total reaction volume	500 µL	1 mL	2.5 mL	5 mL	7.5 mL	15 mL	25 mL

**Note:** Use the Click-iT® reaction cocktail within 15 minutes of preparation.

- 5.3 Add 0.5 mL of Click-iT® reaction cocktail to each tube and mix well.
- 5.4 Incubate the reaction mixture for 30 minutes at room temperature, protected from light.
- 5.5 Wash the cells once with 3 mL of 1X Click-iT® saponin-based permeabilization and wash reagent (prepared in step 1.4), pellet the cells, and remove the supernatant. Dislodge the cell pellet and resuspend the cells in 100 µL of 1X Click-iT® saponin-based permeabilization and wash reagent, if proceeding with intracellular antibody labeling in step 6.1. Otherwise, add 500 µL of 1X Click-iT® saponin-based permeabilization and wash reagent and proceed with step 7.1 for staining the cells for DNA content, or with step 8.1 for analyzing the cells on a flow cytometer.

## Staining Intracellular or Surface Antigens (Optional)

- 6.1 Add antibodies against intracellular antigens or against surface antigens that use RPE, PE-tandem, or Qdot® antibody conjugates. Mix well.
- 6.2 Incubate the tubes for the time and temperature required for antibody staining, protected from light.
- 6.3 Wash each tube with 3 mL of 1X Click-iT® saponin-based permeabilization and wash reagent (prepared in step 1.4), pellet the cells, and remove the supernatant. Dislodge the cell pellet and resuspend the cells in 500 µL of 1X Click-iT® saponin-based permeabilization and wash reagent, and proceed with step 7.1 for staining the cells for DNA content, or with step 8.1 for analyzing the cells on a flow cytometer.

## Staining Cells for DNA Content (Optional)

7.1 If necessary, add Ribonuclease A to each tube and mix (Table 4).

**Table 4** Click-iT® EdU compatibility with DNA content stains

DNA Content Stain	Click-iT® EdU Stain Compatibility			RNase required?
	Pacific Blue™	Alexa Fluor® 647	Alexa Fluor® 488	
FxCycle™ Violet	No	Yes	Yes	No
Propidium iodide (PI)	Yes	Yes	No*	Yes
SYTOX® AADvanced™	Yes	Yes	Yes	Yes
FxCycle™ Far Red	Yes	No	Yes	Yes

\* You can use propidium iodide (PI) for cell cycle analysis with Alexa Fluor® 488 dye; however, this may reduce the Click-iT® Alexa Fluor® 488 EdU signal. Other DNA stains are recommended instead.  
**Note:** The complete excitation and emission spectra for the DNA content stains are available at [www.invitrogen.com/spectraviewer](http://www.invitrogen.com/spectraviewer).

7.2 Add the appropriate DNA stain to each tube, mix well, and incubate as recommended for each DNA stain.

### Analysis by Flow Cytometry

If measuring total DNA content on a traditional flow cytometer using hydrodynamic focusing, use a low flow rate during acquisition. If using the Attune® Acoustic Focusing Cytometer, all collection rates may be used without loss of signal integrity if the event rate is kept below 10,000 events per second. However, for each sample within an experiment, the same collection rate and cell concentration should be used. The fluorescent signal generated by DNA content stains is best detected with linear amplification. The fluorescent signal generated by Click-iT® EdU labeling is best detected with logarithmic amplification.

8.1 Analyze the cells using a flow cytometer.

- For the detection of EdU with Pacific Blue™ azide, use 405 nm excitation with a violet emission filter (450/50 nm or similar).
- For the detection of EdU with Alexa Fluor® 647 azide use 633/635 nm excitation with a red emission filter (660/20 nm or similar).
- For the detection of EdU with Alexa Fluor® 488 azide, use 488 nm excitation with a green emission filter (530/30 nm or similar).




## Dual Pulse Labeling using EdU and BrdU

Follow these guidelines to perform dual labeling of cultured cells by combining EdU with BrdU labeling.

- Use EdU for the first pulse and BrdU for the second pulse.
- Removal of EdU from the cell culture media is not required when BrdU is added as the second label.
- Addition of BrdU to culture media containing EdU results in preferential incorporation of BrdU into the DNA with the exclusion of EdU, while simultaneous addition of EdU with equimolar or half equimolar BrdU to the media results in only BrdU incorporation. This simplifies the dual labeling protocol by eliminating the wash step normally required to remove the first label from the culture media prior to addition of the second label.
- Process the cells after dual pulse labeling using a proven BrdU protocol.
- After the DNA denaturation step in the BrdU protocol, click label the cells first for the detection of EdU, and then follow with an antibody labeling protocol for the detection of BrdU.
- Select a BrdU antibody which does not have cross-reactivity to EdU, such as clone MoBu-1 (Cat. nos. B35129, B35139, B35140, B35141).

**Figure 4** Dual pulse labeling with EdU and BrdU. Jurkat cells were first treated with 20 mM EdU for 1 hour. Without washing or removal of the EdU, BrdU was added at a 10  $\mu$ M concentration, and the cells were incubated for 1 hour. The cells were harvested, washed, fixed with 70% ice-cold ethanol and stored at 4°C for 96 hours, followed by an HCl-based denaturation procedure. EdU was detected with Alexa Fluor® 488 azide using the Click-iT® EdU Flow Cytometry Kit (Cat. no. C10420). BrdU was then detected with anti-BrdU, Alexa Fluor® 647 conjugate (Cat. no. A21305). SYTOX® Blue nucleic acid stain (Cat. no. S11348) with RNase was used to detect DNA content. The labeled cells were analyzed by flow cytometry using 488 nm excitation with a 530/30 nm bandpass to detect EdU, 633 nm excitation with a 660/20 nm bandpass to detect BrdU, and 405 nm excitation with a 450/50 nm bandpass to detect DNA content. Four populations of cells are distinguished in the EdU vs BrdU plot: cells which are positive for both (Q2, upper right), cells which are negative for both (Q3, lower left), EdU positive but BrdU negative (Q4, lower right), and cells which are positive for BrdU but negative for EdU (Q1, upper left).



## References

1. *ChemBiochem* 4, 1147 (2003);
2. *J Am Chem Soc* 125, 3192 (2003);
3. *Angew Chem Int Ed Engl* 41, 2596 (2002);
4. *Angew Chem Int Ed Engl* 40, 2004 (2001);
5. *BioTechniques* 44, 927 (2008);
6. *Curr Protoc Cytom* 55, 7.38.1 (2011).



**Product List** Current prices may be obtained at [www.invitrogen.com](http://www.invitrogen.com) or from our Customer Service Department.

Cat no.	Product Name	Unit Size
C10418	Click-iT® EdU Pacific Blue™ Flow Cytometry Assay Kit *50 assays*	1 kit
C10424	Click-iT® EdU Alexa Fluor® 647 Flow Cytometry Assay Kit *50 assays*	1 kit
C10425	Click-iT® EdU Alexa Fluor® 488 Flow Cytometry Assay Kit *50 assays*	1 kit
<b>Related products</b>		
C10419	Click-iT® EdU Alexa Fluor® 647 Flow Cytometry Assay Kit *100 assays*	1 kit
C10420	Click-iT® EdU Alexa Fluor® 488 Flow Cytometry Assay Kit *100 assays*	1 kit
A10027	Click-iT® EdU Alexa Fluor® 488 High-Throughput Imaging (HCS) Assay *2-plate size*	1 kit
A10028	Click-iT® EdU Alexa Fluor® 488 High-Throughput Imaging (HCS) Assay *10-plate size*	1 kit
A10044	EdU (5-ethynyl-2'-deoxyuridine)	50 mg
A10208	Click-iT® EdU Alexa Fluor® 647 High-Throughput Imaging (HCS) Assay *2-plate size*	1 kit
A10209	Click-iT® EdU Alexa Fluor® 594 High-Throughput Imaging (HCS) Assay *2-plate size*	1 kit
B35129	BrdU mouse monoclonal antibody (Clone MoBU-1) - Pacific Blue™ *for flow cytometry* *100 tests*	1 each
B35139	BrdU mouse monoclonal antibody (Clone MoBU-1) Alexa Fluor® 488 *for flow cytometry* *100 tests*	1 each
B35140	BrdU mouse monoclonal antibody (Clone MoBU-1) Alexa Fluor® 647 *for flow cytometry* *100 tests*	1 each
B35141	BrdU mouse monoclonal antibody (Clone MoBU-1) unconjugated *for flow cytometry* *100 tests*	1 each
F10347	FxCycle™ Violet Stain *for flow cytometry* *500 assays*	1 kit
F10348	FxCycle™ Far Red Stain *for flow cytometry* *500 assays*	1 kit
H3570	Hoechst 33342, trihydrochloride, trihydrate *10 mg/mL solution in water*	10 mL
P3566	Propidium Iodide – 1.0mg/ml solution in water	10 mL
S10349	SYTOX® AADvanced™ dead cell stain *for 488 excitation* *for flow cytometry* *100 tests*	1 kit
V35003	Vybrant® DyeCycle™ Violet stain *5 mM in water* *200 assays*	200 µL
12091-021	RNase A (20 mg/mL)	10 mL
14190-144	Dulbecco's Phosphate Buffered Saline 1X without Calcium Chloride without Magnesium Chloride	500 mL
14190-250	Dulbecco's Phosphate Buffered Saline 1X without Calcium Chloride without Magnesium Chloride	10 x 500 mL

## Contact Information

---

### Corporate Headquarters

5791 Van Allen Way  
Carlsbad, CA 92008  
USA  
Phone: +1 760 603 7200  
Fax: +1 760 602 6500  
Email: techsupport@lifetech.com

### European Headquarters

Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF  
UK  
Phone: +44 141 814 6100  
Toll-Free Phone: 0800 269 210  
Toll-Free Tech: 0800 838 380  
Fax: +44 141 814 6260  
Tech Fax: +44 141 814 6117  
Email: euroinfo@invitrogen.com  
Email Tech: eurotech@invitrogen.com

### Japanese Headquarters

LOOP-X Bldg. 6F  
3-9-15, Kaigan  
Minato-ku, Tokyo 108-0022  
Japan  
Phone: +81 3 5730 6509  
Fax: +81 3 5730 6519  
Email: jpinfo@invitrogen.com

Additional international offices are listed at  
[www.invitrogen.com](http://www.invitrogen.com)

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

### Web Resources

Visit the Invitrogen website at [www.invitrogen.com](http://www.invitrogen.com) for:

- Technical resources, including manuals, vector maps and sequences, application notes, Meds, FAQs, formulations, citations, handbooks, etc.
- Complete technical support contact information
- Access to the Invitrogen Online Catalog
- Additional product information and special offers

### SDS

Safety Data Sheets (SDSs) are available at [www.invitrogen.com/sds](http://www.invitrogen.com/sds).

### Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to [www.invitrogen.com/support](http://www.invitrogen.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the product packaging (tube, pouch, or box).

### Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

### Limited Use Label License No. 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial services of any kind, including, without limitation, reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

©2011 Life Technologies Corporation. All rights reserved.

**For research use only. Not intended for any animal or human therapeutic or diagnostic use.**