



Website: [www.cwbiosciences.com](http://www.cwbiosciences.com)  
Technical Support: [service@cwbiosciences.com](mailto:service@cwbiosciences.com)  
Customer Service: 617-800-6495



# Magbead DNA Purification Kit (for NGS Size Selection)

**Catalog Number:** CW2508S (5 ml)  
CW2508M (50 ml)

**Storage Condition:** store at 2-8°C; shipping at room temperature.

## Kit Components:

| Component    | CW2508S<br>(5 ml) | CW2508M<br>(50 ml) |
|--------------|-------------------|--------------------|
| CMPure beads | 5 ml              | 50 ml              |

## Product Introduction:

This kit provides a simple, fast and efficient nucleic acid purification method. The product can be used for the selective or non-selective recovery of DNA during the construction of NGS library, as well as the purification and recovery of PCR products. After the CMPure beads are mixed with the sample in a certain proportion, the magnetic beads selectively adsorb the nucleic acid. After two steps of washing, the eluted DNA is of high purity. The  $A_{260}/A_{280}$  ratio is between 1.7-1.9, and the  $A_{260}/A_{230}$  ratio is usually above 2.0. The DNA purified by this kit is suitable for PCR, Real-Time PCR, sequencing, southern blotting and other experiments.

## Kit Notes:

| Sample type     | Typical yield |
|-----------------|---------------|
| 5000 bp segment | Up to 90%     |
| 1000 bp segment | Up to 90%     |
| 500 bp segment  | Up to 80%     |
| 200 bp segment  | Up to 70%     |

**Not included in the kit:**

1. Magnetic stand: it is recommended to use DynaMag™-2 (Cat. No. 12321D).
2. 80% ethanol.
3. Elution Buffer: Buffer EB (10 mM Tris-HCl, pH8.0); ddH<sub>2</sub>O (pH between 7.0-8.0).

**Preparation before the experiment and important notes:**

1. Freezing, centrifugation, and ultrasound can cause irreversible damage to the CMPure beads.
2. Magnetic beads in CMPure will aggregate into clusters after being placed for a long time, so that the surface area of magnetic beads will be reduced, and the yield of sample recovery will be reduced. Before use, magnetic beads must be thoroughly mixed by vortexing.
3. Before use, it is recommended that the CMPure beads should be vortexed and aliquoted into 1.5 ml microcentrifuge tubes, each tube containing 1 ml of CMPure beads.
4. This kit is not suitable for the purification of DNA fragments smaller than 100 bp. If the DNA fragments are smaller than 100 bp, it is recommended to increase the amount of CMPure to 4 times of the sample volume.
5. For the selective recovery of DNA, CMPure is more sensitive to the concentration of ions in the DNA solution. Because the

concentration of ions in the adapted DNA solution and the PCR product obtained by the NGS library construction kits from different manufacturers are different, the amount of reagents used varies.

**Protocol:**

1. Vortex CMPure for 20 seconds to thoroughly mix it into a homogeneous solution.
2. Add the purified DNA solution to a 1.5 ml centrifuge tube.
3. Add two times of the sample volume of CMPure to the centrifuge tube in the previous step. Vortex for 5 seconds and allow to stand at room temperature for 5 minutes.
4. Place the centrifuge tube from the previous step on the magnetic stand until the beads are fully adsorbed (approximately 5 minutes).
5. Keep the centrifuge tube on the magnetic stand and discard the solution completely, avoiding touch magnetic beads.
6. Continue to hold the tube on the magnetic stand and add 250  $\mu$ l of freshly prepared 80% ethanol to the tube.
7. Keep the centrifuge tube on the magnetic stand. discard the ethanol completely after the suspended magnetic beads are fully absorbed.
8. Repeat step 6-7 twice.

9. Keep the centrifuge tube on the magnetic stand for 10 minutes to completely evaporate the ethanol.
10. Remove the centrifuge tube from the magnetic stand and add 20-100  $\mu\text{l}$  of EB (self-prepared) or ddH<sub>2</sub>O. After vortexing to completely resuspend the magnetic beads, leave at room temperature for 5 minutes.
11. Place the centrifuge tube from the previous step on the magnetic stand until the beads are fully adsorbed (approximately 5 minutes).
12. Transfer the DNA elution to a new 1.5 ml centrifuge tube. At this point, magnetic beads can be discarded.