

pcDNA™ 3.4-TOPO® TA Cloning Kit

Cat. no. A14697	Amount 20 reactions	Storage See below
Pub. Part no. 100014709	Pub. no. MAN0006708	Rev. 2.0

Kit Contents and Storage

The pcDNA™ 3.4-TOPO® TA Cloning Kit is shipped on dry ice. Each kit contains two boxes; upon receipt, store their contents as detailed below.

	Item	Amount	Storage*
Box 1	pcDNA™ 3.4-TOPO® vector (5–10 ng/μL)	20 μL	-20°C
	10X PCR buffer	100 μL	-20°C
	dNTP Mix (12.5 mM each dATP, dCTP, dGTP, dTTP)	10 μL	-20°C
	Salt solution	50 μL	-20°C
	Sterile water	1 mL	RT
	Control PCR template (50 ng/μL)	10 μL	-20°C
	Control PCR primers (100 ng/μL)	10 μL	-20°C
	CMV forward sequencing primer (100 ng/μL)	20 μL	-20°C
	pcDNA™ 3.4 reverse sequencing primer (2 μg/μL)	20 μL	-20°C
Box 2	One Shot® TOP10 Chemically Competent <i>E. coli</i>	21 × 50 μL	-80°C
	S.O.C. medium	6 mL	4°C or RT
	pUC19 control DNA (10 pg/μL)	50 μL	-20°C

* For convenience, you may store the entire Box 1 at -20°C, and Box 2 at -80°C.
RT: room temperature.

Description

The pcDNA™ 3.4-TOPO® vector allows the rapid TOPO® cloning of a PCR product containing a gene of interest downstream of the CMV promoter for high level expression of the native protein in adherent mammalian cell cultures, or for high level expression of secreted, native protein in suspension mammalian cell cultures following transient transfection.

The instructions provided below are intended for experienced users. For detailed protocols and more information on the pcDNA™ 3.4-TOPO® vector, refer to the pcDNA™ 3.4-TOPO® TA Cloning Kit manual available online at www.lifetechnologies.com/manuals.

Guidelines for Designing PCR Primers

To obtain the pcDNA™ 3.4-TOPO® expression construct containing your gene of interest, your PCR primer design must include:

- a Kozak consensus sequence
- a mammalian secretion signal upstream of your gene of interest (if you wish to produce secreted protein)
- a stop codon at the end of your gene of interest

Amplify Your Gene of Interest

1. Produce PCR products containing your gene of interest using *Taq* polymerase (see **Note** below). Use cycling parameters suitable for your primers and template and include a final extension step of 7–10 minutes to ensure 3' adenylation.
2. Verify there is only one PCR product of the proper size in each reaction by agarose gel electrophoresis.



Note: We recommend using Platinum® *Taq* DNA Polymerase High Fidelity (Cat. no. 11304). You may use an enzyme mixture containing *Taq* polymerase and a proofreading polymerase; however *Taq* must be used in excess of 10:1 to ensure the presence of 3' A-overhangs on the PCR product. If you use polymerase mixtures that do not have enough *Taq* polymerase or use a proofreading polymerase only, you can add 3' A-overhangs after amplification (refer to the online manual for protocol).

Perform TOPO[®] Cloning

1. Set up the following 6 μL TOPO[®] Cloning reaction for each PCR product, adding the reagents in the order shown. For optimal results, use a 0.5:1 to 2:1 molar ratio of PCR product:TOPO[®] vector.

Reagent	Volume
PCR product	0.5–4 μL
Salt solution	1 μL
Sterile water	to 5 μL
pcDNA [™] 3.4-TOPO [®] vector	1 μL

2. Mix gently and incubate for 5 minutes at room temperature.
3. Place tubes on ice and proceed to transformation.

Transform One Shot[®] TOP10 Chemically Competent *E. coli*

1. Thaw One Shot[®] *E. coli* or equivalent competent cells on ice.
2. Add 2 μL of the TOPO[®] Cloning reaction to a vial of One Shot[®] TOP10 Chemically Competent *E. coli* or equivalent, and mix gently.
3. Incubate on ice for 5–30 minutes.
4. Heat-shock the cells for 30 seconds at 42°C without shaking and immediately transfer the cells to ice.
5. Add 250 μL of room temperature S.O.C. Medium.
6. Cap the tube tightly and incubate at 37°C for 1 hour with horizontal shaking at 200 rpm.
7. Spread 10–50 μL from each transformation on a pre-warmed selective LB plate and incubate overnight.

Analyze Positive Clones

1. Pick 10 colonies and culture them overnight in LB medium containing 100 $\mu\text{g}/\text{mL}$ ampicillin.
2. Isolate plasmid DNA using your method of choice.
3. Analyze the plasmid by restriction analysis or by sequencing to confirm the presence and correct orientation of the insert.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Disclaimer

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

NOTICE TO PURCHASER: LIMITED USE LABEL LICENSE: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

NOTICE TO PURCHASER: LIMITED USE LABEL LICENSE No. 308: WPRE Element

Notice to Purchaser: This product contains the Woodchuck Post-transcriptional Regulatory Element ("WPRE") which is the subject of intellectual property owned by The Salk Institute for Biological Studies, and licensed to Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; and/or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. In addition, any use of WPRE outside of this product or the product's authorized use requires a separate license from the Salk Institute. Life Technologies will not assert a claim against the buyer of infringement of patents owned by Life Technologies and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product or for a Commercial Purpose. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008, Phone (760) 603-7200. Fax (760) 602-6500, or The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: (858) 453-4100 extension 1275, Fax: (858) 546-8093.

©2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

For support visit lifetechnologies.com/support or email
techsupport@lifetech.com
lifetechnologies.com

