

# Single Cell-to-CT™ Kit

Catalog Numbers 4458236 and 4458237

Pub. no. 4458354 Rev. E

Component	Quantity		Storage conditions
	50 rxns (Cat. no. 4458237)	400 rxns (Cat. no. 4458236)	
Single Cell DNase I	50 µL	400 µL	-20°C
Single Cell Lysis Solution	0.5 mL	4.0 mL	4°C
Single Cell PreAmp Mix	265 µL	2.1 mL	-20°C
Single Cell Stop Solution	50 µL	400 µL	-20°C
Single Cell SuperScript® RT	75 µL	600 µL	-20°C
Single Cell VILO™ RT Mix	150 µL	1.2 mL	-20°C
TaqMan® Gene Expression Master Mix	5.0 mL	50 mL	4°C

**Note:** For safety and biohazard guidelines, refer to the “Safety” section in the *Single Cell-to-CT™ Kit User Guide* (Pub. no. 4458356). For every chemical, read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

This insert provides abbreviated instructions for using the Single Cell-to-CT™ Kit. For detailed instructions and troubleshooting information, refer to the *Single Cell-to-CT™ Kit User Guide* (Pub. no. 4458356). You can download a PDF version of the protocol from the Life Technologies web site at [www.lifetechnologies.com](http://www.lifetechnologies.com).

## Perform single-cell lysis

1. Add 1 µL Single Cell DNase I to 9 µL Single Cell Lysis Solution.
2. Add 1–10 cells to the 10 µL of Single Cell DNase I/Single Cell Lysis Solution (no mixing required), then incubate at room temperature for 5 minutes.
3. Add 1 µL Single Cell Stop Solution (no mixing required), incubate at room temperature for 2 minutes, then place the sample on ice.

## Perform reverse transcription

1. Prepare sufficient RT master mix for all samples, then add 4.5 µL to each lysed cell sample:

Component	Volume for each reaction
Single Cell VILO™ RT Mix	3 µL
Single Cell SuperScript® RT	1.5 µL
<b>Total volume for each reaction</b>	<b>4.5 µL</b>

2. Perform reverse transcription in a thermal cycler:

Temp	Time
25°C	10 min
42°C	60 min
85°C	5 min

## Perform preamplification

1. Pool the TaqMan® Gene Expression Assays for your targets of interest, then dilute the pooled assays using 1X TE Buffer, pH 8.0 so that each assay is at a final concentration of 0.2X.
2. Prepare sufficient PreAmp reaction mix, then add 11 µL to each reverse-transcribed sample:

Component	Volume for each reaction
Single Cell PreAmp Mix	5 µL
0.2X pooled TaqMan® Gene Expression Assays	6 µL
<b>Total PreAmp reaction mix</b>	<b>11 µL</b>



### 3. Perform preamplification in a thermal cycler:

Stage	Step	Temp	Time
Holding	Enzyme activation	95°C	10 min
Cycling (14 cycles)	Denature	95°C	15 sec
	Anneal/extend	60°C	4 min
Holding	Enzyme deactivation	99°C	10 min

### 4. Place the tubes on ice or store at -20°C.

## Real-time PCR recommendations

### Reaction conditions

**Note:** Dilute the preamplified products with 1X TE Buffer, pH 8.0 before performing real-time PCR.

- Using 96-well or 384-well reaction plates:

Component	Volume for each reaction	
	96-well plate	384-well plate
2X TaqMan® Gene Expression Master Mix	25.0 µL	10.0 µL
Preamplified product diluted 1:20 with 1X TE Buffer, pH 8.0	10.0 µL	4.0 µL
20X TaqMan® Gene Expression Assay	2.5 µL	1.0 µL
Nuclease-free water	12.5 µL	5.0 µL
<b>Total volume for each reaction</b>	<b>50 µL</b>	<b>20 µL</b>

- Using a TaqMan® Array:

Component	Volume for a full array
2X TaqMan® Gene Expression Master Mix	450 µL
Preamplified product diluted 1:40 with 1X TE Buffer, pH 8.0	450 µL
<b>Total volume for the array</b>	<b>900 µL</b>

### Real-time PCR conditions

Stage	Step	Temp	Time
Holding	UDG incubation	50°C	2 min
Holding	Enzyme activation	95°C	10 min
Cycling (40 cycles)	Denature	95°C	5 sec
	Anneal/extend	60°C	1 min

### Analysis conditions

- Use an automatic baseline and set the threshold to 0.2.
- Review the amplification plots and remove outliers.
- Omit samples that are undetectable for all assays tested.
- Use RealTime StatMiner® Software to perform differential expression analysis.

#### For Research Use Only. Not for use in diagnostic procedures.

The information in this guide is subject to change without notice.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

#### Limited Use Label License No: 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com).

#### Limited Use Label License No: 413: PCR Enzymes

Notice to Purchaser: For research purposes only. Diagnostic uses require a separate license from Roche.

© 2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. RealTime StatMiner is a registered trademark of Integromics, S.L.

For support visit [lifetechnologies.com/support](http://lifetechnologies.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[lifetechnologies.com](http://lifetechnologies.com)

30 November 2012

