

Publication No. MAN0016161 Rev B.0

**Package Contents****Catalog Numbers****Amount:**

A33073

4 preps

A31217

10 preps

A31231

25 preps

**Storage Conditions**

- Store all components at room temperature.

**Required Materials**

- Compatible vacuum manifold with vacuum source capable of -600 to -800 mbar (vacuum protocol only)
- 100% isopropanol
- 100% ethanol
- Centrifuges with swinging bucket rotors capable of 1,000 × g at room temperature, and 4,000 × g at 4°C; high speed centrifuge capable of >12,000 × g at 4°C

**Timing**

Bacterial culture: overnight

Purification: 90 minutes

**Selection Guide**

Go online to view related products:

**PureLink™ Nucleic Acid Purification Kits**  
**Expi293™ and ExpiCHO™ Expression Systems****Product Description**

- The PureLink™ Expi Endotoxin-Free Maxi Plasmid Purification Kit enables isolation of highly pure, endotoxin-free (<0.1 EU/μg) advanced transfection grade plasmid DNA for use in your most sensitive downstream applications.
- Isolate up to 1.5 mg of high quality, endotoxin-free ultrapure plasmid DNA from 100–200 mL of bacterial culture.
- Endotoxin-free (<0.1 EU/μg) DNA with ultra-low levels of RNA, gDNA or protein contamination that is ideal for mammalian cell transfection or *in vivo* experiments.
- Lysis Indicator turns blue during lysis and helps to ensure complete lysis and subsequent neutralization.

**Important Guidelines**

- Add Lysis Indicator to Resuspension Buffer (R3) at a 1:1000 ratio if desired.
- If precipitate is observed in the Lysis Buffer (L7), warm the buffer in a 37°C water bath until the solution clears. Swirl contents gently to resuspend.
- Grow transformed *E. coli* in LB medium. Use 100 mL (high copy number plasmid) or 200 mL (low copy number plasmid) of an overnight culture.
- Perform protocol (step 6 onward) using endotoxin-free/pyrogen free components.
- Do not over-dry DNA. If the DNA pellet is difficult to resuspend, allow the pellet to incubate in TE Buffer for a longer period of time.

**Online Resources**

Visit our product page for additional information and protocols. For support, visit [thermofisher.com/support](http://thermofisher.com/support).

## Maxiprep plasmid isolation protocol (centrifuge)

Before first use of the kit, add RNase A to the Resuspension Buffer (R3) according to the label on the bottle (Mark the label to indicate that RNase A is added. Store at 4°C). Prepare endotoxin-free 70% ethanol by adding 100% ethanol to endotoxin-free water.

Steps	Procedure Details
1 Harvest	Sediment cells by centrifugation at $4,000 \times g$ for 15 min at 4°C. Discard all medium.
2 Resuspend	Add 6 mL Resuspension Buffer (R3) with RNase A to the cell pellet and resuspend the pellet until it is homogeneous.
3 Lyse	Add 6 mL Lysis Buffer (L7). Mix gently by inverting the capped tube until the mixture is homogeneous. Do not vortex. Incubate at room temperature for 5 minutes. <b>Note:</b> If using Lysis Indicator, the solution turns blue.
4 Precipitate	Add 6 mL Precipitation Buffer (N3). Mix immediately by inverting the tube until the mixture is homogeneous. Do not vortex. <b>Note:</b> If using Lysis Indicator, the solution turns completely colorless.
5 Clarify	Pour the lysate into a <b>Lysate Clarification Column</b> (pre-inserted in a 50-mL conical tube with white cap). Centrifuge at $1,000 \times g$ for 4 minutes.
6 Remove Endotoxins	Remove and discard the Lysate Clarification Column. Add 2.5 mL Endotoxin Removal Buffer (ER) to the clarified lysate and mix by inverting the 50-mL conical tube (with white cap) 10 times.
7 Bind	Pour the <b>clarified lysate</b> (from step 6) into a <b>DNA-binding Column</b> . Place the column into the 50-mL conical tube (with white cap) from step 6, and Centrifuge in a swinging bucket rotor at $1,000 \times g$ for 1 minute to drain the solution.
8 Wash	Discard flow through, then add 20 mL (one column volume) of Wash Buffer (W8) to the DNA-binding Column and centrifuge in a swinging bucket rotor at $1,000 \times g$ for 1 minute to drain the solution.
9 Elute	Place the DNA-binding Column into an endotoxin-free 50-mL conical tube (with blue cap). Add 15 mL Elution Buffer (E4) to the column and incubate for 1 minute. Centrifuge at $1,000 \times g$ for 1 minute. Remove and discard the column.
10 Precipitate and Wash	Add 10.5 mL (0.7 volume) of isopropanol to the eluate. Mix well. Centrifuge at $>12,000 \times g$ for 30 minutes at 4°C. Remove and discard the supernatant. Wash the DNA pellet with 5 mL 70% ethanol. Centrifuge at $>12,000 \times g$ for 10 minutes at 4°C. Remove the supernatant.
11 Resuspend	Air-dry the pellet for 10 minutes, then resuspend the purified plasmid DNA in TE Buffer (TE). Store plasmid DNA in endotoxin-free tubes at -20°C.

## Maxiprep plasmid isolation protocol (vacuum)

Before first use of the kit, add RNase A to the Resuspension Buffer (R3) according to the label on the bottle (Mark the label to indicate that RNase A is added. Store at 4°C). Prepare endotoxin-free 70% ethanol by adding 100% ethanol to endotoxin-free water.

Steps	Procedure Details
1 Harvest	Sediment cells by centrifugation at $4,000 \times g$ for 15 min at 4°C. Discard all medium.
2 Resuspend	Add 6 mL Resuspension Buffer (R3) with RNase A to the cell pellet and resuspend the pellet until it is homogeneous.
3 Lyse	Add 6 mL Lysis Buffer (L7). Mix gently by inverting the capped tube until the mixture is homogeneous. Do not vortex. Incubate at room temperature for 5 minutes. <b>Note:</b> If using Lysis Indicator, the solution turns blue.
4 Precipitate	Add 6 mL Precipitation Buffer (N3). Mix immediately by inverting the tube until the mixture is homogeneous. Do not vortex. <b>Note:</b> If using Lysis Indicator, the solution turns completely colorless.
5 Clarify	Pour the lysate into a <b>Lysate Clarification Column</b> (pre-inserted in a 50-mL conical tube with white cap). Centrifuge at $1,000 \times g$ for 4 minutes.
6 Remove Endotoxins	Remove and discard the Lysate Clarification Column. Add 2.5 mL Endotoxin Removal Buffer (ER) to the clarified lysate and mix by inverting the 50-mL conical tube (with white cap) 10 times.
7 Bind	Attach the <b>DNA-binding Column</b> to a vacuum manifold and load the <b>clarified lysate</b> (from step 6) onto the column. Apply vacuum to drain the solution.
8 Wash	Add 20 mL (one column volume) of Wash Buffer (W8) to the DNA-binding Column and apply vacuum to drain the solution.
9 Elute	Place the DNA-binding Column into an endotoxin-free 50-mL conical tube (with blue cap). Add 15 mL Elution Buffer (E4) to the column and incubate for 1 minute. Centrifuge at $1,000 \times g$ for 1 minute. Remove and discard the column.
10 Precipitate and Wash	Add 10.5 mL (0.7 volume) of isopropanol to the eluate. Mix well. Centrifuge at $>12,000 \times g$ for 30 minutes at 4°C. Remove and discard the supernatant. Wash the DNA pellet with 5 mL 70% ethanol. Centrifuge at $>12,000 \times g$ for 10 minutes at 4°C. Remove the supernatant.
11 Resuspend	Air-dry the pellet for 10 minutes, then resuspend the purified plasmid DNA in TE Buffer (TE). Store plasmid DNA in endotoxin-free tubes at -20°C.

**Disclaimer:** TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit [thermofisher.com/support](https://www.thermofisher.com/support)

12 February 2017

**ThermoFisher**  
SCIENTIFIC