Conventional in vitro Transcription

This protocol is for the Conventional *in vitro* Transcription

More than 10 µg of RNA transcript can be generated per 1 µg template DNA using the following protocol. The reaction can be scaled up or down. For high yield transcription, generating up to 200 µg RNA, use TranscriptAidTM T7 High Yield Transcription Kit.

- Thaw frozen reagents, mix and centrifuge briefly.
- Keep enzymes and nucleotides on ice.
- Keep the Reaction Buffer at room temperature.
- 1. Prepare the following reaction mixture at room temperature:

5X Transcription buffer	10 µl
ATP/GTP/CTP/UTP Mix, 10 mM each	10 µl (2 mM final concentration)
Linearized template DNA	1 µg
RiboLock™ RNase Inhibitor	1.25 µl (50 u)
T7/T3/SP6 RNA Polymerase	1.5 µl (30 u)
DEPC-treated Water	to 50 μl
Total volume	50 µl

- 2. Incubate at 37°C for 2 hours.
- Optional: To remove template DNA add 2 μl (2 u) of DNase I, RNase-free, mix and incubate at 37°C for 15 min.
- Stop the reaction by addition of 2 μl 0.5 M EDTA, pH 8.0 and incubate at 65°C for 10 min.

Note

RNA hydrolyzes if heated in the absence of a chelating agent.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America Technical Services: techservice.genomics@thermofisher.com Customerservice.genomics@ thermofisher.com Tel 800 235 9880 Fax 800 292 6088 Europe and Asia Technical Services: techservice.emea.genomics@ thermofisher.com Customer Services: customer service.emea.genomics@ thermofisher.com

