

Concentrating RiboMinus™ RNA Using the RiboMinus™ Concentration Module

Pub. Part no. 25-1042

MAN0002509

Rev. Date 3 December 2011

Description

After purifying RiboMinus™ RNA using the RiboMinus™ Eukaryote or Plant Kit for RNA-Seq, you need to concentrate RiboMinus™ RNA for further use in downstream applications. To retain all species of RNA <200 nucleotides (nt), you have two options to concentrate RiboMinus™ RNA:

- Ethanol precipitation, see the manual supplied with RiboMinus™ Eukaryote or Plant Kit for RNA-Seq for the protocol (not recommended for Low-input Protocol)
- RiboMinus™ Concentration Module using silica spin columns (Cat. no. K1550-05), see next page for a **modified** protocol

Generally, RNA species <200 nt are excluded from the standard binding conditions for RNA isolation using silica spin column. Recent studies have shown the importance of these small RNA species which include regulatory RNA molecules such as microRNA (miRNA), short interfering RNA (siRNA), snRNA, and other RNA transcripts of yet unknown function.

To retain all species of RNA <200 nucleotides (nt) during the concentration step of RiboMinus™ RNA isolation, use the modified silica spin column protocol with the RiboMinus Concentration Module wherein the binding of RNA is performed with 60% ethanol to retain all RNA species including <200 nt.

System Overview

The binding conditions of the spin column protocol are optimized for the RiboMinus™ RNA sample with ethanol and buffer. The sample is loaded onto a spin column. The RiboMinus™ RNA binds to the silica-based membrane in the column and impurities are removed by thorough washing with Wash Buffer. The RiboMinus™ RNA is then eluted in RNase-free water.

Important

Do not use the purification protocols included in the RiboMinus™ Concentration Module manuals for concentrating RiboMinus™ RNA as using the standard RNA spin column purification protocols will result in loss of small RNA species. Use the modified spin column purification protocol included in this manual for recovering all small RNA species including <200 nt.

CAUTION! The RiboMinus™ Concentration Module buffers contain guanidine isothiocyanate. Always wear a laboratory coat, disposable gloves, and eye protection when handling buffers. Do not add bleach or acidic solutions directly to solutions containing guanidine isothiocyanate or sample preparation waste as it forms reactive compounds and toxic gases when mixed with bleach or acids. Solutions containing ethanol are considered flammable. Use appropriate precautions when using this chemical. Dispose of the buffers and chemicals in appropriate waste containers.

General Guidelines

- Use disposable, individually wrapped, sterile plasticware and use sterile, new pipette tips and microcentrifuge tubes
- Wear latex gloves while handling reagents and RNA samples to prevent RNase contamination from the skin surface
- Always use proper microbiological aseptic techniques when working with RNA
- Use RNase AWAY® Reagent (Cat. no. 10328-011) to remove RNase contamination from surfaces


Materials Needed

- RiboMinus™ RNA (obtained using RiboMinus™ Eukaryote Kit for RNA-Seq, Cat. no. A10837-02 or A10837-08 and RiboMinus™ Plant Kit for RNA-Seq, Cat. no. A10838-02 or A10838-08)
- RiboMinus™ Concentration Module (Cat. no. K1550-05)
- 96–100% ethanol
- Microcentrifuge capable of centrifuging >12,000 × g

Before Starting

Before using Wash Buffer (W5) for the **first time**, add 6 mL 96–100% ethanol to 1.5 mL Wash Buffer (W5) included with the kit. Check the box on the Wash Buffer label to indicate that ethanol was added. Store Wash Buffer (W5) with ethanol at room temperature.

Workflow



Modified RiboMinus™ Concentration Module Protocol for RiboMinus™ RNA

1. Transfer the RiboMinus™ RNA (~300 µL, Standard/~50 µL, Low-input) to a new tube capable of holding >1.5 mL.
2. Add one sample volume of Binding Buffer L3 (300 µL, Standard/50 µL, Low-input) and 3 volumes of 96–100% ethanol (900 µL, Standard/150 µL, Low-input) for a final ethanol concentration of 60%. Mix by vortexing.
Note: The binding step is performed with a 60% ethanol concentration to ensure recovery of all RNA species, including those of <200 nt, which are typically not retained using standard silica binding conditions.
3. Transfer up to 600 µL of the sample (from step 3) to the spin column (with the collection tube).
4. Centrifuge at 12,000 × g for 1 minute at room temperature. Discard the flow-through, and reinsert the Spin Column into the same Collection Tube.
5. **Repeat** Steps 4–5 until the entire sample is processed.
6. Wash the column by adding 600 µL (Standard) or 200 µL (Low-input) of Wash Buffer (W5) prepared with ethanol. Centrifuge the column at 12,000 × g for 1 minute at room temperature. Discard the flow through.
7. Discard the collection tube and place the column into a clean collection tube, supplied with the kit.
8. Centrifuge the column at maximum speed for 2–3 minutes at room temperature to remove any residual Wash Buffer (W5). Place the column in a clean, 1.5-mL Recovery Tube.
9. Add 10–30 µL of nuclease-free water to the center of the column. Incubate the column at room temperature for 1 minute.
10. Centrifuge the column at maximum speed for 1 minute at room temperature. **The Recovery Tube contains purified RiboMinus™ RNA.**
11. Place RiboMinus™ RNA on ice and proceed to the desired downstream application or store RiboMinus™ RNA at –80°C.

Analyzing RiboMinus™ RNA

The purified RiboMinus™ RNA is easily quantitated using UV absorbance at 260 nm or Quant-iT™ RNA Assay Kit. The RNA isolated is of high-quality and is efficiently depleted in ribosomal RNA (rRNA) species.

To verify rRNA depletion, perform agarose gel electrophoresis of the sample or use a bioanalyzer. Agarose gel electrophoresis analysis shows depletion of 18S and 28S rRNA bands as compared to a control sample. Absence of contaminating DNA and RNA degradation may also be confirmed by agarose gel electrophoresis. The efficiency of rRNA depletion in RiboMinus™ RNA, RNA degradation, and RNA concentration can also be analyzed using a bioanalyzer such as the Agilent 2100 bioanalyzer with Agilent RNA Nano Kit.

Troubleshooting

Observation	Cause	Solution
Low RiboMinus™ RNA yield	Incorrect binding conditions	For efficient binding of RiboMinus™ RNA to the spin column, ensure the sample contains 60% ethanol prior to loading onto the spin column.
	Ethanol not added to Wash Buffer	Be sure to add 96–100% ethanol to Wash Buffer as described above.
	Incorrect elution conditions	Add water to the center of the column and perform incubation for 1 minute with water before centrifugation.
	RNA quantitation performed with water	Be sure the RNA quantitation using UV absorbance is performed with 10 mM Tris-HCl, pH 7.0 to accurately measure the UV absorbance.
RNA degraded	RNase contamination	Follow the guidelines on page 2 to prevent RNase contamination.
Inhibition of downstream enzymatic reactions	Presence of ethanol in purified RNA sample	Remove traces of ethanol from residual Wash Buffer by placing the spin column into a new collection tube after washing the column with Wash Buffer (W5), and centrifuging the column at maximum speed for 2–3 minutes to completely dry the column.

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

©2011 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation or their representative owners. RNase AWAY® is a registered trademark of Molecular Bio-Products, Inc.

For support visit www.lifetechnologies.com/support or contact techsupport@lifetech.com.

www.lifetechnologies.com

Intended Use: For research use only. Not intended for any animal or human therapeutic or diagnostic use.

