

# ChargeSwitch®-Pro PCR Clean-up Kit

Catalog nos. CS32050 and CS32250

Part no. 25-0983 Version A, 15 August 2007

## Centrifugation Protocol

Follow the steps below to purify PCR fragments **using a microcentrifuge**. All steps are performed at room temperature. For more detailed protocols and additional information, refer to the kit manual.

**Note:** After PCR cycling, cool the reaction to room temperature before purification.

### 1. Binding the DNA

1. To the PCR reaction, add an equal volume of ChargeSwitch®-Pro PCR Purification Buffer (e.g., for a 50- $\mu$ l PCR reaction, add 50  $\mu$ l of buffer). Briefly vortex to mix.
2. Transfer the mixture onto the ChargeSwitch®-Pro PCR Clean-up Column inserted in a Collection Tube.
3. Centrifuge the column/tube at 10,000  $\times g$  for 30–60 seconds.
4. Proceed to **Washing the Column**. (**Note:** If the volume of your PCR reaction was >75  $\mu$ l, empty Collection Tube before proceeding to avoid overflow.)

### 2. Washing the Column

1. Add 500  $\mu$ l of ChargeSwitch®-Pro PCR Wash Buffer to the column.
2. Centrifuge the column/tube at 10,000  $\times g$  for 1 minute.
3. Discard the flow-through and the Collection Tube.
4. Insert the column into a new, sterile Elution Tube (provided in the kit). Then proceed to **Eluting the DNA**.

### 3. Eluting the DNA

1. Add 50  $\mu$ l of ChargeSwitch®-Pro PCR Elution Buffer onto the column, and incubate at room temperature for 1 minute.
2. Centrifuge the column/tube at 10,000  $\times g$  for 30–60 seconds. **The flow-through contains the purified DNA.**
3. Store the purified DNA at 4°C for immediate use or at -20°C for long-term storage.

# ChargeSwitch®-Pro PCR Clean-up Kit

Catalog nos. CS32050 and CS32250

Part no. 25-0983 Version A, 15 August 2007

## Vacuum Protocol

Follow the steps below to purify PCR fragments **using a vacuum manifold**. All steps are performed at room temperature. For more detailed protocols and additional information, refer to the kit manual.

**Note:** After PCR cycling, cool the reaction to room temperature before purification.

### 1. Binding the DNA

1. To the PCR reaction, add an equal volume of ChargeSwitch®-Pro PCR Purification Buffer (e.g., for a 50- $\mu$ l PCR reaction, add 50  $\mu$ l of buffer). Briefly vortex to mix.
2. Remove the ChargeSwitch®-Pro PCR Clean-up Column from the Collection Tube and insert it into the luer extension of a vacuum manifold.
3. Transfer the mixture from Step 1 onto the column.
4. Apply vacuum pressure until the liquid has passed through the column, then proceed to **Washing the Column**.

### 2. Washing the Column

1. Add 500  $\mu$ l of ChargeSwitch®-Pro PCR Wash Buffer to the column.
2. Apply vacuum pressure until the liquid has passed through the column.
3. Remove the column from the vacuum manifold and re-insert it into the Collection Tube.
4. Centrifuge the column/tube at 10,000  $\times g$  for 1 minute.
5. Discard the flow-through *and* the Collection Tube.
6. Insert the column into a new, sterile Elution Tube (provided in the kit) and proceed to **Eluting the DNA**.

### 3. Eluting the DNA

1. Add 50  $\mu$ l of ChargeSwitch®-Pro PCR Elution Buffer onto the column, and incubate at room temperature for 1 minute.
2. Centrifuge the column/tube at 10,000  $\times g$  for 30–60 seconds. **The flow-through contains the purified DNA.**
3. Store the purified DNA at 4°C for immediate use or at -20°C for long-term storage.